

Miltex
Redefining Excellence



ORTHODONTIC INSTRUMENTS

Miltex, Inc.
589 Davies Drive
York, PA 17402
toll-free phone 866-854-8300
toll-free fax 866-854-8400
phone 717-840-9335
fax 717-840-9347
email dental@miltex.com
www.miltex.com



INTRODUCTION

Miltex® Orthodontic instruments are crafted from the finest quality German-grade stainless steel. Our product offering includes but is not limited to the following products:

- **Banding Instruments** – include band seating pliers, band contouring pliers, band removing pliers and a band pusher
- **Ortho Cutters** – include distal end cutters, hard wire cutters, wire cutters, pin and ligature cutters
- **Debonding Instruments** – include bracket removers and adhesive removing pliers
- **Utility Pliers** – include How pliers, Weingart pliers, crown crimping pliers and a variety of additional utility pliers
- **Wire Forming Pliers** – include wire bending pliers, loop forming pliers, arch forming pliers and a variety of additional wire forming pliers
- **Wire Instruments and Wire** – include wire twisting instruments and Babcock suture wire
- **Gauges and Dividers** – include boley gauges, crown gauges and spring bow dividers



TABLE OF CONTENTS

Banding Instruments	4
Cutters	5
<i>Distal End</i>	5
<i>Pin and Ligature</i>	6
Debonding Instruments	8
Pliers, Utility	9
Pliers, Wire Forming	10
Wire Instruments and Wire	13
Scissors, Crown and Collar	14
College Pliers	15
Gauges and Dividers	16
Finishing and Polishing	17
Care & Cleaning	22
<i>Cassettes</i>	22
<i>Cleaning Products</i>	29
<i>Color Bands</i>	30
<i>Sharpeners</i>	31
<i>Instrument Guards</i>	31
<i>Sterilizing Items</i>	32
<i>Care and Cleaning Instructions</i>	33
Indexes	35-36

BANDING INSTRUMENTS

REF	Description
74-30	Pliers - Crown and Band Contouring
74-50	Pliers - GORDON, Contouring
74-80	Pliers - BAADE Band and Shell Remover
74-85	Pliers - Posterior Band Removing, Long
74-135	Band Pusher
74-205	Pliers - Band Seating
74-208	Pliers - Posterior Band Removing, Short



**74-80
BAADE**
Used for removing
band and shells



**74-135
Band Pusher**
Used for band pushing,
burnishing & seating bands;
Serrated tip for effective
control of bands



**74-30
Crown and Band Contouring**
Used for fitting
stainless steel crowns



**74-50
GORDON, Contouring**
Used for forming and
contouring bands;
Curved beaks



74-85
Used for posterior
band removal



**74-205
Band Seating**
Serrated tips grip
the bracket better for
placing bands



74-208
Used for band
removal

ORTHO CUTTERS- DISTAL END & HARD WIRE

REF	Description
74-101	Cutter - Universal Cut & Hold Distal End Cutter
74-101L	Cutter - Universal Cut & Hold Distal End Cutter, Long Handle
74-102	Cutter - Distal End Cutter, Flush Cut
74-103	Cutter - Hard Wire, 15°
74-104	Cutter - Hard Wire, Straight



74-101
Universal Cut & Hold
Distal End Cutter
Cuts all types of wires ranging
from .012"/.30mm to .021 x .025"
(.53 mm x .64 mm)



74-102
Distal End Cutter
Universal cutter for all
types of wire up to a maximum
of .021" x .025"
(.53 mm x .64 mm)



74-101L
Universal Cut & Hold
Distal End Cutter
Cuts all types of wires ranging
from .012"/.30 mm to .021 x .025"
(.53 mm x .64 mm)



74-103
Hard Wire 15°
Cuts hard wire up to .022" x .0285"
(.56 mm x .71 mm);
Angled tip for hard to access areas



74-104
Hard Wire, Straight
Cuts hard wire up to .022" x .028"
(.56 mm x .71 mm)

ORTHO CUTTERS- PIN & LIGATURE

REF	Description
74-105	Cutter - Pin & Ligature
74-106	Cutter - Mini Pin & Ligature
74-107	Cutter - Micro-Mini Pin & Ligature
74-107L	Cutter - Micro-Mini Pin & Ligature, Long Handle
74-109	Cutter - Micro-Mini Ligature, 15°



74-107
Micro-Mini Pin & Ligature
Cuts soft wire pins and ligature wire up to .012" (.30 mm);
Slim design for hard to access areas



74-105
Pin & Ligature
Cuts soft wire pins and ligature wire up to .015" (.38 mm)



74-106
Mini Pin & Ligature
Cuts soft wire pins and ligature wire up to .012" (.30 mm)



74-107L
Micro-Mini Pin & Ligature
Cuts soft wire pins and ligature wire up to .012" (.30 mm);
Slim design for hard to access areas



74-109
Micro-Mini Ligature, 15°
Cuts soft wire pins and ligature wire up to .012" (.30 mm);
Slim design and angled tip aid in accessing hard to reach areas

ORTHO CUTTERS- WIRE

REF	Description
9-128	Cutter - Parallel Action Wire Cutter, 4 1/2" Nickel Plated
74-115	Cutter - Diagonal Wire Cutter, 5 1/2"



9-128
Parallel Action Wire Cutter

Cuts and bends large
hard wire up to .090";
Used in cutting and
bending head gear



74-115
Diagonal Wire Cutter

Cuts all types and sizes
of orthodontic wire;
Excellent for cutting heavy
wire with little pressure

DEBONDING INSTRUMENTS

REF	Description
74-206	Pliers - Adhesive Removing
74-219	Bracket Remover
74-220L	Bracket Remover - Angled, Long Handle



74-206
Adhesive Removing
For removing composite material after debonding attachments



74-219
Bracket Remover
Used for removing steel, ceramic and plastic brackets



74-220L
Bracket Remover, Angled
For removing all bracket styles on posterior and anterior teeth

ORTHO UTILITY PLIERS

REF	Description
74-20	Pliers - How, Straight
74-25	Pliers - How, Offset
74-30	Pliers - Crown and Band Contouring
74-82	Pliers - Weingart
74-82L	Pliers - Weingart - Long Handle
74-83	Pliers - Slim Weingart
74-95	Pliers - Crown and Shell Remover
74-100	Pliers - Primary Crown Crimping
74-210	Pliers - Force Module Separating
74-211	Pliers - Coon Style Ligature Tying
74-212	Holder - Direct Bond Bracket



74-20
How, Straight

Utility pliers with serrations, useful for placement and removal of archwires, pins and other auxiliaries



74-25
How, Offset

Tips are angled at 40° for easy access in hard to reach areas



74-30
Crown and Band Contouring

Used for fitting stainless steel crowns



74-82
Weingart
Utility pliers with serrations



74-82L
Weingart - Long Handle
Utility pliers with serrations



74-83
Slim Weingart
Slim taper tips for hard to access areas



74-95
Crown and Shell Remover
Used for placing crown and shells



74-100
Primary Crown Crimping
Used for crimping crowns



74-210
Force Module Separating
For placement of elastic separators



74-211
Coon Style Ligature Tying
For tying stainless steel ligature wire



74-212
Direct Bond Bracket
Used for holding brackets

ORTHO WIRE FORMING PLIERS

REF	Description
74-10	Pliers - Wire Bending
74-35	Pliers - REYNOLDS, Contouring
74-40	Pliers - PEESO, Collar
74-45	Pliers - Lab and Office
74-55	Pliers - Bird Beak
74-75	Pliers - Three Jaw
74-301	Pliers - Jarabak
74-303	Pliers - Omega Loop Forming
74-305	Pliers - Light Wire Bird Beak with Groove



74-10
Wire Bending
Used for bending wire



74-35
REYNOLDS
For contouring wire



74-40
PEESO
Used for forming collars



74-45
Lab and Office
Used for lab and office work



74-55
Bird Beak
Loop forming pliers used for round wire up to .030" (.76 mm)



74-75
Three Jaw
Used for bending and contouring wire up to .030" (.76 mm)



74-301
Jarabak
For light wire bending; Grooves assist with precise bending and closing of loops. Serrated surface grips wire up to .02" (.51 mm)



74-303
Omega Loop Forming
U Graduated cone .045", .060" and .075" (1.14 mm, 1.52 mm, 1.91 mm); used for forming precise omega loops



74-305
Light Wire Bird Beak with Groove
For forming precise loops in light wires up to .020" (.51 mm). The groove in the tip allows for duplicating loops

ORTHO WIRE FORMING PLIERS

REF	Description
74-306	Pliers - Tweed Loop Forming
74-307	Pliers - Arch Bending
74-308	Pliers - Rectangular Arch Bending
74-309	Pliers - Lingual Arch Forming
74-310	Pliers - Hollow Chop Contouring
74-311	Pliers - O'Brien
74-312	Pliers - Mini Three Jaw
74-316	Pliers - Marcotte Looping



74-306

Tweed Loop Forming

Loop forming pliers with .045", .60" and .075" (1.14 mm, 1.52 mm, 1.91 mm); used for creating bull and helical loops in wires up to .022" x .028" (.56 mm x .71 mm)



74-306

Tweed Loop Forming

Loop forming pliers with .045", .60" and .075" (1.14 mm, 1.52 mm, 1.91 mm); used for creating bull and helical loops in wires up to .022" x .028" (.56 mm x .71 mm)



74-307

Arch Bending

For placing first, second and third order bends; .050" (1.27 mm) blade width



74-308

Rectangular Arch Bending

Used for placing first, second and third order bends; .070" (1.78 mm) blade width



74-309

Lingual Arch Forming

To create double back and triple back bends in .030" and .036" lingual arch wires



74-310

Hollow Chop Contouring

For subtle contouring of arch wire of all sizes



74-311

O'Brien

For forming small helical loops and intra-oral tip back bends



74-312

Mini Three Jaw

Rounded tips ease contouring of light wires



74-316

Marcotte Looping

Used for forming helical loops and T loops

ORTHO WIRE FORMING PLIERS

REF	Description
74-317	Pliers - Hollow Chop
74-318	Pliers - Light Wire Bird Beak
74-319	Pliers - Nance
74-320	Pliers - Adams
74-321	Pliers - V-Stop
74-322	Pliers - Three Prong NiTi Bender
74-324	Pliers - Old Style Looping
74-330	Pliers - Mathieu



**74-317
Hollow Chop**
Used for subtle
contouring of arch wire
of all sizes



**74-317
Hollow Chop**
For subtle contouring
of arch wire of all sizes



**74-318
Light Wire Bird Beak**
Used for forming precise
loops in light wires up to
.020" (.51 mm)



**74-319
Nance**
For forming loops in four
lengths - 3, 4, 5 and 6 mm;
Groove at base aids in closing
loops. Wire capacity .022"



**74-320
Adams**
For right angle bends and
appliance adaptation;
Wires up to .040"



**74-321
V-Stop**
For placing V bend in
archwires to prevent
archwire from traveling



**74-322
Three Prong NiTi Bender**
For forming and contouring
archwire, especially nickel
titanium wire up to .02"



**74-324
Old Style Looping**
Used for forming and closing
loops with one cylindrical three-
step beak and one triangular beak



**74-330
Mathieu**
Used for ligature typing
and placing of elastics

WIRE INSTRUMENTS AND WIRE

REF	Description
22-500TC	Wire Twister 6", with Carb-N-Sert Jaws
22-502TC	Wire Twister 8", with Carb-N-Sert Jaws



22-500TC



22-502TC

BABCOCK Suture Wire, Stainless Steel Type AISI 316-L

B&S Gauge	Diameter	1 oz. Spools	4 oz. Spools
18 GA	.040"	9-40-18	9-44-18
20 GA	.032"	9-40-20	9-44-20
22 GA	.025"	9-40-22	9-44-22
24 GA	.020"	9-40-24	9-44-24
25 GA	.018"	9-40-25	9-44-25
26 GA	.016"	9-40-26	9-44-26
28 GA	.012"	9-40-28	9-44-28
30 GA	.010"	9-40-30	9-44-30
32 GA	.008"	9-40-32	9-44-32
34 GA	.006"	9-40-34	9-44-34
36 GA	.004"	9-40-36	9-44-36



9-40-24



9-44-24

CROWN AND COLLAR AND WIRE CUTTING SCISSORS

REF	Description
9D-115	RA Matrix Scissors, 4-1/4" (10.8 cm), one serrated blade
9D-118	Pedo Crown-Festooning Scissors, 3-1/2" (8.9 cm), curved
9-120	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), straight, one serrated blade
9D-121	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), straight
9-122	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), curved, one serrated blade
9D-123	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), curved
9-124	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade
9-124TC	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), angled to side, one serrated blade, CARB-N-SERT blades
9D-132	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), straight, one serrated blade
9D-133	MILTEX Wire Cutting Scissors, 4-3/4" (12.1 cm), straight
9D-134	MILTEX Wire Cutting Scissors, 4-1/4" (10.8 cm), curved, one serrated blade
9D-135	MILTEX Wire Cutting Scissors, 4-1/4" (10.8 cm), curved
V99D-133	VANTAGE Collar Crown Scissors, 4" (10.2 cm), straight
V99D-135	VANTAGE Collar Crown Scissors, 4" (10.2 cm), curved



9D-115
4-1/4"
One Serrated Blade



9D-118
3-1/2" Curved



9-122
4-3/4" Curved, One
Serrated Blade



9D-123
4-3/4" Curved,
Smooth



9-124
4-3/4" Angular

9-124TC
4-3/4" Angular,
One Serrated Blade



9D-132
4-1/4" Straight,
One Serrated Blade

9D-133
5-1/4" Straight, Smooth

V99D-133
Economy



9D-134
4-1/4" Curved,
One Serrated Blade

9D-135
4-1/4" Curved, Smooth

V99D-135
Economy

COLLEGE (COTTON) PLIERS

REF	Description
66-311	311 Dental Cotton (College) Pliers, forked tips, smooth
66-312A	312A Dental Cotton (Perry) Pliers, angled, serrated tips
66-317	317 Dental Cotton (College) Pliers, serrated tips
66-317A	317A Dental Cotton (College) Pliers, extra fine tips, serrated
66-317B	317B College Pliers, X-fine tips, serrated w/long fine points
66-317S	317 Dental Cotton (College) Pliers, smooth tips
66-318	318 MERIAM Pliers, serrated tips
66-318-49	318-49 College Pliers, non-pinch, fine points, serrated
66-49	49 College Pliers, non-pinch, fine points, serrated
66-652	Cotton Pliers with lock, serrated tips
66-656	Cotton Pliers with lock, grooved tips
V966-317	VANTAGE College Pliers, 6" (15.2 cm), with pin, serrated tips
V966-652	VANTAGE Self-Locking College Pliers 6" (15.2 cm), serrated tips



66-311

66-312A



66-318

66-652
V966-652

66-656

66-317
V966-317

66-317A



66-317B



66-317S



66-318-49



66-49

DIVIDERS & GAUGES

REF	Description
68-693	Crown Gauge, stainless steel
68-694	BOLEY Gauge, stainless steel
68-695	Spring Bow Divider, 4-1/2" (11.4 cm) radius
68-696	Spring Bow Divider, 6" (15.2 cm) radius
68-697	Micro BOLEY Gauge, stainless steel



68-697



68-693



68-694



68-695



68-696

LIGHTNING STRIPS

Lightning strips make finishing and polishing filling materials faster and easier. Made of resilient, high-strength materials, Lightning Strips are ideal for sensitive aesthetic finishing, yet tough enough for use on the hardest porcelain or composite materials.

Lightning Strips - Used on Composite Materials

Strips have a safe side and are available in a 10 pack.

	Narrow 2.38 mm	Medium 3.17 mm	Wide 3.97 mm	Assorted
Fine	620-81010	620-81030	620-81050	620-81210
Medium	620-81070	620-81090	620-81110	620-81230
Coarse	620-81130	620-81150	620-81170	620-81250
Assorted	620-81290	620-81270	620-81310	620-81190

Lightning Strips - Used on Natural Dentition

Strips are dual sided and are available in a 10 pack.

	Narrow 2.38 mm	Medium 3.17 mm	Wide 3.97 mm	Assorted
Fine	630-81330	630-81350	630-81370	630-81530
Medium	630-81390	630-81410	630-81430	630-81550
Coarse	630-81450	630-81470	630-81490	630-81570
Assorted	630-81590	630-81610	630-81630	630-81510

FACT:

Lightning Strips are excellent for finishing composite restorations.



620-81190
Composites
Safe-Sided

Lightning Strips - Used for Orthodontics

Strips have a safe side and are available in a 10 pack.

	Narrow 2.38 mm	Medium 3.17 mm	Wide 3.97 mm	X-Wide 6.35 mm	XX-Wide 9.74 mm
Fine	–	–	–	625-82010	625-82030
Medium	–	–	–	625-82050	625-82070
Coarse	–	–	–	625-82090	625-82110
Assorted	–	–	–	625-82130	625-82150

Lightning Strips - Used for Orthodontics

Strips are dual sided and are available in a 10 pack.

	Narrow 2.38 mm	Medium 3.17 mm	Wide 3.97 mm	X-Wide 6.35 mm	XX-Wide 9.74 mm
Fine	–	–	–	635-82510	635-82530
Medium	–	–	–	635-82550	635-82570
Coarse	–	–	–	635-82590	635-82610



630-81510
Natural Teeth
Dual-Sided



625-82010
Orthodontics
Safe-Sided



635-82530
Orthodontics
Dual-Sided

MICRO DIAMONDS

Micro Diamond Strips have a plastic backing, are extra durable and are used to polish and finishing a wide variety of dental materials.

REF	Description
404-40413	Strips, Micro Diamond Strips, Trail Pack, 6 Strips
404-40414	Strips, Micro Diamond Strips, 24 Strips



404-40414

DUAL ALO


#54 & #56 EVENWET™ STRIPS

Dual Alo Evenwet™ Strips are very effective for trimming and polishing composite materials. Evenwet™ Strips are 7" long and have a 1/2" open center that permits safe interproximal entry without destroying the interproximal tooth contact. Evenwet™ Strips are equally effective on wet and dry surfaces.

The #54 strips feature a fine grit on the white end and a medium grit on the gray end. #54 strips are also recommended for hygienist cleaning.

The #56 strips feature a fine grit on the white end and a micro-fine grit on the blue end.

REF	Description
471-87605	Strips, Dual Alo #54 Evenwet Narrow 2.5 mm
471-87610	Strips, Dual Alo #54 Evenwet Medium 3.18 mm
471-87615	Strips, Dual Alo #54 Evenwet Wide 4.00 mm
471-86120	Strips, Dual Alo #56 Evenwet Strips Fine/Microfine




471-87605
471-87610
471-87615

471-86120

SPECTRA-A STRIPS

Spectra-A Strips are made of resilient durable polyester fused with uniform aluminum oxide abrasive grit. The ½" open center of Spectra-A Strips permits safe interproximal entry without destroying the interproximal tooth contact. Spectra-A Strips are designed to be used on the hardest composite materials and are equally effective on wet and dry surfaces.

REF	Description
471-88110	Strips, Spectra-A Fine
471-88140	Strips, Spectra-A Medium
471-88160	Strips, Spectra-A Coarse
471-88100	Strips, Spectra-A Assorted



SPECTRA-F STRIPS

Spectra-F Strips are used for general purpose finishing and polishing. Spectra-F Strips have a ½" open center that permits safe interproximal entry without destroying the interproximal tooth contact. Spectra-F Strips are equally effective on wet and dry surfaces.

REF	Description
471-88020	Strips, Spectra-F Fine
471-88050	Strips, Spectra-F Medium
471-88080	Strips, Spectra-F Coarse
471-88000	Strips, Spectra-F Assorted



STRIP HOLDING FORCEPS

Strip holding forceps are designed to securely hold finishing and polishing strips and matrix bands during procedures, eliminating the need to hold the strip or band by hand. The 1x2 diamond-coated barrel tip gently but firmly clamps the strip, which allows for precise guiding. The forceps are sold as a pair and can be used in finishing adhesive restorations, polishing, positioning matrix bands, orthodontics and prophylaxis procedures.

REF	Description
72-60	Strip Holder Clips - 90 degree R.A.
72-62	Strip Holder Clips - 45 degree
72-70	Strip Holding Forceps, Pair



72-60




72-62



72-70

FEATURES AND BENEFITS of Thompson™ Cassettes

1. Unique flow-through design allows for effective rinsing and cleaning.
2. Medical grade silicone stabilizes the instruments when the cassette is closed.
3. Full length hinge provides stable and secure attachment to cassette.
4. Riveted stainless steel rack provides a secure, non degradable base on which to place instruments.
5. Quick and easy push button release provides a secure lock when closed.
6. Tabs allow for custom coding cassettes.
7. Double-hinged lid folds under for easy storage during use.
8. Self-lubricating plastic hinge pin will not corrode and provides smooth mechanical movement of the lid.
9. Various interior configurations to meet your storage needs.
10. Stainless steel provides aesthetic finish and corrosion protection.



Thompson™ Tuck-Away™ Lid

The Thompson™ Tuck-Away™ lid is a double-hinge lid designed to conveniently fold under the cassette, offering the following advantages:

- Lid cannot be misplaced during procedures.
- Saves valuable space during procedures.
- Provides incline for better instrument viewing and access during procedures.



Lid closed



Lid open



Lid tucked under

FACT:

All Thompson™ Cassettes come with a four year warranty on all moving parts and a **lifetime guarantee against corrosion.**

OPTIONS FOR CUSTOMIZATION

Lid Clips

Lid clips are available on several cassettes by adding "L" to the end of the part number. Please see the following part number charts for proper part number and the respective cassettes with enough room for this option.



Lid Clips

Stone Clips

Stone clips hold instrument sharpening stones in place. Stone clips are available on the same cassettes on which lid clips are optional by adding "S" to the end of the part number. Please see the following part number charts for proper part number and the respective cassettes with enough room for this option.



Stone Clips

Color Coded Tabs

All Thompson™ cassettes feature a tab where silicone tubing can be added. This feature gives your office the option of color coding by procedure, date or dental professional.

*See page 109 for part numbers



Color Coded Tabs

FACT:

To order an eligible cassette with both lid and stone clips, add "A" to the end of the part number. Please see the charts on the following pages for available cassettes

Thompson™ Cassettes

Thompson™ cassettes are manufactured with passivated stainless steel, feature the Thompson™ Tuck-Away™ lid and are available in a variety of sizes and configurations to meet all storage needs.

Thompson™ cassettes are available in several options:

Single Rack Cassettes — Cassette contains a single riveted rack to accommodate various hand instruments. Several sizes are available and hold 5 to 26 instruments.

Double Rack Cassettes — Cassette contains additional removable rack for double stacking instruments.

Utility Cassettes — Cassette does not contain an instrument rack, making it ideal for ring-handled and larger sized instrumentation.



Single Rack cassette



Double Rack cassette





Utility cassette



FACT:

Lid clips are ideal for ring handled instruments such as needleholders.



Thompson™ 5 Cassette

The Thompson™ 5 Cassette is perfect for simple procedures requiring few instruments. Its compact size features a single hinged lid, and it can hold up to 10 hand instruments with the double rack.

REF	Description	Dimensions	
		Inches	Millimeters
3-083005*	Thompson 5 Utility	8 x 3 13/16 x 1 1/8	204 x 81 x 30
3-083105*	Thompson 5 Single Rack	8 x 3 13/16 x 1 1/8	204 x 81 x 30
3-083210*	Thompson 5 Double Rack	8 x 3 13/16 x 1 1/8	204 x 81 x 30

* Thompson™ 5 cassette features single-hinged lid.

Thompson™ 7 Cassette

The Thompson™ 7 Cassette features the double-hinged, Tuck-Away™ lid for easy storage during procedures. It can hold 7 instruments with a single rack, 14 instruments with the double rack, or larger instruments without racks.

REF	Description	Dimensions	
		Inches	Millimeters
3-084007	Thompson 7 Utility	8 x 4 7/8 x 1 1/8	204 x 124 x 28.5
3-084107	Thompson 7 Single Rack	8 x 4 7/8 x 1 1/8	204 x 124 x 28.5
3-084214	Thompson 7 Double Rack	8 x 4 7/8 x 1 1/8	204 x 124 x 28.5

Thompson™ 9 Cassette

The Thompson™ 9 Cassette features the double-hinged, Tuck-Away™ lid for easy storage during procedures. It can hold 9 instruments with a single rack, 18 instruments with the double rack, or larger instruments without racks.

REF	Description	Dimensions	
		Inches	Millimeters
3-085009	Thompson 9 Utility	8 x 6 x 1 3/16	204 x 154 x 30
3-085109	Thompson 9 Single Rack	8 x 6 x 1 3/16	204 x 154 x 30
3-085218	Thompson 9 Double Rack	8 x 6 x 1 3/16	204 x 154 x 30

Thompson™ 10 Cassette

The Thompson™ 10 Cassette can hold 10 instruments and features the double-hinged, Tuck-Away™ lid for easy storage during procedures. This Cassette is also available with lid and/or stone clips to customize to your needs.

REF	Description	Dimensions	
		Inches	Millimeters
3-089110	Thompson 10	8 x 9 x 1 1/8	204 x 229 x 28.5
3-089110L	Thompson 10 w/lid clips	8 X 9 X 1 1/8	204 x 229 x 28.5
3-089110S	Thompson 10 w/ stone clip	8 X 9 X 1 1/8	204 x 229 x 28.5
3-089110A	Thompson 10 w/ lid and stone clips	8 X 9 X 1 1/8	204 x 229 x 28.5



Thompson™ 14 Cassette

The Thompson™ 14 Cassette can hold up to 14 instruments and features the double-hinged, Tuck-Away™ lid for easy storage during procedures. This Cassette is also available with lid and/or stone clips to customize to your needs.

REF	Description	Dimensions	
		Inches	Millimeters
3-118114	Thompson 14	8 x 10 7/8 x 1 3/16	204 x 276.5 x 30
3-118114L	Thompson 14 w/lid clips	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30
3-118114S	Thompson 14 w/ stone clip	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30
3-118114A	Thompson 14 w/ lid and stone clips	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30



Thompson™ 16 Cassette

The Thompson™ 16 Cassette can hold up to 16 instruments and features the double-hinged, Tuck-Away™ lid for easy storage during procedures. This Cassette is also available with lid and/or stone clips to customize to your needs.

REF	Description	Dimensions	
		Inches	Millimeters
3-118116	Thompson 16	8 x 10 7/8 x 1 3/16	204 x 276.5 x 30
3-118116L	Thompson 16 w/lid clip	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30
3-118116S	Thompson 16 w/stone clip	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30
3-118116A	Thompson 16 w/lid and stone clips	8 X 10 7/8 X 1 3/16	204 x 276.5 x 30



SLIMLINE CASSETTES

Thompson™ Slimline Cassettes feature a single-hinged lid and reduced profile. All Slimline cassettes are manufactured from passivated stainless steel that provides corrosion protection and an aesthetic finish.

Slimline 7

The seven instrument capacity of the Slimline 7 Cassette meets small cleaning and equipment requirements. Also available in extended length and deeper versions. The Slimline 7 and Deep Slimline 7 Cassette can be used in Statim™ autoclaves.



REF	Description	Dimensions	
		Inches	Millimeters
3-072007	Slimline 7	7 x 2 5/8 x 11/16	178 x 67 x 17
3-082008	Slimline 7 extended	8 x 2 5/8 x 11/16	203 x 127 x 17
STD5L27	Deep Slimline 7	7 x 2 5/8 x 1	178 x 67 x 25

Slimline 14

The Slimline 14 Cassette is ideal for exam, pedodontic prophylaxis/sealant, small prophylaxis and suture removal. The Slimline 14 Cassette will hold up to 14 instruments, and the extended width allows ring-handled instruments to be stored along with hand instruments.



REF	Description	Dimensions	
		Inches	Millimeters
3-072014	Slimline 14	7 x 4 5/8 x 11/16	178 x 115 x 17

SPECIALTY CASSETTES

In addition to the Slimline and Thompson™ series cassettes, Miltex also offers a variety of specialty cassettes designed for various applications.

Standard Cassette

The Standard Cassette stores up to 14 instruments and also contains an instrument bin for storing large and ring-handled instruments. The standard cassette also features a break-away lid.

REF	Description	Dimensions	
		Inches	Millimeters
STDSTAT	Standard Cassette	7 x 10 5/8 x 1 1/2	178 x 270 x 38



Standard 8 Cassette

The Standard 8 Cassette can store eight instruments and also contains an instrument bin for storing large and ring-handled instruments. The standard 8 cassette features the double-hinged lid.

REF	Description	Dimensions	
		Inches	Millimeters
STDSTAT8	Standard 8 Cassette	6 3/4 x 8 1/8 x 1 1/8	171 x 206 x 29



Standard 814 Cassette

The Standard 814 Cassette can store 14 instruments and also contains an instrument bin for storing large and ring-handled instruments. The standard 814 cassette features the double-hinged lid.

REF	Description	Dimensions	
		Inches	Millimeters
STDSTAT814	Standard 814 Cassette with 14 Rack	6 3/4 x 8 1/8 x 1 1/8	171 x 206 x 29





Twenty-Two Instrument Cassette

This 22 instrument cassette is ideal for storing large quantities of hand instruments and additional larger instruments.

REF	Description	Dimensions	
		Inches	Millimeters
STD222	22 Instrument Cassette	8 x 10 7/8 x 1 3/16	204 x 276.5 x 30

FACT:

Using cassettes for procedural setups is an excellent way to efficiently organize and clean instruments for specific procedures.



Orthodontic Cassette

The upright box design of the Thompson™ Orthodontic Cassette is perfect for storing orthodontic pliers, wire cutters, scissors and wire bending instruments.

REF	Description	Dimensions	
		Inches	Millimeters
STDORTHO	Ortho Cassette	5 x 7 x 2	127 x 178 x 51
STDORTHOS	Small Ortho Cassette	6 1/2 x 2 1/4 x 2	165 x 57 x 51



Handpiece/Accessory Cassette

This cassette is excellent for small specialized setups, handpieces and smaller handled instruments.

REF	Description	Dimensions	
		Inches	Millimeters
3-080205	Handpiece/ Accessory Cassette	8 x 2 3/16 x 1 7/16	204 x 71 x 37



Small Parts Box

Holds small items such as ortho brackets, long shank or lab carbide burs and rubber dam clamps. Features a hinged lid.

REF	Description	Dimensions	
		Inches	Millimeters
STDPB	Small Parts Box	1 2/3 x 1 1/3 x 1	42 x 33 x 26

CASSETTE CODE TUBING

REF	Description	Quantity per Pack
STDCC12BLACK	Cassette Code Tube, Black	12
STDCC12BROWN	Cassette Code Tube, Brown	12
STDCC12DKBLUE	Cassette Code Tube, Blue	12
STDCC12GOLD	Cassette Code Tube, Gold	12
STDCC12GRAY	Cassette Code Tube, Gray	12
STDCC12GREEN	Cassette Code Tube, Green	12
STDCC12LTBLUE	Cassette Code Tube, Light Blue	12
STDCC12PINK	Cassette Code Tube, Pink	12
STDCC12PURPLE	Cassette Code Tube, Purple	12
STDCC12RED	Cassette Code Tube, Red	12
STDCC12WHITE	Cassette Code Tube, White	12
STDCC12YELLOW	Cassette Code Tube, Yellow	12

REF	Description	Quantity per Pack
STDCC24BLACK	Cassette Code Tube, Black	24
STDCC24BROWN	Cassette Code Tube, Brown	24
STDCC24DKBLUE	Cassette Code Tube, Blue	24
STDCC24GOLD	Cassette Code Tube, Gold	24
STDCC24GRAY	Cassette Code Tube, Gray	24
STDCC24GREEN	Cassette Code Tube, Green	24
STDCC24LTBLUE	Cassette Code Tube, Light Blue	24
STDCC24PINK	Cassette Code Tube, Pink	24
STDCC24PURPLE	Cassette Code Tube, Purple	24
STDCC24RED	Cassette Code Tube, Red	24
STDCC24WHITE	Cassette Code Tube, White	24
STDCC24YELLOW	Cassette Code Tube, Yellow	24



Color Coded Tabs

SILITONE™ COLOR BANDS

Silitone™ color bands are made from medical-grade silicone and can withstand temperatures up to 450° F. Silitone™ bands are expandable from 4 mm to 12 mm and are available in a variety of colors. 50 per pack except kits.


REF	Description
VST	Kit 480 pieces, assortment
VSTA	Kit 120 pieces, assortment
VSTSA	GRAY
VSTSB	BLUE
VSTSD	GOLD
VSTSG	GREEN
VSTSK	BLACK
VSTSL	LIGHT BLUE
VSTSN	BROWN
VSTSP	PINK
VSTSR	RED
VSTSU	PURPLE
VSTSW	WHITE
VSTSY	YELLOW



TAPE 'N TELL™ ADHESIVE STRIPS

Ideal for color coding instruments and procedural kits.

REF	Description
017-48220	TAPE'N'TELL WHITE
017-48221	TAPE'N'TELL YELLOW
017-48222	TAPE'N'TELL RED
017-48223	TAPE'N'TELL BLUE
017-48224	TAPE'N'TELL GREEN
017-48225	TAPE'N'TELL BLACK
017-48226	TAPE'N'TELL ORANGE
017-48227	TAPE'N'TELL BROWN
017-48230	TAPE'N'TELL 8 ASST.



MILTEX® INSTRUMENT CARE & CLEANING PRODUCTS

The Miltex® care and cleaning product portfolio is designed to provide optimum cleaning performance and will help to prolong the service life of valuable surgical and dental instrumentation. Many products are available in multiple sizes, including single-dose packets for ease of use, space-saving 8 oz. bottles, and 1-gallon bulk containers with pump dispensers.

REF	Description
3-700	MILTEX Spray Lubricant, 8 fl. oz. spray bottle (0.24 liters) (case of 12)
3-705	MILTEX Lubricant, 1-gallon bottle (3.78 liters).
3-710	MILTEX Surgical Instrument Cleaner/Lubricant, 1-gallon bottle (3.78 liters).
3-720	MILTEX Surgical Instrument Cleaner, 8 fl. oz. bottle (0.24 liters). Contents of 8 oz. bottle makes 32 gallons of solution.
3-725	MILTEX Surgical Instrument Cleaner, 1-gallon bottle (3.78 liters). Contents of 1-gallon bottle makes 512 gallons of solution.
3-726	MILTEX Surgical Instrument Cleaner, Box 80 – 1/4 oz. (7.4 ml) packets. Contents of 1 box makes 80 gallons of solution.
3-740	MILTEX Surgical Instrument Stain Remover, 3 oz. (85g) jar (case of 12)
3-750	MILTEX EZ-Zyme®, All Purpose Enzyme Cleaner, Box 32 – 3/4 oz. (22 ml) packets. Contents of 1 box makes 32 gallons of solution.
3-755	MILTEX EZ-Zyme®, All Purpose Enzyme Cleaner, 1-gallon bottle (3.78 liters). Contents of 1-gallon bottle makes 170 gallons of solution.
3-760	MILTEX Instrument Prep Enzyme Foam, 24 oz. (0.71 liters) spray bottles (case of 12)
3-800	MILTEX Instrument Care System Kit includes one each of Spray Lube, Surgical Instrument Cleaner (8 oz.), Surgical Instrument Stain Remover, Instrument Care Video/DVD, Nylon Brush, Stainless Steel Brush.
3-810	MILTEX Instrument Care DVD, 25-minute instructional DVD designed to help practitioners and staff properly care for instruments
3-1000	MILTEX Nylon Instrument Cleaning Brushes, 7-1/4" (18.4 cm) (set of 3)
3-1001	MILTEX Stainless Steel Instrument Cleaning Brushes, 7-1/4" (18.4 cm) (set of 3)



MILTEX® EDGE SCISSOR SHARPENER

REF	Description
5-5000	MILTEX Edge Scissors Sharpener, not intended for serrated or SuperCut scissors



TIP-IT™ INSTRUMENT GUARDS

Tip-It™ Instrument Guards are used to protect instrument tips and edges from damage during instrument handling. 50 per pack.



REF	REF	REF	REF	Dimension	Size	Color
Clear Vented	Tinted Vented	Solid Vented	Solid Non-Vented		Code	
3-2501C	3-2501T	3-2501V	3-2501	1/16" X 3/4" (1.6mm X 19mm)	1	White*
3-2502C	3-2502T	3-2502V	3-2502	5/64" X 3/4" (2mm X 19mm)	2	Blue*
3-2503C	3-2503T	3-2503V	3-2503	7/64" X 3/4" (2.8mm X 19mm)	3	Green*
3-2504C	3-2504T	3-2504V	3-2504	1/8" X 1" (3.2mm X 25.4mm)	4	Red*
3-2505C	3-2505T	3-2505V	3-2505	3/16" X 1" (4.8mm X 25.4mm)	5	Yellow*
3-2506C	3-2506T	3-2506V	3-2506	13/32" X 3/4" (10.3mm X 19mm)	6	Maroon*
3-2507C	3-2507T	3-2507V	3-2507	1/16" X 3/8" X 1" (1.6mm X 9.5mm X 25.4mm)	7	Orange*
3-2508C	3-2508T	3-2508V	3-2508	5/64" X 5/8" X 1" (2mm X 16mm X 25.4mm)	8	Brown*
3-2509C	3-2509T	3-2509V	3-2509	1/8" X 1" X 1" (3.2mm X 25.4mm X 25.4mm)	9	Black*
3-2510C	3-2510T	3-2510V	3-2510	Assorted, equal number of sizes plus 5	-	-

* Color does not apply to clear vented pattern

STERILIZING INSTRUMENTS

Sterilizing & Grasping Forceps



3-2
3-4

65-25

REF	Description
3-2	DAVIS Sterilizer Forceps, 6" (15.2 cm)
3-4	DAVIS Sterilizer Forceps, 10" (25.4 cm)
3-10	Three prong Sterilizer Forceps, 8" (20.3 cm)
3-20	Utility-Sterilizer Forceps, 8" (20.3 cm), straight
3-22	Utility-Sterilizer Forceps, 8" (20.3 cm), curved
3-24	Utility-Sterilizer Forceps, 11" (27.9 cm), straight
3-26	Utility-Sterilizer Forceps, 11" (27.9 cm), curved
65-25	Grasping Forceps, non-magnetic, gold plated ends



3-10



3-20
3-24



3-22
3-26

Recommended Steps in Cleaning, Sterilization & Maintenance of Miltex® Surgical Instruments

1. HOLDING/PRESOAK

Never hold instruments in a dry container. Doing so allows blood and debris to dry onto instrument surfaces and makes cleaning more difficult. If rinsing and decontamination processes are not immediately available, pre-treat instruments or hold them in a neutral pH holding/presoak enzymatic solution after patient use but before actual cleaning. Miltex® Instrument Prep Enzyme Foam (3-760) is a ready to use foaming spray for pre-cleaning soiled instruments and scopes. Simply spray the foam on instruments until they are ready for processing. Also, Miltex® EZ-Zyme® (3-750 and 3-755) neutral pH all purpose multi-enzyme concentrate is ideal for instrument presoaking and pre-cleaning. As soon as possible, rinse and clean as follows:

2. RINSING

Always wear safety protection gear.

Immediately after surgery, remove organic materials by rinsing instruments under warm (not hot) running water. Rinse should remove most blood, fluids and tissue. Do not process dissimilar metals (stainless, copper, chrome plated, etc.) together.

Note: *Disinfection can be included to protect medical personnel from contamination during cleaning; however, protection also can be accomplished by following standard precautions by wearing appropriate PPE as promulgated by OSHA & AORN.*

3. CLEANING

All blood, dried body fluids and tissue should be completely removed from instruments prior to sterilization. Several methods for removing these materials are available.

- A. **Soak:** An enzymatic cleaner bath (soak) such as Miltex® EZ-Zyme® (3-750 and 3-755) or a solution of water and neutral pH (7) detergent such as Miltex® Surgical Instrument Cleaner (3-720, 3-725 and 3-726) is effective in removing organic material from instruments. Use distilled (demineralized) water, if possible. Instruments should be fully submerged for at least 10 minutes. Do not let "sharps" (scissors, knives, osteotomes, etc.) touch each other. Ensure that dissimilar metal instruments are separated. Rinse instruments under running water to remove solutions and change solutions frequently.
 - B. **Ultrasonic Cleaning:** Most instrument manufacturers recommend ultrasonic cleaning as the most effective way to clean surgical instruments, particularly instruments with hinges, box locks and other moving parts.
 - All instruments must be fully submerged in an open position using distilled (demineralized) water, if possible. To prevent possible surface scratching make sure that "sharps" do not touch other instruments. Also, separate instruments of dissimilar metals.
 - Process instruments for the full recommended ultrasonic cleaning cycle. Change solution frequently or as often as the manufacturer recommends.
 - Rinse instruments with water to remove the cleaning solution. Miltex® EZ-Zyme® (3-750 and 3-755) or Miltex® Surgical Instrument Cleaner (3-720, 3-725 and 3-726) can be safely used with all ultrasonic units.
 - C. **Automatic Washer Sterilizers:** Following manufacturers' recommendations, ensure that instruments are lubricated before the sterilization cycle and after the last rinse cycle.
- CAUTION:** Needle holders and forceps may crack if sterilized with ratchet in closed position.
- D. **Manual Cleaning:** If ultrasonic cleaning is not available, observe the following steps:
 - Use stiff nylon cleaning brushes (3-1000). Do not use steel wool or wire brushes except stainless steel wire brushes (3-1001) specially recommended for instrument serrated areas, bone files, burs or on stained areas of knurled handles.
 - Use only neutral pH (7) detergents such as Miltex® Surgical Instrument Cleaner (3-720, 3-725 and 3-726). If not rinsed off properly, low pH (acidic - less than 6 pH) detergents break down the stainless protective surface resulting in pitting and/or black staining on instruments. High pH detergents (alkaline - more than 8 pH) can cause brown stains (phosphate surface deposit) which can also interfere with the smooth operation of instruments. Most brown stains are not rust and are easily removed with Miltex® Surgical Instrument Stain Remover (3-740).
 - Brush delicate instruments carefully, and if possible, separate them from general instruments.
 - Make sure instrument surfaces are visibly clean and are free from stains and tissue. Miltex® Surgical Instrument Stain Remover (3-740) can help remove residue stains. This is also a good time to inspect each instrument for proper function and condition.
 - Check scissors' blades to ensure proper function. Blades should open and close smoothly. Test cutting performance at 3/4 length of the blade with the materials recommended below. Scissors should cut all the way to the tips. Recommended cutting test materials:
 - Fine/Delicate scissors: Surgical glove
 - Medium scissors: Single layer of stocking/cast netting
 - Large/Utility scissors: Double layer of stocking/cast netting
 - Check forceps (pickups) for proper jaw alignment. Teeth must meet properly without catching.
 - Check hemostats and needle holders to ensure jaw tips close in first ratchet position and that the entire jaw closes in third ratchet position. Check instruments for loose hinges and verify that they lock and unlock easily. Also check instruments for wear on jaw surfaces.

Recommended Steps in Cleaning, Sterilization & Maintenance of Miltex® Surgical Instruments

- Suction tubes should be clean inside.
- Test Biopsy Punches by punching a clean hole in 3-6 mil thick poly-bag material. If poly-bag material is not available, use tissue paper.
- Retractors should function properly.
- Cutting edge instruments and knives should be sharp and free of damage.
- After scrubbing, rinse instruments thoroughly under running water. While rinsing, open and close scissors, hemostats, needle holders and other hinged instruments to ensure that hinged areas are rinsed out and that no debris remains.

4. AFTER CLEANING

Separate dissimilar metals prior to sterilizing/autoclaving. If instruments will be stored, allow them to air-dry and store them in a clean and dry environment.

5. AUTOCLAVING

- A. Lubricate all hinged instruments** which have any "metal to metal" action at the screw or box lock. A non-silicone, water-soluble surgical lubricant such as Miltex® Lube (3-700 or 3-705) is recommended. Do not use industrial oils or lubricants.
- B. Sterilize instruments** either individually or in sets.
- **Individual Instruments:** Disposable paper or plastic pouches are ideal. Make sure to use a wide enough pouch for instruments with ratchet locks so instruments can be sterilized in an open (unlocked) position. Instruments locked during autoclaving can develop cracked hinges (box locks) or other problems because of heat expansion. If wrapping instruments, make sure the towel does not contain detergent residue which can stain instruments.
 - **Instrument Sets:** Unlock all instruments and sterilize in an open position. Place heavy instruments on the bottom of set (when two layers are required). Do not overload the chamber because an air pocket that hinders steam penetration may form.

CAUTION: With most portable tabletop autoclaves, unlock the door and open it no more than a crack, about 1/4" (6.4 mm) (at the end of the autoclave cycle and before the drying cycle). Then run the dry cycle for the period recommended by the autoclave manufacturer. If the autoclave door is fully opened before the drying cycle, cold air will rush into the chamber and will cause condensation on the instruments, potentially resulting in water stains or causing wet packs. Make sure autoclave filters and chambers are cleaned as recommended by the manufacturer. Miltex® Surgical Instrument Stain Remover (3-740) is effective in removing stains and in cleaning the autoclave chamber. If you have any unusual staining on your instruments after sterilization, obtain our Instrument Care DVD (3-810) or contact your local Miltex® representative.

6. CHEMICAL/COLD STERILIZATION

Most chemical/cold sterilization solutions render instruments sterile only after 10-hour immersion. This prolonged chemical action can be more detrimental to instruments than the usual 20-minute autoclave cycle. If the instruments need to be "disinfected" only, a chemical/cold sterilization soak is acceptable, as disinfection will take approximately 10 minutes or more. Check manufacturers' specifications. Also see our warning in using bleach (paragraph 3).

Keep in mind the difference between:

- Sterile - an absolute term (no living organism survives); and
- Disinfected - basically clean. Some organisms may survive. Always use the proper sterilization/cleaning technique to render the instrument in the required condition for use.

CAUTION: For instruments with tungsten carbide insert jaws, the use of chemical/cold sterilization solutions is not recommended as it may deteriorate the instrument's jaw.

We recommend these Miltex® products to properly safeguard and care for your quality instruments.

REF	Description	Use During
3-760	Instrument Prep Enzyme Foam 24 oz. spray	• Holding & presoak
3-750	EZ-Zyme®, All Purpose Enzyme Cleaner Box 32 - 3/4 oz. packets	• Holding & presoak • Soak cleaning
3-755	EZ-Zyme®, All Purpose Enzyme Cleaner 1-gallon pump bottle	• Ultrasonic cleaning • Manual cleaning • Laundry presoak
3-700	Spray Lube 8 oz. spray bottle	• Before autoclaving
3-705	Surgical Instrument Lubricant 1-gallon pump bottle	
3-710	Surgical Instrument Cleaner & Lubricant 1-gallon pump bottle	• Manual cleaning • Washer/decontaminator • Ultrasonic cleaning
3-720	Surgical Instrument Cleaner 8 oz. spray bottle	• Soak cleaning • Ultrasonic cleaning • Manual cleaning
3-725	Surgical Instrument Cleaner 1-gallon pump bottle	
3-726	Surgical Instrument Cleaner Box 80 - 1/4 oz. packets	
3-740	Surgical Instrument Stain Remover 3 oz. jar	• Stain removal
3-1000	Nylon Instrument Cleaning Brushes 3-pack	• Manual cleaning of instruments
3-1001	Stainless Steel Instrument Cleaning Brushes 3-pack	
3-800	Instrument Care System	• Complete instrument care
3-810	Instrument Care DVD	• Training personnel in the care of instruments

Adams Plier	12	Mathieu Plier	12
Adhesive Removing Plier	8	Matrix Scissors.....	14
Arch Bending Plier.....	11	Micro BOLEY Gauge.....	16
BAADE Band and Shell Remover Plier	4	Micro Diamond Strips	19
Babcock Suture Wire	13	Micro-Mini Ligature Cutter, 15°	6
Band Pusher	4	Micro-Mini Pin & Ligature Cutter	6
Band Seating Plier	4	Mini Pin & Ligature Cutter	6
Bird Beak Plier	10	Mini Three Jaw Plier	11
BOLEY Gauge	16	Nance Plier.....	12
Bracket Remover	8	Nylon Instrument Cleaning Brushes	30
Cassette Code Tube.....	29	O'Brien Plier	11
Cassettes	24-28	Old Style Looping Plier.....	12
Color Code Rings	29	Omega Loop Forming Plier	10
Coon Style Ligature Tying Plier.....	9	Parallel Action Wire Cutter, 4 1/2" Nickel Plated	7
Cotton Pliers.....	15	Pedo Crown-Festooning Scissors	14
Crown and Band Contouring Plier	4	PEESO, Collar Plier.....	10
Crown and Band Contouring Plier	9	Pin & Ligature Cutter	6
Crown and Shell Remover	9	Posterior Band Removing, Long Plier.....	4
Crown Gauge.....	16	Posterior Band Removing, Short Plier	4
DAVIS Sterilizer Forceps	32	Primary Crown Crimping Plier.....	9
Diagonal Wire Cutter, 5 1/2"	7	Rectangular Arch Bending Plier	11
Direct Bond Bracket Holder.....	9	REYNOLDS, Contouring Plier.....	10
Distal End Cutter, Flush Cut	5	Slim Weingart Plier	9
Dual Alo #54 Evenwet.....	19	Spectra-A.....	20
Edge Scissors Sharpener	31	Spray Lubricant	30
EZ-Zyme®, All Purpose Enzyme Cleaner.....	30	Spring Bow Divider	16
Force Module Separating Plier.....	9	Stainless Steel Instrument Cleaning Brushes	30
GORDON, Contouring Plier	4	Strip Holder Clips	21
Grasping Forceps.....	32	Strip Holding Forceps.....	21
Hard Wire Cutter, Straight.....	5	Surgical Instrument Cleaner	30
Hard Wire, 15° Cutter.....	5	Surgical Instrument Cleaner/Lubricant	30
Hollow Chop Contouring Plier	11	Surgical Instrument Stain Remover.....	30
Hollow Chop Plier	12	TAPE'N'TELL	30
How, Offset Plier.....	9	Three Jaw Plier.....	10
How, Straight Plier.....	9	Three Prong NiTi Bender	12
Instrument Care DVD	30	Three Prong Sterilizer Forceps	32
Instrument Care System Kit	30	Tip-It Instrument Guards	31
Instrument Prep Enzyme Foam.....	30	Tweed Loop Forming Plier.....	11
Jarabak Plier.....	10	Universal Cut & Hold Distal End Cutter.....	5
Lab and Office Plier	10	V-Stop Plier.....	12
Light Wire Bird Beak Plier	12	Weingart Plier.....	9
Light Wire Bird Beak with Groove.....	10	Wire Bending Plier.....	10
Lightning Strips	17-18	Wire Cutting Scissors.....	14
Lingual Arch Forming Plier	11	Wire Twister 6"	13
Lubricant	30	Wire Twister 8"	13
Marcotte Looping Plier.....	11		

017-48220	30	3-2509	31	9D-134	14	74-324	12	630-81450	17
017-48221	30	3-2510C	31	9D-135	14	74-330	12	630-81470	17
017-48222	30	3-2510T	31	22-500TC	13	66-311	15	630-81590	17
017-48223	30	3-2510V	31	22-502TC	13	66-312A	15	630-81610	17
017-48224	30	3-2510	31	65-25	32	66-317	15	630-81630	17
017-48225	30	3-083005	24	72-60	21	66-317A	15	630-81510	17
017-48226	30	3-083105	24	72-62	21	66-317B	15	635-82510	18
017-48227	30	3-083210	24	72-70	21	66-317S	15	635-82530	18
017-48230	30	3-072007	26	73-2504C	31	66-318	15	635-82550	18
3-2	32	3-072014	26	74-30	4	66-318-49	15	635-82570	18
3-4	32	3-080205	28	74-50	4	66-49	15	635-82590	18
3-10	32	3-082008	26	74-80	4	66-652	15	635-82610	18
3-20	32	3-084007	24	74-85	4	66-656	15	STD222	28
3-22	32	3-084107	24	74-135	4	68-693	16	STDCC12BLACK	29
3-24	32	3-084214	24	74-205	4	68-694	16	STDCC12BROWN	29
3-26	32	3-085009	24	74-208	4	68-695	16	STDCC12DKBLUE	29
3-700	30	3-085109	24	74-101	5	68-696	16	STDCC12GOLD	29
3-705	30	3-085218	24	74-101L	5	68-697	16	STDCC12GRAY	29
3-710	30	3-089110A	25	74-102	5	404-40413	19	STDCC12GREEN	29
3-720	30	3-089110L	25	74-103	5	404-40414	19	STDCC12LTBLUE	29
3-725	30	3-089110S	25	74-104	5	471-87605	19	STDCC12PINK	29
3-726	30	3-089110	25	74-105	6	471-87610	19	STDCC12PURPLE	29
3-740	30	3-118114A	25	74-106	6	471-87615	19	STDCC12RED	29
3-750	30	3-118114L	25	74-107	6	471-86120	19	STDCC12WHITE	29
3-755	30	3-118114S	25	74-107L	6	471-88110	20	STDCC12YELLOW	29
3-760	30	3-118114	25	74-109	6	471-88140	20	STDCC24BLACK	29
3-800	30	3-118116A	25	74-115	7	471-88160	20	STDCC24BROWN	29
3-810	30	3-118116L	25	74-206	8	471-88100	20	STDCC24DKBLUE	29
3-1000	30	3-118116S	25	74-219	8	471-88020	20	STDCC24GOLD	29
3-1001	30	3-118116	25	74-220L	8	471-88050	20	STDCC24GRAY	29
3-2501C	31	5-5000	31	74-20	9	471-88080	20	STDCC24GREEN	29
3-2501T	31	9-40-18	13	74-25	9	471-88000	20	STDCC24LTBLUE	29
3-2501V	31	9-44-18	13	74-30	9	620-81010	17	STDCC24PINK	29
3-2501	31	9-40-20	13	74-82	9	620-81030	17	STDCC24PURPLE	29
3-2502C	31	9-44-20	13	74-82L	9	620-81050	17	STDCC24RED	29
3-2502T	31	9-40-22	13	74-83	9	620-81210	17	STDCC24WHITE	29
3-2502V	31	9-44-22	13	74-95	9	620-81070	17	STDCC24YELLOW	29
3-2502	31	9-40-24	13	74-100	9	620-81090	17	STDORTHO	28
3-2503C	31	9-44-24	13	74-210	9	620-81110	17	STDORTHOS	28
3-2503T	31	9-40-25	13	74-211	9	620-81230	17	STDPB	28
3-2503V	31	9-44-25	13	74-212	9	620-81130	17	STDSL27	26
3-2503	31	9-40-26	13	74-10	10	620-81150	17	STDSTAT8	27
3-2504T	31	9-44-26	13	74-35	10	620-81170	17	STDSTAT814	27
3-2504V	31	9-40-28	13	74-40	10	620-81250	17	STDSTAT	27
3-2504	31	9-44-28	13	74-45	10	620-81290	17	V99D-133	14
3-2505C	31	9-40-30	13	74-55	10	620-81270	17	V99D-135	14
3-2505T	31	9-44-30	13	74-75	10	620-81310	17	V966-317	15
3-2505V	31	9-40-32	13	74-301	10	620-81190	17	V966-652	15
3-2505	31	9-44-32	13	74-303	10	625-82010	18	VST	29
3-2506C	31	9-40-34	13	74-305	10	625-82030	18	VSTA	29
3-2506T	31	9-44-34	13	74-306	11	625-82050	18	VSTSA	29
3-2506V	31	9-40-36	13	74-307	11	625-82070	18	VSTSB	29
3-2506	31	9-44-36	13	74-308	11	625-82090	18	VSTSD	29
3-2507C	31	9-120	14	74-309	11	625-82110	18	VSTSG	29
3-2507T	31	9-122	14	74-310	11	625-82130	18	VSTSK	29
3-2507V	31	9-124	14	74-311	11	625-82150	18	VSTSL	29
3-2507	31	9-124TC	14	74-312	11	630-81330	17	VSTSN	29
3-2508C	31	9-128	7	74-316	11	630-81350	17	VSTSP	29
3-2508T	31	9D-115	14	74-317	12	630-81370	17	VSTSR	29
3-2508V	31	9D-118	14	74-318	12	630-81530	17	VSTSU	29
3-2508	31	9D-121	14	74-319	12	630-81390	17	VSTSW	29
3-2509C	31	9D-123	14	74-320	12	630-81410	17	VSTSY	29
3-2509T	31	9D-132	14	74-321	12	630-81430	17		
3-2509V	31	9D-133	14	74-322	12	630-81550	17		



MILTEX WARRANTY

INSTRUMENTS AND PRODUCTS

Miltex hereby warrants all Miltex® products purchased directly from Miltex or an Authorized Dealer in accordance with the following terms and conditions:

All Miltex® instruments and products are guaranteed to be free from defects in material and workmanship when used for their intended surgical or dental purpose. Any Miltex® product that proves defective in workmanship or material will either be repaired or replaced, at Miltex's discretion, without charge. Products not purchased from an authorized source may not be authentic or meet Miltex's quality standards and may have been tampered with.

Other Instrument and Product Warranties:

- . Union Broach Chiropractic Mallet – 90 days
- . Union Broach Electrical Devices – 1 year except for Sprint TCM – 2 years
- . Cassettes - 4-year warranty on all moving parts and lifetime guarantee against corrosion

REPAIRS

Instruments and products repaired by the Miltex Instrument Repair Service are guaranteed to be free from defects in material and workmanship for 90 days when used for their intended surgical or dental purpose. Any Miltex repair that proves defective in workmanship or material within this 90-day period will either be repaired again or replaced, at Miltex's discretion, without charge. Instruments must be cleaned and sterilized prior to returning to Miltex for repair.

WARRANTY EXCEPTIONS

This warranty is void for Miltex® instruments and products serviced by any person or facility other than the Miltex Instrument Repair Services. Any custom made, custom etched or neutralized products are not returnable. Warranties are not valid for Miltex® instruments and products that prove defective as a result of improper care and cleaning or misuse. Instruments or products that are damaged in fire, flood or other acts or disasters are also excluded from warranty. Except as set forth in these Terms and Conditions, Miltex makes no other representations and warranties regarding its instruments or products, either expressed or implied. In addition, in no event will Miltex be liable for any incidental, special or exemplary damages or loss of profits in connection with the use of Miltex® instruments and products.



Thompson
Dental Manufacturing Co., Inc.

Moyco Union Broach.

www.miltex.com
toll-free phone 866-854-8300
toll-free fax 866-854-8400

**Please call Customer
Service today for additional
information, an in-office
demonstration, to place an
order, or to receive other
Miltex catalogs.**



**311 Enterprise Drive
Plainsboro, New Jersey 08536**

Miltex, Miltex Redefining Excellence and EZ-Zyme are registered trademarks of Miltex, Inc. a subsidiary of Integra LifeSciences Corporation.

EvenWet, Thompson, Tuck-Away, Tape 'N Tell, Tip It, Silitone and Moyco Union Broach are trade marks of Miltex, Inc.

The Integra wave logo is a trademark of Integra LifeSciences Corporation.

© Copyright 2009 Integra LifeSciences Corporation

