



Integra[®] Jarit[®]
Laparoscopic Instruments

Guarantee/Warranty.....	2
One-Piece Laparoscopic Instrumentation.....	3
Grasping Forceps.....	4
Dissecting Forceps.....	16
Biopsy Forceps.....	19
Scissors	20
Needle Holders, Receivers, and Ligation	22
Suction/Irrigation	25
Retracting Instruments	26
Electrosurgical Coagulation	28
Specialization	30
Take-Apart Laparoscopic Instruments	33
Take-Apart Laparoscopic Grasping Forceps	36
Take-Apart Laparoscopic Dissecting Forceps.....	45
Take-Apart Laparoscopic Biopsy Forceps.....	48
Take-Apart Laparoscopic Scissors	50
Jarit® Detach® System	53
Trocar/Cannula Systems	55
Storage and Accessories	57
Laparoscopes.....	58
3.5mm Laparoscopic Instruments	59
Replacement Parts Listing	63
Tip Comparison Chart.....	68
Numerical Index.....	72
Alphabetical Index.....	74

All Integra Jarit Instruments are unconditionally guaranteed against defects in material and workmanship. Integra Jarit will, at its option and without charge, either repair or replace any instrument which Integra Jarit determines to be defective in material or workmanship when used for its intended surgical purpose.

All instruments sold by Integra Jarit are warranted only to the original purchaser and only against defects in workmanship or materials which under the intended use render such instruments inoperable.

All Integra Jarit Carb-Bite and Integra Jarit Carb-Edge instruments are further warranted for a period of three (3) years under normal use. Should any adjustment become necessary during this warranty period, CARB-EDGE Scissors will be resharpened, and Carb-Bite Needle Holder jaw inserts will be replaced without charge.

All Integra Jarit Supercut Scissors are further warranted for a period of one year under normal use. Should any resharpening become necessary during this warranty period, Integra Jarit will provide this service without charge. A nominal charge will be made for sharpening after expiration of this warranty. Integra Jarit Supercut Scissors will be date-stamped for control.

The autoclavable feature of the Autoclavable Laparoscope is covered under a special one (1) year limited warranty.

This Guarantee/Warranty shall be void in the event any Integra Jarit instrument is used for purposes other than as originally intended, is not maintained or cleaned properly, is damaged as a result of misuse or accident or is repaired or altered by unauthorized persons.

The above is a limited warranty, and it is the only warranty made by Integra Jarit. Integra Jarit makes no other warranty, express or implied, including any warranty of merchantability or fitness for a particular purpose, as well as any express or implied warranty to patients. Integra Jarit shall have no liability for consequential, exemplary or incidental damages even if it is aware of the possibility of such damages. No warranty or guarantee may be created by any act or statement nor may this warranty be modified in any way, except as a result of a writing signed by an officer of Jarit. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct. The stated express warranty is in lieu of all liabilities or obligations of Integra Jarit arising out of, or in connection with, the delivery, use or performance of any Integra Jarit instrument. Integra Jarit neither assumes nor authorizes any representative or other person to assume for it any other liability in connection with Integra Jarit's products.

All products are subject to slight changes in dimensions and finishes due to changes in manufacturing standards.

Ordering Procedure

A written purchase order on customer's form may be requested for purchase of Integra products.

Acceptance of Order

Orders are accepted upon approval by Integra.

Returned Goods Policy

Products must be returned in unopened packages with manufacturer's seals intact to be accepted for replacement or credit unless returned due to a complaint of product defect. Determination of a product defect will be made by Integra. Products will not be accepted for replacement if they have been in the possession of the customer for more than 90 days.

Repairs and Maintenance

Should your instruments require repair or maintenance, contact Integra for return authorization and address. Instruments returned to Integra for repair must have a statement testifying that each instrument has been thoroughly cleaned and sterilized. Failure to supply evidence of cleaning and disinfection will result in a cleaning charge and delayed processing of your instrument repair.

Product Information Disclosure

Integra and its subsidiaries ("Integra") and manufacturer exclude all warranties, except Integra's applicable standard warranty whether expressed or implied, including but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Neither Integra nor manufacturer shall be liable for any incidental or consequential loss, damage, or expense, directly or indirectly arising from use of these products. Neither Integra nor manufacturer assume nor authorize any person to assume for them any other or additional liability or responsibility in connection with these products.

Care and Handling/IFU

For Care and Handling information and Instructions for use, please visit integralife.com/jarit.



Integra York PA, Inc.
589 Davies Dr.
York, PA 17402 USA
integralife.com/jarit




Miltex GmbH
Muehlstrasse 6
78604 Rietheim-Weilheim, Germany




Integra, the Integra logo, Jarit, Jarit Detach, and Jarit Stardust are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Carb-Bite, Carb-Edge, Jarit Articulator, Jarit Rotator, Jarit Supercut, P.E.E.R., Roto-Cam, and Roto-Lok are trademarks of Integra LifeSciences Corporation or its subsidiaries. Lap-Band is a registered trademark of Allergan, Inc.

One-Piece



One-Piece




Stringent quality standards are the most important aspect in the development of Integra Jarit instruments. Laparoscopic instruments are designed, manufactured and tested to provide surgeons with superior performance and dependability. High grade surgical stainless steel, used in the manufacture of all Integra Jarit laparoscopic instruments, offers high tensile strength and durability. For effective and efficient processing, instruments are made with internal seals. To confirm electrosurgical safety, electrosurgical instruments are tested per IEC International Electromechanical Commission standards. Integra Jarit offers laparoscopic instruments with a choice of features to meet both the preference of the surgeon and the economic demands of the hospital.


Rotation

The Rotator design provides comfortable, single-handed, bi-directional rotation of the distal tip. The Rotator Instruments offer exclusive features to meet the technical requirements of blunt dissection, sharp dissection and grasping. Precise German engineering, combined with handcrafted workmanship, has resulted in the features of the Roto-Lok design for grasping forceps and the Roto-Cam design for dissectors and scissors.

Roto-Lok: When the traction/countertraction technique is applied during dissection, it is important to maintain an effective grasp on tissue. If the instrument rotates freely, the distal tip could slip, making dissection difficult. With Integra Jarit's unique Roto-Lok feature, rotation locks into place once the tissue is grasped. The result is a secure grasp that will not slip or spin.



Roto-Cam: The Roto-Cam design for dissectors and scissors allows the surgeon to easily turn the rotating wheel, yet provides enough tension to prevent unnecessary movement. The design enables rotation in both open and closed positions to facilitate precise dissection.



Locking Mechanisms

Computer-generated production techniques employed in the precisely calibrated ratchet, hinge and locking mechanisms result in smooth, secure jaw closure.

Ratchet: The advanced design of the half-moon ratchet offers fine-teeth and a sturdy leaf-spring that is specially manufactured to maintain a firm hold.



Locking Lever: As an alternative to the ratchet, the sophisticated Integra Jarit locking lever engages and releases easily, yet maintains a firm hold. The locking lever offers the surgeon the ability to open and close the jaws freely, locking only when the surgeon chooses to engage the lever. The unique internal design adds built-in strength for a long, useful life.



Slide Lock: Featuring an ergonomic handle design, the slide locking mechanism is a different style of ratchet. It allows surgeons to operate the jaws freely until the sliding ratchet mechanism is engaged.



Integra Jarit grasping forceps feature a hinge mechanism to assure strength and durability. Precisely calibrated ratchets provide smooth, gradual closure of jaws. The easy-release locking lever offers the surgeon the ability to open and close the jaws freely, locking only when the lever is engaged.



Grasping Forceps


Straight, serrated jaws with proximal recess, double action.



600-100
Ratchet,
insulated shaft,
5mm, 332.5mm
12.5mm jaw



625-100
Roto-Lok,
ratchet,
insulated shaft,
5mm, 334mm
14mm jaw



600-105
Monopolar,
5mm, 332.5mm
12.5mm jaw

600-106
Monopolar,
5mm, 462.5mm
12.5mm jaw



610-100
Roto-Lok
Slide lock,
insulated shaft,
5mm, 321mm
14mm jaw




605-100
Slide lock,
insulated shaft,
5mm, 322mm
14mm jaw



625-100LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 334mm
14mm jaw



625-105
Roto-Cam,
monopolar,
5mm, 342.5mm
12.5mm jaw



TC Grasping Forceps


Uniquely designed grasping forceps feature double action Carb-Bite tips.



625-102
Roto-Lok,
ratchet,
insulated shaft,
5mm, 343mm




625-102LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 343mm



Strong TC Grasping Forceps

Heavy Carb-Bite inserts, double action.



625-101

Roto-Lok, ratchet, insulated shaft, 5mm, 326mm


625-103

Roto-Lok, ratchet, insulated shaft, 5mm, 377mm



625-103LL

Roto-Lok, locking lever, insulated shaft, 5mm, 377mm



Grabber Forceps

Long, fenestrated jaws with interlocking wave-like serrations, double action.




610-104

Roto-Lok, slide lock, insulated shaft, 5mm, 380mm



600-104IL

In-line, ratchet, insulated shaft, 5mm, 378mm



(Locking Lever Shown)

625-104


Roto-Lok, ratchet, insulated shaft, 5mm, 382mm

625-104LL



Roto-Lok, locking lever, insulated shaft, 5mm, 382mm

625-106LL

Roto-Lok, locking lever, insulated shaft, 5mm, 462mm



Fundus Grasping Forceps
Long, beveled-edge jaws.




(Locking Lever Shown)

625-107DA
Double action,
Roto-Lok,
locking lever,
insulated shaft,
5mm, 384mm

625-107RD
Double action,
Roto-Lok,
ratchet,
insulated shaft,
5mm, 384mm


625-107D5
Double action,
Roto-Lok,
locking lever,
insulated shaft,
5mm, 464mm



Dual-Cup Grasping Forceps
Double cup design
increases surface area
contact, double action.



610-107
Roto-Lok,
slide-lock,
insulated shaft,
5mm, 381mm



610-107DA
Double action,
Roto-Lok,
slide lock,
insulated shaft,
5mm, 382mm



625-107
Roto-Lok,
ratchet,
insulated shaft,
5mm, 383mm


625-107LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 383mm



600-107IL
Double action,
in-line, ratchet,
insulated shaft,
5mm, 380mm



610-108
Roto-Lok,
slide lock,
insulated shaft,
5mm, 326mm




625-108
Roto-Lok,
ratchet,
insulated shaft,
5mm, 338mm




625-108LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 338mm




600-108IL
In-line, ratchet,
insulated shaft,
5mm, 374mm


Retraction Grasping Forceps


Long-ridged jaws with 2x3 blunt interdigitating teeth.


Atrau-Allis Grasping Forceps


Unique jaw design offers strong atraumatic grasp, double action.



Atrau-Allis Long Jaw Grasping Forceps

Unique jaw design offers strong atraumatic grasp, double action.



Allis Grasping Forceps

Interlocking distal teeth with fine-serrated inner jaws, double action.






Retrograde Grasping Forceps
Retrograde, deep-serrated teeth, double action.



600-116
Roto-Lok, ratchet, insulated shaft, 5mm, 331mm 15mm jaw


600-118
Roto-Lok, ratchet, insulated shaft, 5mm, 462mm 15mm jaw




Grasping and Dissecting Forceps
Cup-shaped, toothed jaws, double action.




600-120
Monopolar, 5mm, 328mm



Winer Grasping Forceps
Cross-serrated, disc-shaped jaws with wide opening between the jaws, in-line ratcheted handle.




600-123
In-line, ratchet, insulated shaft, 5mm, 344mm



Grasping and Extracting Forceps
Retrograde, deep-serrated teeth.



600-130
Spring-loaded, 10mm, 325mm




Claw Forceps

2x3 interlocking teeth for penetrating hold and extraction of tissue.



600-135

Spring-loaded,
10mm, 330mm




Claw Forceps

2x3 interlocking teeth for penetrating hold and extracting tissue.




600-137

Spring-loaded,
insulated shaft,
5mm, 325mm



DJB Grasping Forceps

Wide fenestration and straight fine serrations of single action jaws.



(Locking Lever Shown)

625-139

Roto-Lok,
ratchet,
insulated shaft,
5mm, 379mm

625-139LL

Roto-Lok,
locking lever,
insulated shaft,
5mm, 379mm



DeBakey Forceps

Fenestrated 40mm long jaws with DeBakey-type teeth, double action.



610-140

Roto-Lok, slide lock,
insulated shaft,
5mm, 398mm

610-140L

Roto-Lok, slide lock,
insulated shaft,
5mm, 474mm




600-140

Ratchet,
insulated shaft,
5mm, 384mm



625-140LL

Roto-Lok,
locking lever,
insulated shaft,
5mm, 480mm




DeBakey Forceps, Curved

48mm long jaws with DeBakey-type teeth, double action.



625-141

Roto-Lok,
ratchet,
10mm, 373mm




Strong Hold Babcock Forceps

Traditional design, double action.



610-149
Roto-Lok, slide lock, insulated shaft, 5mm, 379mm




Paddle Babcock Forceps, 10mm

Flat paddles with fine, straight serrations, double action.




625-142
Roto-Lok, ratchet, 10mm, 370mm



Paddle Babcock Forceps, 5mm

Flat paddles with fine, straight serrations, double action.




625-143
Roto-Lok, ratchet, insulated shaft, 5mm, 393mm

625-143LL
Roto-Lok, locking lever, insulated shaft, 5mm, 393mm

625-146LL
Roto-Lok, locking lever, insulated shaft, 5mm, 473mm



600-143IL
In-line, insulated shaft, ratchet, 5mm, 389mm




Babcock Forceps

Traditional style Babcock jaws designed with six rows of fine serrations, double action.




625-144
Roto-Lok, ratchet, 10mm, 370mm



Babcock Forceps, 5mm

Traditional design, double action.

600-147
Ratchet,
insulated shaft,
5mm, 350mm




610-147
Roto-Lok,
slide lock,
insulated shaft,
5mm, 388mm



625-147
Roto-Lok,
ratchet,
insulated shaft,
5mm, 390mm

625-147LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 390mm

(Ratchet Shown)


Babcock Forceps, Round Jaw


Rounded jaws with 2x3 rows of DeBakey-type teeth, double action.



600-145
Ratchet,
10mm, 346mm




625-145
Roto-Lok,
ratchet,
10mm, 330mm


Swanström Grasping Forceps

Unique design of 1x2 DeBakey-type teeth with triangular fenestration, double action.



610-148
Roto-Lok,
slide lock,
insulated shaft,
5mm, 376mm




625-148
Roto-Lok,
ratchet,
insulated shaft,
5mm, 378mm

625-148LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 378mm

(Ratchet Shown)



600-148IL
In-line,
insulated shaft,
ratchet,
5mm, 374mm




Duckbill Grasping Forceps

Wide fenestration and straight serrations of jaws.



625-118
Roto-Cam,
monopolar,
5mm, 350mm



Duval Grasping Forceps


Triangular jaws with deep-serrated surfaces, double action.



600-150
Ratchet,
10mm, 331mm



625-150
Roto-Lok,
ratchet,
10mm, 335mm




Russian Grasping Forceps

40mm long jaws with heavy serrations, double action.



600-154
Ratchet,
10mm, 331mm




Anvil Grasping Forceps

Jaws specially shaped to grasp and hold the anvil of a circular stapler.



625-154
Roto-Lok,
ratchet,
10mm, 370mm




**Foerster
Sponge Forceps**

45mm jaws,
double action.



600-158
Ratchet,
10mm, 335mm




Hunter Clamp

Right-angle shaft with
Atrau-Allis teeth,
double action.




600-160
Ratchet,
5mm, 321mm




Hunter Bowel Grasper

Fenestrated jaws with
DeBakey-type teeth,
double action.




610-163
Roto-Lok,
slide lock,
insulated shaft,
5mm, 330mm

610-164
Roto-Lok,
slide lock,
insulated shaft,
5mm, 381mm



600-164IL
In line,
non-ratcheted
insulated shaft,
5mm, 379mm

610-165
Roto-Lok,
slide lock,
insulated shaft,
5mm, 457mm



(Locking Lever Shown)

625-163
Roto-Lok,
non-locking,
insulated shaft,
5mm, 322mm


625-164
Roto-Lok,
non-locking,
insulated shaft,
5mm, 383mm

625-164LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 383mm

625-165LL
Roto-Lok,
locking lever,
insulated shaft,
5mm, 463mm



610-166
Roto-Lok,
slide lock,
insulated shaft,
5mm, 403mm,
45mm jaw length



Tenaculum Forceps

Single-toothed
long jaws,
double action.




605-152
Slide Lock,
insulated shaft,
5mm 339mm
31mm jaw




615-150
Ratchet,
insulated shaft
10mm, 363mm
37mm jaw

615-152
Ratchet
insulated shaft
5mm, 350mm
31mm jaw



Strong Forceps

Fenestrated, serrated,
oval-shaped jaws that
converge at distal end,
double action.




615-110
Monopolar,
5mm, 330mm
12mm jaw

615-111
Monopolar,
5mm, 462mm
12mm jaw




615-109
Roto-Lok,
5mm, 329mm
13.5mm jaw




Stone Extractor

To facilitate stone
removal once neck of
gallbladder has been
externalized.




600-520
7 1/2" (193mm)




Lobe Spleen Extractor

Open cupped jaws
for extraction of the
spleen during
laparoscopic
splenectomy.



600-522
8 3/4" (221mm)


Strong Atraumatic Grasping Forceps


Fenestrated, serrated, curved jaws with split-steel construction.


615-113

Spring handle,
5mm, 340mm

615-114

Spring handle,
5mm, 480mm


Grasping Forceps, Fallopian Tube


Curved jaws with split-steel construction.


615-116

Spring handle,
5mm, 340mm

615-118

Spring handle,
5mm, 460mm


Bowel Forceps

Cross-serrated, disc-shaped jaws with split-steel construction.



615-143

Spring handle,
5mm, 336mm

615-144

Spring handle,
5mm, 456mm

Integra Jarit dissecting forceps are specially designed with the varied needs and preferences of the surgeon in mind — from fine, precise dissection around delicate tissue to gross dissection of large structures.




Micro Dissecting Forceps


Delicate, fenestrated, serrated jaws, double action.



600-110
Monopolar,
5mm, 332mm



625-110
Roto-Cam, monopolar,
5mm, 337mm



Needle Nose Forceps


Slender, tapered and fine-serrated jaws.



600-112
Monopolar,
5mm, 330mm
18.5mm jaw



625-112
Roto-Cam,
monopolar,
5mm, 345mm
19mm jaw



Tapered Dissecting Forceps


Tips of serrated, tapered jaws meet at fine point, double action.



600-115
Monopolar,
5mm, 332.5mm



625-115
Roto-Cam,
monopolar,
5mm, 342.5mm




Micro Grasping Forceps

Slender, tapered jaws for dissection of small structures.



600-117
Monopolar,
5mm, 330mm




Bullet-Nose Dissecting Forceps

Bullet-shaped, tapered, serrated jaws for dissection, double action.



625-117
 Roto-Cam,
 monopolar,
 5mm, 341mm




Tapered Maryland Dissector

Curved, finely-tapered, serrated jaws with wide spread, double action.



625-124
 Roto-Cam,
 monopolar,
 5mm, 345.5mm




Maryland Dissector

Slender, full-curved, long, serrated jaws with wide spread, double action.



625-125
 Roto-Cam,
 monopolar,
 10mm, 339mm



Maryland Dissectors

Curved, tapered, serrated jaws with wide spread, double action.



600-126
 Monopolar,
 5mm, 326mm




600-125
 Monopolar,
 10mm, 325mm



625-126
 Roto-Cam,
 monopolar,
 5mm, 346mm

625-127
 Roto-Cam,
 monopolar,
 5mm, 450mm




Mixer Dissecting Forceps

Slender, full-curved, long, serrated jaws with wide spread, double action.



600-127
10mm, 335mm



Mixer Spreader-Dissectors

Jaw design offers a true right-angle for dissection around delicate structures, double action.




625-128
Roto-Cam, monopolar, 10mm, 320mm 26mm jaw



600-128
10mm, 322mm 27mm jaw



625-128RH
Roto-Cam, monopolar, ring-handle, 10mm, 319mm 26mm jaw



Mixer Spreader-Dissectors

Jaw design offers a true right-angle for dissection around delicate structures; double action.




625-129
Roto-Cam, monopolar, 5mm, 320mm



625-129RH
Roto-Cam, monopolar, ring-handle, 5mm, 318mm

625-130RH
Roto-Cam, monopolar, ring-handle, 5mm, 463mm




Right Angle Dissector

Cross-serrated jaws with full right angle; double action



625-132
Roto-Cam, monopolar, 10mm, 360mm




Appel Dissecting Hook

Right-angle tip allows blunt dissection and isolation of structures.




600-500
10mm, 406.5mm



Biopsy Forceps

Two sharp teeth, one on each jaw, for biopsy or hold of thick-walled tissue.




615-160

1x1 teeth,
Spring-loaded,
5mm, 330mm
11.5mm jaw


615-162

1x1 teeth,
Spring-loaded,
5mm, 450mm
10mm jaw



Biopsy Punches

Serrated lower jaw.




615-164

Monopolar,
5mm, 330mm


615-166

Monopolar,
5mm, 450mm



Hook Biopsy Punches

Frangenheim-type electro-surgical biopsy punch with hook blade.




615-168

Monopolar,
5mm, 320mm


615-169

Monopolar,
5mm, 450mm



Biopsy Punches

Single tooth in upper jaw.




615-170

Monopolar,
5mm, 330mm
9.5mm jaw

615-172

Monopolar,
5mm, 450mm
12mm jaw



Biopsy Forceps

Oval-shaped fenestrated jaws.




615-174

Monopolar,
5mm, 325mm

615-176

Monopolar,
5mm, 445mm



Micro Biopsy Forceps

For biopsy of small, delicate tissues.




615-178

Spring-loaded,
5mm, 330mm

The Integra Jarit scissors line was developed by utilizing advanced metallurgy, as well as innovative methods of grinding and sharpening. The result — scissors that cut cleanly and stay sharp.

The exclusive Integra Jarit Supercut blade is razor sharp and glides effortlessly through tissue with a clean, precise cut. The finely-serrated lower blade helps prevent tissue slippage.



Supercut Scissors


Razor-sharp upper blade with finely-serrated lower blade.



600-200 5mm, 330mm
600-201 5mm, 450mm



625-200
Roto-Cam,
monopolar,
5mm, 343mm




Supercut Endo Scissors, Straight

Electrosurgical instrument with straight, Integra Jarit Supercut blades and precise cutting action offers mechanical cutting, electrosurgical cutting, blunt dissection, electrosurgical coagulation—all in one reusable, cost-effective instrument; double action.



625-201
Roto-Cam
monopolar,
5mm, 343.5mm




Supercut Endo Scissors, Curved


Electrosurgical instrument with curved, Integra Jarit Supercut blades and precise cutting action offers mechanical cutting, electrosurgical cutting, blunt dissection, electrosurgical coagulation—all in one reusable, cost-effective instrument; double action.



600-202 Monopolar, 5mm, 330mm
600-203 Monopolar, 5mm, 465mm




625-202 Roto-Cam, monopolar, 5mm, 344mm
625-203 Roto-Cam, monopolar, 5mm, 454mm



Hook Scissors


Blades designed to lift, isolate and transect tissue.



600-205 Monopolar, 5mm, 320mm
600-206 Monopolar, 5mm, 450mm
600-207 Non-insulated, Non-monopolar 5mm, 325mm




625-205
Roto-Cam,
monopolar,
5mm, 330mm

**Micro Scissors**


Fine, delicate, 3mm blades for small incision.

**600-210**

Curved Left,
5mm, 333mm

**600-212**


Straight,
5mm, 333mm

**Mayo Operating Scissors**

Large, sturdy blades for incising heavy tissue.

**600-220**


10mm, 330mm

**Metzenbaum Scissors**

Long, curved, tapered blades with wide spread, double action.

**625-230**

Roto-Cam,
monopolar,
5mm, 338mm

**Mini Scissors**

Curved, tapered, delicate blades, double action.

**625-240**



Roto-Cam,
monopolar,
5mm, 335mm

The difficult task of laparoscopic suturing is more manageable with Integra Jarit's Ligation Line of unique instruments for both intracorporeal and extracorporeal suturing. The line features the Appel Needle Holders, which give the surgeon improved control of curved needles. These high quality laparoscopic instruments offer an ergonomically designed handle that places the surgeon's hand in line with the needle, resulting in a natural suturing motion, along with a leaf-spring mechanism that is engaged with speed and ease. The curved, tapered jaws permit fine control of the needle together with improved visibility at the suture site. Integra Jarit's Carb-Bite jaws allow for a firm, secure grip on the needle and long, useful life of the instrument. Alternative needle holder options include the In-line ring-handle Needle Holder and the Self-Righting Romeo Needle Holder.


Also featured are the Lyons Needle Forceps and Juliet Receiving Forceps. These unique instruments are specially designed to work in conjunction with the needle holder to achieve an efficient laparoscopic suturing technique. Functioning as a tissue forceps, the tapered jaws grasp and retract tissue for needle placement without obstructing visualization. Functioning as a receiver, the driven needle is grasped at the tip. The design then allows the needle to be easily repositioned within the jaw for re-grasping by the needle holder at the appropriate angle.

Appel Needle Holders


Tapered jaws stabilize the needle while allowing better visibility at the suture site; Carb-Bite jaws are designed for a firm, secure grip on the needle and long life of the instrument.




600-250
Left Curve,
5mm, 330mm




600-252
Right Curve,
5mm, 330mm



600-255
Up-Angle,
10mm, 330mm



600-260
Straight,
5mm, 330mm





600-264
Micro, Left Curve,
5mm, 330mm


600-249
Left Curve, 5mm, 450mm

In-line Carb-Bite Needle Holders


In-line ring handle offers a handle design similar to that used in open surgery with the same superior hold of needles.



600-251
Left Curve,
5mm, 330mm




600-253
Right Curve,
5mm, 330mm




600-261
Straight,
5mm, 330mm

600-254
Left Curve, 5mm, 450mm




Romeo Self-Righting Needle Holder

Up-curved, Jarit Stardust jaws self-right curved needles; needles held very securely at any angle within the jaws.




600-266
Up-curve,
5mm, 325mm

600-267
Up-curve,
5mm, 450mm




Lyons Needle Forceps

Tapered, cross-serrated jaws grasp and retract tissue for needle placement without obstructing visualization; allows the needle to be easily repositioned within the jaw for re-grasping by the needle holder at the appropriate angle; in-line handle places hand in same line as that with needle holder for more efficient suturing process, double action.



600-268
5mm, 330mm




Juliet Receiving Forceps

Right curved jaws with contained hinges facilitate intracorporeal knot tying; cross-serrated jaws for tissue control during suturing.




600-269
Curved right,
5mm, 330mm

600-271
Curved right,
5mm, 450mm



Needle Holders

Straight serrations on jaws designed to grasp straight needle securely.



600-273
Spring-loaded,
3mm, 330mm
7mm jaw



600-275
Spring-loaded,
5mm, 330mm
9mm jaw

Suturing Set

The Integra Jarit Suturing Set offers instruments required for endoscopic suturing. The standard set, provided in a wooden case includes:

Appel Needle Holder (600-250)
Lyons Needle Forceps (600-268)
Knot Pusher (600-289).

600-270



Appel Knot Pushers


Unique design facilitates the guiding of extracorporeal knots to suture site; knots can be tied and passed without repositioning the knot pusher; the twisting of suture that is experienced with standard knot pushers is minimized.

600-285

10mm,
410mm total length

600-287

5mm,
410mm total length



Knot Pusher

Classic design for guiding of extracorporeal knots to suture site.

600-289

4mm,
350mm total length




Endoloop Applicator

Allows placement of endoloop suture intraperitoneally.

600-240

5mm to 3mm, 175mm






Suction/Irrigation Cannulas

Dual trumpet valves for suction and irrigation; laser guide (fiber insert) and three additional seals included.




600-320 5mm, 350mm
600-322 5mm, 280mm




Suction/Irrigation Set

Set includes handpiece, 10mm cannula for aspiration of hematomas, gallstones, and viscous drainage, and 5mm cannula for standard suction/irrigation; suction connection straight off end minimizes clogging of cannula and sticking of valves and allows for stronger vacuum pressures; 3.5mm cannulas available.



600-340
Suction/irrigation cannula
10mm, 330mm color coded
Includes 5mm and 10mm cannulas



Puncture Needles with Irrigation

Large bore needle (19 Gauge to 17 Gauge) with Luer-Lok fitting for irrigation or aspiration.




600-490 5mm, 330mm
600-491 5mm, 450mm

Integra Jarit Articulators offer unique articulating instruments for retraction and dissection of veins, nerves, esophagus, and other structures. With the open/close action of the handle, the distal arm smoothly angles from 0° to 90° improving the speed and ease of articulation.

Articulators


Articulate rapidly from 0° to 90°; bi-directional ROTO-LOK wheel.



Herron Articulator Bar Retractor

Straight arm for dissection and retraction; keyhole for holding suture (e.g., anvil placement).


600-580
Roto-Lok,
locking lever,
5mm, 405mm



Vein/Nerve Retractor

U-shaped distal tip for atraumatic retraction of veins, nerves, etc.


600-582
Roto-Lok,
locking lever,
10mm, 365mm



Articulator Retractor

S-shaped distal tip for atraumatic retraction of the esophagus; keyhole for maneuvering penrose drain or LAP-BAND® System.


600-584
Roto-Lok,
locking lever,
10mm, 410mm




Articulator Retractor

S-shaped distal tip for atraumatic retraction of the esophagus; keyhole for maneuvering penrose drain or LAP-BAND® System.

600-585
Roto-Lok,
locking lever,
10mm, 420mm


P.E.E.R. Retractors

Atraumatic, telescoping blades allow for elevation, retraction and mobilization of organs and tissues, while permitting optimum visualization of surgical field; fully rotatable shaft and locking mechanism combine for versatile instrument placement with stationary positioning; bar ratchet with lock secures jaws in open position; black ring at distal shaft provides identifiable mark to assure jaws have cleared cannula.


600-550

Roto-Lok,
ratchet,
5mm, 320mm
48mm retractor jaw length


600-560

Roto-Lok,
ratchet,
10mm, 320mm
68mm retractor jaw length


600-540
Fan Retractor


Working length 380mm
Five Finger Liver Elevator
with up to 60° spread and
80° Lift Functions, 10mm
58mm fan blade length, in
closed position.

For surgeons who prefer electrosurgical dissection, Integra Jarit has designed electrosurgical instruments to meet their specific needs. On the monopolar electrodes, the distal tip is partially insulated to assure that electricity is delivered to the tissue the surgeon has targeted. The Integra Jarit coagulator handle offers a large, comfortable surface for the surgeon's hand. The channel and trumpet valve allow quick and easy aspiration of smoke. For the surgeon who does not require smoke aspiration, cost-effective, fully insulated electrodes are available.

Bipolar instruments are designed with micro tips for effective coagulation of fine, superficial vessels, and with grasping action for electrosurgical coagulation of blood vessels.

(Note: The bipolar grasping forceps are not for use with tubal ligation.)

Suction Coagulators




Spatula Tip Suction Coagulator

Smooth, arched, partially insulated spatula tip.

600-300

5mm, 330mm




L-shaped Hook Suction Coagulator

Partially insulated right angle tip.

600-305

5mm, 330mm




Hook End Suction Coagulator

Partially insulated hook tip.

600-310

5mm, 330mm




Blunt Suction Coagulator

Blunt, open-ended distal tip.

600-315

5mm, 330mm

Electrodes



Spatula Electrodes


Smooth, arched, partially insulated spatula tip.

600-317

5mm, 330mm

600-317L

5mm, 450mm



L-shaped Hook Electrodes


Partially insulated right angle tip.

600-318

5mm, 330mm

600-318L

5mm, 450mm



Hook Electrodes


Partially insulated hook tip.

600-319

5mm, 330mm

600-319L

5mm, 450mm




Ball Electrode

Ball tip.

600-316

5mm, 370mm

Bipolar Forceps



Micro Bipolar Forceps


Fine, pointed distal tip for coagulation of surface blood vessels.

615-330

5mm, 336mm

615-331

5mm, 446mm



Large Bipolar Forceps


Flat 3mm wide jaws for bipolar coagulation of blood vessels.

615-334

5mm, 347mm

615-335

5mm, 457mm



Grasping Bipolar Forceps

Wave-like jaw for bipolar coagulation of blood vessels.

615-338


5mm, 343mm

615-339

5mm, 453mm

Electrosurgical Accessories

IMPORTANT: When disconnecting electrosurgical cables from the generator or instrument, always grasp the cables by the connector only. Never pull the cable by the cord. Do not use instruments or cables if the insulation is not fully intact.



Monopolar Cable

600-290


10' (3m)



Monopolar Adapter

600-291

male to female



Bipolar Cable

615-295

10' (3m)




Bipolar Cable, Wolf Generator

615-296B

10' (3m)


Integra Jarit offers a diversified selection of specialized instrumentation for cholangiography, morcellation, and tissue manipulation.



Spoon Forceps
Wide cup-shaped jaws
for tissue or stone
removal.




600-138
Spring-loaded,
10mm, 335mm



Cholangiogram Clamp
Channel for placement of
5Fr or smaller cholangiogram
catheter.




600-400
Ratchet,
5mm, 330mm



Cholangiogram Guide
Curved tip for simple placement
of cholangiogram catheter.




600-410
5mm, 320mm total length


Dilation Set

For dilation of 5mm port to 10mm diameter.



600-475
5mm to 10mm



Veress Needles

Internal safety shield allows initial puncture and introduction of CO₂ for pneumoperitoneum, 14 gauge needle.



600-480
120mm

600-481
150mm



Palpation Probes

For atraumatic manipulation of tissue.



600-495
5mm, 300mm

600-496
5mm, 450mm


Myoma Screws

For manipulation and securing of thick-walled myoma.




615-210
5mm, 330mm

615-212
10mm, 330mm



Take-Aparts



Take-Aparts

Take-Aparts

Integra Jarit Take-Apart Laparoscopic Instruments were designed with the surgeon's flexibility needs in mind. Handles are available with and without ratchet mechanisms and with and without rotation. The shafts come in 350mm and 425mm lengths. Take-Apart instruments may be ordered separately by component as well as in the assembled configuration. These features, along with many grasping, dissecting, cutting, and biopsy tip options allows the surgeon to select the instrument that will meet their needs.

Handles



Passive Ratchet Mechanism

Precise teeth easily engage and disengage for surgical flexibility.

Monopolar Connection

All handles have monopolar connection for eletrosurgical flexibility.



642-200
w/ Rotation, w/ Ratchet



642-225
w/o Rotation, w/ Ratchet



642-250
w/ Rotation, w/o Ratchet



642-275
w/o Rotation, w/o Ratchet

Shafts



643-350
5mm, 350mm



643-425
5mm, 425mm



643-650
10mm, 350mm



643-725
10mm, 425mm

Cleaning Brush

643-100
For all shafts


Directions for Disassembly and Re-Assembly

Integra Jarit Take-Apart Laparoscopic Instruments are reusable, non-sterile, fully-assembled, monopolar electrosurgical instruments comprised of three individual parts – handle, shaft and stainless steel insert. Disassembly is completed in a few short steps; however, inserts must always be disengaged prior to removing the shaft from the handle or damage may result. For re-assembly, the shaft must be connected to the handle for the insert to be seated fully within the handle assembly.

A. DISASSEMBLY

Step 1:

For handles with ratchets: Push gold-plated lever down to detach/disassemble ratchet connection. For handles without ratchets: Go directly to Step 2.



Step 2:


Note: Insert should be disengaged from the handle prior to unscrewing the shaft or damage to the handle may result.

Press and hold rounded button on the top leg of the handle. Open handle by lifting until an audible pop is heard and push button is visibly tucked underneath the handle assembly.




Step 3:

Insert is now disengaged from handle and may be removed by grasping distal end and pulling from shaft.



Step 4:


For handles with rotation: Hold rotating wheel in place, loosen screw connection nut, then separate shaft from handle. For handles without rotation: Hold handle in place, loosen screw connection nut, then separate shaft from handle.



B. RE-ASSEMBLY

Step 1:

For handles with rotation: Screw shaft with connection nut to handle's rotating wheel while holding rotating wheel in a fixed position. For handles without rotation: Screw shaft with connection nut to handle while holding handle in a fixed position.





Step 2:

For handles with ratchets: Push gold-plated lever down to detach disassemble ratchet connection. For handles without ratchets: Go directly to Step 3.




Step 3:

Open the handle per Disassembly: Step 2. Handle must remain in the fully open position to properly seat the insert into the handle assembly. Slide insert into shaft with jaws in the closed position. Align the raised pins below the insert tip with the slots on the shaft. Push the raised pins into the slots. When an audible pop is heard, observe the push button on the handle is visible.




Step 4:

Close handle against resistance until insert engages into shaft and an audible **snap** is heard. Insert will be seated flush to the shaft and the handle is easily manipulated.



Fundus Grasping Forceps - Single Action


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-508	645-508	-	-	-
425mm	644-509	645-509	-	-	-

Fundus Grasping Forceps - Double Action


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-510	645-510	-	-	-
425mm	644-511	645-511	-	-	-

Retraction Grasping Forceps - Single Action


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-514	645-514	-	-	-
425mm	644-515	645-515	-	-	-

Retraction Grasping Forceps - Double Action


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-516	645-516	-	-	-
425mm	644-517	645-517	-	-	-

Allis Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-524	645-524	646-524	-	-
425mm	644-525	645-525	646-525	-	-

Atrau-Allis Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-518	645-518	646-518	-	-
425mm	644-519	645-519	646-519	-	-

Atrau-Allis Long Jaw Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-520	645-520	-	-	-
425mm	644-521	645-521	-	-	-

Retrograde Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-522	645-522	646-522	-	-
425mm	644-523	645-523	646-523	-	-

DJB Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-526	645-526	-	-	-
425mm	644-527	645-527	-	-	-

DeBakey Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-528	645-528	646-528	-	-
425mm	644-529	645-529	646-529	-	-

DeBakey Forceps, Curved


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-530	645-530	-	-	-
425mm	644-531	645-531	-	-	-

Babcock Forceps - 2mm


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-540	645-540	646-540	-	-
425mm	644-541	645-541	646-541	-	-

Babcock Forceps - 3mm


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-536	645-536	-	-	-
425mm	644-537	645-537	-	-	-

Babcock Forceps, Round Jaw - 5mm


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-538	645-538	646-538	-	-
425mm	644-539	645-539	646-539	-	-

Paddle Babcock Forceps - 8mm


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-532	645-532	-	-	-
425mm	644-533	645-533	-	-	-

Paddle Babcock Forceps - 5mm


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-534	645-534	-	-	-
425mm	644-535	645-535	-	-	-

Claw Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-557	-	646-557	-	-
425mm	644-558	-	646-558	-	-

Claw Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-559	-	646-559	-	-
425mm	644-560	-	646-560	-	-

Tenaculum Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-569		646-569	-	-
425mm	644-570		646-570	-	-

Tenaculum Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-571	-	646-571	-	-
425mm	644-572	-	646-572	-	-

Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-500	645-500	646-500	647-500	648-500
425mm	644-501	645-501	646-501	647-501	648-501

TC Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-502	645-502	-	-	-
425mm	644-503	645-503	-	-	-

Strong TC Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-504	645-504	-	-	-
425mm	644-505	645-505	-	-	-

Grabber Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-506	645-506	-	-	-
425mm	644-507	645-507	-	-	-

Foerster Sponge Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-565	-	646-565	-	-
425mm	644-566	-	646-566	-	-

Strong Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-567	-	-	647-567	648-567
425mm	644-568	-	-	647-568	648-568

Duckbill Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-561	-	-	647-561	-
425mm	644-562	-	-	647-562	-

Russian Grasping Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-563	-	646-563	-	-
425mm	644-564	-	646-564	-	-

Dual-Cup Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-512	645-512	-	-	-
425mm	644-513	645-513	-	-	-

Swanstrom Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-542	645-542	-	-	-
425mm	644-543	645-543	-	-	-

Strong Hold Babcock Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-544	645-544	-	-	-
425mm	644-545	645-545	-	-	-

Duval Grasping Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-546	645-546	646-546	-	-
425mm	644-547	645-547	646-547	-	-

Anvil Grasping Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-548	645-548	-	-	-
425mm	644-549	645-549	-	-	-

Hunter Bowel Grasper


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-550	645-550	-	647-550	-
425mm	644-551	645-551	-	647-551	-

Grasping and Dissecting Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-552	-	-	-	648-552
425mm	644-553	-	-	-	648-553

Grasping and Extracting Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-554	-	646-554	-	-
425mm	644-555	-	646-555	-	-

Tapered Maryland Dissectors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-610	647-610	-
425mm	644-611	647-611	-

Maryland Dissectors - 30mm


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-612	647-612	-
425mm	644-613	647-613	-

Maryland Dissectors - 16.5mm


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-614	647-614	648-614
425mm	644-615	647-615	648-615

Maryland Dissectors - 16.5mm


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-616	-	648-616
425mm	644-617	-	648-617

Mixer Spreader-Dissectors


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-620	647-620	-
425mm	644-621	647-621	-

Mixer Spreader-Dissectors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-622	647-622	-
425mm	644-623	647-623	-

Mixer Dissecting Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-618	-	648-618
425mm	644-619	-	648-619

Right Angle Dissectors


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-624	647-624	-
425mm	644-625	647-625	-

Micro Dissecting Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-600	647-600	648-600
425mm	644-601	647-601	648-601

Needle Nose Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-602	647-602	648-602
425mm	644-603	647-603	648-603

Tapered Dissecting Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-604	647-604	648-604
425mm	644-605	647-605	648-605

Micro Grasping Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-606	-	648-606
425mm	644-607	-	648-607

Bullet-Nose Dissecting Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-608	647-608	-
425mm	644-609	647-609	-

Biopsy Punches


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-802	-	-	647-802	648-802
425mm	644-803	-	-	647-803	648-803

Biopsy Punches, 1 Tooth


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-806	-	-	-	648-806
425mm	644-807	-	-	-	648-807

Biopsy Forceps, 2 Teeth


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-800	-	646-800	647-800	-
425mm	644-801	-	646-801	647-801	-

Biopsy Forceps, Oval Fenestrated


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/Rotation	W/O Ratchet W/O Rotation
350mm	644-808	-	-	647-808	648-808
425mm	644-809	-	-	647-809	648-809

Hook Biopsy Punches


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-804	-	-	-	648-804
425mm	644-805	-	-	-	648-805

Micro Biopsy Forceps


Diameter 5mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-810	-	646-810	-	-
425mm	644-811	-	646-811	-	-

Spoon Forceps


Diameter 10mm



Shaft Length	Insert	Full Instrument			
		W/ Ratchet W/ Rotation	W/ Ratchet W/O Rotation	W/O Ratchet W/ Rotation	W/O Ratchet W/O Rotation
350mm	644-850	-	646-850	-	-
425mm	644-851	-	646-851	-	-

Supercut Scissors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-700	647-700	648-700
425mm	644-701	647-701	648-701

Supercut Endo Scissors, Straight


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-702	647-702	-
425mm	644-703	647-703	-

Supercut Endo Scissors, Curved


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-704	647-704	648-704
425mm	644-705	647-705	648-705

Micro Scissors, Straight


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-709	-	648-709
425mm	644-710	-	648-710

Micro Scissors, Curved Left


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-711	-	648-711
425mm	644-712	647-712	-

Metzenbaum Scissors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-715	647-715	-
425mm	644-716	647-716	-

Baby Metzenbaum Scissors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-719	647-719	-
425mm	644-720	647-720	-

Hook Scissors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-706	647-706	648-706
425mm	644-707	647-707	648-707

Mini Scissors


Diameter 5mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-717	647-717	-
425mm	644-718	647-718	-

Mayo Operating Scissors


Diameter 10mm



Shaft Length	Insert	Full Instrument (no ratchet)	
		W/ Rotation	W/O Rotation
350mm	644-713	-	648-713
425mm	644-714	-	648-714



Jarrit Detach



Jarrit Detach

Jarrit Detach

Integra's two-part Jarit Detach® System easily separates into a handle/shaft component and a scissors insert. Using backup inserts, which represent only a fraction of the cost of a complete Jarit Detach instrument, substantial cost savings are achieved. Each reusable endoscopic scissors offers electrosurgical coagulation and ergonomic rotation — all in one high-quality, cost-effective instrument. For easy cleaning, a cleaning brush is included with each instrument.

Jarit Detach Scissors Insert


625-0903

Cleaning brushes, pkg/6

Jarit Detach Handle/Shaft


625-090H

Roto-Cam,
monopolar,
5mm, 307mm



Jarit Detach Supercut Scissors

Unique Integra Jarit Supercut blade has a finely-serrated lower blade that helps prevent tissue slippage; curved blades meet the contour of structures while the blunt-shaped distal tip allows tissue manipulation and cauterization with blades closed; double action.



625-202E


Roto-Cam,
monopolar,
5mm

625-202N

Insert,
5mm


625-2023

Insert,
5mm, pkg/3



Jarit Detach Hook Scissors

Classic hook blades designed to lift, isolate and transect tissue.



625-205E


Roto-Cam,
monopolar,
5mm

625-205N

Insert,
5mm


625-2053

Insert,
5mm, pkg/3



Jarit Detach Micro Scissors Curved

Delicate, sharp-sharp,
curved blades



625-210E


Roto-Cam,
monopolar,
5mm

625-210N

Insert,
5mm


625-2103

Insert,
5mm, pkg/3



Jarit Detach Micro Scissors Straight

Delicate, sharp-sharp,
straight blades



625-212E


Roto-Cam,
monopolar,
5mm

625-212N

Insert,
5mm


625-2123

Insert,
5mm, pkg/3




Jarrit Detach Metzenbaum Scissors

Long, curved, tapered blades with wide spread, double action.



625-230E Roto-Cam, monopolar, 5mm	625-230N Insert, 5mm	625-2303 Insert, 5mm, pkg/3
---	-----------------------------------	--




Jarrit Detach Baby Metzenbaum Scissors

Shorter, curved, very tapered blades reach more narrow and delicate dissection areas, double action.




625-235E Roto-Cam, monopolar, 5mm	625-235N Insert, 5mm	625-2353 Insert, 5mm, pkg/3
---	-----------------------------------	--




Jarrit Detach Mini Scissors

Curved, tapered, delicate blades, double action.




625-240E Roto-Cam, monopolar, 5mm	625-240N Insert, 5mm	625-2403 Insert, 5mm, pkg/3
---	-----------------------------------	--




Jarrit Detach Biopsy Forceps

Two sharp teeth, one on each jaw, for biopsy or secure hold of thick-walled tissue




615-160E 1 x 1 Teeth, Roto-Cam, monopolar, 5mm	615-160N Insert, 1 x 1 teeth, 5mm, pkg/1	615-1603 Insert, 1 x 1 teeth, 5mm, pkg/3
---	---	--




Jarrit Detach Biopsy Punch

Punch with serrated lower jaw.




615-164E Roto-Cam, monopolar, 5mm	615-164N Insert, 5mm, pkg/1	615-1643 Insert, 5mm, pkg/3
---	---	--



Jarrit Detach Biopsy Forceps

Oval-shaped fenestrated jaws.




615-174E Roto-Cam, monopolar, 5mm	615-174N Insert, 5mm, pkg/1	615-1743 Insert, 5mm, pkg/3
---	---	--

As ORs continue to look for ways to control spending, many have looked to reduce the dollars spent on disposable cannula systems. This demand is met by Integra Jarit's reusable auto-valve cannulas which also allow for the convenience of hands-free operation. The built-in fascial threads prevent the cannula from slipping, and the lightweight construction makes it easier to insert instruments.

Auto-Valve Cannulas and Trocars, Non-threaded

Hands-free flapper valve; color-coded seals with size and color reference stamped on cannula; arrow on cannula indicates location of valve opening; stopcock for insufflation tubing; pyramidal trocar.


600-423

10/11mm, 101mm

600-428


5mm, 104mm

600-444

12mm, 98mm

Auto-Valve Cannulas and Trocars, Threaded

Fascial threads minimize in/out movement of cannula; hands-free flapper valve; color-coded seals with size and color reference stamped on cannula; arrow on cannula indicates location of valve opening; stopcock for insufflation tubing; pyramidal trocar.


600-424


10/11mm, 106mm

600-426

5mm, 107.5mm

Trocars

Pyramidal trocar available separately for backup.


600-425

 Auto-Valve,
11mm

600-431

 Piston-Valve,
5mm

600-427

 Auto-Valve,
5mm

600-441


 Piston-Valve,
11mm

600-447

 Auto-Valve,
12mm

Auto-Valve Hasson Cannulas

For open technique of cannula placement; Auto-valve allows hands-free insertion and removal of instruments; with suture ties on cone, cannula can be repositioned without retying suture; blunt obturator.




600-452
10/11mm, 101mm

630-452
5mm, 105mm

S-Retractors (set of 2)


For retraction of tissue during open technique placement of Hasson cannula.



600-455
10mm

Thoracic Cannula, Threaded


Fascial threads minimize in/out movement of cannula; blunt obturator.



600-448
12mm

Reducer Sleeves

Reduces piston valve cannulas and facilitates extraction of tissue.




600-460
10/11mm to 5mm

600-468
12mm to 5mm

600-465
5mm to 3mm

600-469
12mm to 10mm



600-461
Reducer "Flip Top"
for auto-valve
12.5/11/10mm to 5.5mm



Accessories

Accessories

Accessories

Laparoscopic Case**600-618**

Length 57.9 cm

22-3/4" x 11" x 8"


(579mm x 282mm x 203mm)

Laparoscopic Bariatric Case**600-620**

Length 68.1 cm

26-3/4" x 11" x 8"

(681mm x 282mm x 203mm)

**Instrument Holder**

The articulated Integra Jarit Endoscope & Instrument Holder System is a robust and simple system that will provide added stability and precision to the OR setting. The Integra Jarit Endoscopic & Instrument Holder System can be used with instruments with 5mm or 10mm diameter.

620-800 Endoscope & Instrument Holder

Complete set includes the following components:

- Articulated arm with Endoscope & Instrument Clamp 5/10mm (one piece), 85cm total length
- Rail Clamping Base, double clamp (**620-804**)

620-805 Bariatric Endoscope & Instrument Holder

Complete set includes the following components:


- Articulated arm with Endoscope & Instrument Clamp 5/10mm (one piece), 100cm total length
- Rail Clamping Base, double clamp (**620-804**)

**620-800**

Precision engineering with a rod lens optical system allows these high-quality scopes to deliver exceptional brilliance and clarity. Each scope transmits crisp, clear images with more natural color and enhanced contrast.

Autoclavable Laparoscopes


Gravity and prevacuum steam autoclavable; panoramic view with high-quality optics.



600-700
10mm, 0°
323mm

600-750
5mm, 0°
306mm

630-700
2.7mm, 0°
175mm




600-730
10mm, 30°
325mm

600-780
5mm, 30°
306mm


630-730
2.7mm, 30°
175mm

Operating Laparoscopes

Offset laparoscope with instrument channel which accommodates 5mm, 45cm instruments or laser; high-quality optics with panoramic view.




615-700
Autoclavable,
10mm, 0°



615-700

Light Post Adapters

- | | |
|----------------|---------------------------------|
| 600-292 | Laparoscope Light Post, Storz |
| 600-293 | Laparoscope Light Post, Wolf |
| 600-294 | Laparoscope Light Post, Olympus |





3.5mm

3.5mm

3.5mm

Roto-Lok, ratchet, insulated shaft




Grasping Forceps

630-100

3.5mm, 240mm

632-100

3.5mm, 320mm




Grabber Forceps

630-104

3.5mm, 240mm

632-104

3.5mm, 320mm




Dual-cup Grasping Forceps

630-108

3.5mm, 240mm

632-108

3.5mm, 320mm




Atrau-Allis Grasping Forceps

630-113

3.5mm, 240mm

632-113

3.5mm, 320mm



Babcock Forceps


For atraumatic, yet secure grasp of delicate tissue.

630-147

3.5mm, 240mm

632-147

3.5mm, 320mm



Hunter Bowel Grasper


Fenestrated, 14mm jaws for atraumatic grasping.

630-164

Non-ratcheted, 3.5mm, 240mm

632-164

Non-ratcheted, 3.5mm, 320mm



Fundus Grasping Forceps


Beveled-edge jaws.

630-107

3.5mm, 240mm

632-107

3.5mm, 320mm




Fundus Grasping Forceps

630-107IL

Ratchet, 3.5mm, 240mm

632-107IL

Ratchet, 3.5mm, 320mm



Hunter Bowel Grasper

630-164IL

Non-ratcheted, 3.5mm, 240mm

632-164IL

Non-ratcheted, 3.5mm, 320mm



In-Line,
insulated shaft


Roto-Cam, monopolar



Grasping Forceps

630-105
Monopolar,
3.5mm, 240mm


632-105
Monopolar,
3.5mm, 320mm



Micro Dissecting Forceps

630-110
Monopolar,
3.5mm, 240mm


632-110
Monopolar,
3.5mm, 320mm



Strong Forceps


630-111
Monopolar,
3.5mm, 240mm

632-111
Monopolar,
3.5mm, 320mm



Pyloromyotomy Dissector


630-116
Monopolar,
3.5mm, 140mm



Maryland Dissectors

630-126
Monopolar,
3.5mm, 240mm


632-126
Monopolar,
3.5mm, 320mm



Mixer Spreader-Dissectors

630-129RH
Monopolar, ring-
handle, 3.5mm, 240mm


632-129RH
Monopolar, ring-
handle, 3.5mm, 320mm



Endo Supercut Scissors, Curved

630-202
Monopolar,
3.5mm, 240mm


632-202
Monopolar,
3.5mm, 320mm



Hook Endo Scissors


630-205
Monopolar,
3.0mm, 240mm

632-205
Monopolar,
3.0mm, 320mm




Biopsy Forceps

630-178
Monopolar,
3.5mm, 240mm




Mixer Spreader-Dissectors




630-129
Roto-Cam,
monopolar,
3.5mm, 240mm

632-129
Roto-Cam,
monopolar,
3.5mm, 320mm



Appel Needle Holders

For use with 5-0 suture and smaller.




630-250


Left Curve,
3.5mm, 240mm

632-250

Left Curve,
3.5mm, 320mm



In-line Carb-Bite Needle Holder




630-251

Left Curve,
3.5mm, 240mm

632-251

Left Curve,
3.5mm, 330mm



Romeo Self-Righting Needle Holder




630-266


3.5mm, 241mm

632-266

3.5mm, 325mm



Lyons Needle Forceps




630-268


2.7mm, 246mm

632-268

2.7mm, 326mm



Juliet Receiving Forceps




630-269

2.7mm, 245mm

632-269

2.7mm, 320mm




Appel Knot Pusher



630-287

2.7mm, 340mm overall length




L-shaped Hook Electrodes

Partially insulated,
right angle tip.



630-318
3.5mm, 243mm

632-318
3.5mm, 330mm




Suction/Irrigation Cannula



630-340
3.5mm, 240mm

Auto-Valve Cannula and Trocar

Auto-valve cannula and pyramidal trocar 3.5mm with
silicone valve and luer-lock for insufflation tubing



630-427
Auto-valve cannula and pyramidal trocar
3.5mm




630-429
Pyramidal trocar for
auto-valve 3.5mm



630-4300
Blunt trocar for
auto-valve 3.5mm



630-4512
Nose Cone with wings,
for auto-valve cannula
3.5mm




630-4350
Silicone valve for auto-
valve cannula 3.5mm,
drop-shaped (package of 5)
for single use only!

Reducer Sleeve


Reduces piston valve cannulas and facilitates
extraction of tissue.

600-465
5mm to 3mm




Autoclavable Laparoscope


Panoramic view with high-quality optics.



630-700
2.7mm, 0°, 175mm



630-730
2.7mm, 30°, 175mm



Catalog #	Description
600-300/305/310/315	Suction Coagulators
600-TVA	T-V Replacement Parts Kit
600-3003	T-V Guide Pin, Pkg/6 (Tappet/Cap/Valve)
600-3004	T-V Spring, Pkg/6
600-3005	Trumpet Valve Case
600-3006	T-V Bottom Cap, Pkg/6
600-3007	T-V Set Screw, Pkg/12
600-320/322	Suction/Irrigation Cannulas
600-3010	Rubber Cap W/O Hole, Pkg/12 (Suction Cannulas)
600-3021	Trumpet Valve (320/322), Pkg/6
600-3022	Trumpet Valve Top Cap (320/322), Pkg/6
600-3023	Trumpet Valve Pin Cap (320/322), Pkg/6
600-3027	Trumpet Valve Set Screw (320/322), Pkg/6
600-3029	Metal Base Cap (320/322), Pkg/12
600-CTVA	Replacement Parts Kit (Includes 600-3024, 3033 And 3036)
600-3024	T-V Spring (320/322), Pkg/6
600-3033	Ctv Assembly (320/322), Pkg/3
600-3036	Ctv Bottom Cap (320/322), Pkg/6
X600-320LASER	Laser Channel Only For 600-320
X600-322LASER	Laser Channel Only For 600-322
600-330	Suction/Irrigation Set, 10mm, 33cm **Discontinued
630-330	Suction/Irrigation Cannula, 3.5mm **Discontinued
600-3301	Irrigation Valve (330)
600-3302	Suction Valve (330)
600-3303	Valve Top Cap (330), Pkg/6
600-3304	Valve Spring (330/340), Pkg/6
600-3305	Valve Bottom Cap (330), Pkg/6
600-3306	Replacement Cannula 5mm (600-330/600-340)
600-3307	Replacement Cannula 10mm (600-330/600-340)
630-3308	Replacement Cannula 3.5mm (630-330)
600-330HOUS	Main Housing For (330)
600-3309	O-Ring (330/340), Pkg/6
600-340	Suction/Irrigation Set, 10mm
630-340	Suction/Irrigation Cannula, 3.5mm
600-3401	Irrigation Valve (340) Blue
600-3402	Suction Valve (340) Green
600-3304	Valve Spring (330/340), Pkg/6
600-3405	Valve Bottom Cap (340), Pkg/6
600-3306	Replacement Cannula 5mm 600-330/600-340
600-3307	Replacement Cannula 10mm 600-330/600-340
600-340HOUS	Main Housing For (340)
600-3309	O-Ring (330/340), Pkg/6
630-3308	Replacement Cannula 3.5mm 240mm
632-3308	Replacement Cannula 3.5mm 340mm
600-423	Auto-Valve Cannula and Trocar, Non-threaded, 10/11mm
600-424	Auto-Valve Cannula and Trocar, Threaded, 10/11mm
600-4240	Inner Seal (Yellow) (423/424/444/452), Pkg/6
600-4241	Knurled Ring (423/424/444/452), Pkg/6
600-4242	O-Ring (423/424/444/452), Pkg/6
600-4420	Outer Seal (Blue) (423/424/440/450/452/471), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6

Catalog #	Description
600-426	Auto-Valve Cannula and Trocar, Non-threaded, 5mm
600-428	Auto-Valve Cannula and Trocar, Threaded, 5mm
600-4260	Inner Seal (Tan) (426/428), Pkg/6
600-4261	Knurled Ring (426/428), Pkg/6
600-4262	O-Ring (426/428), Pkg/6
600-4320	Outer Seal (Red) (426/428/430/462/470), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-444	Auto-Valve Cannulas and Trocars, Non-threaded, 12mm
600-4240	Inner Seal (Yellow) (423/424/444/452), Pkg/6
600-4241	Knurled Ring (423/424/444/452), Pkg/6
600-4242	O-Ring (423/424/444/452), Pkg/6
600-4470	Outer Seal (Green) (444/445), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-430	Cannula & Trocar Pyramid Tip 5mm **Discontinued
600-3030	Valve Spring (430), Pkg/6
600-3031	K-Cap W/Hole (430), Pkg/6
600-3032	K-Cap W/O Hole (430), Pkg/6
600-4320	Outer Seal (Red) (426/428/430/462/470), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-440	Cannula & Trocar Pyramid Tip 10/11mm **Discontinued
600-3040	Valve Spring (440/450), Pkg/6
600-3041	K-Cap W/Hole (440/450), Pkg/6
600-3042	K-Cap W/O Hole (440/450), Pkg/6
600-4420	Outer Seal (Blue) (423/424/440/450/452/471), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-445	Cannula & Trocar Pyramid Tip 12mm **Discontinued
600-3050	Valve Spring (445), Pkg/6
600-3051	K-Cap W/Hole (445), Pkg/6
600-3052	K-Cap W/O Hole (445), Pkg/6
600-4470	Outer Seal (Green) (444/445), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-448	Thoracic Cannula Threaded 12mm
600-4481	Stabilizing Cylinder (448), Pkg/2
600-450	Hasson Cannula 10/11mm, Blunt **Discontinued
600-3040	Valve Spring (440/450), Pkg/6
600-3041	K-Cap W/Hole (440/450), Pkg/6
600-3042	K-Cap W/O Hole (440/450), Pkg/6
600-4420	Outer Seal (Blue) (423/424/440/450/452/471), Pkg/6
600-4501	Blunt Trocar (450/452)
600-4502	Nose Cone (450/452)
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-4504	Suture Attachment (450/452)
600-4520	Cone Seal (450/452), Pkg/6

Catalog #	Description
600-452	Auto-Valve Hasson Cannula, 10/11mm
600-4240	Inner Seal (Yellow) (423/424/444/452), Pkg/6
600-4241	Knurled Ring (423/424/444/452), Pkg/6
600-4242	O-Ring (423/424/444/452), Pkg/6
600-4420	Outer Seal (Blue) (423/424/440/450/452/471), Pkg/6
600-4501	Blunt Trocar (450/452)
600-4502	Nose Cone (450/452)
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
600-4520	Cone Seal (450/452), Pkg/6
630-452	Auto-Valve Hasson Cannula, 5mm
600-4260	Replace Inner Seal(426/428), Pkg/6
600-4261	Knurled Ring (426/428), Pkg/6
600-4262	O-Ring (426/428), Pkg/6
600-4320	Outer Seal (Red) (426/428/430/462/470), Pkg/6
600-4503	Stopcock Cap (All Trocars/Cannulas), Pkg/6
630-4501	Blunt Trocar (452)
630-4502	Nose Cone (452)
630-4520	Cone Seal (452), Pkg/6
630-427	Auto-Valve Cannula and Trocar, 3.5mm
630-4300	Blunt Trocar
630-4320	Outer Seal (Red), Pkg/6
630-4350	Silicone Valve For Auto-Valve, Pkg/5
600-240	Endoloop Applicator
600-460/465/468/469	Reducer Sleeves
600-4620	Replacement Seal (460), Pkg/6
600-4670	Replacement Seal Diam. 2mm (240/465), Pkg/6
600-4680	Replacement Seal (468), Pkg/6
600-4690	Replacement Seal (469), Pkg/6
600-461	Reducer Flip-Top
600-4320	Outer Seal (Red) (426/428/430/462/470), Pkg/6

Catalog #	Description
600-550/560 600-SPS-550	P.E.E.R. Retractors Spring And Screw (550/560)
600-580/582/584/585 X600-580CPCAP	Articulators Plastic Cap For Cleaning Port
615-330/334/338 615-331/335/339 615-334HC 615-EC-330	Bipolar Forceps 32cm Bipolar Forceps 45cm Plastic Brushing For Bipolar Connector (White), Pkg/3 End Caps For Bipolars, Pkg/3
600-410 600-3011	Cholangiogram Guide Rubber Cap W/Hole, Pkg/12 (410)
600-700/730/750/780 615-700 630-700/730 615-3011	5/10mm Laparoscopes 10mm Offset Laparoscope 2.7mm Laparoscopes Rubber Cap For Light Cable Connectors, Pkg/6
600-130/135/137/138/273/275 615-160/162/178 600-SCR-138 600-3017 X600-135SPRING	Spring-Loaded Instruments Screw Only For Spring-Loaded Handles Spring And Screw For Spring-Loaded Handles, Pkg/6 Spring Only For Spring-Loaded Handles
620-800/805 620-804 620-EC-800 620-WHS-800	Endoscope and Instrument Holders Rails Clamping Base, Double Clamp End Screw (800/805) Wing Handle Screw (800/805)



Index

Index

Index

3-Piece Take-Apart Quick Reference

Handles



Passive Ratchet Mechanism

Precise teeth easily engage and disengage for surgical flexibility.

Monopolar Connection

All handles have monopolar connection for eletrosurgical flexibility.



642-200
w/ Rotation, w/ Ratchet



642-225
w/o Rotation, w/ Ratchet



642-250
w/ Rotation, w/o Ratchet



642-275
w/o Rotation, w/o Ratchet

Shafts



643-350
5mm, 350mm



643-425
5mm, 425mm












643-650
10mm, 350mm























643-725
10mm, 425mm





Cleaning Brush

643-100
For all shafts

Tip Style		Dia	One-Piece		3-Piece Take-Apart		
			SKU's	Page	Inserts Only	Assembled Instrument	
Graspers					SKU's	SKU's	Page
	Allis Grasping Forceps	5mm	600-122 625-122	7	644-524 644-525	645-524, 645-525 646-524, 646-525	37
	Anvil Grasping Forceps	10mm	625-154	12	644-548 644-549	645-548, 645-549 646-518, 646-519	44
	Atrau-Allis Grasping Forceps	5mm	600-113, 600-114 625-113, 625-113LL	7	644-518 644-519	645-518, 645-519 646-518, 646-519	37
	Atrau-Allis Long Jaw Grasping Forceps	5mm	625-114 625-114LL	7	644-520 644-521	645-520 645-521	37
	Babcock Forceps, 5mm	5mm	600-147, 610-147 625-147, 625-147LL	11	644-540 644-541	645-540, 645-541 646-540, 646-541	38
	Babcock Forceps, 10mm	10mm	625-144	10	644-536 644-537	645-536 645-537	39
	Babcock Forceps, Round Jaw	10mm	600-145 625-145	11	644-538 644-539	645-538, 645-539 646-538, 646-539	39
	Claw Forceps	5mm	600-137	9	644-559 644-560	646-559 646-560	40
	Claw Forceps	10mm	600-135	9	644-557 644-558	646-557 646-558	40
	DeBakey Forceps	5mm	600-140, 610-140 610-140L, 625-140LL	9	644-528 644-529	645-528, 645-529 646-528, 646-529	38
	DeBakey Forceps, Curved	10mm	625-141	9	644-530 644-531	645-530 645-531	38
	DJB Grasping Forceps	5mm	625-139 625-139LL	9	644-526 644-527	645-526 645-527	38
	Dual-Cup Grasping Forceps	5mm	600-108IL, 610-108 625-108, 625-108LL	6	644-512 644-513	645-512 645-513	43
	Duckbill Grasping Forceps	5mm	625-118	12	644-561 644-562	647-561 647-562	42
	Duval Grasping Forceps	10mm	600-150 625-150	12	644-546 644-547	645-546, 645-547 646-546, 646-547	43
	Foerster Sponge Forceps	10mm	600-158	13	644-565 644-566	646-565 646-566	42
	Fundus Grasping Forceps - Double Action	5mm	600-107IL, 610-107DA 625-107D5, 625-107DA 625-107RD	6	644-510 644-511	645-510 645-511	36
	Fundus Grasping Forceps - Single Action	5mm	610-107, 625-107 625-107LL	6	644-508 644-509	645-508 645-509	36
	Grabber Forceps	5mm	600-104IL, 610-104 625-104, 625-104LL 625-106LL	5	644-506 644-507	645-506 645-507	41

Tip Style		Dia	One-Piece		3-Piece Take-Apart		
			SKU's	Page	Inserts Only	Assembled Instrument	
Graspers					SKU's	SKU's	Page
	Grasping and Dissecting Forceps	5mm	600-120	8	644-552 644-553	648-552 648-553	44
	Grasping and Extracting Forceps	10mm	600-130	8	644-554 644-555	646-554 646-555	44
	Grasping Forceps	5mm	600-100, 600-105 600-106, 605-100 610-100, 625-100 625-100LL, 625-105	4	644-500 644-501	645-500, 645-501 646-500, 646-501 647-500, 647-501 648-500, 648-501	41
	Hunter Bowel Grasper	5mm	600-164IL, 610-163 610-164, 610-165 625-163, 625-164 625-164LL, 625-165LL	13	644-550 644-551	645-550, 645-551 647-550, 647-551	44
	Paddle Babcock Forceps	5mm	600-143IL, 625-143 625-143LL, 625-146LL	10	644-534 644-535	645-534 645-535	39
	Paddle Babcock Forceps	10mm	625-142	10	644-532 644-533	645-532 645-533	39
	Retraction Grasping Forceps - Double Action	5mm	625-109DA, 625-109IL 652-109RD	7	644-516 644-517	645-516 645-517	36
	Retraction Grasping Forceps - Single Action	5mm	625-109 625-109LL	7	644-514 644-515	645-514 645-515	36
	Retrograde Grasping Forceps	5mm	600-116, 600-118 600-119	8	644-522 644-523	645-522, 645-523 646-522, 646-523	37
	Russian Grasping Forceps	10mm	600-154	12	644-563 644-564	646-563 646-564	42
	Strong Forceps	5mm	615-109, 615-110 615-111	14	644-567 644-568	647-567, 647-568 648-567, 648-568	42
	Strong Hold Babcock Forceps	5mm	610-149	10	644-544 644-545	645-544 645-545	43
	Strong TC Grasping Forceps	5mm	625-101, 625-103 625-103LL	5	644-504 644-505	645-504 645-505	41
	Swanstrom Grasping Forceps	5mm	600-148IL, 610-148 625-148, 625-148LL	11	644-542 644-543	645-542 645-543	43
	TC Grasping Forceps	5mm	625-102 625-102LL	4	644-502 644-503	645-502 645-503	41
	Tenaculum Forceps	5mm	605-152	14	644-569 644-570	646-569 646-570	40
	Tenaculum Forceps	10mm	615-150 615-152	14	644-571 644-572	646-571 646-572	40

Tip Style		Dia	One-Piece		3-Piece Take-Apart		
			SKU's	Page	Inserts Only	Assembled Instrument	Page
Dissectors					SKU's	SKU's	
	Bullet-Nose Dissecting Forceps	5mm	625-117	17	644-608 644-609	647-608 647-609	47
	Maryland Dissector - Long Jaw	10mm	625-125	17	644-612 644-613	647-612 647-613	45
	Maryland Dissector - Standard	5mm	600-126, 625-126 625-127	17	644-614 644-615	647-614, 647-615 648-614, 648-615	45
	Maryland Dissector - Standard	10mm	600-125	17	644-616 644-617	648-616 648-617	45
	Micro Dissecting Forceps	5mm	600-110 625-110	16	644-600 644-601	647-600, 647-601 648-600, 648-601	47
	Micro Grasping Forceps	5mm	600-117	16	644-606 644-607	648-606 648-607	47
	Mixer Dissecting Forceps	10mm	600-127	18	644-618 644-619	648-618 648-619	46
	Mixer Spreader-Dissector	5mm	625-129, 625-129RH 625-130RH	18	644-622 644-623	647-622 647-623	46
	Mixer Spreader-Dissector	10mm	600-128, 625-128 625-128RH	18	644-620 644-621	647-620 647-621	46
	Needle Nose Forceps	5mm	600-112 625-112	16	644-602 644-603	647-602, 647-603 648-602, 648-603	47
	Right Angle Dissector	10mm	625-132	18	644-624 644-625	647-624 647-625	46
	Tapered Dissecting Forceps	5mm	600-115 625-115	16	644-604 644-605	647-604, 647-605 648-604, 648-605	47
	Tapered Maryland Dissector	5mm	625-124	17	644-610 644-611	647-610 647-611	45
Biopsy Forceps							
	Biopsy Forceps, 2 Teeth	5mm	615-160 615-162	19	644-800 644-801	646-800 646-801	48
	Biopsy Forceps, Oval Fenestrated	5mm	615-174 615-176	19	644-808 644-809	647-808, 647-809 648-808, 648-809	48
	Biopsy Punch	5mm	615-164 615-166	19	644-802 644-803	647-802, 647-803 648-802, 648-803	48
	Biopsy Punch, 1 Tooth	5mm	615-170 615-172	19	644-806 644-807	648-806 648-807	48
	Hook Biopsy Punch	5mm	615-168 615-169	19	644-804 644-805	648-804 648-805	49
	Micro Biopsy Forceps	5mm	615-178	19	644-810 644-811	646-810 646-811	49
	Spoon Forceps	10mm	600-138	30	644-850 644-851	646-850 646-851	49

Tip Style		Dia	One-Piece		3-Piece Take-Apart		
			SKU's	Page	Inserts Only	Assembled Instrument	Page
Scissors					SKU's	SKU's	
	Baby Metzenbaum	5mm	-	-	644-719 644-720	647-719 647-720	51
	Hook Scissors	5mm	600-205 600-206 600-207 625-205	20	644-706 644-707	647-706 647-707 648-706 648-707	51
	Mayo Operating Scissors	10mm	600-220	21	644-713 644-714	648-713 648-714	51
	Metzenbaum Scissors	5mm	625-230	21	644-715 644-716	647-715 647-716	51
	Micro Scissors, Curved Left	5mm	600-212	21	644-711 644-712	648-711 648-712	50
	Micro Scissors, Straight	5mm	600-210	21	644-709 644-710	648-709 648-710	50
	Mini Scissors	5mm	625-240	21	644-717 644-718	647-717 647-718	51
	Supercut Endo Scissors, Curved	5mm	600-202 600-203	20	644-704 644-705	647-704 647-705 648-704 648-705	50
	Supercut Endo Scissors, Straight	5mm	625-201	20	644-702 644-703	647-702 647-703	50
	Supercut Scissors	5mm	600-200 600-201	20	644-700 644-701	647-700 647-701 648-700 648-701	50

600-100	4, 69	600-3003	63	600-4470	64	615-1603	54	625-127	17, 70
600-104IL	5, 68	600-3004	63	600-448	65, 64	615-160E	54	625-128	18, 70
600-105	4, 69	600-3005	63	600-4481	64	615-160N	54	625-128RH	18, 70
600-106	4, 69	600-3006	63	600-450	64	615-162	19, 66, 70	625-129	18, 70
600-107IL	6, 68	600-3007	63	600-4501	64, 65	615-164	19, 70	625-129RH	18, 70
600-108IL	6, 68	600-3010	63	600-4502	64, 65	615-1643	54	625-130RH	18, 70
600-109IL	7	600-3011	66	600-4503	64, 65	615-164E	54	625-132	18, 70
600-110	16, 70	600-3017	66	600-4504	64	615-164N	54	625-139	9, 68
600-112	16, 70	600-3021	63	600-452	56, 65	615-166	19, 70	625-139LL	9, 68
600-113	7, 68	600-3022	63	600-4520	64, 65	615-168	19, 70	625-140LL	9, 68
600-114	7, 68	600-3023	63	600-455	56	615-169	19, 70	625-141	9, 68
600-115	16, 70	600-3024	63	600-460	56, 65	615-170	19, 70	625-142	10, 69
600-116	8, 69	600-3027	63	600-461	56, 65	615-172	19, 70	625-143	10, 69
600-117	16, 70	600-3029	63	600-4620	65	615-174	19, 70	625-143LL	10, 69
600-118	8, 69	600-3030	64	600-465	56, 62, 65	615-1743	54	625-144	10, 68
600-119	8, 69	600-3031	64	600-4670	65	615-174E	54	625-145	11, 68
600-120	8, 69	600-3032	64	600-468	56, 65	615-174N	54	625-146LL	10, 69
600-122	7, 68	600-3033	65	600-4680	65	615-176	19, 70	625-147	11, 68
600-123	8	600-3036	63	600-469	56, 65	615-178	19, 70	625-147LL	11, 68
600-125	17, 70	600-3040	64	600-4690	65	615-210	31	625-148	11, 69
600-126	17, 70	600-3041	64	600-475	31	615-212	31	625-148LL	11, 69
600-127	18, 70	600-3042	64	600-480	31	615-295	29	625-150	12, 68
600-128	18, 70	600-305	28, 63	600-481	31	615-296B	29	625-154	12, 68
600-130	8, 66, 69	600-3050	64	600-490	25	615-3011	66	625-163	13, 69
600-135	9, 68	600-3051	64	600-491	25	615-330	29, 66	625-164	13, 69
600-137	9, 68	600-3052	64	600-495	31	615-331	29, 66	625-164LL	13, 69
600-138	30, 70	600-310	28, 63	600-496	31	615-334	29, 66	625-165LL	13, 69
600-140	9, 68	600-315	28, 63	600-500	18	615-334HC	66	625-200	20
600-143IL	10, 69	600-316	28	600-520	14	615-335	29, 66	625-201	20, 71
600-145	11, 68	600-317	28	600-522	14	615-338	29, 66	625-202	20
600-147	11, 68	600-317L	28	600-540	27	615-339	29, 66	625-2023	53
600-148IL	11, 69	600-318	28	600-550	27, 66	615-700	55, 66	625-202E	53
600-150	12, 68	600-318L	28	600-560	27, 66	615-EC-330	66	625-202N	53
600-154	12, 69	600-319	28	600-580	26, 66	620-800	57, 66	625-203	20
600-158	13, 68	600-319L	28	600-582	26, 66	620-804	57, 66	625-205	20, 71
600-160	13	600-320	25, 63	600-584	26, 66	620-805	54, 66	625-2053	53
600-164IL	13, 69	600-322	25, 60	600-585	26, 66	620-EC-800	66	625-205E	53
600-200	20, 71	600-330	63	600-618	57	620-WHS-800	66	625-205N	53
600-201	20, 71	600-3301	63	600-620	57	625-0903	53	625-2103	53
600-202	20, 71	600-3302	63	600-700	58, 66	625-090H	53	625-210E	53
600-203	20, 71	600-3303	63	600-730	58, 66	625-100	4, 69	625-210N	53
600-205	20, 71	600-3304	63	600-750	58, 66	625-100LL	4, 69	625-2123	53
600-206	20, 71	600-3305	63	600-780	58, 66	625-101	5, 69	625-212E	53
600-207	20, 71	600-3306	63	600-CTVA	63	625-102	4, 69	625-212N	53
600-210	21, 71	600-3307	63	600-SCR-138	66	625-102LL	4, 69	625-230	21, 71
600-212	21, 71	600-3309	63	600-SPS-550	66	625-103	5, 69	625-2303	54
600-220	21, 71	600-330HOUS	63	600-TVA	63	625-103LL	5, 69	625-230E	54
600-240	24, 65	600-340	25, 63	605-100	4, 69	625-104	5, 68	625-230N	54
600-249	22	600-340HOUS	63	605-152	14, 69	625-104LL	5, 68	625-2353	54
600-250	22, 24	600-3401	63	610-100	4, 69	625-105	4, 69	625-235E	54
600-251	22	600-3402	63	610-104	5, 68	625-106LL	5, 68	625-235N	54
600-252	22	600-3405	63	610-107	6, 68	625-107	6, 68	625-240	21, 71
600-253	22	600-400	30	610-107DA	6, 68	625-107D5	6, 68	625-2403	54
600-254	22	600-410	30, 66	610-108	6, 68	625-107DA	6, 68	625-240E	54
600-255	22	600-423	55, 63	610-140	9, 68	625-107LL	6, 68	625-240N	54
600-260	22	600-424	55, 63	610-140L	9, 68	625-107RD	6, 68	630-100	59
600-261	22	600-4240	63-65	610-147	11, 68	625-108	6, 68	630-104	59
600-264	22	600-4241	63-65	610-148	11, 69	625-108LL	6, 68	630-105	60
600-266	23	600-4242	63-65	610-149	10, 69	625-109	7, 69	630-107	59
600-267	23	600-425	55	610-163	13, 69	625-109DA	7, 69	630-107IL	59
600-268	23, 24	600-426	55, 64	610-164	13, 69	625-109LL	7, 69	630-108	59
600-269	23	600-4260	64, 65	610-165	13, 69	625-109RD	7, 69	630-110	60
600-270	24	600-4261	64, 65	610-166	13	625-110	16, 70	630-111	60
600-271	23	600-4262	64, 65	615-109	14, 69	625-112	16, 70	630-113	59
600-273	23	600-427	55	615-110	14, 69	625-113	7, 68	630-116	60
600-275	23	600-428	55, 64	615-111	14, 69	625-113LL	7, 68	630-126	60
600-285	24	600-430	64	615-113	15	625-114	7, 68	630-129	60
600-287	24	600-431	55	615-114	15	625-114LL	7, 68	630-129RH	60
600-289	24	600-4320	64, 65	615-116	15	625-115	16, 70	630-147	59
600-290	29	600-440	64	615-118	15	625-117	17, 70	630-164	59
600-291	29	600-441	55	615-143	15	625-118	12, 68	630-164IL	59
600-292	58	600-4420	63-65	615-144	15	625-122	7, 68	630-178	60
600-293	58	600-444	55, 64	615-150	14, 69	625-124	17, 70	630-202	60
600-294	58	600-445	64	615-152	14, 69	625-125	17, 70	630-205	60
600-300	28, 63	600-447	55	615-160	19, 66, 70	625-126	17, 70	630-250	61

630-251.....	61	644-521.....	37, 68	644-624.....	46, 70	645-539.....	39, 68	647-624.....	46, 70
630-266.....	61	644-522.....	37, 69	644-625.....	46, 70	645-540.....	38, 68	647-625.....	46, 70
630-268.....	61	644-523.....	37, 69	644-700.....	50, 71	645-541.....	38, 68	647-700.....	50, 71
630-269.....	61	644-524.....	37, 68	644-701.....	50, 71	645-542.....	43, 69	647-701.....	50, 71
630-287.....	61	644-525.....	37, 68	644-702.....	50, 71	645-543.....	43, 69	647-702.....	50, 71
630-318.....	62	644-526.....	38, 68	644-703.....	50, 71	645-544.....	43, 69	647-703.....	50, 71
630-330.....	63	644-527.....	38, 68	644-704.....	50, 71	645-545.....	43, 69	647-704.....	50, 71
630-3308.....	63	644-528.....	38, 68	644-705.....	50, 71	645-546.....	43, 68	647-705.....	50, 71
630-340.....	62, 63	644-529.....	38, 68	644-706.....	51, 71	645-547.....	43, 68	647-706.....	51, 71
630-427.....	62, 65	644-530.....	38, 68	644-707.....	51, 71	645-548.....	44, 68	647-707.....	51, 71
630-429.....	62	644-531.....	38, 68	644-709.....	50, 71	645-549.....	44, 68	647-712.....	50
630-4300.....	62, 65	644-532.....	39, 69	644-710.....	50, 71	645-550.....	44, 69	647-715.....	51, 71
630-4320.....	65	644-533.....	39, 69	644-711.....	50, 71	645-551.....	44, 69	647-716.....	51, 71
630-4350.....	62, 65	644-534.....	39, 69	644-712.....	50, 71	646-500.....	41, 69	647-717.....	51, 71
630-4501.....	65	644-535.....	39, 69	644-713.....	51, 71	646-501.....	41, 69	647-718.....	51, 71
630-4502.....	65	644-536.....	39, 68	644-714.....	51, 71	646-518.....	37, 68	647-719.....	51, 71
630-4512.....	62	644-537.....	39, 68	644-715.....	51, 71	646-519.....	37, 68	647-720.....	51, 71
630-452.....	56, 65	644-538.....	39, 68	644-716.....	51, 71	646-522.....	37, 69	647-800.....	48
630-4520.....	65	644-539.....	39, 68	644-717.....	51, 71	646-523.....	37, 69	647-801.....	48
630-700.....	58, 62, 66	644-540.....	38, 68	644-718.....	51, 71	646-524.....	37, 68	647-802.....	48, 70
630-730.....	58, 62, 66	644-541.....	38, 68	644-719.....	51, 71	646-525.....	37, 68	647-803.....	48, 70
632-100.....	59	644-542.....	43, 69	644-720.....	51, 71	646-528.....	38, 68	647-808.....	48, 70
632-104.....	59	644-543.....	43, 69	644-800.....	48, 70	646-529.....	38, 68	647-809.....	48, 70
632-105.....	60	644-544.....	43, 69	644-801.....	48, 70	646-538.....	39, 68	648-500.....	41, 69
632-107.....	59	644-545.....	43, 69	644-802.....	48, 70	646-539.....	39, 68	648-501.....	41, 69
632-107IL.....	59	644-546.....	43, 68	644-803.....	48, 70	646-540.....	38, 68	648-552.....	44, 69
632-108.....	59	644-547.....	43, 68	644-804.....	49, 70	646-541.....	38, 68	648-553.....	44, 69
632-110.....	60	644-548.....	44, 68	644-805.....	49, 70	646-546.....	43, 68	648-567.....	42, 69
632-111.....	60	644-549.....	44, 68	644-806.....	48, 70	646-547.....	43, 68	648-568.....	42, 69
632-113.....	59	644-550.....	44, 69	644-807.....	48, 70	646-554.....	44, 69	648-600.....	47, 70
632-126.....	60	644-551.....	44, 69	644-808.....	48, 70	646-555.....	44, 69	648-601.....	47, 70
632-129.....	60	644-552.....	44, 69	644-809.....	48, 70	646-557.....	40, 68	648-602.....	47, 70
632-129RH.....	60	644-553.....	44, 69	644-810.....	49, 70	646-558.....	40, 68	648-603.....	47, 70
632-147.....	59	644-554.....	44, 69	644-811.....	49, 70	646-559.....	40, 68	648-604.....	47, 70
632-164.....	59	644-555.....	44, 69	644-850.....	49, 70	646-560.....	40, 68	648-605.....	47, 70
632-164IL.....	59	644-557.....	40, 68	644-851.....	49, 70	646-563.....	42, 69	648-606.....	47, 70
632-202.....	60	644-558.....	40, 68	645-500.....	41, 69	646-564.....	42, 69	648-607.....	47, 70
632-205.....	60	644-559.....	40, 68	645-501.....	41, 69	646-565.....	42, 68	648-614.....	45, 70
632-250.....	61	644-560.....	40, 68	645-502.....	41, 69	646-566.....	42, 68	648-615.....	45, 70
632-251.....	61	644-561.....	42, 68	645-503.....	41, 69	646-569.....	40, 69	648-616.....	45, 70
632-266.....	61	644-562.....	42, 68	645-504.....	41, 69	646-570.....	40, 69	648-617.....	45, 70
632-268.....	61	644-563.....	42, 69	645-505.....	41, 69	646-571.....	40, 69	648-618.....	46, 70
632-269.....	61	644-564.....	42, 69	645-506.....	41, 68	646-572.....	40, 69	648-619.....	46, 70
632-318.....	62	644-565.....	42, 68	645-507.....	41, 68	646-800.....	48, 70	648-700.....	50, 71
632-3308.....	63	644-566.....	42, 68	645-508.....	36, 68	646-801.....	48, 70	648-701.....	50, 71
642-200.....	33, 67	644-567.....	42, 69	645-509.....	36, 68	646-810.....	49, 70	648-704.....	50, 71
642-225.....	33, 67	644-568.....	42, 69	645-510.....	36, 68	646-811.....	49, 70	648-705.....	50, 71
642-250.....	33, 67	644-569.....	40, 69	645-511.....	36, 68	646-850.....	49, 70	648-706.....	51, 71
642-275.....	33, 67	644-570.....	40, 69	645-512.....	43, 68	646-851.....	49, 70	648-707.....	51, 71
643-100.....	33, 67	644-571.....	40, 69	645-513.....	43, 68	647-500.....	41, 69	648-709.....	50, 71
643-350.....	33, 67	644-572.....	40, 69	645-514.....	36, 69	647-501.....	41, 69	648-710.....	50, 71
643-425.....	33, 67	644-600.....	47, 70	645-515.....	36, 69	647-550.....	44, 69	648-711.....	50, 71
643-650.....	33, 67	644-601.....	47, 70	645-516.....	36, 69	647-551.....	44, 69	648-713.....	51, 71
643-725.....	33, 67	644-602.....	47, 70	645-517.....	36, 69	647-561.....	42, 68	648-714.....	51, 71
644-500.....	41, 69	644-603.....	47, 70	645-518.....	37, 68	647-562.....	42, 68	648-802.....	48, 70
644-501.....	41, 69	644-604.....	47, 70	645-519.....	37, 68	647-567.....	42, 69	648-803.....	48, 70
644-502.....	41, 69	644-605.....	47, 70	645-520.....	37, 68	647-568.....	42, 69	648-804.....	49, 70
644-503.....	41, 69	644-606.....	47, 70	645-521.....	37, 68	647-600.....	47, 70	648-805.....	49, 70
644-504.....	41, 69	644-607.....	47, 70	645-522.....	37, 69	647-601.....	47, 70	648-806.....	48, 70
644-505.....	41, 69	644-608.....	47, 70	645-523.....	37, 69	647-602.....	47, 70	648-807.....	48, 70
644-506.....	41, 69	644-609.....	47, 70	645-524.....	37, 68	647-603.....	47, 70	648-808.....	48, 70
644-507.....	41, 69	644-610.....	45, 70	645-525.....	37, 68	647-604.....	47, 70	648-809.....	48, 70
644-508.....	36, 68	644-611.....	45, 70	645-526.....	38, 68	647-605.....	47, 70	X600-135SPRING.....	66
644-509.....	36, 68	644-612.....	45, 70	645-527.....	38, 68	647-608.....	47, 70	X600-320LASER.....	63
644-510.....	36, 68	644-613.....	45, 70	645-528.....	38, 68	647-609.....	47, 70	X600-322LASER.....	63
644-511.....	36, 68	644-614.....	45, 70	645-529.....	38, 68	647-610.....	45, 70	X600-580CPCAP.....	66
644-512.....	43, 68	644-615.....	45, 70	645-530.....	38, 68	647-611.....	45, 70		
644-513.....	43, 68	644-616.....	45, 70	645-531.....	38, 68	647-612.....	45, 70		
644-514.....	36, 69	644-617.....	45, 70	645-532.....	39, 69	647-613.....	45, 70		
644-515.....	36, 69	644-618.....	46, 70	645-533.....	39, 69	647-614.....	45, 70		
644-516.....	36, 69	644-619.....	46, 70	645-534.....	39, 69	647-615.....	45, 70		
644-517.....	36, 69	644-620.....	46, 70	645-535.....	39, 69	647-620.....	46, 70		
644-518.....	37, 68	644-621.....	46, 70	645-536.....	39, 68	647-621.....	46, 70		
644-519.....	37, 68	644-622.....	46, 70	645-537.....	39, 68	647-622.....	46, 70		
644-520.....	37, 68	644-623.....	46, 70	645-538.....	39, 68	647-623.....	46, 70		

3.5mm Biopsy Forceps	60	Hasson Cannulas	56, 64-65
3.5mm Dissecting Forceps	60	Herron Articulator Bar Retractor	26
3.5mm Grasping Forceps	59-60	Hook Biopsy Punches	19, 49
3.5mm Laparoscopic Instruments.....	59-62	Hook Electrodes	28
3.5mm Ligation.....	61	Hook End Suction Coagulator	28
3.5mm Needle Holders.....	61	Hook Endo Scissors	60
3.5mm Receivers	61	Hook Scissors.....	20, 51, 53
3.5mm Scissors	60	Hook, Appel Dissecting	18
Allis Grasping Forceps.....	7, 37	Hunter Bowel Graspers.....	13, 44, 59
Anvil Grasping Forceps	12, 44	Hunter Clamp.....	13
Appel Dissecting Hook	18	In-line Carb-Bite Needle Holders.....	22, 61
Appel Knot Pushers.....	24, 61	Instrument Holder.....	57, 66
Appel Needle Holders.....	22, 61	Jarit Detach Baby Metzenbaum Scissors	54
Articulator Retractors	26, 68	Jarit Detach Biopsy Forceps.....	54
Articulators.....	26, 68	Jarit Detach Biopsy Punch	54
Atrau-Allis Grasping Forceps.....	7, 37, 59	Jarit Detach Hook Scissors.....	53
Atrau-Allis Long Jaw Grasping Forceps.....	7, 37, 58	Jarit Detach Metzenbaum Scissors.....	54
Autoclavable Laparoscopes	58, 62, 68	Jarit Detach Micro Scissors	53
Auto-Valve Cannula and Trocar	66, 68	Jarit Detach Mini Scissors	54
Auto-Valve Cannulas and Trocars.....	57, 67, 70	Jarit Detach Scissors.....	53-54
Auto-Valve Hasson Cannulas.....	56, 67	Jarit Detach Supercut Scissors.....	53
Babcock Forceps	10-11, 38-39, 59	Jarit Detach System	53-54
Babcock Forceps, Round Jaw	11, 39	Juliet Receiving Forceps.....	23, 61
Babcock Forceps, Strong Hold.....	10, 43	Knot Pushers	24, 61
Baby Metzenbaum Scissors	51, 54	Laparoscopes	58, 62
Ball Electrode	28	Laparoscopic Case	57
Bariatric Case	57	Large Bipolar Forceps.....	29
Bariatric Endoscope & Instrument Holder	57	Ligation	22-24
Biopsy Forceps	19, 48-49, 55, 60	Light Post Adapters	58
Biopsy Instruments	19, 48-49, 55, 60	Lobe Spleen Extractor.....	14
Biopsy Punch.....	55	L-shaped Hook Electrodes	28, 62
Biopsy Punches.....	19, 48	L-shaped Hook Suction Coagulator	28
Bipolar Cable, Wolf Generator	29	Lyons Needle Forceps	23, 61
Bipolar Cables	29	Maryland Dissectors	17, 45, 60
Bipolar Forceps.....	29, 66	Mayo Operating Scissors.....	21, 51
Blunt Suction Coagulator.....	28	Metzenbaum Scissors.....	21, 51, 54
Bowel Forceps	15	Micro Biopsy Forceps	19, 49
Bullet-Nose Dissecting Forceps	17, 47	Micro Bipolar Forceps	29
Cannula & Trocar Pyramid Tip.....	64	Micro Dissecting Forceps	16, 47, 60
Carb-Bite Needle Holders	22, 61	Micro Grasping Forceps.....	16, 47
Cholangiogram Clamp.....	30	Micro Scissors	21, 50, 53
Cholangiogram Guide.....	30, 66	Mini Scissors	21, 51, 54
Claw Forceps	9, 40	Mixer Dissecting Forceps.....	18, 46
DeBakey Forceps.....	9, 38	Mixer Spreader-Dissectors	18, 46, 60
Dilation Set.....	31	Monopolar Adapter.....	29
Dissecting Forceps	16-18, 45-47, 60	Monopolar Cable.....	29
DJB Grasping Forceps.....	9, 38	Myoma Screws.....	31
Dual-Cup Grasping Forceps.....	6, 43, 59	Needle Holders	22-24, 61
Duckbill Grasping Forceps.....	12, 42	Needle Nose Forceps	16, 47
Duval Grasping Forceps.....	12, 43	Offset Laparoscope	58, 66
Electrodes	28	One-Piece Biopsy Instruments.....	19
Electrosurgical Accessories.....	29	One-Piece Dissecting Forceps.....	16-18
Electrosurgical Coagulation.....	28-29	One-Piece Electrosurgical Coagulation	28-29
Endoloop Applicator	24, 65	One-Piece Grasping Forceps.....	4-15
Endoscope & Instrument Holder	57, 66	One-Piece Laparoscopic Instrumentation.....	3-31
Fallopian Tube, Grasping Forceps.....	15	One-Piece Ligation	22-24
Fan Retractor	27	One-Piece Needle Holders	22-24
Foerster Sponge Forceps	13, 42	One-Piece Receivers	22-24
Fundus Grasping Forceps	6, 36, 59	One-Piece Retracting Instruments	26-27
Grabber Forceps.....	5, 41, 59	One-Piece Retractors.....	26-27
Grasping and Dissecting Forceps.....	8, 44	One-Piece Scissors.....	20-21
Grasping and Extracting Forceps	8, 44	One-Piece Specialization	30-31
Grasping Bipolar Forceps.....	29	Operating Laparoscopes	58
Grasping Forceps	4-15, 34-44, 59-60	P.E.E.R. Retractors	27, 66
Grasping Forceps, Fallopian Tube.....	15	Paddle Babcock Forceps	10, 39

Palpation Probes	31
Puncture Needles with Irrigation	25
Pyloromyotomy Dissector	60
Receivers	22-24
Reducer Flip Top	65
Reducer Sleeves	56, 62, 65
Replacement Parts	63-66
Retracting Instruments	26-27, 56
Retraction Grasping Forceps	7, 36
Retractors	26-27, 56
Retrograde Grasping Forceps	8, 37
Right Angle Dissectors	18, 46
Romeo Self-Righting Needle Holders	23, 61
Round Jaw Babcock Forceps	11, 39
Russian Grasping Forceps	12, 42
Scissors	20-21, 50-51, 53-54, 60
Spare Parts	63-66
Spatula Electrodes	28
Spatula Tip Suction Coagulator	28
Specialization	30-31
Spoon Forceps	30, 49
S-Retractors	56
Stone Extractor	14
Storage and Accessories	57
Strong Atraumatic Grasping Forceps	15
Strong Forceps	14, 42, 60
Strong Hold Babcock Forceps	10, 43
Strong TC Grasping Forceps	5, 41
Suction Coagulators	28, 63
Suction/Irrigation Cannulas	25, 62, 63
Suction/Irrigation Set	25, 63
Supercut Endo Scissors	20, 50
Supercut Scissors	20, 50, 53, 60
Suturing Set	24
Swanstrom Grasping Forceps	11, 43
Take-Apart Biopsy Instruments	48-49
Take-Apart Dissecting Forceps	45-47
Take-Apart Grasping Forceps	36-44
Take-Apart Laparoscopic Instruments	33-51
Take-Apart Scissors	50-51
Tapered Dissecting Forceps	16, 47
Tapered Maryland Dissectors	17, 45
TC Grasping Forceps	4, 41
Tenaculum Forceps	14, 40
Thoracic Cannula	56, 64
Trocar/Cannula System	55-56, 62, 63-65
Trocars	55, 62
Vein/Nerve Retractor	26
Veress Needles	31
Winer Grasping Forceps	8
Wolf Generator Bipolar Cable	29



For more information or to place an order, please contact:
Integra ■ 589 Davies Drive, York, PA 17402
877-468-5572 USA ■ 800-275-7192 fax
914-789-7094 outside USA ■ 609-750-4218 outside USA fax
integralife.com

