




Integra®
Forefoot
Mid & HindFoot
Products Catalogue

Table of Contents

Warranty Information	5
Integra® Forefoot Solutions	7
Forefoot Set	8
Bold® Compression Screw dia. 2.5 mm / 3.0 mm / 3.7 mm	14
QWIX® Fixation Screw dia. 3.0 mm	14
SPIN® Screw dia. 2.0 mm / 2.7 mm	15
UNI-CLIP® Compressive Staple	16
SOLUSTAPLE® Standard Staple	16
Osteoguides	17
K-FIX® - K-wire End Protector	17
Weil Clamp	17
B-BOP® LOCK plate	18
IPP-ON® Interphalangeal Arthrodesis Implant	19
METIS® Metatarso-phalangeal prosthesis	20
K2™ Hemi Toe Implant	22
Movement Great Toe system	23
HALLU®-FIX MTP Arthrodesis System	24
HALLU® Plates	24
HALLU® SNAP-OFF® Screws	25
HALLU®-FIX Instruments	25
HALLU®-LOCK Arthrodesis System	26
HALLU®-REAM set	28
DPR System Minimal Invasive Foot Surgery	29
Post-Operative Shoes	30

Integra® Mid & HindFoot Solutions	33
Subtalar Arthrodesis	34
QWIX® Fixation screw	34
Large QWIX® Fixation and Positioning Screws	35
I.CO.S.® Ideal Compression Screw	36
Flatfoot implants	39
KALIX® II Flat Foot Implant	39
MBA™ Titanium Subtalar Implant	40
Dorsal Lisfranc plate & Medial Lisfranc plate	41
ADVANSYS® Lisfranc Plate	41
UNI-CP™ Compression Plate	43
Ankle Solutions	44
Total Ankle Replacement Hintegra® prosthesis	44
Ankle Fusion	54
PANTA® Nail	54
ADVANSYS® Plate	59
TIBIAXYS® Plating System	61
Minimally Invasive Achilles Tendon Ruptures	64
ACHILLON® System	64
Surgical Instruments	65
Hintermann® Distractors	65
Osteoguides	66
K-FIX® Protector	66
Weil Clamp	66
General Optional Instruments	67
Post-Operatives Shoes	68
Product References	71



Pricing:

Prices stated in specific written quotations are firm for thirty days from the date given, and are otherwise subject to change without prior notice. Pricing terms stated on written agreements are governed by such agreements.

Minimum Order Requirements:

Minimum order requirement is €250 (a €20 minimum order requirement charge will be added to any order under €250).

Ordering Procedure:

A written purchase order on the customer's form may be requested for purchases of Integra LifeSciences products.

Acceptance of Orders:

Orders are accepted upon approval by Integra Customer Service.

Return Policy:

- Authorization from Customer Service must be obtained prior to returning product.
- Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.
- Non-sterile product must be returned in unused saleable condition in original package.
- Custom or special order products will not be accepted for credit.
- Credit will be issued for goods returned prior to ninety days from ship date with a 20% restocking charge. This assumes that the product returned is not damaged and can be verified to have not been used or opened.

Integra Limited Warranty – Integra products

INTEGRA LIFESCIENCES CORPORATION and its wholly owned subsidiaries ("INTEGRA") warrant to INTEGRA authorized distributors and the original purchaser only that each new INTEGRA Surfix[®] product is free from manufacturing defects in material and workmanship when used for its intended surgical purpose. For purposes of products sold by INTEGRA through an authorized distributor of INTEGRA, "original purchaser" shall include the purchaser of INTEGRA products to whom the distributor first sells the product.

If any covered defect occurs, the purchaser or distributor should communicate directly with INTEGRA. If purchaser or distributor seeks to invoke the terms of this warranty, the product must be returned to INTEGRA. The defective product should be returned promptly, properly packaged and postage prepaid. Loss or damage in return shipment to INTEGRA shall be at sender's risk. INTEGRA's sole responsibility under this warranty shall be repair or replacement, at INTEGRA's sole discretion at INTEGRA's expense, subject to the terms of this warranty and applicable agreements.

INTEGRA provides for no commercial or additional warranty whatsoever, except as may be otherwise expressly agreed in a separate commercial agreement or undertaking duly signed by an INTEGRA authorised representative.

IN NO EVENT SHALL INTEGRA BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA shall not be responsible for, any loss arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA.

Customer Service Contact Information:

International

+33 (0)4 37 47 59 50
+33 (0)4 37 47 59 25 (Fax)
csemea@Integralife.com

Benelux

+32 (0)2 257 4130
+32 (0)2 253 2466 (Fax)
custsvcbenelux@Integralife.com

France

+33 (0)4 37 47 59 10
+33 (0)4 37 47 59 29 (Fax)
cs-ortho@Integralife.com

Switzerland

+41 (0)2 27 21 23 30
+41 (0)2 27 21 23 99 (Fax)
custsvcsuisse@Integralife.com

United Kingdom

+44 (0)1 264 345 780
+44 (0)1 264 363 782 (Fax)
custsvcs.uk@integralife.com

THIS INTEGRA LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART OR THE PART OF ITS DISTRIBUTORS, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS.

INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS. No warranty or guarantee may be created by any act or statement nor may this Standard Warranty be modified in any way, except as a result of a writing signed by an officer of INTEGRA. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct.

IN NO EVENT SHALL INTEGRA AUTHORIZED DISTRIBUTORS BE LIABLE TOWARDS THE ORIGINAL PURCHASER FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA authorized distributors shall not be responsible towards the original purchaser for, any loss, arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA. THIS INTEGRA DISTRIBUTOR LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES TOWARDS THE ORIGINAL PURCHASER ON INTEGRA AUTHORIZED DISTRIBUTOR'S PART.

INTEGRA AUTHORISED DISTRIBUTORS DISCLAIM ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS.

All the references numbers mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED".






Integra® Forefoot Solutions

Forefoot Set

A new modular and scalable forefoot instruments set:

- Bold® 2.5 – 3.0 – 3.7 Compression screw.
- QWIX® 3.0 Fixation screw.
- Spin® 2.0 – 2.7 Spin break off screws.
- Solustaple® Standard staple.
- Uni-clip® Compression Staple.



Screw instrumentation module
(details page 8)



Staple instrumentation module
(details page 9)



Generic module
(details page 10)



Implant rack module
(details page 10)



Frame

229996	Frame - 2 floors
229990	Frame - 3 floors
229970	Forefoot II- lid- generic



- Staple instrumentation
- Screw instrumentation




- Non sterile Implants
- Staple instrumentation
- Screw instrumentation



- Non sterile Implants
- Generic module
- Screw instrumentation

Forefoot Set – Screw Module

Reference	Description
119105	Scarf clamp
229002	Screwdriver - Hexa 2.0
229005	AO Screwdriver - Hexa 2.0 - AO Attachment - Cannulated
229003	Screwdriver - T7 Star - Cannulated
229004	Screwdriver - AO - T7 Star - Cannulated
997301	Depth Gauge
997311	AO Straight Handle
229011	Longitudinal Cutting Guide (Optional)
229012	Transversal Cutting Guide (Optional)
229101	Ruler
229951	Sterilization Screws Container including:
229961	Tray
229970	Lid
229971	Mat



Drill Reference	Corresponding Product	2 in 1	Diameter	Cannulated	Length	AO attachment
159 027 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	12 mm	■
159 025 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	22 mm	■
159 023 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	32 mm	■
119 227 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		12 mm	
119 225 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		22 mm	
119 223 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		32 mm	
119 027 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	12 mm	
119 025 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	22 mm	
119 023 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	32 mm	
159 227 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		12 mm	■
159 225 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		22 mm	■
159 223 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		32 mm	■
229 027 S	Bold® 2.5	■	1.9 - 2.5 mm	■	12 mm	■
229 037 S	Bold® 2.5	■	1.9 - 2.5 mm	■	22 mm	■
229 227 S	Bold® 2.5	■	1.9 - 2.5 mm	■	12 mm	
229 237 S	Bold® 2.5	■	1.9 - 2.5 mm	■	22 mm	
229 028 S	Bold® 2.5	■	1.9 - 2.5 mm		12 mm	■
229 038 S	Bold® 2.5	■	1.9 - 2.5 mm		22 mm	■
229 228 S	Bold® 2.5	■	1.9 - 2.5 mm		12 mm	
229 238 S	Bold® 2.5	■	1.9 - 2.5 mm		22 mm	

K-wire Reference	Corresponding Product	Diameter	Length
115 070 (S)	Bold® 3.0 / Bold® 3.7	1.0 mm	70 mm
115 100 (S)	Bold® 3.0 / Bold® 3.7	1.0 mm	100 mm
115 008 S	Bold® 2.5	0.8 mm	70 mm

All drills and K-wires are delivered sterile or non sterile except Bold® 2.5 K-wire which is only delivered in sterile.

(S) means available in sterile and non sterile

Forefoot Set – Staple Module

UNI-CLIP® Associated Instrument


229201	Uni-Clip® Drilling guide
229202	Uni-Clip® impactor
119311	Uni-Clip® spreading forceps

Solustaple® Associated Instrument

229102	Solustaple® impactor 90° & 26°
229101	Solustaple® drilling guide
997301	Depth gauge - dia 2.0 X L. 10 - 60 mm

Sterilization Container

229952	Staples container includes:	
229962	Tray	
229970	lid	



Drill Reference	Corresponding Product	Diameter	Cannulated	Length	AO attachment
119 004 (S)	Uni-Clip®	2.2 mm	■	34 mm	
119 006 (S)	Uni-Clip®	2.2 mm		34 mm	
159 004 (S)	Uni-Clip®	2.2 mm	■	34 mm	■
119 016 (S)	Uni-Clip®	2.2 mm		34 mm	■


K-wire Reference	Corresponding Product	Diameter	Length
115 070 (S)	Uni-Clip®	1.0 mm	70 mm
115 100 (S)	Uni-Clip®	1.0 mm	100 mm

All drills and K-wires are delivered sterile or non sterile except Bold® 2.5 K-wire which is only delivered in sterile.

(S) means available in sterile and non sterile

Sterilization Container


229965	Forefoot II - tray - non sterile implants
229975	Forefoot II - rack - non sterile implants



Forefoot Set – Generic Module

Sterilization Container

229966	Forefoot II - tray - generic
229976	Forefoot II - mat - generic



Bold® Compression Screw Dia. 2.5 mm / 3.0 mm / 3.7 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Totally intra-osseous implants.
- Large range of sizes:
 - > 10 to 30 mm for Bold® 2.5 screws by 2 mm increment.
 - > 10 to 34 mm for Bold® 3.0 screws by 2 mm increment.
 - > 14 to 34 mm for Bold® 3.7 screws by 2 mm increment.
- K-Wire guided drilling and insertion.
- Self tapping.
- True & controlled compression.
- Dual threads.
- Material: Titanium alloy TA6V.



Bold® dia 2.5 mm

Reference	Length
200010(S)	L. 10 mm
200012(S)	L. 12 mm
200014(S)	L. 14 mm
200016(S)	L. 16 mm
200018(S)	L. 18 mm
200020(S)	L. 20 mm
200022(S)	L. 22 mm
200024(S)	L. 24 mm
200026(S)	L. 26 mm
200028(S)	L. 28 mm
200030(S)	L. 30 mm



Bold® dia 3.0 mm

Reference	Length
111010(S)	L. 10 mm
111012(S)	L. 12 mm
111014(S)	L. 14 mm
111016(S)	L. 16 mm
111018(S)	L. 18 mm
111020(S)	L. 20 mm
111022(S)	L. 22 mm
111024(S)	L. 24 mm
111026(S)	L. 26 mm
111028(S)	L. 28 mm
111030(S)	L. 30 mm
111032(S)	L. 32 mm
111034(S)	L. 34 mm



Bold® dia 3.7 mm

Reference	Length
131014(S)	L. 14 mm
131016(S)	L. 16 mm
131018(S)	L. 18 mm
131020(S)	L. 20 mm
131022(S)	L. 22 mm
131024(S)	L. 24 mm
131026(S)	L. 26 mm
131028(S)	L. 28 mm
131030(S)	L. 30 mm
131032(S)	L. 32 mm
131034(S)	L. 34 mm

QWIX® Fixation Screw, Dia. 3.0 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Self-drilling and self-tapping.
- Totally intra-osseous fixation.
- K-Wire guided drilling and insertion.
- Material: Titanium alloy TA6V.



QWIX® Dia. 3.0 mm

Reference	Length
111312(S)	L. 12 mm
111314(S)	L. 14 mm
111316(S)	L. 16 mm
111318(S)	L. 18 mm
111320(S)	L. 20 mm
111322(S)	L. 22 mm
111324(S)	L. 24 mm
111326(S)	L. 26 mm
111328(S)	L. 28 mm
111330(S)	L. 30 mm
111332(S)	L. 32 mm
111334(S)	L. 34 mm

(S) means available in sterile and non sterile

SPIN[®] Break Off Screw, Dia. 2.0 mm (for associated instrumentation, see p. 8)

For bone fixation.

- Break-off barrel.
- Self-drilling and self-tapping.
- Material: Titanium alloy TA6V.

Reference	Description
112011(S)	Spin [®] Dia. 2.0 x L. 11 mm
112012(S)	Spin [®] Dia. 2.0 x L. 12 mm
112013(S)	Spin [®] Dia. 2.0 x L. 13 mm
112014(S)	Spin [®] Dia. 2.0 x L. 14 mm



SPIN[®] Break Off Screw, Dia. 2.7 mm (for associated instrumentation, see p. 8)

Spin[®] 2.7 screw is recommended for soft or osteoporotic bones fixation.

- Break-off barrel.
- Self-drilling and self-tapping.
- Material: Titanium alloy TA6V.

Reference	Length
112111(S)	Spin [®] Dia. 2.7 x L. 11 mm
112112(S)	Spin [®] Dia. 2.7 x L. 12 mm
112113(S)	Spin [®] Dia. 2.7 x L. 13 mm
112114(S)	Spin [®] Dia. 2.7 x L. 14 mm



(S) means available in sterile and non sterile

UNI-CLIP® Notched Staples (for associated instrumentation, see p. 9)

For fixation of bone fractures or for bone reconstruction.

- Range sizes:
 - > Interaxis: 11, 12, 13, 15 and 20 mm.
 - > Leg lengths: 12, 13, 14, 15, 16, 17 and 20 mm.
- Material: Stainless steel 316L.



Reference	Description
213113(S)	Uni-Clip® interaxis 11.0 x L. 13 mm
213114(S)	Uni-Clip® interaxis 11.0 x L. 14 mm
213115(S)	Uni-Clip® interaxis 11.0 x L. 15 mm
213116(S)	Uni-Clip® interaxis 11.0 x L. 16 mm
213117(S)	Uni-Clip® interaxis 11.0 x L. 17 mm
213213(S)	Uni-Clip® interaxis 12.0 x L. 13 mm
213214(S)	Uni-Clip® interaxis 12.0 x L. 14 mm
213215(S)	Uni-Clip® interaxis 12.0 x L. 15 mm
213216(S)	Uni-Clip® interaxis 12.0 x L. 16 mm
213217(S)	Uni-Clip® interaxis 12.0 x L. 17 mm
213313(S)	Uni-Clip® interaxis 13.0 x L. 13 mm
213314(S)	Uni-Clip® interaxis 13.0 x L. 14 mm
213315(S)	Uni-Clip® interaxis 13.0 x L. 15 mm
213316(S)	Uni-Clip® interaxis 13.0 x L. 16 mm
213317(S)	Uni-Clip® interaxis 13.0 x L. 17 mm
213512(S)	Uni-Clip® interaxis 15.0 x L. 12 mm
213820(S)	Uni-Clip® interaxis 20.0 x L. 20 mm

SOLUSTAPLE® Standard Staple (for associated instrumentation, see p. 9)

Static staple for fixation of mono cortical phalangeal osteotomies.

- Two shapes of staples (26° and 90°) are available.
- Self-perforating.
- Material: stainless steel 316L.

Reference	Description
114002(S)	Solustaple® 90° interaxis 8 mm
114004(S)	Solustaple® 90° interaxis 10 mm
114023(S)	Solustaple® 26° interaxis 8 mm
114025(S)	Solustaple® 26° interaxis 10 mm



(S) means available in sterile and non sterile

Osteoguides

- Suggested use:
 - > Scarf osteotomy.
 - > Chevron osteotomy.
 - > Basal osteotomy.
 - > Arthrodesis.
 - > Wedge osteotomies (Akin...).
- Details:
 - > Hole for K-Wire: Diameter 1.0 to 1.6 mm.
 - > Knurled handle for optimal grip.
 - > Blade thickness maxi 0.8 mm.
 - > Positioning tip.

Reference	Description
229011	Longitudinal
229012	Transversal



K-FIX[®] – K-Wire End Protector

The K-Fix[®] pin protector is a simple system to protect patients or hospital staff from K-Wires sharp ends.

- Delivered sterile.

Reference	Length
119050	K-FIX [®] for K-WIRE - Dia. 0.8 mm to 2.0 mm (5 units per pack)



Weil Clamp

Reference	Length
119106	Weil Clamp

B-BOP® LOCK Plate

B-BOP®-LOCK plate is indicated for fixation of osteotomy of the basis of the first metatarsal. Examples include:

- Moderate to severe Hallux valgus.
- Hallux varus.


Implants

Reference	Description
320 100S	Basal plantar plate B-Bop® Lock right
320 200S	Basal plantar plate B-Bop® Lock left
295 110S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 10 mm
295 112S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 12 mm
295 114S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 14 mm
295 116S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 16 mm
295 118S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 18 mm
295 120S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 20 mm
295 122S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 22 mm
295 124S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 24 mm
295 126S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 26 mm
295 128S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 28 mm
295 130S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 30 mm
295 132S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 32 mm
295 134S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 34 mm
195 100S	Lock screw B-Bop® Lock



Instrumentation

Reference	Description
219 027	Drilling guide Surfix Alpha Dia 2.0 mm
219 127	Screwdriver Torx T7
219 227	Screwdriver AO blade Torx T7
299 050	Drill AO Dia 2.0 mm
299 005	Drill Dia 2.0 mm
329 500	Bending forceps
329 002	Forceps
115 070	K-wire Dia 1.0 mm x L 70 mm
115 116	K-wire Dia 1.6 mm x L 150 mm
299 020	Depth gauge B-Bop Lock®
189 210	Wedge 10°
189 215	Wedge 15°
189 220	Wedge 20°



Container

Reference	Description
329 900	Container
including:	
329 901	Basis
329 902	Silicone mat
996200	Lid

IPP-ON[®] Interphalangeal Arthrodesis Implant

IPP-ON[®] interphalangeal implant is an implant which respects anatomical flexion to achieve interphalangeal arthrodesis of the lesser toes in case of:

- Fixed or semi-fixed hammer toe.
- Revision in case of arthrodesis or arthroplasty failure.
- Shortening in case of length excess of the second toe.
- Material: Stainless steel: 316 LVM.
- 2 sizes.




Reference	Description
230001S	IPP-ON [®] Interphalangeal Implant Size 1
230002S	IPP-ON [®] Interphalangeal Implant Size 2

IPP-ON[®] Instruments

Reference	Description
239011	Distal Hand Drill Size 1
239012	Proximal Hand Drill Size 1
239021	Distal Hand Drill Size 2
239022	Proximal Hand Drill Size 2
239500	Clamp
239030	Implant Holder

Sterilization Container	
239000	IPP-ON [®] sterilization container including:
239001	Base
996100	Lid
278902	Mat
119909	Silicone wedge



METIS® Metatarso-Phalangeal Prosthesis

The three-component METIS® prosthesis replaces the articular surfaces of the great toe joint.

- Metatarsal component:
 - > 7 sizes (0 to 6).
 - > Material: Cobalt-Chromium alloy (CoCr).
 - > Coating: Porous Titanium & Hydroxyapatite (Ti and HAP).
- PE inlay:
 - > Ultra High Molecular Weight Polyethylene (UHMWPE).
 - > 5 sizes (1 to 5), 3 thicknesses (1 to 3).
- Phalangeal implant:
 - > 5 sizes (1 to 5).
 - > Material: Titanium alloy TA6V.
 - > Coating: Porous Titanium (Ti).



Reference	Description
250100	Metatarsal component size 0
250101	Metatarsal component size 1
250102	Metatarsal component size 2
250103	Metatarsal component size 3
250104	Metatarsal component size 4
250105	Metatarsal component size 5
250106	Metatarsal component size 6
250201	Phalangeal component size 1
250202	Phalangeal component size 2
250203	Phalangeal component size 3

Reference	Description
250204	Phalangeal component size 4
250205	Phalangeal component size 5
250311	METIS® Polyethylene Inlay - Size 1 - Thickness 1
250312	METIS® Polyethylene Inlay - Size 1 - Thickness 2
250313	METIS® Polyethylene Inlay - Size 1 - Thickness 3
250321	METIS® Polyethylene Inlay - Size 2 - Thickness 1
250322	METIS® Polyethylene Inlay - Size 2 - Thickness 2
250323	METIS® Polyethylene Inlay - Size 2 - Thickness 3
250331	METIS® Polyethylene Inlay - Size 3 - Thickness 1
250332	METIS® Polyethylene Inlay - Size 3 - Thickness 2
250333	METIS® Polyethylene Inlay - Size 3 - Thickness 3
250341	METIS® Polyethylene Inlay - Size 4 - Thickness 1
250342	METIS® Polyethylene Inlay - Size 4 - Thickness 2
250343	METIS® Polyethylene Inlay - Size 4 - Thickness 3
250351	METIS® Polyethylene Inlay - Size 5 - Thickness 1
250352	METIS® Polyethylene Inlay - Size 5 - Thickness 2
250353	METIS® Polyethylene Inlay - Size 5 - Thickness 3


METIS® Instruments

METIS® Instruments

Reference	Description
139107	Metatarsal Osteotome
139120	Implant holder
139131	Locking screw
139210	Phalangeal reamer starter size
139211–215	Phalangeal reamer size 1 to 5
139140	Metatarsal sizer size 0
139121–126	Metatarsal sizer size 1 to 6
139150	Metatarsal cutting guide size 0
139111–116	Metatarsal cutting guide size 1 to 6
139216	Metatarsal impactor
139128	Visualization guide
139129	Visualization guide wire
119645	Screwdriver hexa 3.5 mm
119 601	Screwdriver hexa 2.5 mm
115216	K-Wire Dia. 1.6 x L. 150 mm
119611	Drill Dia. 2.5 x L. 150 mm
139940	METIS® sterilization container including:
139941	Metatarsal base
139942	Phalangeal base
139943	Base
139953	Lid
119909	Silicone wedge

METIS® Trial Implants

Reference	Description
139100–106	Metatarsal trial implant size 0 to 6
139311–313	Trial inlay size 1 thickness 1 to 3
139321–323	Trial inlay size 2 thickness 1 to 3
139331–333	Trial inlay size 3 thickness 1 to 3
139341–343	Trial inlay size 4 thickness 1 to 3
139351–353	Trial inlay size 5 thickness 1 to 3
139201–205	Phalangeal trial implant size 1 to 5



K2™ Hemi Toe Implant


The K2™ Hemi Toe Implant System is an anatomically designed, metallic, one-piece surface replacement system for the base of the hallux. 4 sizes are available.

K2™ Implants

Reference	Description
201010	Phalangeal, Size 1
201020	Phalangeal, Size 2
201030	Phalangeal, Size 3
201040	Phalangeal, Size 4

K2™ Instruments

Reference	Description
202010	Size 1 Trial
202020	Size 2 Trial
202030	Size 3 Trial
202040	Size 4 Trial
203010	Size 1 Shaper
203020	Size 2 Shaper
203030	Size 3 Shaper
203040	Size 4 Shaper
203100	Cutting Guide
203300	Sizer/Drill Guide
203500	Broach Handle
203600	Impactor
205000	Sterilization tray
203200	Reamer
203415	1.5 mm 80 mm drill bit



Movement Great Toe System

- Hemi Proximal Phalanx Implant.
- Hemi Metatarsal Implant.
- Total Great Toe Implant.

Reference	Description
MGT-890-10PPH	Hemi Proximal Phalanx, Size 10
MGT-890-20PPH	Hemi Proximal Phalanx, Size 20
MGT-890-30PPH	Hemi Proximal Phalanx, Size 30
MGT-890-40PPH	Hemi Proximal Phalanx, Size 40
MGT-890-10MT	Hemi Metatarsal, Size 10
MGT-890-20MT	Hemi Metatarsal, Size 20
MGT-890-30MT	Hemi Metatarsal, Size 30
MGT-890-40MT	Hemi Metatarsal, Size 40
MGT-890-10PPT	Phalangeal Total, Size 10
MGT-890-20PPT	Phalangeal Total, Size 20
MGT-890-30PPT	Phalangeal Total, Size 30
MGT-890-40PPT	Phalangeal Total, Size 40

Reference	Description
INS-890-00	Instrument Set
MIS-890-00	Metatarsal Implant Sizer
PIS-890-00	Proximal Phalanx Implant Sizer
GDW-890-00	Guide Pin, 2.0 mm X 125 mm
MSR-890-10/20	Metatarsal Surface Reamer, Size 10/20
MSR-890-30/40	Metatarsal Surface Reamer, Size 30/40
PSR-890-10/20	Proximal Phalanx Surface Reamer, Size 10/20
PSR-890-30/40	Proximal Phalanx Surface Reamer, Size 30/40
DRL-890-00	Cannulated Drill, 4.5 mm
DCG-890-00	Metatarsal Dorsal Cutting Guide
TRL-890-10MT	Metatarsal Trial, Size 10
TRL-890-20MT	Metatarsal Trial, Size 20
TRL-890-30MT	Metatarsal Trial, Size 30
TRL-890-40MT	Metatarsal Trial, Size 40
TRL-890-10PPH	Proximal Phalanx Hemi Trial, Size 10
TRL-890-20PPH	Proximal Phalanx Hemi Trial, Size 20
TRL-890-30PPH	Proximal Phalanx Hemi Trial, Size 30
TRL-890-40PPH	Proximal Phalanx Hemi Trial, Size 40
TRL-890-10PPT	Phalangeal Total Trial, Size 10
TRL-890-20PPT	Phalangeal Total Trial, Size 20
TRL-890-30PPT	Phalangeal Total Trial, Size 30
TRL-890-40PPT	Phalangeal Total Trial, Size 40
IMP-890-00MT	Metatarsal Impactor
IMP-890-00PP	Proximal Phalanx Impactor
ALG-890-00	Sizer Alignment Guide



HALLU®-FIX MTP Arthrodesis System

The HALLU®-FIX solution is a system dedicated to first MTP arthrodesis.

- Pre-bent plate.
- SNAP-OFF® screw fixation.
- Color coded.
- 3 sizes (3 left / 3 right).
- Material: Titanium alloy TA6V.

HALLU®-C Plate

HALLU®-C Plates

Reference	Description
117340(S)	HALLU®-C Plate Right-Size 1: L. 40 mm - 4 holes
117345(S)	HALLU®-C Plate Right-Size 2: L. 45 mm - 5 holes
117350(S)	HALLU®-C Plate Right-Size 3: L. 50 mm - 6 holes
117440(S)	HALLU®-C Plate Left-Size 1: L. 40 mm - 4 holes
117445(S)	HALLU®-C Plate Left-Size 2: L. 45 mm - 5 holes
117450(S)	HALLU®-C Plate Left-Size 3: L. 50 mm - 6 holes



HALLU®-S Plate

HALLU®-S Plates

Reference	Description
117150(S)	HALLU®-S Plate Right-Size 1: L. 45 mm - 3 holes
117155(S)	HALLU®-S Plate Right-Size 2: L. 50 mm - 4 holes
117160(S)	HALLU®-S Plate Right-Size 3: L. 55 mm - 5 holes
117250(S)	HALLU®-S Plate Left-Size 1: L. 45 mm - 3 holes
117255(S)	HALLU®-S Plate Left-Size 2: L. 50 mm - 4 holes
117260(S)	HALLU®-S Plate Left-Size 3: L. 55 mm - 5 holes



HALLU® SNAP-OFF® Screws

- 2 diameters: 2.7 mm & 3.0 mm (for porous bone).
- Self-tapping.
- Color code for length identification.
- Material: Titanium alloy TA6V.

HALLU® Snap-Off® screws Dia. 3.0 mm

Reference	Length
117110(S)	10 mm
117112(S)	12 mm
117114(S)	14 mm
117116(S)	16 mm
117118(S)	18 mm

HALLU® Snap-Off® screws Dia. 2.7 mm

Reference	Length
117010(S)	10 mm
117012(S)	12 mm
117014(S)	14 mm
117016(S)	16 mm
117018(S)	18 mm
117020(S)	20 mm
117022(S)	22 mm
117024(S)	24 mm
117026(S)	26 mm
117028(S)	28 mm
117030(S)	30 mm
117032(S)	32 mm
117034(S)	34 mm



HALLU®-FIX Instruments



HALLU®-Ream Metatarsal

Reference	Diameter
129714	Dia. 14 mm
129716	Dia. 16 mm
129718	Dia. 18 mm
129720	Dia. 20 mm
129722	Dia. 22 mm

HALLU®-Ream Phalangeal

Reference	Diameter
129724	Dia. 14 mm
129726	Dia. 16 mm
129728	Dia. 18 mm
129730	Dia. 20 mm
129752	Dia. 22 mm

HALLU® Associated instruments

Reference	Description
129710	Quick coupling
129733	SNAP-OFF® Cannulated screwdriver - Dia. 2.7 mm and 3.0 mm
129731	Plate Bender Left
129732	Plate Bender Right
129734	Drilling Guide Dia. 1.9 mm
129735	Screwdriver tip
129736	Depth gauge
119618	Drill Dia. 1.9 x L. 100 mm
115216	K-Wire Dia. 1.6 x L. 150 mm
115100	K-Wire Dia. 1.0 x L. 100 mm



Sterilization container

129930	HALLU®-FIX Sterilization Container including:
129931	HALLU®-FIX Basis
129933	HALLU®-S Plate Insert
129932	HALLU®-C Plate Insert
996200	Lid
129936	HALLU®-FIX screws rack
119909	Silicone wedge (Optional)

(S) means available in sterile and non sterile

HALLU®-LOCK MTP Arthrodesis System

The HALLU®-LOCK solution is a locked system dedicated to first MTP arthrodesis:

Pre-bent plate.

- Color coded.
- 3 sizes (3 left / 3 right).
- Material: Titanium alloy TA 6V.
- Monobloc fixation with Surfix® locking system.

Surfix® locking screw

- 2 diameters: 2.7 mm & 3.0 mm (for porous bone).
- Standard and orientable screws.
- Color code for length identification.
- Material: Titanium alloy TA 6V.
- Sterily Packaged.



Surfix® standard screw
Dia. 2.7 mm + lock-screw



Reference	Length
285210S	10 mm
285212S	12 mm
285214S	14 mm
285216S	16 mm
285218S	18 mm
285220S	20 mm
285222S	22 mm
285224S	24 mm
285226S	26 mm
185200S	lock-screw

Surfix® standard screw
Dia. 3.0 mm + lock-screw



Reference	Length
285110S	10 mm
285112S	12 mm
285114S	14 mm
285116S	16 mm
285118S	18 mm
285120S	20 mm
285122S	22 mm
285124S	24 mm
285126S	26 mm
185100S	lock-screw

Surfix® Alpha
variable angle screw
Dia. 2.7 mm + lock-screw



Reference	Length
295210S	10 mm
295212S	12 mm
295214S	14 mm
295216S	16 mm
295218S	18 mm
295220S	20 mm
295222S	22 mm
295224S	24 mm
295226S	26 mm
195200S	lock-screw

Surfix® Alpha
variable angle screw
Dia. 3.0 mm + lock-screw



Reference	Length
295110S	10 mm
295112S	12 mm
295114S	14 mm
295116S	16 mm
295118S	18 mm
295120S	20 mm
295122S	22 mm
295124S	24 mm
295126S	26 mm
195100S	lock-screw

Hallu®-Lock C plate 10°

Reference	Description
290340S	Right-Size 1: L.40 mm – 4 holes
290345S	Right-Size 2: L.45 mm – 5 holes
290350S	Right-Size 3: L.50 mm – 6 holes
290440S	Left-Size 1: L.40 mm – 4 holes
290445S	Left-Size 2: L.45 mm – 5 holes
290450S	Left-Size 3: L.50 mm – 6 holes

Hallu®-Lock S plate 10°
Number of holes excludes the 4 extremity holes

Reference	Description
290150S	Right-Size 1: L.45 mm – 3 holes
290155S	Right-Size 2: L.50 mm – 3 holes
290160S	Right-Size 3: L.55 mm – 3 holes
290250S	Left-Size 1: L.45 mm – 3 holes
290255S	Left-Size 2: L.50 mm – 3 holes
290260S	Left-Size 3: L.55 mm – 3 holes

QWIX® transarticular
compressive screw Dia.
3.0 mm

Reference	Length
111328S	28 mm
111330S	30 mm
111332S	32 mm
111334S	34 mm





HALLU®-LOCK Instruments

HALLU®-REAM System

129714	Metatarsal Reamer Dia. 14 mm
129716	Metatarsal Reamer Dia. 16 mm
129718	Metatarsal Reamer Dia. 18 mm
129720	Metatarsal Reamer Dia. 20 mm
129722	Metatarsal Reamer Dia. 22 mm
129724	Phalangeal Reamer Dia. 14 mm
129726	Phalangeal Reamer Dia. 16 mm
129728	Phalangeal Reamer Dia. 18 mm
129730	Phalangeal Reamer Dia. 20 mm
129752	Phalangeal Reamer Dia. 22 mm

HALLU®-Lock associated instruments

129710	Quick coupling
997100	Quick-Coupling AO
299070	Left bender
299080	Right bender
219027	SURFIX® Alpha drilling guide Dia. 2.0 mm
219627	SURFIX® drilling guide Dia. 2.0 mm
219127	Screwdriver torx T7
219227	AO Screwdriver tip torx T7
299020	Depth Gauge
119133	Handle with AO Attachment Dia. 3.0 mm
119135	AO Screwdriver tip dia 3.0 mm
299010	Trial plate holder
299150	Trial plate-S Right - Size 1
299155	Trial plate-S Right - Size 2
299160	Trial plate-S Right - Size 3
299250	Trial plate-S Left - Size 1
299255	Trial plate-S Left - Size 2
299260	Trial plate-S Left - Size 3
299340	Trial plate-C Right - Size 1
299345	Trial plate-C Right - Size 2
299350	Trial plate-C Right - Size 3
299440	Trial plate-C Left - Size 1
299445	Trial plate-C Left - Size 2
299450	Trial plate-C Left - Size 3
299005	Drill Dia. 2.0 x L. 100 mm
299050	AO Drill AO Dia. 2.0 x L. 100 mm
159027S	2 in 1 cannulated drill AO dia 2.2/3.0 x L. 12 mm sterile
115070	K-Wire Dia. 1.0 x L. 70 mm
115100	K-wire Dia. 1.0 x L. 100 mm
115116	K-Wire Dia. 1.6 x L. 150 mm



HALLU®-LOCK sterilization container

299900	HALLU®-Lock sterilization container
299901	Base
299903	Trial implants tray
299904	K-Wire tray
996200	Lid
119909	Silicone wedge

HALLU®-REAM Set

Reference	Description
129714	Metatarsal Reamer Dia. 14 mm
129716	Metatarsal Reamer Dia. 16 mm
129718	Metatarsal Reamer Dia. 18 mm
129720	Metatarsal Reamer Dia. 20 mm
129722	Metatarsal Reamer Dia. 22 mm

129724	Phalangeal Reamer Dia. 14 mm
129726	Phalangeal Reamer Dia. 16 mm
129728	Phalangeal Reamer Dia. 18 mm
129730	Phalangeal Reamer Dia. 20 mm
129752	Phalangeal Reamer Dia. 22 mm

Instruments

129710	Quick coupling
115216	K-Wire Dia. 1.6 x L. 150 mm

Sterilization container

129935	HALLU®-Ream sterilization container including:
129938	HALLU®-Ream Basis
996100	Lid



DPR System Minimal Invasive Foot Surgery

DPR System minimal invasive foot surgery was born from our collaboration with minimal invasive foot surgeons.

- Single use: Sterile packaging.
- Safety: Protection by silicone ring.
- Identification: Color coded.
- Material: Stainless Steel 316L.



DPR System Minimal Invasive Foot Surgery burs

278002S	Burr straight Dia. 2.0 x L. 15 mm
278007S	Burr conical Dia. 3.1/1.0 x L. 12 mm
278005S	Burr straight Dia. 1.9 x L. 10 mm
278008S	Burr conical Dia. 4.1/1.0 x L. 12 mm
278009S	Burr cylindrical Dia. 4.1 x L. 15 mm
278010S	Burr cylindrical Dia. 3.1 x L. 15 mm

DPR System Minimal Invasive Foot Surgery curettes

278121	Curette Medium - Down
278122	Curette Medium - Up
278123	Curette Large - Down
278124	Curette Large - Up
278100	Sleeve protection Cap
278110	Sleeve protection Dia. 3.1 mm
278109	Sleeve protection Dia. 4.1 mm

Sterilization container

278900	DPR System Minimal Invasive Foot Surgery sterilization container
278901	Base
996100	Lid
278902	Mat
119909	Silicone wedge



Please contact us if you need a MIS motor offer.

Post-Operative Shoes

Reference Description

Flat Post Operative Shoe


999001	Small
999002	Medium
999003	Large
999004	X-Large

Reference Description

Wedged Heel Post Operative Shoe

999100	Size 4 to 6.5
999101	Size 7 to 8.5
999102	Size 9 to 10.5





Integra® Mid & HindFoot Solutions

QWIX® Fixation Screw, dia. 4.3 mm

For the fixation of bone fractures or for bone reconstruction.

- Self-drilling and self-tapping.
- Intra-osseous fixation.
- Variable lag length for extremity fixation.
- Material: Titanium alloy TA6V.

QWIX® dia. 4.3 mm

Reference	Length
111422S	L. 22 mm
111424S	L. 24 mm
111426S	L. 26 mm
111428S	L. 28 mm
111430S	L. 30 mm
111432S	L. 32 mm
111434S	L. 34 mm
111436S	L. 36 mm
111438S	L. 38 mm
111440S	L. 40 mm
111445S	L. 45 mm
111450S	L. 50 mm
111455S	L. 55 mm
111460S	L. 60 mm




QWIX® dia. 4.3 mm Instruments

Reference	Description
119104	Screwdriver tip
119139	Short Drill dia. 2.6/5.0 x L. 100 mm
119138	Quick coupling handle
119516	Long cannulated drill dia. 2.6 x L. 140 mm
115116	K-Wire dia. 1.6 x L. 150 mm

Sterilization container

Reference	Description
119935	QWIX® 4.3 sterilization container including:
119937	Cylinder
119936	Base
119938	Mat
996400	Lid



Large QWIX® Fixation and Positioning Screws dia. 5.5 mm / 7.5 mm

- Fixation and positioning screws dia. 5.5 and 7.5 mm.
- Self-drilling and self-tapping.
- Intra-osseous fixation.
- Material: Titanium alloy TA6V.
- Sterile packaged.

Large QWIX® Fixation screws
dia. 7.5 mm



Reference	Length
111740s	L. 40 mm
111745s	L. 45 mm
111750s	L. 50 mm
111755s	L. 55 mm
111760s	L. 60 mm
111765s	L. 65 mm
111770s	L. 70 mm
111775s	L. 75 mm
111780s	L. 80 mm
111785s	L. 85 mm
111790s	L. 90 mm
111795s	L. 95 mm
111800s	L. 100 mm
111805s	L. 105 mm
111810s	L. 110 mm
111815s	L. 115 mm
111820s	L. 120 mm



Large QWIX® Fixation screws
dia. 5.5 mm

Reference	Length
111530s	L. 30 mm
111535s	L. 35 mm
111540s	L. 40 mm
111545s	L. 45 mm
111550s	L. 50 mm
111555s	L. 55 mm
111560s	L. 60 mm
111565s	L. 65 mm
111570s	L. 70 mm
111575s	L. 75 mm
111580s	L. 80 mm

Large QWIX® Positioning screws
dia. 7.5 mm



Reference	Length
121740s	L. 40 mm
121745s	L. 45 mm
121750s	L. 50 mm
121755s	L. 55 mm
121760s	L. 60 mm
121765s	L. 65 mm
121770s	L. 70 mm
121775s	L. 75 mm
121780s	L. 80 mm
121785s	L. 85 mm
121790s	L. 90 mm
121795s	L. 95 mm
121800s	L. 100 mm
121805s	L. 105 mm
121810s	L. 110 mm
121815s	L. 115 mm
121820s	L. 120 mm




Large QWIX® Positioning screws
dia. 5.5 mm

Reference	Length
121530s	L. 30 mm
121535s	L. 35 mm
121540s	L. 40 mm
121545s	L. 45 mm
121550s	L. 50 mm
121555s	L. 55 mm
121560s	L. 60 mm
121565s	L. 65 mm
121570s	L. 70 mm
121575s	L. 75 mm
121580s	L. 80 mm

Large QWIX® Instruments

Reference	Length
119150	Measurer 200 and 250 mm
119151	Internal protection sleeve dia. 1.8 mm
119152	Drilling guide dia. 4.5 mm
119153	External protection sleeve dia. 7.3 mm
119154	Cannulated drill dia. 1.8/4.2 x L. 180 mm
119155	Cannulated 2 in 1 drill dia. 4.2/5.8 x L. 113 mm
119156	Hexa screwdriver Tip 3 mm
119171	Internal protection sleeve dia. 2.7 mm
119172	Drilling guide dia. 5.8 mm
119173	External protection sleeve dia. 9.3 mm
119174	Cannulated drill dia. 2.7/5.5 x L. 210 mm
119175	Cannulated 2 in 1 drill dia. 5.5/7.0 x L. 128 mm
119176	Hexa screwdriver Tip 4 mm
119177	T handle quick coupling
115516	K-Wire dia. 1.6 x L. 200 mm
529061	K-Wire dia. 2.5 x L. 250 mm



Sterilization container

Reference	Description
121950	LARGE QWIX® sterilization container including:
121951	Base
996200	Lid
119170	Cylinder
119909	Silicone wedge

I.CO.S.® Ideal Compression Screw

Cannulated Titanium 4.0 mm and 6.5 mm diameter compression screw are dedicated to fixation of osteotomies or for fractures stabilization.

- K-Wire guidance for drilling, tapping and screw setting (dia. 1.6 or 2.5 mm).
- Adjustable and controlled static compression.
- Dynamic compression by additional head screwing.
- Intra osseous implant.
- Large range of sizes :
 - > I.CO.S.® dia. 4.0 mm:
 - » L. 26 to 40 mm by 2 mm
 - » L. 40 to 60 mm by 5 mm
 - > I.CO.S.® dia. 6.5 mm:
 - » L. 40 to 90 mm by 5 mm
- Material: Titanium alloy TA6V.
- Cortico-spongyous thread.
- Long non threaded shaft (lag screws) allowing compression.



I.CO.S.® dia. 4.0 mm

Reference	Length
105426	L. 26 mm
105428	L. 28 mm
105430	L. 30 mm
105432	L. 32 mm
105434	L. 34 mm
105436	L. 36 mm
105438	L. 38 mm
105440	L. 40 mm
105445	L. 45 mm
105450	L. 50 mm
105455	L. 55 mm
105460	L. 60 mm



I.CO.S.® dia. 6.5 mm


Reference	Length
105640	L. 40 mm
105645	L. 45 mm
105650	L. 50 mm
105655	L. 55 mm
105660	L. 60 mm
105665	L. 65 mm
105670	L. 70 mm
105675	L. 75 mm
105680	L. 80 mm
105685	L. 85 mm
105690	L. 90 mm

I.CO.S.® dia. 4.0 mm Instruments

Reference	Description
119501	Protection Sleeve dia. 6 mm
119506	Drilling Guide dia. 1.8 mm
119510	Measurer
119516	Long cannulated drill dia. 2.6 x L. 140 mm
119522	Short cannulated Drill dia. 5.0 x L. 110 mm
129710	Quick-coupling
119525	Tap
119531	Internal screwdriver dia. 4.0 x L. 100 mm
119532	Screwdriver tip dia. 4.0 x L. 100 mm
119536	External screwdriver dia. 6.4 x L. 170 mm

I.CO.S.® 4.0 sterilization Container

Reference	Description
119974	I.CO.S.® 4.0 sterilization container including:
119956	Base
119951	Screw rack
119949	Metallic lid for screws (Optional)
996200	Lid
119909	Silicone wedge



I.CO.S.® dia. 6.5 mm Instruments

Reference	Description
119541	Drilling Guide dia. 6.5 x L. 100 mm
119551	Internal protection sleeve dia. 6.5 mm
119552	External protection sleeve dia. 6.5 mm
119576	Long cannulated drill dia. 4.6 x L. 180 mm
119546	Short cannulated drill dia. 6.5 x L. 118 mm
119561	Internal screwdriver dia. 6.5 x L. 170 mm
119562	Screwdriver tip dia. 6.5 x L. 110 mm
119566	External screwdriver dia. 10.1 x L. 170 mm
119555	Tap
119510	Measurer
129710	Quick coupling
115225	K-Wire dia. 2.5 x L. 200 mm

I.CO.S.® 6.5 sterilization Container

Reference	Description
119976	I.CO.S.® 6.5 sterilization Container including:
119958	Base
119953	Screws rack
119949	Metallic lid for screws (Optional)
119957	Lid
119909	Silicone wedge



I.CO.S.® Removal Set

This removal set contains instrumentation required to remove I.CO.S.® screws dia. 4.0 and 6.5 mm

Removal set Instruments (for all diameters)

Reference	Description
Instruments for I.CO.S.® dia. 4.0 mm	
115116	K-Wire dia. 1.6 x L. 150 mm
119531	Internal screwdriver dia. 4.0 x L. 100 mm
119580	Trephine dia. 4.0 x L. 70 mm
119582	Internal Extractor dia. 4.0 x L. 60 mm (single use)
Instruments for I.CO.S.® dia. 6.5 mm	
Reference	Description
115225	K-Wire dia. 2.5 x L. 200 mm
119561	Internal screwdriver dia. 6.5 x L. 170 mm
119590	Trephine dia. 6.5 x L. 80 mm
119592	Internal Extractor dia. 6.5 x L. 70 mm (single use)

Sterilization container

Reference	Description
119975	I.CO.S.® Removal Sterilization container including
119959	Lid
119955	Base
119977	K-Wire cylinder



KALIX® II Flat Foot Implant

Expansion endorthesis for flat foot correction in child or adult.

Progressive expansion.

- Material:
 - > Titanium alloy TA6V
 - > Polyethylene (UHMWPE)
- Color coded.
- Sterily packaged.

Reference	Ext. diameter
141009	dia. 9 mm
141010	dia. 10 mm
141011	dia. 11 mm
141012	dia. 12 mm
141013	dia. 13 mm
141014	dia. 14 mm
141015	dia. 15 mm
141017	dia. 17 mm



KALIX® II Flat Foot Implant Instruments

Reference	Description
119825	Viladot Lever
119850	Impactor
119831	Trial implant holder
119859	KALIX® II Trial implant dia. 9.0 mm
119860	KALIX® II Trial implant dia. 10.0 mm
119861	KALIX® II Trial implant dia. 11 mm
119862	KALIX® II Trial implant dia. 12 mm
119863	KALIX® II Trial implant dia. 13 mm
119864	KALIX® II Trial implant dia. 14 mm
119865	KALIX® II Trial implant dia. 15 mm
119867	KALIX® II Trial implant dia. 17 mm



KALIX® II sterilization container

Reference	Description
119880	KALIX® II sterilization container including
119882	Lid
119881	Base

MBA™ Titanium Subtalar Implant

The Subtalar MBA™ System is indicated for use in treatment of the hyperpronated foot and stabilization of the subtalar joint. It is designed to block the posterior and inferior displacement of the talus, thus allowing normal subtalar joint motion while blocking excessive pronation and the resulting sequela.

- Severely pronated foot.
- Walking intemperance.
- Calcaneal stance position greater than 5°.
- Manually correctable deformities.
- Mid-tarsal breech (arch pain).
- Forefoot varus greater than 10°.



Removal set Instruments (for all diameters)

Reference	Description
Instruments tray	
05-0016	Guide pin insert holder
05-0206	Trial 6 mm
05-0208	Trial 8 mm
05-0209	Trial 9 mm
05-0210	Trial 10 mm
05-0212	Trial 12 mm
05-1002	Yellow probe
05-1006	6 mm Sizer
05-1008	8 mm Sizer
05-1009	9 mm Sizer
05-1010	10 mm Sizer
05-1012	12 mm Sizer
05-1020	Insertion Device
05-0017	Guide pin 2 mm
05-5001	Sterilization tray
Implants	
05-0106	MBA™ 6 mm implant
05-0108	MBA™ 8 mm implant
05-0109	MBA™ 9 mm implant
05-0110	MBA™ 10 mm implant
05-0112	MBA™ 12 mm implant



ADVANSYS® Dorsal Lisfranc Plate

The Dorsal Lisfranc Plates are intended for fractures, fusions, osteotomies and replantations of small bones at the tarso-metatarsal joints (Lisfranc joints).

- 3 sizes of implants:
3 right & 3 left plates (small, medium and large).
- 3.5 mm diameter screws available in lengths from 8 mm to 34 mm in 2 mm increments.
- 3.5 mm lock-screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.



ADVANSYS® Dorsal Lisfranc Plates

Reference	Description
181021S	ADVANSYS® DLP small-right
181022S	ADVANSYS® DLP medium-right
181023S	ADVANSYS® DLP large-right
181031S	ADVANSYS® DLP small-left
181032S	ADVANSYS® DLP medium-left
181033S	ADVANSYS® DLP large-left



Surfix® Locking screw Sterile Titanium dia.3.5 mm + lock-screw

Reference	Description
285310S	L. 10 mm
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
185300S	Lock-screw dia.3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

ADVANSYS® Medial Lisfranc Plate

The Medial Lisfranc Plates are intended to be used for bone fixation such as:

- Lisfranc arthrodesis.
- Mono or bi-cortical osteotomies or fractures near the 1st metatarso-cuneiform joint.
- 2 sizes of implants:
2 right & 2 left plates (small and large sizes).
- 3.5 mm diameter screws available in lengths from 8 mm to 34 mm in 2 mm increments.
- 3.5 mm lock-screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.
- Sterily packaged.



ADVANSYS® Medial Lisfranc Plate

Reference	Description
181041S	ADVANSYS® MLP small-right
181042S	ADVANSYS® MLP large-right
181051S	ADVANSYS® MLP small-left
181052S	ADVANSYS® MLP large-left



Surfix® Locking screw Sterile Titanium dia. 3.5 mm + lock-screw

Reference	Description
285310S	L. 10 mm
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
185300S	Lock-screw dia.3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

ADVANSYS® Medial Lisfranc Plate Instruments

Reference	Description
115100	K-Wire dia. 1.0 x L. 100 mm
219335	Depth gauge hole dia. 3.5 mm
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219535	AO drill dia. 2.7 x L. 125 mm
219635	Drilling guide dia. 2.7 mm
219735	Bending forceps L. 171 hole - dia 3.5 mm
219835	Screwdriver hexa 2.0 x L. 180 mm
169041	MLP trial implant right-small
169042	MLP trial implant right-large
169051	MLP trial implant left-small
169052	MLP trial implant left-large

ADVANSYS® Midfoot sterilization container


Reference	Description
169102	ADVANSYS® Midfoot sterilization container includes:
169112	Base
996200	Lid
169106	Cylinder
169107	DLP/MLP Trial implant module
119909	Silicone wedge

DLP trial plates

Reference	Description
169021	Dorsal trial plate, right, small
169022	Dorsal trial plate, right, medium
169023	Dorsal trial plate, right, large
169031	Dorsal trial plate, left, small
169032	Dorsal trial plate, left, medium
169033	Dorsal trial plate, left, large

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm



UNI-CP™ Compression plate

The UNI-CP™ compression plate is indicated for fixation of bone fractures or for bone reconstruction:

- 4 plate configurations:
 - > 2 hole design (17, 20, 25, and 30 mm interaxis)
 - > 4 hole design (20, 25, and 30 mm interaxis)
 - > 4 hole T-shape design (20 mm interaxis)
 - > 4 hole U-shape design (17x20, 19x22.5, 21x25 mm interaxis)
- Range of screw lengths: from 12 mm to 40 mm in 2 mm increments.
- Surfix® Locking Technology.
- Material
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.

UNI-CP™ Compression Plate

Reference	Description
330430S	UNI-CP™ 4 holes - 30 mm
330425S	UNI-CP™ 4 holes - 25 mm
330420S	UNI-CP™ 4 holes - 20 mm
330230S	UNI-CP™ 2 holes - 30 mm
330225S	UNI-CP™ 2 holes - 25 mm
330220S	UNI-CP™ 2 holes - 20 mm
330217S	UNI-CP™ 2 holes - 17 mm
330030S	UNI-CP™ T shaped plate
330021S	4 holes - 17x20 mm U shaped plate
330023S	4 holes - 19x22.5 mm U shaped plate
330025S	4 holes - 21x25 mm U shaped plate



Surfix® Sterile Titanium screw
dia. 3.5 mm + lock-screw

Reference	Description
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
285336S	L. 36 mm
285338S	L. 38 mm
285340S	L. 40 mm
185300S	Surfix® lock screw dia. 3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm



UNI-CP™ plate 4 holes


UNI-CP™ plate 2 holes

CP T-shape plate

UNI-CP™ U-shape plate

UNI-CP™ Compression plate Instruments

Reference	Description
219335	Depth gauge hole dia. 3.5 mm
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219535	AO drill dia. 2.7 x L. 125 mm
219635	Drilling guide dia. 2.7 mm
219735	Bending forceps L. 171 hole - dia 3.5 mm
219835	Screwdriver hexa 2.0 x L. 180 mm
339001	Compression forceps
339004	Trial implant for T and U plates
339005	Trial implant for 2 and 4 holes plates
339003	Holder
115101	K-Wire dia. 1.0 x L. 100 mm



Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

UNI-CP™ sterilization container

Reference	Description
339950	UNI-CP™ sterilization container including:
996200	Lid
339920	Base
169106	Cylinder
119909	Silicone wedge

HINTEGRA® Total Ankle Prosthesis – Primary Surgery

The HINTEGRA® Total Ankle Prosthesis is a three components design anatomical prosthesis.

- Materials:
 - > Talar & tibial components:
 - Cobalt Chromium alloy (CoCr), with double coating:
 - » Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
 - > PE inlay: UHMWPE.
 - > Fixation screws (optional): Stainless steel M30NW.



HINTEGRA® Total Ankle Prosthesis – Tibial Component – Primary Surgery

Minimum contact stress to the bone.

- Materials:
 - > Cobalt Chromium alloy (CoCr) with double coating:
 - Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- 7 sizes right & 7 sizes left.

Reference	Description
HINTEGRA® - Tibial Component - Right	
301200	HINTEGRA® tibial component-right - size 0
301201	HINTEGRA® tibial component-right - size 1
301202	HINTEGRA® tibial component-right - size 2
301203	HINTEGRA® tibial component-right - size 3
301204	HINTEGRA® tibial component-right - size 4
301205	HINTEGRA® tibial component-right - size 5
301206	HINTEGRA® tibial component-right - size 6
HINTEGRA® - Tibial Component - Left	
302200	HINTEGRA® tibial component-left - size 0
302201	HINTEGRA® tibial component-left - size 1
302202	HINTEGRA® tibial component-left - size 2
302203	HINTEGRA® tibial component-left - size 3
302204	HINTEGRA® tibial component-left - size 4
302205	HINTEGRA® tibial component-left - size 5
302206	HINTEGRA® tibial component-left - size 6



HINTEGRA® Total Ankle Prosthesis – PE Inlay – Primary Surgery

- Ultra High Molecular Weight Polyethylene (UHMWPE).
- X-rays markers (Titanium TA6V rods).
- Same size as talar component: 6 sizes with 4 thicknesses (5, 6, 7, 9 mm).

HINTEGRA® - PE inlay - size 0

300005	HINTEGRA® PE inlay - size 0 - thickness 5 mm
300006	HINTEGRA® PE inlay - size 0 - thickness 6 mm
300007	HINTEGRA® PE inlay - size 0 - thickness 7 mm
300009	HINTEGRA® PE inlay - size 0 - thickness 9 mm

HINTEGRA® - PE inlay - size 1

300105	HINTEGRA® PE inlay - size 1 - thickness 5 mm
300106	HINTEGRA® PE inlay - size 1 - thickness 6 mm
300107	HINTEGRA® PE inlay - size 1 - thickness 7 mm
300109	HINTEGRA® PE inlay - size 1 - thickness 9 mm

HINTEGRA® - PE inlay - size 2

300205	HINTEGRA® PE inlay - size 2 - thickness 5 mm
300206	HINTEGRA® PE inlay - size 2 - thickness 6 mm
300207	HINTEGRA® PE inlay - size 2 - thickness 7 mm
300209	HINTEGRA® PE inlay - size 2 - thickness 9 mm

HINTEGRA® - PE inlay - size 3

300305	HINTEGRA® PE inlay - size 3 - thickness 5 mm
300306	HINTEGRA® PE inlay - size 3 - thickness 6 mm
300307	HINTEGRA® PE inlay - size 3 - thickness 7 mm
300309	HINTEGRA® PE inlay - size 3 - thickness 9 mm

HINTEGRA® - PE inlay - size 4

300405	HINTEGRA® PE inlay - size 4 - thickness 5 mm
300406	HINTEGRA® PE inlay - size 4 - thickness 6 mm
300407	HINTEGRA® PE inlay - size 4 - thickness 7 mm
300409	HINTEGRA® PE inlay - size 4 - thickness 9 mm

HINTEGRA® - PE inlay - size 5

300505	HINTEGRA® PE inlay - size 5 - thickness 5 mm
300506	HINTEGRA® PE inlay - size 5 - thickness 6 mm
300507	HINTEGRA® PE inlay - size 5 - thickness 7 mm
300509	HINTEGRA® PE inlay - size 5 - thickness 9 mm

HINTEGRA® - PE inlay - size 6

300605	HINTEGRA® PE inlay - size 6 - thickness 5 mm
300606	HINTEGRA® PE inlay - size 6 - thickness 6 mm
300607	HINTEGRA® PE inlay - size 6 - thickness 7 mm
300609	HINTEGRA® PE inlay - size 6 - thickness 9 mm



HINTEGRA® Total Ankle Prosthesis – Talar component with Pegs – Primary Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- Cobalt Chromium alloy (CoCr).
- Anatomical conical shape.
- 7 sizes right and 7 sizes left components.



Reference Description

HINTEGRA® - Talar Component with Pegs - Right

301110	HINTEGRA® talar component with pegs-right - size 0
301111	HINTEGRA® talar component with pegs-right - size 1
301112	HINTEGRA® talar component with pegs-right - size 2
301113	HINTEGRA® talar component with pegs-right - size 3
301114	HINTEGRA® talar component with pegs-right - size 4
301115	HINTEGRA® talar component with pegs-right - size 5
301116	HINTEGRA® talar component with pegs-right - size 6

HINTEGRA® - Talar Component with Pegs - Left

302110	HINTEGRA® talar component with pegs-left - size 0
302111	HINTEGRA® talar component with pegs-left - size 1
302112	HINTEGRA® talar component with pegs-left - size 2
302113	HINTEGRA® talar component with pegs-left - size 3
302114	HINTEGRA® talar component with pegs-left - size 4
302115	HINTEGRA® talar component with pegs-left - size 5
302116	HINTEGRA® talar component with pegs-left - size 6

HINTEGRA® Total Ankle Prosthesis – (Optional) Tibial cortical screw – Primary Surgery


- Stainless steel M30NW.
- Optionally to be put only with tibial components.



Reference Description

303328S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 28 mm - tibia - size 0
303330S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 30 mm - size 1
303334S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 34 mm - size 2
303338S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 38 mm - size 3
303340S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 40 mm - size 4
303342S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 42 mm - size 5 & 6


HINTEGRA® Total Ankle Prosthesis Associated Instruments for Primary Surgery



Generic base
(upper tray)



Trial tray



Talar base

HINTEGRA® Total Ankle Prosthesis Associated Instruments for Primary Surgery

Instruments for LOWER TRAY

Talar Base	
309300	Drilling guide talus - Size 0
309301	Drilling guide talus - Size 1
309302	Drilling guide talus - Size 2
309303	Drilling guide talus - Size 3
309304	Drilling guide talus - Size 4
309305	Drilling guide talus - Size 5
309306	Drilling guide talus - Size 6
309360	Tal. cut. guide-Right - Size 0
309361	Tal. cut. guide-Right - Size 1
309362	Tal. cut. guide-Right - Size 2
309363	Tal. cut. guide-Right - Size 3
309364	Tal. cut. guide-Right - Size 4
309365	Tal. cut. guide-Right - Size 5
309366	Tal. cut. guide-Right - Size 6
309370	Tal. cut. guide-Left - Size 0
309371	Tal. cut. guide-Left - Size 1
309372	Tal. cut. guide-Left - Size 2
309373	Tal. cut. guide-Left - Size 3
309374	Tal. cut. guide-Left - Size 4
309375	Tal. cut. guide-Left - Size 5
309376	Tal. cut. guide-Left - Size 6
309601	Screwdriver tibial screws
309608	Spacer h 12 mm
119609	Talar impactor
119641	Cutting guide tibial revision
309200	Reamer talar dia. 6 mm
309309	Peg drill talus - dia. 4.5 mm
309697	Holder trial talar implant
309699	Impactor trial talar implant

Trial Tray

HINTEGRA® Tibial Trial-Right	
119646	Tibial Trial-Right - Size 0
119690	Tibial Trial-Right - Size 1
119691	Tibial Trial-Right - Size 2
119692	Tibial Trial-Right - Size 3
119693	Tibial Trial-Right - Size 4
119694	Tibial Trial-Right - Size 5
119778	Tibial Trial-Right - Size 6
HINTEGRA® Tibial Trial-Left	
119647	Tibial Trial-Left - Size 0
119695	Tibial Trial-Left - Size 1
119696	Tibial Trial-Left - Size 2
119697	Tibial Trial-Left - Size 3
119698	Tibial Trial-Left - Size 4
119699	Tibial Trial-Left - Size 5
119779	Tibial Trial-Left - Size 6
HINTEGRA® Talar Trial-Right	
309680	Talar Trial-Right - Size 0
309681	Talar Trial-Right - Size 1
309682	Talar Trial-Right - Size 2
309683	Talar Trial-Right - Size 3
309684	Talar Trial-Right - Size 4
309685	Talar Trial-Right - Size 5
309686	Talar Trial-Right - Size 6
HINTEGRA® Talar Trial-Left	
309690	Talar Trial-Left - Size 0
309691	Talar Trial-Left - Size 1
309692	Talar Trial-Left - Size 2
309693	Talar Trial-Left - Size 3
309694	Talar Trial-Left - Size 4
309695	Talar Trial-Left - Size 5
309696	Talar Trial-Left - Size 6
HINTEGRA® Trial inlay	
119665	Trial inlay - thickness 5mm
119666	Trial inlay - thickness 6mm
119667	Trial inlay - thickness 7mm
119669	Trial inlay - thickness 9 mm

Sterilization container

309900	HINTEGRA® Container (Upper + Lower) including:		
309902	Stainless steel base upper, including:	309904	Stainless steel base lower, including:
309920	Generic base	309960	Talar base
309940	Pins insert	309980	Thermo insert trial implants
309942	Lid pins insert	309970	Lid trial implants
		996350	Lid
		119909	Silicone wedge

Instruments for UPPER TRAY

309625	Tibial v positioning
309630	V translation block
119664	Distractor Hintermann® big model
119662	Tibial - Talar Holder
119751	Tibial impactor
309604	Pins L. 30 mm
309605	Pins L. 70 mm
309606	Pins L. 45 mm
309607	Depth gauge tibial
159127	Drilling guide dia. 2.7mm
309615	Alignment rod tibial
309620	Connector v - tibial rod
309622	Saw blade Aesculap attachment
309623	Saw blade AO Synthes attachment
309624	Saw blade Stryker 6 attachment
309626	Saw blade Stryker attachment
309627	Saw blade Conmed Linvatec attachment
309635	Tib. cut. guide large
309636	Tib. cut. guide medium / standard
309637	Tib. cut. guide small
309645	Screwdriver Hex dia. 3.5 mm
309652	Imp / Extract pins
309656	Cutting block-talar cut / Medium - Standard
309657	Cutting block-talar cut - Small
309661	Holder PE inlays
309753	Cutting block proximal - tibial
309773	Cutting block distal - tibial
997034	Curved chisel 20mm r40mm
997035	Straight chisel 10mm
275560	Parallel pliers 7-1/4"
119611	Drill dia. 2.5 mm x L. 150 mm
115116	K-Wire dia. 1.6 x L. 150 mm
115225	K-Wire dia. 2.5 x L. 200 mm

HINTEGRA® Total Ankle Prosthesis – Flat Cut – Revision Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).



HINTEGRA® Total Ankle Prosthesis Flat Cut – Tibial Component – Revision Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- Two thicknesses: 8 and 12 mm.



Reference Description

HINTEGRA® - revision tibial component - right - 8 mm

301222	HINTEGRA® revision tibial component-right-8 mm - size 2
301223	HINTEGRA® revision tibial component-right-8 mm - size 3
301224	HINTEGRA® revision tibial component-right-8 mm - size 4
301225	HINTEGRA® revision tibial component-right-8 mm - size 5

HINTEGRA® - revision tibial component - left - 8 mm

302222	HINTEGRA® revision tibial component-left-8 mm - size 2
302223	HINTEGRA® revision tibial component-left-8 mm - size 3
302224	HINTEGRA® revision tibial component-left-8 mm - size 4
302225	HINTEGRA® revision tibial component-left-8 mm - size 5

Reference Description

HINTEGRA® - revision tibial component - right - 12 mm

301232	HINTEGRA® revision tibial component-right-12 mm - size 2
301233	HINTEGRA® revision tibial component-right-12 mm - size 3
301234	HINTEGRA® revision tibial component-right-12 mm - size 4
301235	HINTEGRA® revision tibial component-right-12 mm - size 5

HINTEGRA® - revision tibial component - left - 12 mm

302232	HINTEGRA® revision tibial component-left-12 mm - size 2
302233	HINTEGRA® revision tibial component-left-12 mm - size 3
302234	HINTEGRA® revision tibial component-left-12 mm - size 4
302235	HINTEGRA® revision tibial component-left-12 mm - size 5

HINTEGRA® Total Ankle Prosthesis Flat Cut – Talar Component – Revision Surgery / Flat Talus

- For flat talus and revision surgeries.
- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).



Reference Description

HINTEGRA® - flat cut talar component - right

301121	HINTEGRA® revision talar component-right - size 1
301122	HINTEGRA® revision talar component-right - size 2
301123	HINTEGRA® revision talar component-right - size 3
301124	HINTEGRA® revision talar component-right - size 4
301125	HINTEGRA® revision talar component-right - size 5

Reference Description

HINTEGRA® - flat cut talar component - left

302121	HINTEGRA® revision talar component-left - size 1
302122	HINTEGRA® revision talar component-left - size 2
302123	HINTEGRA® revision talar component-left - size 3
302124	HINTEGRA® revision talar component-left - size 4
302125	HINTEGRA® revision talar component-left - size 5

HINTEGRA® Total Ankle Prosthesis – Cortical screw

- Stainless steel: M30NW.
- Sterile.



Reference Description

303312S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 12 mm - talus
303316S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 16 mm - talus
303320S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 20 mm - talus
303334S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 34 mm - tibia - size 2
303338S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 38 mm - tibia - size 3
303340S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 40 mm - tibia - size 4
303342S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 42 mm - tibia - size 5 & 6

HINTEGRA® Total Ankle Prosthesis – PE Inlay – Revision Surgery (See page 19 for product references.)

Size of the PE inlay has to be the same as the talar component.




HINTEGRA® Total Ankle Prosthesis – Flat Cut Instruments for Revision Surgery

Instruments

119790	FC Talar Trial for Drilling-Right - Size 1
119791	FC Talar Trial for Drilling-Right - Size 2
119792	FC Talar Trial for Drilling-Right - Size 3
119793	FC Talar Trial for Drilling-Right - Size 4
119794	FC Talar Trial for Drilling-Right - Size 5
119657	Drill dia. 4.5 mm x L. 22 mm
119659	FC Talar trial for drilling Impactor
119795	FC Talar Trial for Drilling-Left - Size 1
119796	FC Talar Trial for Drilling-Left - Size 2
119797	FC Talar Trial for Drilling-Left - Size 3
119798	FC Talar Trial for Drilling-Left - Size 4
119799	FC Talar Trial for Drilling-Left - Size 5
119631	FC Spacer 4 mm for revision Tibial Component 8 mm
119632	FC Spacer 8 mm for revision Tibial Component 12 mm
119757	FC spacer
119760	FC Talar Cutting Guide-Right - Size 1
119761	FC Talar Cutting Guide-Right - Size 2
119762	FC Talar Cutting Guide-Right - Size 3
119763	FC Talar Cutting Guide-Right - Size 4
119764	FC Talar Cutting Guide-Right - Size 5
119765	FC Talar Cutting Guide-Left - Size 1
119766	FC Talar Cutting Guide-Left - Size 2
119767	FC Talar Cutting Guide-Left - Size 3
119768	FC Talar Cutting Guide-Left - Size 4
119769	FC Talar Cutting Guide-Left - Size 5
309 708	Tibial depth gauge for tibial component 8 mm
309712	Tibial depth gauge for tibial component 12 mm

Sterilization container

119965	HINTEGRA® FC revision container includes:
119978	Container pocket
119966	HINTEGRA® FC container base
119967	Lid for HINTEGRA® FC container base



HINTEGRA® Sensitive Total Ankle Prosthesis – Primary Surgery

The HINTEGRA® Sensitive Total Ankle Prosthesis is also a three component design specially made for allergic patients.

- Materials:
 - > Talar & Tibial components:
 - » Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
 - > PE inlay:
 - » Polyethylene (UHMWPE).



HINTEGRA® Sensitive Total Ankle Prosthesis – Tibial Component – Primary Surgery

- Anatomical shape for minimal contact stress to the bone.
- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 6 right and 6 left components.

Reference	Description
HINTEGRA® sensitive - tibial component - right	
351251	HINTEGRA® sensitive - tibia-right - size 1
351252	HINTEGRA® sensitive - tibia-right - size 2
351253	HINTEGRA® sensitive - tibia-right - size 3
351254	HINTEGRA® sensitive - tibia-right - size 4
351255	HINTEGRA® sensitive - tibia-right - size 5
351256	HINTEGRA® sensitive - tibia-right - size 6
HINTEGRA® sensitive - tibial component - left	
352251	HINTEGRA® sensitive - tibia-left - size 1
352252	HINTEGRA® sensitive - tibia-left - size 2
352253	HINTEGRA® sensitive - tibia-left - size 3
352254	HINTEGRA® sensitive - tibia-left - size 4
352255	HINTEGRA® sensitive - tibia-left - size 5
352256	HINTEGRA® sensitive - tibia-left - size 6



HINTEGRA® Sensitive Total Ankle Prosthesis – Talar component with Pegs – Primary Surgery

- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- Anatomical conical shape.
- 5 sizes and 5 left components.

Reference	Description
HINTEGRA® - talar component with pegs - right	
351152	HINTEGRA® sensitive - talus -right - size 2
351153	HINTEGRA® sensitive - talus -right - size 3
351154	HINTEGRA® sensitive - talus -right - size 4
351155	HINTEGRA® sensitive - talus -right - size 5
351156	HINTEGRA® sensitive - talus -right - size 6
HINTEGRA® - talar component with pegs - left	
352152	HINTEGRA® sensitive - talus-left - size 2
352153	HINTEGRA® sensitive - talus-left - size 3
352154	HINTEGRA® sensitive - talus-left - size 4
352155	HINTEGRA® sensitive - talus-left - size 5
352156	HINTEGRA® sensitive - talus-left - size 6



HINTEGRA® Sensitive Total Ankle Prosthesis – PE Inlay (See page 45 for product references.)

- HINTEGRA® PE inlays are “sensitive compatible” since metallic rods are made of titanium Alloy (TA6V).
- Size of the PE inlay has to be the same as the talar component.



HINTEGRA® Sensitive Total Ankle Prosthesis – (Optional) Tibial cortical screw

- Titanium alloy (TA6V).
- Optionally to be put only with tibial components.

Reference	Description
303530S	Cortical Screw Ti - Dia. 3.5 mm x L. 30 mm - Tibia - Size 1
303534S	Cortical Screw Ti - Dia. 3.5 mm x L. 34 mm - Tibia - Size 2
303538S	Cortical Screw Ti - Dia. 3.5 mm x L. 38 mm - Tibia - Size 3
303540S	Cortical Screw Ti - Dia. 3.5 mm x L. 40 mm - Tibia - Size 4
303542S	Cortical Screw Ti - Dia. 3.5 mm x L. 42 mm - Tibia - Size 5 & 6



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Revision Surgery

The HINTEGRA® Sensitive Total Ankle Prosthesis is also a three component design specially made for allergic patients.

- Materials:
 - > Talar & Tibial components:
 - » Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
 - > PE inlay:
 - » Polyethylene (UHMWPE).



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Tibial Component – Revision Surgery

- Anatomical shape for minimal contact stress to the bone.
- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 4 right and 4 left components.

Reference	Description	Reference	Description
HINTEGRA® - Revision Tibial Component - Right - 8 mm		HINTEGRA® - Revision Tibial Component - Right - 12 mm	
351522	HINTEGRA® sensitive FC - tibia-right 8 mm - size 2	351532	HINTEGRA® sensitive FC - tibia-right 12 mm - size 2
351523	HINTEGRA® sensitive FC - tibia-right 8 mm - size 3	351533	HINTEGRA® sensitive FC - tibia-right 12 mm - size 3
351524	HINTEGRA® sensitive FC - tibia-right 8 mm - size 4	351534	HINTEGRA® sensitive FC - tibia-right 12 mm - size 4
351525	HINTEGRA® sensitive FC - tibia-right 8 mm - size 5	351535	HINTEGRA® sensitive FC - tibia-right 12 mm - size 5
HINTEGRA® - Revision Tibial Component - Left - 8 mm		HINTEGRA® - Revision Tibial Component - Left - 12 mm	
352522	HINTEGRA® sensitive FC - tibia-left 8 mm - size 2	352532	HINTEGRA® sensitive FC - tibia-left 12 mm - size 2
352523	HINTEGRA® sensitive FC - tibia-left 8 mm - size 3	352533	HINTEGRA® sensitive FC - tibia-left 12 mm - size 3
352524	HINTEGRA® sensitive FC - tibia-left 8 mm - size 4	352534	HINTEGRA® sensitive FC - tibia-left 12 mm - size 4
352525	HINTEGRA® sensitive FC - tibia-left 8 mm - size 5	352535	HINTEGRA® sensitive FC - tibia-left 12 mm - size 5



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut Talar component with Pegs – Revision Surgery / Flat talus

- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 4 right and 4 left components.

Reference	Description
HINTEGRA® - Talar Component with Pegs - Right	
351512	HINTEGRA® sensitive FC - talus-right - size 2
351513	HINTEGRA® sensitive FC - talus-right - size 3
351514	HINTEGRA® sensitive FC - talus-right - size 4
351515	HINTEGRA® sensitive FC - talus-right - size 5
HINTEGRA® - Talar Component with Pegs - Left	
352512	HINTEGRA® sensitive FC - talus-left - size 2
352513	HINTEGRA® sensitive FC - talus-left - size 3
352514	HINTEGRA® sensitive FC - talus-left - size 4
352515	HINTEGRA® sensitive FC - talus-left - size 5



HINTEGRA® Sensitive Total Ankle Prosthesis – PE Inlay (See page 45 for product references.)

- HINTEGRA® PE inlays are “sensitive compatible” since metallic rods are made of titanium Alloy (TA6V).
- Size of the PE inlay has to be the same as the talar component.



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Cortical screw

Titanium alloy (TA6V).

Reference	Description
303512S	Cortical Screw Ti - dia.3.5 x L. 12 mm - Talus
303516S	Cortical Screw Ti - dia.3.5 x L. 16 mm - Talus
303520S	Cortical Screw Ti - dia.3.5 x L. 20 mm - Talus
303534S	Cortical Screw Ti - dia. 3.5 x L. 34 mm - Tibia - Size 2
303538S	Cortical Screw Ti - dia. 3.5 x L. 38 mm - Tibia - Size 3
303540S	Cortical Screw Ti - dia. 3.5 x L. 40 mm - Tibia - Size 4
303542S	Cortical Screw Ti - dia. 3.5 x L. 42 mm - Tibia - Size 5



PANTA® / PANTA® XL Arthrodesis Nail

The PANTA® Nail is intended for use in tibiototalcalcaneal arthrodesis and treatment of trauma of the hindfoot and distal tibia.

- Cannulated Titanium Nail.
- Radio-transparent instrumentation system to apply compression.
- Multi-axial fixation:
 - > Oblique antero-posterior calcaneal screws.
 - > Medio-lateral tibial screws.
- 2 kind of screws:
 - > Oblique antero-posterior calcaneal screws.
 - » fully threaded
 - » partially threaded
- 5 mm diameter Titanium screw.
- Screw lengths: 20 to 110 mm.
- Material of the nail, the screws and the end cap: Titanium alloy TA6V.

PANTA® Nails (sterile)

Reference Dia. - length

500050	Dia. 10 mm - L. 150 mm
500080	Dia. 10 mm - L. 180 mm
500150	Dia. 11 mm - L. 150 mm
500180	Dia. 11 mm - L. 180 mm
500250	Dia. 12 mm - L. 150 mm
500280	Dia. 12 mm - L. 180 mm
500350	Dia. 13 mm - L. 150 mm
500380	Dia. 13 mm - L. 180 mm

PANTA® XL Nails (sterile)

To be used only with Partially threaded Screws

Reference Dia. - length

510111	Dia. 11mm - L 210 mm
510141	Dia. 11mm - L 240 mm
510211	Dia. 12 mm - L 210 mm
510241	Dia. 12 mm - L 240 mm
510311	Dia. 13mm - L 210 mm
510341	Dia. 13mm - L 240 mm

End cap (sterile)

Reference Description

500001	End cap
510005	Locking end cap



PANTA® Arthrodesis Nail – Fully threaded screws dia. 5 mm (not for use with PANTA® XL Arthrodesis Nail)

Reference	Length
501020	L. 20 mm
501022	L. 22.5 mm
501025	L. 25 mm
501027	L. 27.5 mm
501030	L. 30 mm
501035	L. 35 mm
501040	L. 40 mm
501045	L. 45 mm
501050	L. 50 mm
501055	L. 55 mm
501060	L. 60 mm
501065	L. 65 mm
501070	L. 70 mm
501075	L. 75 mm
501080	L. 80 mm
501085	L. 85 mm
501090	L. 90 mm
501095	L. 95 mm
501100	L. 100 mm
501105	L. 105 mm
501110	L. 110 mm


PANTA® / PANTA® XL Arthrodesis Nail – Partially threaded lag screws dia. 5 mm

Reference	Length
511020	L. 20 mm
511022	L. 22.5 mm
511025	L. 25 mm
511027	L. 27.5 mm
511030	L. 30 mm
511035	L. 35 mm
511040	L. 40 mm
511045	L. 45 mm
511050	L. 50 mm
511055	L. 55 mm
511060	L. 60 mm
511065	L. 65 mm
511070	L. 70 mm
511075	L. 75 mm
511080	L. 80 mm
511085	L. 85 mm
511090	L. 90 mm
511095	L. 95 mm
511100	L. 100 mm
511105	L. 105 mm
511110	L. 110 mm



PANTA® Arthrodesis Nail – Instruments

Reference	Description
519003	Drill dia. 4.3 mm for Calcaneum Screws
519004	Drill dia. 4.3 mm for Tibial Screws
519005	Drill dia. 5.0 mm – AO Attachment
519006	Long drill for calcaneus screws dia. 4.3 mm
519007	Cannulated Drill dia. 7.0 mm
519009	Cannulated Drill dia. 9.0 mm
519011	Reamer dia. 11 mm
519012	Reamer dia. 12 mm
519013	Reamer dia. 13 mm
519014	Reamer dia. 10.5 mm
519015	Reamer dia. 11.5 mm
519016	Reamer dia. 12.5 mm
519017	Reamer dia. 13.5 mm
519021	T Handle AO Attachment
519028	Internal Protecting Sleeve dia. 3.2 mm
519029	Central Internal Protecting Sleeve dia. 9.0 mm
519030	External Protecting Sleeve
519032	Pin dia. 3.2 mm x L. 400 mm
519040	Trocard Awl
519110	Support Device
519120	Nail Fixation Axis
519121	Toothed Wheel
519131	Threaded Axis
519130	Compression Device
519133	Teflon Ring
519135	Compression Wheel
519160	Depth Gauge
519175	Compression Rod
519179	Short Drill Guide dia. 4.3 mm L. 70 mm
519180	Long Drill Guide dia. 4.3 mm L. 120 mm
519181	Drill Guide dia. 5.0 mm
519184	Short Drill Guide dia. 7.0 mm L. 70 mm
519185	Drill Guide dia. 7.0 mm L. 120 mm
519290	Hexagonal screwdriver tip dia. 3.5 mm
519295	Hexagonal screwdriver dia. 3.5 mm
539015	Tap dia. 6.5 mm

PANTA® – Sterilization Container

519920	PANTA® sterilization container including:
519910	PANTA® base
519911	PANTA® Insert
519920	Lid
119909	Silicone wedge

PANTA® – Optional Instruments

519010	Reamer dia. 10 mm
519020	Quick Coupling AO
519002	Short drill for calcaneus screws used with medium drill guide dia. 4.3 mm
519008	Long drill for calcaneus screws used with medium drill guide dia. 4.3 mm
519034(S)	Guide wire dia. 3.2 mm, L. 600 mm (sterile upon request)
519178	Medium drill guide dia. 4.3 mm
519183	Medium soft tissue protector dia. 7 mm





PANTA® XL – Associated Instruments (optional)

539020	Conical reamer dia. 10.5 mm
539021	Conical reamer dia. 11.5 mm
539022	Conical reamer dia. 12.5 mm
539023	Conical reamer dia. 13.5 mm
519912	Reamer holder



(S) means available in sterile and non sterile

PANTA[®] Arthrodesis Nail - Associated Instruments



PANTA® Removal Set

Reference	Description
519209	Sliding hammer
519210	Handle extractor
309645	Screwdriver Hex dia.3.5 mm
119552	External protection sleeve dia.6.5 mm
Sterilization container	
519950	PANTA® Removal sterilization container
119909	Silicone wedge
519951	Basis
996100	Lid



ADVANSYS® TTC Fusion Plate

The TTC Plates are intended for arthrodesis of the ankle joint and calcanea, osteotomies, fusions and replantations of small bones in the ankle and distal tibia.

- 4 sizes of implants
 - > 4 right & 4 left plates (in 6,7,8,9 holes configurations).
- Screw dia.4.5 mm for fixation of the tibia and talus – Length 14 mm to 65 mm.
- Screw dia.6.5 mm for fixation of the calcaneus – Length 20 mm to 65 mm.
- Lock-screws for 4.5 mm dia. screws.
- Lock-screws for 6.5 mm dia. screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.



Advansys® TTC Fusion Plate

Reference	Description
181001S	ADVANSYS® TTCP Right-6 holes
181002S	ADVANSYS® TTCP Right-8 holes
181003S	ADVANSYS® TTCP Right-7 holes
181004S	ADVANSYS® TTCP Right-9 holes
181011S	ADVANSYS® TTCP Left-6 holes
181012S	ADVANSYS® TTCP Left-8 holes
181013S	ADVANSYS® TTCP Left-7 holes
181014S	ADVANSYS® TTCP Left-9 holes



Surfix® Locking screw
Sterile Titanium dia. 4.5 mm + lock-screw

Reference	Length
285414s	L. 14 mm
285416s	L. 16 mm
285418s	L. 18 mm
285420s	L. 20 mm
285422s	L. 22 mm
285424s	L. 24 mm
285426s	L. 26 mm
285428s	L. 28 mm
285430s	L. 30 mm
285435s	L. 35 mm
285440s	L. 40 mm
285445s	L. 45 mm
285450s	L. 50 mm
285455s	L. 55 mm
285460s	L. 60 mm
285465s	L. 65 mm
185400s	SURFIX® Lock Screw dia. 4.5 mm



Surfix® Locking screw
Sterile Titanium dia. 6.5 mm + lock-screw

Reference	Length
285620s	L. 20 mm
285625s	L. 25 mm
285630s	L. 30 mm
285635s	L. 35 mm
285640s	L. 40 mm
285645s	L. 45 mm
285650s	L. 50 mm
285655s	L. 55 mm
285660s	L. 60 mm
285665s	L. 65 mm
185600s	SURFIX® lock Screw dia. 6.5 mm


Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

Reference	Description
219845	Screwdriver hexa 2.5 x L. 191 mm
219865	Screwdriver hexa dia.3.5 x L. 230 mm
219445	AO Screwdriver tip hexa dia.2.5 x L. 76 mm
219465	AO Screwdriver tip hexa dia.3.5 x L. 76 mm
219645	Drilling guide dia.3.0 mm
219665	Drilling guide dia.3.5 mm
219545	AO Drill dia.3.0 x L. 190 mm
219565	Drill dia.3.5 x L. 200 mm
219900	Compression forceps L. 260 mm
219765	Bending forceps L. 193 - hole dia.6.5 mm
219345	Depth gauge hole dia.4.5 mm
219365	Depth gauge hole dia.6.5 mm
115116	K-Wire dia.1.6 x L. 150 mm
115425	K-Wire dia.2.5 x L. 200 mm
169001	TTCP trial implant 6 holes
169002	TTCP trial implant 8 holes
169003	TTCP trial implant 7 holes
169004	TTCP trial implant 9 holes

169101	ADVANSYS® Hindfoot sterilization container including:
169111	Base
996200	Lid
169104	Cylinder
169105	TTCP trial implant module
119909	Silicone wedge

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm





TIBIAXYS® Ankle Arthrodesis Plating System

TIBIAXYS® Distal Tibia and Fibula Osteotomy Plating Systems

The TIBIAXYS® Plates are intended for fixation of arthrodesis osteotomies and fractures of ankle joint, distal tibia and fibula.

- Each plate is anatomically contoured and available in 1 size.
- 2 systems of 2 anterior plates for ankle fusion (green left and blue right systems, including 1 antero-medial and 1 antero-lateral plates).
- 3 different plates for distal tibia osteotomies (medial plate, left lateral / right medial anterior plate, left medial/ right lateral anterior plate).
- 2 different plates for fibula osteotomies (4 and 6 holes).
- dia. 3.5 mm Surfix® or Surfix®-Alpha screws and lock-screws.
- dia. 4.0 mm cortical screws (only for 1 hole in each plate for the ankle fusion).
- Material:
 - > Plates: Titanium alloy TA6V.
 - > Surfix® and Surfix®-Alpha screws and lock-screws: Titanium alloy TA6V.
 - > Cortical screws: Titanium alloy TA6V.



Reference	Description
150010S	TIBIAXYS® Arthrodesis Plate / Left antero-medial
150020S	TIBIAXYS® Arthrodesis Plate / Left antero-lateral
150110S	TIBIAXYS® Arthrodesis Plate / right antero-medial
150120S	TIBIAXYS® Arthrodesis Plate / right antero-lateral

Reference	Description
150030S	TIBIAXYS® Osteotomy Plate / Lateral-right / Medial-left
150130S	TIBIAXYS® Osteotomy Plate / Medial-right / Lateral-left
150040S	TIBIAXYS® Osteotomy Plate / Medial

Reference	Description
150514S	TIBIAXYS® Fibula Osteotomy Plate 4 holes
150516S	TIBIAXYS® Fibula Osteotomy Plate 6 holes

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

TIBIAXYS® Ankle Arthrodesis Plating System TIBIAXYS® Distal Tibia and Fibula Osteotomy Plating Systems



Surfix®-Alpha variable angle locking screws – Sterile Titanium
dia. 3.5 mm + lock-screw

Reference	Length
295310s	L. 10 mm
295312s	L. 12 mm
295314s	L. 14 mm
295316s	L. 16 mm
295318s	L. 18 mm
295320s	L. 20 mm
295322s	L. 22 mm
295324s	L. 24 mm
295326s	L. 26 mm
295328s	L. 28 mm
295330s	L. 30 mm
295332s	L. 32 mm
295334s	L. 34 mm
295336s	L. 36 mm
295338s	L. 38 mm
295340s	L. 40 mm
295344s	L. 44 mm
295348s	L. 48 mm
295350s	L. 50 mm
195300s	Lock-screw dia. 3.5 mm

Surfix® Locking screws
Sterile Titanium
dia. 3.5 mm + lock-screw

Reference	Length
285310s	L. 10 mm
285312s	L. 12 mm
285314s	L. 14 mm
285316s	L. 16 mm
285318s	L. 18 mm
285320s	L. 20 mm
285322s	L. 22 mm
285324s	L. 24 mm
285326s	L. 26 mm
285328s	L. 28 mm
285330s	L. 30 mm
285332s	L. 32 mm
285334s	L. 34 mm
285336s	L. 36 mm
285338s	L. 38 mm
285340s	L. 40 mm
285344s	L. 44 mm
285348s	L. 48 mm
285350s	L. 50 mm
185300s	Lock-screw dia. 3.5 mm

Sterile Titanium cortical screw
dia. 4.0 mm

Reference	Length
150240s	L. 40 mm
150242s	L. 42 mm
150246s	L. 46 mm
150250s	L. 50 mm
150255s	L. 55 mm
150260s	L. 60 mm
150265s	L. 65 mm
150270s	L. 70 mm
150275s	L. 75 mm
150280s	L. 80 mm
150285s	L. 85 mm
150290s	L. 90 mm
150200s	L. 100 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

TIBIAXYS® Instruments

Reference	Description
115116	K-Wire dia.1.6 x L. 150 mm
115225	K-Wire dia.2.5 x L. 200 mm
159103	Wedge thickness 3 mm
159106	Wedge thickness 6 mm
159109	Wedge thickness 9 mm
159130	Cortical screw drilling guide dia.3.0 mm
159400	Depth gauge for dia.4.0 mm screws
159635	Compression forceps guide
159740	Compression forceps screw dia.4.0 x L. 40 mm
159755	Compression forceps screw dia.4.0 x L. 55 mm
159760	Compression forceps screw dia.4.0 x L. 60 mm
219035	Surfix® Alpha drilling guide dia.2.7 mm
219135	Screwdriver hexalobular T10
219335	Depth gauge for dia. 3.5 mm screws
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219445	AO Screwdriver tip hexa dia.2.5 x L. 76 mm
219535	AO drill dia.2.7 x L. 125 mm
219545	AO Drill dia.3.0 x L. 190 mm
219635	Drilling guide dia.2.7 mm
219835	Screwdriver hexa dia. 2.0 x L. 180 mm
219845	Screwdriver hexa dia. 2.5 x L. 191 mm
219960	Compression forceps L. 260 mm
159010	Trial plate Left Antero-Medial
159020	Trial plate Left Antero-Lateral
159110	Trial plate Right Antero-Medial (optional)
159120	Trial plate Right Antero-Lateral (optional)

Sterilization container

159991	Tibixys® sterilization container including
159960	Base
996200	Lid
119909	Silicone wedge

Trial container


159010	Left medial anterior trial plate
159020	Left lateral anterior trial plate
159110	Right medial anterior trial plate
159020	Right lateral anterior trial plate

Container

159071	Base
309942	Lid

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm




ACHILLON® Minimal Invasive Achilles Tendon Suture System – Single Use

The ACHILLON® System is a method to treat acute Achilles tendon ruptures.

- Needle driver.
- Suture size: USP (3.5 metric).
- Materials:
 - > Polycarbonate.
 - > Stainless steel 316LVM.

Reference	Description
119700	ACHILLON® + 2 Needles and 1 Needle Driver



Hintermann® Distractors

Range of 4 spreading forceps.

For each of the 4 models, two pins (2.5 mm for the large model, and 1.6 or 2.5 mm diameter for the small one) are slid into the appropriate holes of the spreading forceps, and are then introduced at appropriate localisation into each of the two bone fragments to be outspread. The spreader is then distracted, and locked when the distance between the two fragments is reached.

The flexible positioning of these spreading forceps enables them not to interfere with the main procedure. These spreaders are available in two different sizes (15 cm long, or 21 cm long), and with closed or outspread arms, making them fitting many possible indication.

Hintermann® Small Distractor – 15 cm

Reference	Description
119653	Hintermann® Distractor Small open handle
119654	Hintermann® Distractor Small closed handle

Hintermann® Large Distractor – 21 cm

Reference	Description
119664	Hintermann® Distractor Big model
119651	Hintermann® Distractor Large closed handle

Hintermann® Distractor open handle



Hintermann® Distractor closed handle



Osteoguides

- Suggested use:
 - > Scarf osteotomy.
 - > Chevron osteotomy.
 - > Basal osteotomy.
 - > Arthrodesis.
 - > Wedge osteotomies (Akin...).
- Details:
 - > Hole for K-Wire: Diameter 1.0 to 1.6 mm.
 - > Knurled handle for optimal grip.
 - > Blade thickness maxi 0.8 mm.
 - > Positioning tip.

Reference	Description
229011	Longitudinal
229012	Transversal



K-FIX® – K-Wire End Protector

The K-Fix® pin protector is a simple system to protect patients or hospital staff from K-Wires sharp ends.

- Delivered sterile.

Reference	Length
119050	K-FIX® for K-WIRE - Dia. 0.8 mm to 2.0 mm (5 units per pack)



Weil Clamp

Reference	Length
119106	Weil Clamp

General Optional Instruments

Reference	Description
998001	Scalpel Handle N° 3
998002	Scalpel Handle N° 4
998003	Beyer Bone Rongeur Forceps
998004	Backhaus Towel Forceps
998005	Beckmann Mini Retractor
998006	Museux Forceps (14 cm)
998007	Lange- Hohmann Retractor
998008	Liston Bone Cutting Forceps
998010	Ombredanne Square Points
998011	Wolkman Retractor
998012	Wolkman Curette
998013	Lambotte Raspatory 12 mm
998014	Guillaume Gouge 8 mm
998015	Mac Even Hammer
998016	Redon Needle
998017	Finochietto Needle Holder 13 cm
998018	Pauwels Osteotomes 15 mm
998019	Farabeuf Retractor 12 cm
998020	Dissection Forceps (GR) 14 cm
998021	Dissection Forceps Curved 14 cm
998022	Metzenbaum Scissors
998023	Straight Scissors
998024	Mayo Scissors
998025	Leriche Forceps
998026	Verbrugge Bone Holding Forceps
998027	Poirier Gouge (14 mm)
998031	Metatarsal Head Extractor
MA3511	Surfix® screw 3.5 mm extractor
119651	Hintermann® Distractor / Large Closed Jaws
119653	Hintermann® Distractor / Small Open Jaws
119654	Hintermann® Distractor / Small Closed Jaws
119106	Weil clamp
998028	Inge Retractor
998029	Banaleck Forceps Right
998030	Banaleck Forceps Left
119909	Silicone wedge
119989	Instrument Tray



MA3511

Surfix® screw extractor 3.5 mm

Post-Operative Shoes

Reference	Description
-----------	-------------

Flat Post Operative Shoe

999001	Small
999002	Medium
999003	Large
999004	X-Large



Reference	Description
-----------	-------------

Wedged Heel Post Operative Shoe

999100	Size 4 to 6.5
999101	Size 7 to 8.5
999102	Size 9 to 10.5





Integra[®] References

Index

05-0016.....	40	195 100S	18	105675	36
05-0017.....	40	219 027	18	105680.....	36
05-0106	40	219 127	18	105685	36
05-0108.....	40	219 227	18	105690	36
05-0109	40	229 027 S	10	111010.....	14
05-0110	40	229 028 S	10	111012.....	14
05-0112.....	40	229 037 S	10	111014.....	14
05-0206.....	40	229 038 S	10	111016.....	14
05-0208.....	40	229 227 S	10	111018.....	14
05-0209.....	40	229 228 S	10	111020.....	14
05-0210.....	40	229 237 S	10	111022	14
05-0212	40	229 238 S	10	111024	14
05-1002.....	40	295 110S.....	18	111026	14
05-1006.....	40	295 112S.....	18	111028	14
05-1008.....	40	295 114S.....	18	111312	14
05-1009.....	40	295 116S.....	18	111314.....	14
05-1010	40	295 118S.....	18	111316	14
05-1012.....	40	295 120S.....	18	111318.....	14
05-1020	40	295 122S.....	18	111320.....	14
05-5001.....	40	295 124S.....	18	111322	14
115 008 S	10	295 126S.....	18	111324	14
115 070.....	18	295 128S	18	111326.....	14
115 070 (S).....	10, 11	295 130S.....	18	111328	14
115 100 (S).....	10, 11	295 132S	18	111328S.....	26
115 116	18	295 134S	18	111330.....	14
119 004 (S).....	11	299 005.....	18	111330S.....	26
119 006 (S)	11	299 020.....	18	111332	14
119 016 (S).....	11	299 050.....	18	111332S.....	26
119 023 (S)	10	309 708	49	111334	14
119 025 (S).....	10	320 100S	18	111334S.....	26
119 027 (S)	10	320 200S.....	18	111424.....	34
119 223 (S).....	10	329 002	18	111426.....	34
119 225 (S)	10	329 500	18	111428.....	34
119 227 (S).....	10	329 900	18	111430.....	34
119 501	37	329 901.....	18	111432.....	34
119 506	37	329 902	18	111434.....	34
119 510	37	105426	36	111436.....	34
119 516	37	105428	36	111440.....	34
119 522.....	37	105430.....	36	111445.....	34
119 525.....	37	105432	36	111450.....	34
119 531.....	37	105434.....	36	111455.....	34
119 532	37	105436.....	36	111460.....	34
119 536	37	105438	36	111530S	35
129 710.....	37	105440.....	36	111535S.....	35
159 004 (S)	11	105445	36	111540S	35
159 023 (S).....	10	105450.....	36	111545S.....	35
159 025 (S).....	10	105455	36	111550S	35
159 027 (S).....	10	105460.....	36	111555S.....	35
159 223 (S).....	10	105640.....	36	111560S.....	35
159 225 (S).....	10	105645.....	36	111565S.....	35
159 227 (S).....	10	105650.....	36	111570S.....	35
189 210.....	18	105655.....	36	111575S.....	35
189 215.....	18	105660	36	111580S	35
189 220.....	18	105665.....	36	111740S.....	35
		105670.....	36	111745S	35

111750S.....	35	117155.....	24	119632.....	49
111755S.....	35	117160.....	24	119641.....	47
111760S.....	35	117250.....	24	119645.....	21
111765S.....	35	117255.....	24	119646.....	47
111770S.....	35	117260.....	24	119647.....	47
111775S.....	35	117340.....	24	119651.....	65, 67
111780S.....	35	117345.....	24	119653.....	65, 67
111785S.....	35	117350.....	24	119654.....	65, 67
111790S.....	35	117440.....	24	119657.....	49
111795S.....	35	117445.....	24	119659.....	49
111800S.....	35	117450.....	24	119662.....	47
111805S.....	35	119050.....	17, 66	119664.....	47, 65
111810S.....	35	119104.....	34	119665.....	47
111815S.....	35	119105.....	10	119666.....	47
111820S.....	35	119106.....	17, 66, 67	119667.....	47
112011.....	15, 16	119133.....	27	119669.....	47
112012.....	15, 16	119135.....	27	119690.....	47
112013.....	15, 16	119138.....	34	119691.....	47
112014.....	15, 16	119139.....	34	119692.....	47
112111.....	15	119150.....	35	119693.....	47
112111S.....	15	119151.....	35	119694.....	47
112112.....	15	119152.....	35	119695.....	47
112112S.....	15	119153.....	35	119696.....	47
112113.....	15	119154.....	35	119697.....	47
112113S.....	15	119155.....	35	119698.....	47
112114(s).....	15	119156.....	35	119699.....	47
112114S.....	15	119170.....	35	119700.....	64
114438.....	34	119171.....	35	119751.....	47
115070.....	27	119172.....	35	119757.....	49
115100.....	25, 27, 42	119173.....	35	119760.....	49
115101.....	43	119174.....	35	119761.....	49
115116.....	27, 34, 38, 47, 63	119175.....	35	119762.....	49
115216.....	21, 25, 28, 60	119176.....	35	119763.....	49
115225.....	37, 38, 47, 63	119177.....	35	119764.....	49
115425.....	60	119311.....	11	119765.....	49
115516.....	35	119510.....	37	119766.....	49
117010.....	25	119516.....	34	119767.....	49
117012.....	25	119531.....	38	119768.....	49
117014.....	25	119541.....	37	119769.....	49
117016.....	25	119546.....	37	119778.....	47
117018.....	25	119551.....	37	119779.....	47
117020.....	25	119552.....	37, 58	119790.....	49
117022.....	25	119555.....	37	119791.....	49
117024.....	25	119561.....	37, 38	119792.....	49
117026.....	25	119562.....	37	119793.....	49
117028.....	25	119566.....	37	119794.....	49
117030.....	25	119580.....	38	119795.....	49
117032.....	25	119582.....	38	119796.....	49
117034.....	25	119590.....	38	119797.....	49
117110.....	25	119592.....	38	119798.....	49
117112.....	25	119601.....	21	119799.....	49
117114.....	25	119609.....	47	119825.....	39
117116.....	25	119611.....	21, 47	119831.....	39
117118.....	25	119618.....	25	119850.....	39
117150.....	24	119631.....	49	119859.....	39

119860	39	121800S	35	139211-215	21
119861	39	121805S	35	139216	21
119862	39	121810S	35	139311-313	21
119863	39	121815S	35	139321-323	21
119864	39	121820S	35	139331-333	21
119865	39	121950	35	139341-343	21
119867	39	121951	35	139351-353	21
119880	39	129710	25, 27, 28, 37	139940	21
119881	39	129714	25, 27, 28	139941	21
119882	39	129716	25, 27, 28	139942	21
119909	19, 21, 25, 27, 29, 35, 37, 42, 43, 47, 56, 58, 60, 63, 67	129718	25, 27, 28	139943	21
119935	34	129720	25, 27, 28	139953	21
119936	34	129722	27, 28	141009	39
119937	34	129724	25, 27, 28	141010	39
119938	34	129726	25, 27, 28	141011	39
119949	37	129728	25, 27, 28	141012	39
119951	37	129730	25, 27, 28	141013	39
119953	37	129731	25	141014	39
119955	38	129732	25	141015	39
119956	37	129733	25	141017	39
119957	37	129734	25	150010S	61
119958	37	129735	25	150020S	61
119959	38	129736	25	150030S	61
119965	49	129752	27, 28	150040S	61
119966	49	129930	25	150110S	61
119967	49	129931	25	150120S	61
119974	37	129932	25	150130S	61
119975	38	129933	25	150200S	62
119976	37	129935	28	150240S	62
119977	38	129936	25	150242S	62
119978	49	129938	28	150246S	62
119989	67	131014S	14	150250S	62
121530S	35	131016S	14	150255S	62
121535S	35	131018S	14	150260S	62
121540S	35	131020S	14	150265S	62
121545S	35	131022S	14	150270S	62
121550S	35	131024S	14	150275S	62
121555S	35	131026S	14	150280S	62
121560S	35	131028S	14	150285S	62
121565S	35	131030S	14	150290S	62
121570S	35	131032S	14	150514S	61
121575S	35	131034S	14	150516S	61
121580S	35	139100-106	21	159010	63
121740S	35	139107	21	159020	63
121745S	35	139111-116	21	159027S	27
121750S	35	139120	21	159071	63
121755S	35	139121	21	159103	63
121760S	35	139121-126	21	159106	63
121765S	35	139128	21	159109	63
121770S	35	139129	21	159110	63
121775S	35	139140	21	159120	63
121780S	35	139150	21	159127	47
121785S	35	139201-205	21	159130	63
121790S	35	139210	21	159400	63
121795S	35	139211	21	159635	63

159740.....	63	200016S.....	14	219535.....	42, 43, 63
159755.....	63	200018S.....	14	219545.....	60, 63
159760.....	63	200020S.....	14	219565.....	60
159960.....	63	200022S.....	14	219627.....	27
159991.....	63	200024S.....	14	219635.....	42, 43, 63
169001.....	60	200026S.....	14	219645.....	60
169002.....	60	200028S.....	14	219665.....	60
169003.....	60	200030S.....	14	219735.....	42, 43
169004.....	60	201010.....	22	219765.....	60
169021.....	42	201020.....	22	219835.....	42, 43, 63
169022.....	42	201030.....	22	219845.....	60, 63
169023.....	42	201040.....	22	219865.....	60
169031.....	42	202010.....	22	219900.....	60
169032.....	42	202020.....	22	219960.....	63
169033.....	42	202030.....	22	229002.....	10
169041.....	42	202040.....	22	229003.....	10
169042.....	42	203010.....	22	229004.....	10
169051.....	42	203020.....	22	229005.....	10
169052.....	42	203030.....	22	229011.....	10, 17, 66
169101.....	60	203040.....	22	229012.....	10, 17, 66
169102.....	42	203100.....	22	229101.....	10, 11
169104.....	60	203200.....	22	229102.....	11
169105.....	60	203300.....	22	229201.....	11
169106.....	42, 43	203415.....	22	229202.....	11
169107.....	42	203500.....	22	229951.....	10
169111.....	60	203600.....	22	229952.....	11
169112.....	42	205000.....	22	229961.....	10
181001S.....	59	213113.....	16	229962.....	11
181002S.....	59	213114.....	16	229965.....	12
181003S.....	59	213115.....	16	229966.....	13
181004S.....	59	213116.....	16	229970.....	8, 10, 11
181011S.....	59	213117.....	16	229971.....	10
181012S.....	59	213213.....	16	229975.....	12
181013S.....	59	213214.....	16	229976.....	13
181014S.....	59	213215.....	16	229990.....	8
181021S.....	41	213216.....	16	229996.....	8
181022S.....	41	213217.....	16	230001S.....	19
181023S.....	41	213313.....	16	230002S.....	19
181031S.....	41	213314.....	16	239000.....	19
181032S.....	41	213315.....	16	239001.....	19
181033S.....	41	213316.....	16	239011.....	19
181041S.....	41	213317.....	16	239012.....	19
181042S.....	41	213512.....	16	239021.....	19
181051S.....	41	213820.....	16	239022.....	19
181052S.....	41	219027.....	27	239030.....	19
185100S.....	26	219035.....	63	239500.....	19
185200S.....	26	219127.....	27	250100.....	20
185300S.....	41, 43, 62	219135.....	63	250101.....	20
185400S.....	59	219227.....	27	250102S.....	20
185600S.....	59	219335.....	42, 43, 63	250103.....	20
195100S.....	26	219345.....	60	250104.....	20
195200S.....	26	219365.....	60	250105.....	20
200010(S).....	14	219435.....	42, 43, 63	250106.....	20
200012S.....	14	219445.....	60, 63	250201.....	20
200014S.....	14	219465.....	60	250202.....	20

250203.....	20	285316S.....	41, 43, 62	290445S.....	26
250204.....	20	285318S.....	41, 43, 62	290450S.....	26
250205.....	20	285320S.....	41, 43, 62	295110S.....	26
250311.....	20	285322S.....	41, 43, 62	295112S.....	26
250312.....	20	285324S.....	41, 43, 62	295114S.....	26
250321.....	20	285326S.....	41, 43, 62	295116S.....	26
250322.....	20	285328S.....	41, 43, 62	295118S.....	26
250323.....	20	285330S.....	41, 43	295120S.....	26
250331.....	20	285330S.....	62	295122S.....	26
250332.....	20	285332S.....	41, 43, 62	295124S.....	26
250333.....	20	285334S.....	41, 43, 62	295126S.....	26
250341.....	20	285336S.....	43	295210S.....	26
250342.....	20	285336S.....	62	295212S.....	26
250343.....	20	285338S.....	43	295214S.....	26
250351.....	20	285338S.....	62	295216S.....	26
250352.....	20	285340S.....	43	295218S.....	26
250353.....	20	285340S.....	62	295220S.....	26
275560.....	47	285344S.....	62	295222S.....	26
278002S.....	29	285348S.....	62	295224S.....	26
278005S.....	29	285350S.....	62	295226S.....	26
278007S.....	29	285414S.....	59	295310S.....	62
278008S.....	29	285416S.....	59	295310S.....	62
278009S.....	29	285418S.....	59	295312S.....	62
278010S.....	29	285420S.....	59	295314S.....	62
278100.....	29	285422S.....	59	295316S.....	62
278110.....	29	285424S.....	59	295318S.....	62
278119.....	29	285426S.....	59	295320S.....	62
278121.....	29	285428S.....	59	295322S.....	62
278122.....	29	285430S.....	59	295324S.....	62
278123.....	29	285435S.....	59	295326S.....	62
278124.....	29	285440S.....	59	295328S.....	62
278900.....	29	285445S.....	59	295330S.....	62
278901.....	29	285450S.....	59	295332S.....	62
278902.....	19, 29	285455S.....	59	295334S.....	62
285110S.....	26	285460S.....	59	295336S.....	62
285112S.....	26	285465S.....	59	295338S.....	62
285114S.....	26	285625S.....	59	295340S.....	62
285116S.....	26	285630S.....	59	295344S.....	62
285118S.....	26	285635S.....	59	295348S.....	62
285120S.....	26	285640S.....	59	295350S.....	62
285122S.....	26	285645S.....	59	299005.....	27
285124S.....	26	285650S.....	59	299010.....	27
285126S.....	26	285655S.....	59	299020.....	27
285210S.....	26	285660S.....	59	299050.....	27
285212S.....	26	285665S.....	59	299070.....	27
285214S.....	26	290150S.....	26	299080.....	27
285216S.....	26	290155S.....	26	299150.....	27
285218S.....	26	290160S.....	26	299155.....	27
285220S.....	26	290250S.....	26	299160.....	27
285222S.....	26	290255S.....	26	299250.....	27
285224S.....	26	290260S.....	26	299255.....	27
285226S.....	26	290340S.....	26	299260.....	27
285310S.....	41, 62	290345S.....	26	299340.....	27
285312S.....	41, 43, 62	290350S.....	26	299345.....	27
285314S.....	41, 43, 62	290440S.....	26	299350.....	27

299440	27	301224	48	309305	47
299445	27	301225	48	309306	47
299450	27	301232	48	309309	47
299900	27	301233	48	309360	47
299901	27	301234	48	309361	47
299903	27	301235	48	309362	47
299904	27	302110	45	309363	47
300005	45	302111	45	309364	47
300006	45	302112	45	309365	47
300007	45	302113	45	309366	47
300009	45	302114	45	309370	47
300105	45	302115	45	309371	47
300106	45	302116	45	309372	47
300107	45	302121	48	309373	47
300109	45	302122	48	309374	47
300205	45	302123	48	309375	47
300206	45	302124	48	309376	47
300207	45	302125	48	309601	47
300209	45	302201	44	309604	47
300305	45	302202	44	309605	47
300306	45	302203	44	309606	47
300307	45	302204	44	309607	47
300309	45	302205	44	309608	47
300405	45	302206	44	309615	47
300406	45	302222	48	309620	47
300407	45	302223	48	309622	47
300409	45	302224	48	309623	47
300505	45	302225	48	309624	47
300506	45	302232	48	309625	47
300507	45	302233	48	309626	47
300509	45	302234	48	309627	47
300605	45	302235	48	309630	47
300606	45	303312S	48	309635	47
300607	45	303316S	48	309636	47
300609	45	303320S	48	309637	47
301110	45	303328S	45	309645	47, 58
301111	45	303330S	45	309652	47
301112	45	303334S	45, 48	309656	47
301113	45	303338S	45, 48	309657	47
301114	45	303340S	45, 48	309661	47
301115	45	303342S	45, 48	309680	47
301116	45	303512S	53	309681	47
301121	48	303516S	53	309682	47
301122	48	303520S	53	309683	47
301123	48	303530S	51	309684	47
301124	48	303534S	51, 53	309685	47
301125	48	303538S	51, 53	309686	47
301201	44	303540S	51, 53	309690	47
301202	44	303542S	51, 53	309691	47
301203	44	309200	47	309692	47
301204	44	309300	47	309693	47
301205	44	309301	47	309694	47
301206	44	309302	47	309695	47
301222	48	309303	47	309696	47
301223	48	309304	47	309697	47

309699	47	352251	50	511020	55
309712	49	352252	50	511022	55
309753	47	352253	50	511025	55
309773	47	352254	50	511027	55
309900	47	352255	50	511030	55
309902	47	352256	50	511035	55
309904	47	352512	53	511040	55
309920	47	352513	53	511045	55
309940	47	352514	53	511055	55
309942	47	352515	53	511060	55
309960	47	352522	52	511065	55
309970	47	352523	52	511070	55
309980	47	352524	52	511075	55
330021S	43	352525	52	511080	55
330023S	43	352532	52	511085	55
330025S	43	352533	52	511090	55
330217S	43	352534	52	511095	55
330220S	43	352535	52	511100	55
330225S	43	500001	54	511105	55
330230S	43	500050	54	511110	55
330420S	43	500080	54	519002	56
330425S	43	500150	54	519003	56
330430S	43	500180	54	519004	56
339001	43	500250	54	519005	56
339003	43	500280	54	519006	56
339004	43	500350	54	519007	56
339005	43	500380	54	519008	56
339920	43	501020	55	519009	56
339950	43	501022	55	519010	56
351152	51	501025	55	519011	56
351153	51	501027	55	519012	56
351154	51	501030	55	519013	56
351155	51	501035	55	519014	56
351251	50	501040	55	519015	56
351252	50	501045	55	519016	56
351253	50	501050	55	519017	56
351254	50	501055	55	519020	56
351255	50	501060	55	519021	56
351256	50	501065	55	519028	56
351512	53	501070	55	519029	56
351513	53	501075	55	519030	56
351514	53	501080	55	519032	56
351515	53	501085	55	519034	56
351522	52	501090	55	519040	56
351523	52	501095	55	519110	56
351524	52	501100	55	519120	56
351525	52	501105	55	519121	56
351532	52	501110	55	519130	56
351533	52	510005	54	519131	56
351534	52	510111	54	519133	56
351535	52	510141	54	519135	56
352152	51	510211	54	519160	56
352153	51	510241	54	519175	56
352154	51	510311	54	519178	56
352155	51	510341	54	519179	56

250
240
230
220
210
200
190
180
170
160
150
140
130
120
110
100
90
80
70
60
50
40
30
20
10
MILLIMETERS



Distributed by

Integra LifeSciences Services (France) SAS
 Sales & Marketing EMEA
 Immeuble Séquoia 2 • 97 allée Alexandre Borodine
 Parc technologique de la Porte des Alpes
 69800 Saint Priest • FRANCE
 Phone: +33 (0)4 37 47 59 00 • Fax: +33 (0)4 37 47 59 99
 emea.info@integralife.com • integralife.eu

Customer Services
 International: +33 (0)4 37 47 59 50 • +33 (0)4 37 47 59 25 (Fax) • csmea@integralife.com
 United Kingdom: +44 (0)1 264 345 781 • +44 (0)1 264 363 782 (Fax) • csuk.ortho@integralife.com
 France: +33 (0)4 37 47 59 10 • +33 (0)4 37 47 59 29 (Fax) • cs-ortho@integralife.com
 Benelux: +32 (0)2 257 4130 • +32 (0)2 253 2466 (Fax) • custsvcbenelux@integralife.com
 Switzerland: +41 (0)2 27 21 23 30 • +41 (0)2 27 21 23 99 (Fax) • custsvcsuisse@integralife.com

©2013 Integra LifeSciences Corporation. All rights reserved. ILS 02-02-004-01-12 / ILS 02-02-006-01-12
 PRODUCTS FOR SALE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.