



Integra®
Forefoot
Mid & HindFoot
Upper Extremities
Large Bones
Products Catalogue

Products for sale in Europe, Middle-East and Africa only
Catalogue des ventes

Table of Contents

Warranty Information	5
Integra® Forefoot Solutions	7
Forefoot Set	8
Bold® Compression Screw dia. 2.5 mm / 3.0 mm / 3.7 mm	14
QWIX® Fixation Screw dia. 3.0 mm	14
SPIN® Screw dia. 2.0 mm / 2.7 mm	15
UNI-CLIP® Compressive Staple	16
SOLUSTAPLE® Standard Staple	16
Osteoguides	17
K-FIX® - K-wire End Protector	17
Weil Clamp	17
B-BOP® LOCK plate	18
IPP-ON® Interphalangeal Arthrodesis Implant	19
METIS® Metatarso-phalangeal prosthesis	20
K2™ Hemi Toe Implant	22
Movement Great Toe system	23
HALLU®-FIX MTP Arthrodesis System	24
HALLU® Plates	24
HALLU® SNAP-OFF® Screws	25
HALLU®-FIX Instruments	25
HALLU®-LOCK Arthrodesis System	26
HALLU®-REAM set	28
DPR System Minimal Invasive Foot Surgery	29
Post-Operative Shoes	30
Integra® Mid & HindFoot Solutions	33
Subtalar Arthrodesis	34
QWIX® Fixation screw	34
Large QWIX® Fixation and Positioning Screws	35
I.CO.S.® Ideal Compression Screw	36
Flatfoot implants	39
KALIX® II Flat Foot Implant	39
MBA™ Titanium Subtalar Implant	40
Dorsal Lisfranc plate & Medial Lisfranc plate	41

ADVANSYS® Lisfranc Plate	41
UNI-CP™ Compression Plate	43
Ankle Solutions	44
Total Ankle Replacement Hintegra® prosthesis	44
Ankle Fusion	54
PANTA® Nail	54
ADVANSYS® Plate	59
TIBIAXYS® Plating System	61
Minimally Invasive Achilles Tendon Ruptures	64
ACHILLON® System	64
Surgical Instruments	65
Hintermann® Distractors	65
Osteoguides	66
K-FIX® Protector	66
Weil Clamp	66
General Optional Instruments	67
Post-Operatives Shoes	68
Integra® Upper Extremities Solutions	71
Hand Deformity & Fracture Repair	72
Pyrocarbon PIP Implant	73
Silicon PIP Implant	74
Pyrocarbon MCP Implant	76
Silicon MCP Implant	78
The Thumb carpometacarpal (CMC) joint	80
Pyrodisk™	80
NuGrip™ hemi arthroplasty implant	81
Saddle™ CMC Implant	82
PCS™ spacer	83
Hand Fracture Repair	84
Safeguard™ Mini-open carpal tunnel release system	84
Kompressor™, screw system	85
Hand Fractures and Deformity repair set	86
Upper modules	87
Other possible drills	88
QWIX® Fixation Screw	89
BOLD® Compression Screw	89
UNI-CLIP® Notched & Standard Staples	90
SPIN® SNAP-OFF® screw	91
Wrist Arthritis & Deformity Repair	92
Universal2™ total wrist implant system	92
Integra Total Wrist Fusion Plating System - ITWFP	94

Spider™ – Mini Spider™ / limited wrist fusion system	96
PyroCarbon Lunate Implant	98
First Choice DRUJ system™	99
Wifix® plate – Surfix® volar plate for distal radius	102
Elbow pathologies solutions	103
Elbow pathologies solutions	103
Pycarbon Radial Head implants (CMRH)	104
Metal radial head	105
RADFx™ radius fixation system	106
Katalyst™ telescoping bipolar radial head system	108
EndoRelease™ endoscopic cubital tunnel release system	109
Shoulder repair	110
Titan™ Humeral Resurfacing Arthroplasty	110
Titan™ Modular Shoulder system	112
PyroTitan™ Humeral Resurfacing Arthroplasty	115
PSH plate – Surfix® Proximal humeral plate	117
Surfix® osteotomy solutions	120
Knee osteotomy plates	121
POF plate	122
PTO plate	123
POM plate	124
PTC plate	125
PTF plate	126
ATAO® Bone Graft Substitute	127
Set osteotomy	127
Pediatric osteotomy plates	128
PSFE plates	128
PDTL plates	130
Pediatric set	132
Product References	135

Pricing:

Prices stated in specific written quotations are firm for thirty days from the date given, and are otherwise subject to change without prior notice. Pricing terms stated on written agreements are governed by such agreements.

Minimum Order Requirements:

Minimum order requirement is €250 (a €20 minimum order requirement charge will be added to any order under €250).

Ordering Procedure:

A written purchase order on the customer's form may be requested for purchases of Integra LifeSciences products.

Acceptance of Orders:

Orders are accepted upon approval by Integra Customer Service.

Return Policy:

- Authorization from Customer Service must be obtained prior to returning product.
- Sterile product must be returned in unopened, undamaged cartons, packed to prevent damage.
- Non-sterile product must be returned in unused saleable condition in original package.
- Custom or special order products will not be accepted for credit.
- Credit will be issued for goods returned prior to ninety days from ship date with a 20% restocking charge. This assumes that the product returned is not damaged and can be verified to have not been used or opened.

Integra Limited Warranty – Integra products

INTEGRA LIFESCIENCES CORPORATION and its wholly owned subsidiaries ("INTEGRA") warrant to INTEGRA authorized distributors and the original purchaser only that each new INTEGRA Surfix[®] product is free from manufacturing defects in material and workmanship when used for its intended surgical purpose. For purposes of products sold by INTEGRA through an authorized distributor of INTEGRA, "original purchaser" shall include the purchaser of INTEGRA products to whom the distributor first sells the product.

If any covered defect occurs, the purchaser or distributor should communicate directly with INTEGRA. If purchaser or distributor seeks to invoke the terms of this warranty, the product must be returned to INTEGRA. The defective product should be returned promptly, properly packaged and postage prepaid. Loss or damage in return shipment to INTEGRA shall be at sender's risk. INTEGRA's sole responsibility under this warranty shall be repair or replacement, at INTEGRA's sole discretion at INTEGRA's expense, subject to the terms of this warranty and applicable agreements.

INTEGRA provides for no commercial or additional warranty whatsoever, except as may be otherwise expressly agreed in a separate commercial agreement or undertaking duly signed by an INTEGRA authorised representative.

IN NO EVENT SHALL INTEGRA BE LIABLE FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA shall not be responsible for, any loss arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as, in INTEGRA's judgment, to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA.

Customer Service Contact Information:

International

+33 (0)4 37 47 59 50
+33 (0)4 37 47 59 25 (Fax)
csemea@Integralife.com

Benelux

+32 (0)2 257 4130
+32 (0)2 253 2466 (Fax)
custsvcbenelux@Integralife.com

France

+33 (0)4 37 47 59 10
+33 (0)4 37 47 59 29 (Fax)
cs-ortho@Integralife.com

Switzerland

+41 (0)2 27 21 23 30
+41 (0)2 27 21 23 99 (Fax)
custsvcsuisse@Integralife.com

United Kingdom

+44 (0)1 264 345 780
+44 (0)1 264 363 782 (Fax)
custsvcs.uk@integralife.com

THIS INTEGRA LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON INTEGRA'S PART OR THE PART OF ITS DISTRIBUTORS, AND INTEGRA NEITHER ASSUMES NOR AUTHORIZES ANY REPRESENTATIVE OR OTHER PERSON TO ASSUME FOR IT ANY OTHER LIABILITY IN CONNECTION WITH INTEGRA'S PRODUCTS.

INTEGRA DISCLAIMS ALL OTHER WARRANTIES, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS. No warranty or guarantee may be created by any act or statement nor may this Standard Warranty be modified in any way, except as a result of a writing signed by an officer of INTEGRA. These limitations on the creation or modification of this warranty may not be waived or modified orally or by any conduct.

IN NO EVENT SHALL INTEGRA AUTHORIZED DISTRIBUTORS BE LIABLE TOWARDS THE ORIGINAL PURCHASER FOR ANY INCIDENTAL, INDIRECT, CONSEQUENTIAL OR PUNITIVE DAMAGES IN CONNECTION WITH THE ACQUISITION OR USE OF ANY INTEGRA PRODUCT. Further, this warranty shall not apply to, and INTEGRA authorized distributors shall not be responsible towards the original purchaser for, any loss, arising in connection with the purchase or use of any INTEGRA product that has been repaired by anyone other than an authorized INTEGRA service representative or altered in any way so as to affect its stability or reliability, or which has been subject to misuse, negligence or accident, or which has been used otherwise than in accordance with the instructions furnished by INTEGRA. THIS INTEGRA DISTRIBUTOR LIMITED WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES TOWARDS THE ORIGINAL PURCHASER ON INTEGRA AUTHORIZED DISTRIBUTOR'S PART.

INTEGRA AUTHORISED DISTRIBUTORS DISCLAIM ALL OTHER WARRANTIES TOWARDS THE ORIGINAL PURCHASER, EXPRESS OR IMPLIED INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE OR APPLICATION OR WARRANTY OF QUALITY AS WELL AS ANY EXPRESS OR IMPLIED WARRANTY TO PATIENTS.

All the references numbers mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED".






Integra® Forefoot Solutions

Forefoot Set

A new modular and scalable forefoot instruments set:

- Bold® 2.5 – 3.0 – 3.7 Compression screw.
- QWIX® 3.0 Fixation screw.
- Spin® 2.0 – 2.7 Spin break off screws.
- Solustaple® Standard staple.
- Uni-clip® Compression Staple.



Screw instrumentation module
(details page 8)



Staple instrumentation module
(details page 9)



Generic module
(details page 10)



Implant rack module
(details page 10)



Frame

229996	Frame - 2 floors
229990	Frame - 3 floors
229970	Forefoot II- lid- generic



- Staple instrumentation
- Screw instrumentation




- Non sterile Implants
- Staple instrumentation
- Screw instrumentation



- Non sterile Implants
- Generic module
- Screw instrumentation

Forefoot Set – Screw Module

Reference	Description
119105	Scarf clamp
229002	Screwdriver - Hexa 2.0
229005	AO Screwdriver - Hexa 2.0 - AO Attachment - Cannulated
229003	Screwdriver - T7 Star - Cannulated
229004	Screwdriver - AO - T7 Star - Cannulated
997301	Depth Gauge
997311	AO Straight Handle
229011	Longitudinal Cutting Guide (Optional)
229012	Transversal Cutting Guide (Optional)
229101	Ruler
229951	Sterilization Screws Container including:
229961	Tray
229970	Lid
229971	Mat



Drill Reference	Corresponding Product	2 in 1	Diameter	Cannulated	Length	AO attachment
159 027 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	12 mm	■
159 025 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	22 mm	■
159 023 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	32 mm	■
119 227 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		12 mm	
119 225 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		22 mm	
119 223 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		32 mm	
119 027 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	12 mm	
119 025 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	22 mm	
119 023 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm	■	32 mm	
159 227 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		12 mm	■
159 225 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		22 mm	■
159 223 (S)	Bold® 3.0 / Bold® 3.7	■	2.2 - 3.0 mm		32 mm	■
229 027 S	Bold® 2.5	■	1.9 - 2.5 mm	■	12 mm	■
229 037 S	Bold® 2.5	■	1.9 - 2.5 mm	■	22 mm	■
229 227 S	Bold® 2.5	■	1.9 - 2.5 mm	■	12 mm	
229 237 S	Bold® 2.5	■	1.9 - 2.5 mm	■	22 mm	
229 028 S	Bold® 2.5	■	1.9 - 2.5 mm		12 mm	■
229 038 S	Bold® 2.5	■	1.9 - 2.5 mm		22 mm	■
229 228 S	Bold® 2.5	■	1.9 - 2.5 mm		12 mm	
229 238 S	Bold® 2.5	■	1.9 - 2.5 mm		22 mm	

K-wire Reference	Corresponding Product	Diameter	Length
115 070 (S)	Bold® 3.0 / Bold® 3.7	1.0 mm	70 mm
115 100 (S)	Bold® 3.0 / Bold® 3.7	1.0 mm	100 mm
115 008 S	Bold® 2.5	0.8 mm	70 mm

All drills and K-wires are delivered sterile or non sterile except Bold® 2.5 K-wire which is only delivered in sterile.

(S) means available in sterile and non sterile

Forefoot Set – Staple Module

UNI-CLIP® Associated Instrument


229201	Uni-Clip® Drilling guide
229202	Uni-Clip® impactor
119311	Uni-Clip® spreading forceps

Solustaple® Associated Instrument

229102	Solustaple® impactor 90° & 26°
229101	Solustaple® drilling guide
997301	Depth gauge - dia 2.0 X L. 10 - 60 mm

Sterilization Container

229952	Staples container includes:	
229962	Tray	
229970	lid	



Drill Reference	Corresponding Product	Diameter	Cannulated	Length	AO attachment
119 004 (S)	Uni-Clip®	2.2 mm	■	34 mm	
119 006 (S)	Uni-Clip®	2.2 mm		34 mm	
159 004 (S)	Uni-Clip®	2.2 mm	■	34 mm	■
119 016 (S)	Uni-Clip®	2.2 mm		34 mm	■

K-wire Reference	Corresponding Product	Diameter	Length
115 070 (S)	Uni-Clip®	1.0 mm	70 mm
115 100 (S)	Uni-Clip®	1.0 mm	100 mm


All drills and K-wires are delivered sterile or non sterile except Bold® 2.5 K-wire which is only delivered in sterile.

(S) means available in sterile and non sterile

Forefoot Set – Implant Module

Sterilization Container


229965	Forefoot II - tray - non sterile implants
229975	Forefoot II - rack - non sterile implants



Forefoot Set – Generic Module

Sterilization Container

229966	Forefoot II - tray - generic
229976	Forefoot II - mat - generic



Bold® Compression Screw Dia. 2.5 mm / 3.0 mm / 3.7 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Totally intra-osseous implants.
- Large range of sizes:
 - > 10 to 30 mm for Bold® 2.5 screws by 2 mm increment.
 - > 10 to 34 mm for Bold® 3.0 screws by 2 mm increment.
 - > 14 to 34 mm for Bold® 3.7 screws by 2 mm increment.
- K-Wire guided drilling and insertion.
- Self tapping.
- True & controlled compression.
- Dual threads.
- Material: Titanium alloy TA6V.



Bold® dia 2.5 mm

Reference	Length
200010(S)	L. 10 mm
200012(S)	L. 12 mm
200014(S)	L. 14 mm
200016(S)	L. 16 mm
200018(S)	L. 18 mm
200020(S)	L. 20 mm
200022(S)	L. 22 mm
200024(S)	L. 24 mm
200026(S)	L. 26 mm
200028(S)	L. 28 mm
200030(S)	L. 30 mm



Bold® dia 3.0 mm

Reference	Length
111010(S)	L. 10 mm
111012(S)	L. 12 mm
111014(S)	L. 14 mm
111016(S)	L. 16 mm
111018(S)	L. 18 mm
111020(S)	L. 20 mm
111022(S)	L. 22 mm
111024(S)	L. 24 mm
111026(S)	L. 26 mm
111028(S)	L. 28 mm
111030(S)	L. 30 mm
111032(S)	L. 32 mm
111034(S)	L. 34 mm



Bold® dia 3.7 mm

Reference	Length
131014(S)	L. 14 mm
131016(S)	L. 16 mm
131018(S)	L. 18 mm
131020(S)	L. 20 mm
131022(S)	L. 22 mm
131024(S)	L. 24 mm
131026(S)	L. 26 mm
131028(S)	L. 28 mm
131030(S)	L. 30 mm
131032(S)	L. 32 mm
131034(S)	L. 34 mm

QWIX® Fixation Screw, Dia. 3.0 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction.

- Self-drilling and self-tapping.
- Totally intra-osseous fixation.
- K-Wire guided drilling and insertion.
- Material: Titanium alloy TA6V.



QWIX® Dia. 3.0 mm

Reference	Length
111312(S)	L. 12 mm
111314(S)	L. 14 mm
111316(S)	L. 16 mm
111318(S)	L. 18 mm
111320(S)	L. 20 mm
111322(S)	L. 22 mm
111324(S)	L. 24 mm
111326(S)	L. 26 mm
111328(S)	L. 28 mm
111330(S)	L. 30 mm
111332(S)	L. 32 mm
111334(S)	L. 34 mm

(S) means available in sterile and non sterile

SPIN[®] Break Off Screw, Dia. 2.0 mm (for associated instrumentation, see p. 8)

For bone fixation.

- Break-off barrel.
- Self-drilling and self-tapping.
- Material: Titanium alloy TA6V.

Reference	Description
112011(S)	Spin [®] Dia. 2.0 x L. 11 mm
112012(S)	Spin [®] Dia. 2.0 x L. 12 mm
112013(S)	Spin [®] Dia. 2.0 x L. 13 mm
112014(S)	Spin [®] Dia. 2.0 x L. 14 mm



SPIN[®] Break Off Screw, Dia. 2.7 mm (for associated instrumentation, see p. 8)

Spin[®] 2.7 screw is recommended for soft or osteoporotic bones fixation.

- Break-off barrel.
- Self-drilling and self-tapping.
- Material: Titanium alloy TA6V.

Reference	Length
112111(S)	Spin [®] Dia. 2.7 x L. 11 mm
112112(S)	Spin [®] Dia. 2.7 x L. 12 mm
112113(S)	Spin [®] Dia. 2.7 x L. 13 mm
112114(S)	Spin [®] Dia. 2.7 x L. 14 mm



(S) means available in sterile and non sterile

UNI-CLIP® Notched Staples (for associated instrumentation, see p. 9)

For fixation of bone fractures or for bone reconstruction.

- Range sizes:
 - > Interaxis: 11, 12, 13, 15 and 20 mm.
 - > Leg lengths: 12, 13, 14, 15, 16, 17 and 20 mm.
- Material: Stainless steel 316L.



Reference	Description
213113(S)	Uni-Clip® interaxis 11.0 x L. 13 mm
213114(S)	Uni-Clip® interaxis 11.0 x L. 14 mm
213115(S)	Uni-Clip® interaxis 11.0 x L. 15 mm
213116(S)	Uni-Clip® interaxis 11.0 x L. 16 mm
213117(S)	Uni-Clip® interaxis 11.0 x L. 17 mm
213213(S)	Uni-Clip® interaxis 12.0 x L. 13 mm
213214(S)	Uni-Clip® interaxis 12.0 x L. 14 mm
213215(S)	Uni-Clip® interaxis 12.0 x L. 15 mm
213216(S)	Uni-Clip® interaxis 12.0 x L. 16 mm
213217(S)	Uni-Clip® interaxis 12.0 x L. 17 mm
213313(S)	Uni-Clip® interaxis 13.0 x L. 13 mm
213314(S)	Uni-Clip® interaxis 13.0 x L. 14 mm
213315(S)	Uni-Clip® interaxis 13.0 x L. 15 mm
213316(S)	Uni-Clip® interaxis 13.0 x L. 16 mm
213317(S)	Uni-Clip® interaxis 13.0 x L. 17 mm
213512(S)	Uni-Clip® interaxis 15.0 x L. 12 mm
213820(S)	Uni-Clip® interaxis 20.0 x L. 20 mm

SOLUSTAPLE® Standard Staple (for associated instrumentation, see p. 9)

Static staple for fixation of mono cortical phalangeal osteotomies.

- Two shapes of staples (26° and 90°) are available.
- Self-perforating.
- Material: stainless steel 316L.

Reference	Description
114002(S)	Solustaple® 90° interaxis 8 mm
114004(S)	Solustaple® 90° interaxis 10 mm
114023(S)	Solustaple® 26° interaxis 8 mm
114025(S)	Solustaple® 26° interaxis 10 mm



(S) means available in sterile and non sterile

Osteoguides

- Suggested use:
 - > Scarf osteotomy.
 - > Chevron osteotomy.
 - > Basal osteotomy.
 - > Arthrodesis.
 - > Wedge osteotomies (Akin...).
- Details:
 - > Hole for K-Wire: Diameter 1.0 to 1.6 mm.
 - > Knurled handle for optimal grip.
 - > Blade thickness maxi 0.8 mm.
 - > Positioning tip.

Reference	Description
229011	Longitudinal
229012	Transversal



K-FIX[®] – K-Wire End Protector

The K-Fix[®] pin protector is a simple system to protect patients or hospital staff from K-Wires sharp ends.

- Delivered sterile.

Reference	Length
119050	K-FIX [®] for K-WIRE - Dia. 0.8 mm to 2.0 mm (5 units per pack)



Weil Clamp

Reference	Length
119106	Weil Clamp

B-BOP® LOCK Plate

B-BOP®-LOCK plate is indicated for fixation of osteotomy of the basis of the first metatarsal. Examples include:

- Moderate to severe Hallux valgus.
- Hallux varus.


Implants

Reference	Description
320 100S	Basal plantar plate B-Bop® Lock right
320 200S	Basal plantar plate B-Bop® Lock left
295 110S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 10 mm
295 112S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 12 mm
295 114S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 14 mm
295 116S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 16 mm
295 118S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 18 mm
295 120S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 20 mm
295 122S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 22 mm
295 124S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 24 mm
295 126S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 26 mm
295 128S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 28 mm
295 130S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 30 mm
295 132S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 32 mm
295 134S	Screw and Lock screw titanium B-Bop® Lock Dia 3.0 x L 34 mm
195 100S	Lock screw B-Bop® Lock



Instrumentation

Reference	Description
219 027	Drilling guide Surfix Alpha Dia 2.0 mm
219 127	Screwdriver Torx T7
219 227	Screwdriver AO blade Torx T7
299 050	Drill AO Dia 2.0 mm
299 005	Drill Dia 2.0 mm
329 500	Bending forceps
329 002	Forceps
115 070	K-wire Dia 1.0 mm x L 70 mm
115 116	K-wire Dia 1.6 mm x L 150 mm
299 020	Depth gauge B-Bop Lock®
189 210	Wedge 10°
189 215	Wedge 15°
189 220	Wedge 20°



Container

Reference	Description
329 900	Container
including:	
329 901	Basis
329 902	Silicone mat
996200	Lid

IPP-ON[®] Interphalangeal Arthrodesis Implant

IPP-ON[®] interphalangeal implant is an implant which respects anatomical flexion to achieve interphalangeal arthrodesis of the lesser toes in case of:

- Fixed or semi-fixed hammer toe.
- Revision in case of arthrodesis or arthroplasty failure.
- Shortening in case of length excess of the second toe.
- Material: Stainless steel: 316 LVM.
- 2 sizes.




Reference	Description
230001S	IPP-ON [®] Interphalangeal Implant Size 1
230002S	IPP-ON [®] Interphalangeal Implant Size 2

IPP-ON[®] Instruments

Reference	Description
239011	Distal Hand Drill Size 1
239012	Proximal Hand Drill Size 1
239021	Distal Hand Drill Size 2
239022	Proximal Hand Drill Size 2
239500	Clamp
239030	Implant Holder

Sterilization Container	
239000	IPP-ON [®] sterilization container including:
239001	Base
996100	Lid
278902	Mat
119909	Silicone wedge



METIS® Metatarso-Phalangeal Prosthesis

The three-component METIS® prosthesis replaces the articular surfaces of the great toe joint.

- Metatarsal component:
 - > 7 sizes (0 to 6).
 - > Material: Cobalt-Chromium alloy (CoCr).
 - > Coating: Porous Titanium & Hydroxyapatite (Ti and HAP).
- PE inlay:
 - > Ultra High Molecular Weight Polyethylene (UHMWPE).
 - > 5 sizes (1 to 5), 3 thicknesses (1 to 3).
- Phalangeal implant:
 - > 5 sizes (1 to 5).
 - > Material: Titanium alloy TA6V.
 - > Coating: Porous Titanium (Ti).



Reference	Description
250100	Metatarsal component size 0
250101	Metatarsal component size 1
250102	Metatarsal component size 2
250103	Metatarsal component size 3
250104	Metatarsal component size 4
250105	Metatarsal component size 5
250106	Metatarsal component size 6
250201	Phalangeal component size 1
250202	Phalangeal component size 2
250203	Phalangeal component size 3

Reference	Description
250204	Phalangeal component size 4
250205	Phalangeal component size 5
250311	METIS® Polyethylene Inlay - Size 1 - Thickness 1
250312	METIS® Polyethylene Inlay - Size 1 - Thickness 2
250313	METIS® Polyethylene Inlay - Size 1 - Thickness 3
250321	METIS® Polyethylene Inlay - Size 2 - Thickness 1
250322	METIS® Polyethylene Inlay - Size 2 - Thickness 2
250323	METIS® Polyethylene Inlay - Size 2 - Thickness 3
250331	METIS® Polyethylene Inlay - Size 3 - Thickness 1
250332	METIS® Polyethylene Inlay - Size 3 - Thickness 2
250333	METIS® Polyethylene Inlay - Size 3 - Thickness 3
250341	METIS® Polyethylene Inlay - Size 4 - Thickness 1
250342	METIS® Polyethylene Inlay - Size 4 - Thickness 2
250343	METIS® Polyethylene Inlay - Size 4 - Thickness 3
250351	METIS® Polyethylene Inlay - Size 5 - Thickness 1
250352	METIS® Polyethylene Inlay - Size 5 - Thickness 2
250353	METIS® Polyethylene Inlay - Size 5 - Thickness 3


METIS® Instruments

METIS® Instruments

Reference	Description
139107	Metatarsal Osteotome
139120	Implant holder
139131	Locking screw
139210	Phalangeal reamer starter size
139211–215	Phalangeal reamer size 1 to 5
139140	Metatarsal sizer size 0
139121–126	Metatarsal sizer size 1 to 6
139150	Metatarsal cutting guide size 0
139111–116	Metatarsal cutting guide size 1 to 6
139216	Metatarsal impactor
139128	Visualization guide
139129	Visualization guide wire
119645	Screwdriver hexa 3.5 mm
119 601	Screwdriver hexa 2.5 mm
115216	K-Wire Dia. 1.6 x L. 150 mm
119611	Drill Dia. 2.5 x L. 150 mm
139940	METIS® sterilization container including:
139941	Metatarsal base
139942	Phalangeal base
139943	Base
139953	Lid
119909	Silicone wedge

METIS® Trial Implants

Reference	Description
139100–106	Metatarsal trial implant size 0 to 6
139311–313	Trial inlay size 1 thickness 1 to 3
139321–323	Trial inlay size 2 thickness 1 to 3
139331–333	Trial inlay size 3 thickness 1 to 3
139341–343	Trial inlay size 4 thickness 1 to 3
139351–353	Trial inlay size 5 thickness 1 to 3
139201–205	Phalangeal trial implant size 1 to 5



K2™ Hemi Toe Implant


The K2™ Hemi Toe Implant System is an anatomically designed, metallic, one-piece surface replacement system for the base of the hallux. 4 sizes are available.

K2™ Implants

Reference	Description
201010	Phalangeal, Size 1
201020	Phalangeal, Size 2
201030	Phalangeal, Size 3
201040	Phalangeal, Size 4

K2™ Instruments

Reference	Description
202010	Size 1 Trial
202020	Size 2 Trial
202030	Size 3 Trial
202040	Size 4 Trial
203010	Size 1 Shaper
203020	Size 2 Shaper
203030	Size 3 Shaper
203040	Size 4 Shaper
203100	Cutting Guide
203300	Sizer/Drill Guide
203500	Broach Handle
203600	Impactor
205000	Sterilization tray
203200	Reamer
203415	1.5 mm 80 mm drill bit



Movement Great Toe System

- Hemi Proximal Phalanx Implant.
- Hemi Metatarsal Implant.
- Total Great Toe Implant.

Reference	Description
MGT-890-10PPH	Hemi Proximal Phalanx, Size 10
MGT-890-20PPH	Hemi Proximal Phalanx, Size 20
MGT-890-30PPH	Hemi Proximal Phalanx, Size 30
MGT-890-40PPH	Hemi Proximal Phalanx, Size 40
MGT-890-10MT	Hemi Metatarsal, Size 10
MGT-890-20MT	Hemi Metatarsal, Size 20
MGT-890-30MT	Hemi Metatarsal, Size 30
MGT-890-40MT	Hemi Metatarsal, Size 40
MGT-890-10PPT	Phalangeal Total, Size 10
MGT-890-20PPT	Phalangeal Total, Size 20
MGT-890-30PPT	Phalangeal Total, Size 30
MGT-890-40PPT	Phalangeal Total, Size 40

Reference	Description
INS-890-00	Instrument Set
MIS-890-00	Metatarsal Implant Sizer
PIS-890-00	Proximal Phalanx Implant Sizer
GDW-890-00	Guide Pin, 2.0 mm X 125 mm
MSR-890-10/20	Metatarsal Surface Reamer, Size 10/20
MSR-890-30/40	Metatarsal Surface Reamer, Size 30/40
PSR-890-10/20	Proximal Phalanx Surface Reamer, Size 10/20
PSR-890-30/40	Proximal Phalanx Surface Reamer, Size 30/40
DRL-890-00	Cannulated Drill, 4.5 mm
DCG-890-00	Metatarsal Dorsal Cutting Guide
TRL-890-10MT	Metatarsal Trial, Size 10
TRL-890-20MT	Metatarsal Trial, Size 20
TRL-890-30MT	Metatarsal Trial, Size 30
TRL-890-40MT	Metatarsal Trial, Size 40
TRL-890-10PPH	Proximal Phalanx Hemi Trial, Size 10
TRL-890-20PPH	Proximal Phalanx Hemi Trial, Size 20
TRL-890-30PPH	Proximal Phalanx Hemi Trial, Size 30
TRL-890-40PPH	Proximal Phalanx Hemi Trial, Size 40
TRL-890-10PPT	Phalangeal Total Trial, Size 10
TRL-890-20PPT	Phalangeal Total Trial, Size 20
TRL-890-30PPT	Phalangeal Total Trial, Size 30
TRL-890-40PPT	Phalangeal Total Trial, Size 40
IMP-890-00MT	Metatarsal Impactor
IMP-890-00PP	Proximal Phalanx Impactor
ALG-890-00	Sizer Alignment Guide



HALLU®-FIX MTP Arthrodesis System

The HALLU®-FIX solution is a system dedicated to first MTP arthrodesis.

- Pre-bent plate.
- SNAP-OFF® screw fixation.
- Color coded.
- 3 sizes (3 left / 3 right).
- Material: Titanium alloy TA6V.

HALLU®-C Plate

HALLU®-C Plates

Reference	Description
117340(S)	HALLU®-C Plate Right-Size 1: L. 40 mm - 4 holes
117345(S)	HALLU®-C Plate Right-Size 2: L. 45 mm - 5 holes
117350(S)	HALLU®-C Plate Right-Size 3: L. 50 mm - 6 holes
117440(S)	HALLU®-C Plate Left-Size 1: L. 40 mm - 4 holes
117445(S)	HALLU®-C Plate Left-Size 2: L. 45 mm - 5 holes
117450(S)	HALLU®-C Plate Left-Size 3: L. 50 mm - 6 holes



HALLU®-S Plate

HALLU®-S Plates

Reference	Description
117150(S)	HALLU®-S Plate Right-Size 1: L. 45 mm - 3 holes
117155(S)	HALLU®-S Plate Right-Size 2: L. 50 mm - 4 holes
117160(S)	HALLU®-S Plate Right-Size 3: L. 55 mm - 5 holes
117250(S)	HALLU®-S Plate Left-Size 1: L. 45 mm - 3 holes
117255(S)	HALLU®-S Plate Left-Size 2: L. 50 mm - 4 holes
117260(S)	HALLU®-S Plate Left-Size 3: L. 55 mm - 5 holes



HALLU® SNAP-OFF® Screws

- 2 diameters: 2.7 mm & 3.0 mm (for porous bone).
- Self-tapping.
- Color code for length identification.
- Material: Titanium alloy TA6V.

HALLU® Snap-Off® screws Dia. 3.0 mm

Reference	Length
117110(S)	10 mm
117112(S)	12 mm
117114(S)	14 mm
117116(S)	16 mm
117118(S)	18 mm

HALLU® Snap-Off® screws Dia. 2.7 mm

Reference	Length
117010(S)	10 mm
117012(S)	12 mm
117014(S)	14 mm
117016(S)	16 mm
117018(S)	18 mm
117020(S)	20 mm
117022(S)	22 mm
117024(S)	24 mm
117026(S)	26 mm
117028(S)	28 mm
117030(S)	30 mm
117032(S)	32 mm
117034(S)	34 mm



HALLU®-FIX Instruments



HALLU®-Ream Metatarsal

Reference	Diameter
129714	Dia. 14 mm
129716	Dia. 16 mm
129718	Dia. 18 mm
129720	Dia. 20 mm
129722	Dia. 22 mm

HALLU®-Ream Phalangeal

Reference	Diameter
129724	Dia. 14 mm
129726	Dia. 16 mm
129728	Dia. 18 mm
129730	Dia. 20 mm
129752	Dia. 22 mm

HALLU® Associated instruments

Reference	Description
129710	Quick coupling
129733	SNAP-OFF® Cannulated screwdriver - Dia. 2.7 mm and 3.0 mm
129731	Plate Bender Left
129732	Plate Bender Right
129734	Drilling Guide Dia. 1.9 mm
129735	Screwdriver tip
129736	Depth gauge
119618	Drill Dia. 1.9 x L. 100 mm
115216	K-Wire Dia. 1.6 x L. 150 mm
115100	K-Wire Dia. 1.0 x L. 100 mm



Sterilization container

129930	HALLU®-FIX Sterilization Container including:	
129931	HALLU®-FIX Basis	
129933	HALLU®-S Plate Insert	
129932	HALLU®-C Plate Insert	
996200	Lid	
129936	HALLU®-FIX screws rack	
119909	Silicone wedge (Optional)	

(S) means available in sterile and non sterile

HALLU®-LOCK MTP Arthrodesis System

The HALLU®-LOCK solution is a locked system dedicated to first MTP arthrodesis:

Pre-bent plate.

- Color coded.
- 3 sizes (3 left / 3 right).
- Material: Titanium alloy TA 6V.
- Monobloc fixation with Surfix® locking system.

Surfix® locking screw

- 2 diameters: 2.7 mm & 3.0 mm (for porous bone).
- Standard and orientable screws.
- Color code for length identification.
- Material: Titanium alloy TA 6V.
- Sterily Packaged.



Surfix® standard screw
Dia. 2.7 mm + lock-screw



Reference	Length
285210S	10 mm
285212S	12 mm
285214S	14 mm
285216S	16 mm
285218S	18 mm
285220S	20 mm
285222S	22 mm
285224S	24 mm
285226S	26 mm
185200S	lock-screw

Surfix® standard screw
Dia. 3.0 mm + lock-screw



Reference	Length
285110S	10 mm
285112S	12 mm
285114S	14 mm
285116S	16 mm
285118S	18 mm
285120S	20 mm
285122S	22 mm
285124S	24 mm
285126S	26 mm
185100S	lock-screw

Surfix® Alpha
variable angle screw
Dia. 2.7 mm + lock-screw



Reference	Length
295210S	10 mm
295212S	12 mm
295214S	14 mm
295216S	16 mm
295218S	18 mm
295220S	20 mm
295222S	22 mm
295224S	24 mm
295226S	26 mm
195200S	lock-screw

Surfix® Alpha
variable angle screw
Dia. 3.0 mm + lock-screw



Reference	Length
295110S	10 mm
295112S	12 mm
295114S	14 mm
295116S	16 mm
295118S	18 mm
295120S	20 mm
295122S	22 mm
295124S	24 mm
295126S	26 mm
195100S	lock-screw

Hallu®-Lock C plate 10°

Reference	Description
290340S	Right-Size 1: L.40 mm – 4 holes
290345S	Right-Size 2: L.45 mm – 5 holes
290350S	Right-Size 3: L.50 mm – 6 holes
290440S	Left-Size 1: L.40 mm – 4 holes
290445S	Left-Size 2: L.45 mm – 5 holes
290450S	Left-Size 3: L.50 mm – 6 holes



Hallu®-Lock S plate 10°
Number of holes excludes the 4 extremity holes

Reference	Description
290150S	Right-Size 1: L.45 mm – 3 holes
290155S	Right-Size 2: L.50 mm – 3 holes
290160S	Right-Size 3: L.55 mm – 3 holes
290250S	Left-Size 1: L.45 mm – 3 holes
290255S	Left-Size 2: L.50 mm – 3 holes
290260S	Left-Size 3: L.55 mm – 3 holes



QWIX® transarticular
compressive screw Dia.
3.0 mm

Reference	Length
111328S	28 mm
111330S	30 mm
111332S	32 mm
111334S	34 mm





HALLU®-LOCK Instruments

HALLU®-REAM System

129714	Metatarsal Reamer Dia. 14 mm
129716	Metatarsal Reamer Dia. 16 mm
129718	Metatarsal Reamer Dia. 18 mm
129720	Metatarsal Reamer Dia. 20 mm
129722	Metatarsal Reamer Dia. 22 mm
129724	Phalangeal Reamer Dia. 14 mm
129726	Phalangeal Reamer Dia. 16 mm
129728	Phalangeal Reamer Dia. 18 mm
129730	Phalangeal Reamer Dia. 20 mm
129752	Phalangeal Reamer Dia. 22 mm

HALLU®-Lock associated instruments

129710	Quick coupling
997100	Quick-Coupling AO
299070	Left bender
299080	Right bender
219027	SURFIX® Alpha drilling guide Dia. 2.0 mm
219627	SURFIX® drilling guide Dia. 2.0 mm
219127	Screwdriver torx T7
219227	AO Screwdriver tip torx T7
299020	Depth Gauge
119133	Handle with AO Attachment Dia. 3.0 mm
119135	AO Screwdriver tip dia 3.0 mm
299010	Trial plate holder
299150	Trial plate-S Right - Size 1
299155	Trial plate-S Right - Size 2
299160	Trial plate-S Right - Size 3
299250	Trial plate-S Left - Size 1
299255	Trial plate-S Left - Size 2
299260	Trial plate-S Left - Size 3
299340	Trial plate-C Right - Size 1
299345	Trial plate-C Right - Size 2
299350	Trial plate-C Right - Size 3
299440	Trial plate-C Left - Size 1
299445	Trial plate-C Left - Size 2
299450	Trial plate-C Left - Size 3
299005	Drill Dia. 2.0 x L. 100 mm
299050	AO Drill AO Dia. 2.0 x L. 100 mm
159027S	2 in 1 cannulated drill AO dia 2.2/3.0 x L. 12 mm sterile
115070	K-Wire Dia. 1.0 x L. 70 mm
115100	K-wire Dia. 1.0 x L. 100 mm
115116	K-Wire Dia. 1.6 x L. 150 mm



HALLU®-LOCK sterilization container

299900	HALLU®-Lock sterilization container
299901	Base
299903	Trial implants tray
299904	K-Wire tray
996200	Lid
119909	Silicone wedge

HALLU®-REAM Set

Reference	Description
129714	Metatarsal Reamer Dia. 14 mm
129716	Metatarsal Reamer Dia. 16 mm
129718	Metatarsal Reamer Dia. 18 mm
129720	Metatarsal Reamer Dia. 20 mm
129722	Metatarsal Reamer Dia. 22 mm

129724	Phalangeal Reamer Dia. 14 mm
129726	Phalangeal Reamer Dia. 16 mm
129728	Phalangeal Reamer Dia. 18 mm
129730	Phalangeal Reamer Dia. 20 mm
129752	Phalangeal Reamer Dia. 22 mm

Instruments

129710	Quick coupling
115216	K-Wire Dia. 1.6 x L. 150 mm

Sterilization container

129935	HALLU®-Ream sterilization container including:
129938	HALLU®-Ream Basis
996100	Lid



DPR System Minimal Invasive Foot Surgery

DPR System minimal invasive foot surgery was born from our collaboration with minimal invasive foot surgeons.

- Single use: Sterile packaging.
- Safety: Protection by silicone ring.
- Identification: Color coded.
- Material: Stainless Steel 316L.



DPR System Minimal Invasive Foot Surgery burs

278002S	Burr straight Dia. 2.0 x L. 15 mm
278007S	Burr conical Dia. 3.1/1.0 x L. 12 mm
278005S	Burr straight Dia. 1.9 x L. 10 mm
278008S	Burr conical Dia. 4.1/1.0 x L. 12 mm
278009S	Burr cylindrical Dia. 4.1 x L. 15 mm
278010S	Burr cylindrical Dia. 3.1 x L. 15 mm

DPR System Minimal Invasive Foot Surgery curettes

278121	Curette Medium - Down
278122	Curette Medium - Up
278123	Curette Large - Down
278124	Curette Large - Up
278100	Sleeve protection Cap
278110	Sleeve protection Dia. 3.1 mm
278109	Sleeve protection Dia. 4.1 mm

Sterilization container

278900	DPR System Minimal Invasive Foot Surgery sterilization container
278901	Base
996100	Lid
278902	Mat
119909	Silicone wedge



Please contact us if you need a MIS motor offer.


Post-Operative Shoes

Reference	Description
Flat Post Operative Shoe	
999001	Small
999002	Medium
999003	Large
999004	X-Large



Reference	Description
Wedged Heel Post Operative Shoe	
999100	Size 4 to 6.5
999101	Size 7 to 8.5
999102	Size 9 to 10.5





Integra® Mid & HindFoot Solutions

QWIX® Fixation Screw, dia. 4.3 mm

For the fixation of bone fractures or for bone reconstruction.

- Self-drilling and self-tapping.
- Intra-osseous fixation.
- Variable lag length for extremity fixation.
- Material: Titanium alloy TA6V.

QWIX® dia. 4.3 mm

Reference	Length
111422S	L. 22 mm
111424S	L. 24 mm
111426S	L. 26 mm
111428S	L. 28 mm
111430S	L. 30 mm
111432S	L. 32 mm
111434S	L. 34 mm
111436S	L. 36 mm
111438S	L. 38 mm
111440S	L. 40 mm
111445S	L. 45 mm
111450S	L. 50 mm
111455S	L. 55 mm
111460S	L. 60 mm




QWIX® dia. 4.3 mm Instruments

Reference	Description
119104	Screwdriver tip
119139	Short Drill dia. 2.6/5.0 x L. 100 mm
119138	Quick coupling handle
119516	Long cannulated drill dia. 2.6 x L. 140 mm
115116	K-Wire dia. 1.6 x L. 150 mm

Sterilization container

Reference	Description
119935	QWIX® 4.3 sterilization container including:
119937	Cylinder
119936	Base
119938	Mat
996400	Lid



Large QWIX® Fixation and Positioning Screws dia. 5.5 mm / 7.5 mm

- Fixation and positioning screws dia. 5.5 and 7.5 mm.
- Self-drilling and self-tapping.
- Intra-osseous fixation.
- Material: Titanium alloy TA6V.
- Sterile packaged.

Large QWIX® Fixation screws
dia. 7.5 mm



Reference	Length
111740s	L. 40 mm
111745s	L. 45 mm
111750s	L. 50 mm
111755s	L. 55 mm
111760s	L. 60 mm
111765s	L. 65 mm
111770s	L. 70 mm
111775s	L. 75 mm
111780s	L. 80 mm
111785s	L. 85 mm
111790s	L. 90 mm
111795s	L. 95 mm
111800s	L. 100 mm
111805s	L. 105 mm
111810s	L. 110 mm
111815s	L. 115 mm
111820s	L. 120 mm



Large QWIX® Fixation screws
dia. 5.5 mm

Reference	Length
111530s	L. 30 mm
111535s	L. 35 mm
111540s	L. 40 mm
111545s	L. 45 mm
111550s	L. 50 mm
111555s	L. 55 mm
111560s	L. 60 mm
111565s	L. 65 mm
111570s	L. 70 mm
111575s	L. 75 mm
111580s	L. 80 mm

Large QWIX® Positioning screws
dia. 7.5 mm



Reference	Length
121740s	L. 40 mm
121745s	L. 45 mm
121750s	L. 50 mm
121755s	L. 55 mm
121760s	L. 60 mm
121765s	L. 65 mm
121770s	L. 70 mm
121775s	L. 75 mm
121780s	L. 80 mm
121785s	L. 85 mm
121790s	L. 90 mm
121795s	L. 95 mm
121800s	L. 100 mm
121805s	L. 105 mm
121810s	L. 110 mm
121815s	L. 115 mm
121820s	L. 120 mm




Large QWIX® Positioning screws
dia. 5.5 mm

Reference	Length
121530s	L. 30 mm
121535s	L. 35 mm
121540s	L. 40 mm
121545s	L. 45 mm
121550s	L. 50 mm
121555s	L. 55 mm
121560s	L. 60 mm
121565s	L. 65 mm
121570s	L. 70 mm
121575s	L. 75 mm
121580s	L. 80 mm

Large QWIX® Instruments

Reference	Length
119150	Measurer 200 and 250 mm
119151	Internal protection sleeve dia. 1.8 mm
119152	Drilling guide dia. 4.5 mm
119153	External protection sleeve dia. 7.3 mm
119154	Cannulated drill dia. 1.8/4.2 x L. 180 mm
119155	Cannulated 2 in 1 drill dia. 4.2/5.8 x L. 113 mm
119156	Hexa screwdriver Tip 3 mm
119171	Internal protection sleeve dia. 2.7 mm
119172	Drilling guide dia. 5.8 mm
119173	External protection sleeve dia. 9.3 mm
119174	Cannulated drill dia. 2.7/5.5 x L. 210 mm
119175	Cannulated 2 in 1 drill dia. 5.5/7.0 x L. 128 mm
119176	Hexa screwdriver Tip 4 mm
119177	T handle quick coupling
115516	K-Wire dia. 1.6 x L. 200 mm
529061	K-Wire dia. 2.5 x L. 250 mm



Sterilization container

Reference	Description
121950	LARGE QWIX® sterilization container including:
121951	Base
996200	Lid
119170	Cylinder
119909	Silicone wedge

I.CO.S.® Ideal Compression Screw

Cannulated Titanium 4.0 mm and 6.5 mm diameter compression screw are dedicated to fixation of osteotomies or for fractures stabilization.

- K-Wire guidance for drilling, tapping and screw setting (dia. 1.6 or 2.5 mm).
- Adjustable and controlled static compression.
- Dynamic compression by additional head screwing.
- Intra osseous implant.
- Large range of sizes :
 - > I.CO.S.® dia. 4.0 mm:
 - » L. 26 to 40 mm by 2 mm
 - » L. 40 to 60 mm by 5 mm
 - > I.CO.S.® dia. 6.5 mm:
 - » L. 40 to 90 mm by 5 mm
- Material: Titanium alloy TA6V.
- Cortico-spongyous thread.
- Long non threaded shaft (lag screws) allowing compression.



I.CO.S.® dia. 4.0 mm

Reference	Length
105426	L. 26 mm
105428	L. 28 mm
105430	L. 30 mm
105432	L. 32 mm
105434	L. 34 mm
105436	L. 36 mm
105438	L. 38 mm
105440	L. 40 mm
105445	L. 45 mm
105450	L. 50 mm
105455	L. 55 mm
105460	L. 60 mm



I.CO.S.® dia. 6.5 mm


Reference	Length
105640	L. 40 mm
105645	L. 45 mm
105650	L. 50 mm
105655	L. 55 mm
105660	L. 60 mm
105665	L. 65 mm
105670	L. 70 mm
105675	L. 75 mm
105680	L. 80 mm
105685	L. 85 mm
105690	L. 90 mm

I.CO.S.® dia. 4.0 mm Instruments

Reference	Description
119501	Protection Sleeve dia. 6 mm
119506	Drilling Guide dia. 1.8 mm
119510	Measurer
119516	Long cannulated drill dia. 2.6 x L. 140 mm
119522	Short cannulated Drill dia. 5.0 x L. 110 mm
129710	Quick-coupling
119525	Tap
119531	Internal screwdriver dia. 4.0 x L. 100 mm
119532	Screwdriver tip dia. 4.0 x L. 100 mm
119536	External screwdriver dia. 6.4 x L. 170 mm

I.CO.S.® 4.0 sterilization Container

Reference	Description
119974	I.CO.S.® 4.0 sterilization container including:
119956	Base
119951	Screw rack
119949	Metallic lid for screws (Optional)
996200	Lid
119909	Silicone wedge



I.CO.S.® dia. 6.5 mm Instruments

Reference	Description
119541	Drilling Guide dia. 6.5 x L. 100 mm
119551	Internal protection sleeve dia. 6.5 mm
119552	External protection sleeve dia. 6.5 mm
119576	Long cannulated drill dia. 4.6 x L. 180 mm
119546	Short cannulated drill dia. 6.5 x L. 118 mm
119561	Internal screwdriver dia. 6.5 x L. 170 mm
119562	Screwdriver tip dia. 6.5 x L. 110 mm
119566	External screwdriver dia. 10.1 x L. 170 mm
119555	Tap
119510	Measurer
129710	Quick coupling
115225	K-Wire dia. 2.5 x L. 200 mm

I.CO.S.® 6.5 sterilization Container

Reference	Description
119976	I.CO.S.® 6.5 sterilization Container including:
119958	Base
119953	Screws rack
119949	Metallic lid for screws (Optional)
119957	Lid
119909	Silicone wedge



I.CO.S.[®] Removal Set

This removal set contains instrumentation required to remove I.CO.S.[®] screws dia. 4.0 and 6.5 mm

Removal set Instruments (for all diameters)

Reference	Description
Instruments for I.CO.S.[®] dia. 4.0 mm	
115116	K-Wire dia. 1.6 x L. 150 mm
119531	Internal screwdriver dia. 4.0 x L. 100 mm
119580	Trephine dia. 4.0 x L. 70 mm
119582	Internal Extractor dia. 4.0 x L. 60 mm (single use)
Instruments for I.CO.S.[®] dia. 6.5 mm	
Reference	Description
115225	K-Wire dia. 2.5 x L. 200 mm
119561	Internal screwdriver dia. 6.5 x L. 170 mm
119590	Trephine dia. 6.5 x L. 80 mm
119592	Internal Extractor dia. 6.5 x L. 70 mm (single use)

Sterilization container

Reference	Description
119975	I.CO.S. [®] Removal Sterilization container including
119959	Lid
119955	Base
119977	K-Wire cylinder



KALIX® II Flat Foot Implant

Expansion endorthesis for flat foot correction in child or adult.

Progressive expansion.

- Material:
 - > Titanium alloy TA6V
 - > Polyethylene (UHMWPE)
- Color coded.
- Sterily packaged.

Reference	Ext. diameter
141009	dia. 9 mm
141010	dia. 10 mm
141011	dia. 11 mm
141012	dia. 12 mm
141013	dia. 13 mm
141014	dia. 14 mm
141015	dia. 15 mm
141017	dia. 17 mm



KALIX® II Flat Foot Implant Instruments

Reference	Description
119825	Viladot Lever
119850	Impactor
119831	Trial implant holder
119859	KALIX® II Trial implant dia. 9.0 mm
119860	KALIX® II Trial implant dia. 10.0 mm
119861	KALIX® II Trial implant dia. 11 mm
119862	KALIX® II Trial implant dia. 12 mm
119863	KALIX® II Trial implant dia. 13 mm
119864	KALIX® II Trial implant dia. 14 mm
119865	KALIX® II Trial implant dia. 15 mm
119867	KALIX® II Trial implant dia. 17 mm



KALIX® II sterilization container

Reference	Description
119880	KALIX® II sterilization container including
119882	Lid
119881	Base

MBA™ Titanium Subtalar Implant

The Subtalar MBA™ System is indicated for use in treatment of the hyperpronated foot and stabilization of the subtalar joint. It is designed to block the posterior and inferior displacement of the talus, thus allowing normal subtalar joint motion while blocking excessive pronation and the resulting sequela.

- Severely pronated foot.
- Walking intemperance.
- Calcaneal stance position greater than 5°.
- Manually correctable deformities.
- Mid-tarsal breech (arch pain).
- Forefoot varus greater than 10°.



Removal set Instruments (for all diameters)

Reference	Description
Instruments tray	
05-0016	Guide pin insert holder
05-0206	Trial 6 mm
05-0208	Trial 8 mm
05-0209	Trial 9 mm
05-0210	Trial 10 mm
05-0212	Trial 12 mm
05-1002	Yellow probe
05-1006	6 mm Sizer
05-1008	8 mm Sizer
05-1009	9 mm Sizer
05-1010	10 mm Sizer
05-1012	12 mm Sizer
05-1020	Insertion Device
05-0017	Guide pin 2 mm
05-5001	Sterilization tray
Implants	
05-0106	MBA™ 6 mm implant
05-0108	MBA™ 8 mm implant
05-0109	MBA™ 9 mm implant
05-0110	MBA™ 10 mm implant
05-0112	MBA™ 12 mm implant



ADVANSYS® Dorsal Lisfranc Plate

The Dorsal Lisfranc Plates are intended for fractures, fusions, osteotomies and replantations of small bones at the tarso-metatarsal joints (Lisfranc joints).

- 3 sizes of implants:
3 right & 3 left plates (small, medium and large).
- 3.5 mm diameter screws available in lengths from 8 mm to 34 mm in 2 mm increments.
- 3.5 mm lock-screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.



ADVANSYS® Dorsal Lisfranc Plates

Reference	Description
181021S	ADVANSYS® DLP small-right
181022S	ADVANSYS® DLP medium-right
181023S	ADVANSYS® DLP large-right
181031S	ADVANSYS® DLP small-left
181032S	ADVANSYS® DLP medium-left
181033S	ADVANSYS® DLP large-left



Surfix® Locking screw Sterile Titanium dia.3.5 mm + lock-screw

Reference	Description
285310S	L. 10 mm
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
185300S	Lock-screw dia.3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

ADVANSYS® Medial Lisfranc Plate

The Medial Lisfranc Plates are intended to be used for bone fixation such as:

- Lisfranc arthrodesis.
- Mono or bi-cortical osteotomies or fractures near the 1st metatarso-cuneiform joint.
- 2 sizes of implants:
2 right & 2 left plates (small and large sizes).
- 3.5 mm diameter screws available in lengths from 8 mm to 34 mm in 2 mm increments.
- 3.5 mm lock-screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.
- Sterily packaged.



ADVANSYS® Medial Lisfranc Plate

Reference	Description
181041S	ADVANSYS® MLP small-right
181042S	ADVANSYS® MLP large-right
181051S	ADVANSYS® MLP small-left
181052S	ADVANSYS® MLP large-left



Surfix® Locking screw Sterile Titanium dia. 3.5 mm + lock-screw

Reference	Description
285310S	L. 10 mm
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
185300S	Lock-screw dia.3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

ADVANSYS® Medial Lisfranc Plate Instruments

Reference	Description
115100	K-Wire dia. 1.0 x L. 100 mm
219335	Depth gauge hole dia. 3.5 mm
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219535	AO drill dia. 2.7 x L. 125 mm
219635	Drilling guide dia. 2.7 mm
219735	Bending forceps L. 171 hole - dia 3.5 mm
219835	Screwdriver hexa 2.0 x L. 180 mm
169041	MLP trial implant right-small
169042	MLP trial implant right-large
169051	MLP trial implant left-small
169052	MLP trial implant left-large

ADVANSYS® Midfoot sterilization container


Reference	Description
169102	ADVANSYS® Midfoot sterilization container includes:
169112	Base
996200	Lid
169106	Cylinder
169107	DLP/MLP Trial implant module
119909	Silicone wedge

DLP trial plates

Reference	Description
169021	Dorsal trial plate, right, small
169022	Dorsal trial plate, right, medium
169023	Dorsal trial plate, right, large
169031	Dorsal trial plate, left, small
169032	Dorsal trial plate, left, medium
169033	Dorsal trial plate, left, large

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm



UNI-CP™ Compression plate

The UNI-CP™ compression plate is indicated for fixation of bone fractures or for bone reconstruction:

- 4 plate configurations:
 - > 2 hole design (17, 20, 25, and 30 mm interaxis)
 - > 4 hole design (20, 25, and 30 mm interaxis)
 - > 4 hole T-shape design (20 mm interaxis)
 - > 4 hole U-shape design (17x20, 19x22.5, 21x25 mm interaxis)
- Range of screw lengths: from 12 mm to 40 mm in 2 mm increments.
- Surfix® Locking Technology.
- Material
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.

UNI-CP™ Compression Plate

Reference	Description
330430S	UNI-CP™ 4 holes - 30 mm
330425S	UNI-CP™ 4 holes - 25 mm
330420S	UNI-CP™ 4 holes - 20 mm
330230S	UNI-CP™ 2 holes - 30 mm
330225S	UNI-CP™ 2 holes - 25 mm
330220S	UNI-CP™ 2 holes - 20 mm
330217S	UNI-CP™ 2 holes - 17 mm
330030S	UNI-CP™ T shaped plate
330021S	4 holes - 17x20 mm U shaped plate
330023S	4 holes - 19x22.5 mm U shaped plate
330025S	4 holes - 21x25 mm U shaped plate



Surfix® Sterile Titanium screw
dia. 3.5 mm + lock-screw

Reference	Description
285312S	L. 12 mm
285314S	L. 14 mm
285316S	L. 16 mm
285318S	L. 18 mm
285320S	L. 20 mm
285322S	L. 22 mm
285324S	L. 24 mm
285326S	L. 26 mm
285328S	L. 28 mm
285330S	L. 30 mm
285332S	L. 32 mm
285334S	L. 34 mm
285336S	L. 36 mm
285338S	L. 38 mm
285340S	L. 40 mm
185300S	Surfix® lock screw dia. 3.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm



UNI-CP™ plate 4 holes


UNI-CP™ plate 2 holes

CP T-shape plate

UNI-CP™ U-shape plate

UNI-CP™ Compression plate Instruments

Reference	Description
219335	Depth gauge hole dia. 3.5 mm
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219535	AO drill dia. 2.7 x L. 125 mm
219635	Drilling guide dia. 2.7 mm
219735	Bending forceps L. 171 hole - dia 3.5 mm
219835	Screwdriver hexa 2.0 x L. 180 mm
339001	Compression forceps
339004	Trial implant for T and U plates
339005	Trial implant for 2 and 4 holes plates
339003	Holder
115101	K-Wire dia. 1.0 x L. 100 mm



Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

UNI-CP™ sterilization container

Reference	Description
339950	UNI-CP™ sterilization container including:
996200	Lid
339920	Base
169106	Cylinder
119909	Silicone wedge

HINTEGRA® Total Ankle Prosthesis – Primary Surgery

The HINTEGRA® Total Ankle Prosthesis is a three components design anatomical prosthesis.

- Materials:
 - > Talar & tibial components:
 - Cobalt Chromium alloy (CoCr), with double coating:
 - » Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
 - > PE inlay: UHMWPE.
 - > Fixation screws (optional): Stainless steel M30NW.



HINTEGRA® Total Ankle Prosthesis – Tibial Component – Primary Surgery

Minimum contact stress to the bone.

- Materials:
 - > Cobalt Chromium alloy (CoCr) with double coating:
 - Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- 7 sizes right & 7 sizes left.

Reference	Description
HINTEGRA® - Tibial Component - Right	
301200	HINTEGRA® tibial component-right - size 0
301201	HINTEGRA® tibial component-right - size 1
301202	HINTEGRA® tibial component-right - size 2
301203	HINTEGRA® tibial component-right - size 3
301204	HINTEGRA® tibial component-right - size 4
301205	HINTEGRA® tibial component-right - size 5
301206	HINTEGRA® tibial component-right - size 6
HINTEGRA® - Tibial Component - Left	
302200	HINTEGRA® tibial component-left - size 0
302201	HINTEGRA® tibial component-left - size 1
302202	HINTEGRA® tibial component-left - size 2
302203	HINTEGRA® tibial component-left - size 3
302204	HINTEGRA® tibial component-left - size 4
302205	HINTEGRA® tibial component-left - size 5
302206	HINTEGRA® tibial component-left - size 6



HINTEGRA® Total Ankle Prosthesis – PE Inlay – Primary Surgery

- Ultra High Molecular Weight Polyethylene (UHMWPE).
- X-rays markers (Titanium TA6V rods).
- Same size as talar component: 6 sizes with 4 thicknesses (5, 6, 7, 9 mm).

HINTEGRA® - PE inlay - size 0

300005	HINTEGRA® PE inlay - size 0 - thickness 5 mm
300006	HINTEGRA® PE inlay - size 0 - thickness 6 mm
300007	HINTEGRA® PE inlay - size 0 - thickness 7 mm
300009	HINTEGRA® PE inlay - size 0 - thickness 9 mm

HINTEGRA® - PE inlay - size 1

300105	HINTEGRA® PE inlay - size 1 - thickness 5 mm
300106	HINTEGRA® PE inlay - size 1 - thickness 6 mm
300107	HINTEGRA® PE inlay - size 1 - thickness 7 mm
300109	HINTEGRA® PE inlay - size 1 - thickness 9 mm

HINTEGRA® - PE inlay - size 2

300205	HINTEGRA® PE inlay - size 2 - thickness 5 mm
300206	HINTEGRA® PE inlay - size 2 - thickness 6 mm
300207	HINTEGRA® PE inlay - size 2 - thickness 7 mm
300209	HINTEGRA® PE inlay - size 2 - thickness 9 mm

HINTEGRA® - PE inlay - size 3

300305	HINTEGRA® PE inlay - size 3 - thickness 5 mm
300306	HINTEGRA® PE inlay - size 3 - thickness 6 mm
300307	HINTEGRA® PE inlay - size 3 - thickness 7 mm
300309	HINTEGRA® PE inlay - size 3 - thickness 9 mm

HINTEGRA® - PE inlay - size 4

300405	HINTEGRA® PE inlay - size 4 - thickness 5 mm
300406	HINTEGRA® PE inlay - size 4 - thickness 6 mm
300407	HINTEGRA® PE inlay - size 4 - thickness 7 mm
300409	HINTEGRA® PE inlay - size 4 - thickness 9 mm

HINTEGRA® - PE inlay - size 5

300505	HINTEGRA® PE inlay - size 5 - thickness 5 mm
300506	HINTEGRA® PE inlay - size 5 - thickness 6 mm
300507	HINTEGRA® PE inlay - size 5 - thickness 7 mm
300509	HINTEGRA® PE inlay - size 5 - thickness 9 mm

HINTEGRA® - PE inlay - size 6

300605	HINTEGRA® PE inlay - size 6 - thickness 5 mm
300606	HINTEGRA® PE inlay - size 6 - thickness 6 mm
300607	HINTEGRA® PE inlay - size 6 - thickness 7 mm
300609	HINTEGRA® PE inlay - size 6 - thickness 9 mm



HINTEGRA® Total Ankle Prosthesis – Talar component with Pegs – Primary Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- Cobalt Chromium alloy (CoCr).
- Anatomical conical shape.
- 7 sizes right and 7 sizes left components.



Reference Description

HINTEGRA® - Talar Component with Pegs - Right

301110	HINTEGRA® talar component with pegs-right - size 0
301111	HINTEGRA® talar component with pegs-right - size 1
301112	HINTEGRA® talar component with pegs-right - size 2
301113	HINTEGRA® talar component with pegs-right - size 3
301114	HINTEGRA® talar component with pegs-right - size 4
301115	HINTEGRA® talar component with pegs-right - size 5
301116	HINTEGRA® talar component with pegs-right - size 6

HINTEGRA® - Talar Component with Pegs - Left

302110	HINTEGRA® talar component with pegs-left - size 0
302111	HINTEGRA® talar component with pegs-left - size 1
302112	HINTEGRA® talar component with pegs-left - size 2
302113	HINTEGRA® talar component with pegs-left - size 3
302114	HINTEGRA® talar component with pegs-left - size 4
302115	HINTEGRA® talar component with pegs-left - size 5
302116	HINTEGRA® talar component with pegs-left - size 6

HINTEGRA® Total Ankle Prosthesis – (Optional) Tibial cortical screw – Primary Surgery


- Stainless steel M30NW.
- Optionally to be put only with tibial components.



Reference Description

303328S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 28 mm - tibia - size 0
303330S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 30 mm - size 1
303334S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 34 mm - size 2
303338S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 38 mm - size 3
303340S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 40 mm - size 4
303342S	HINTEGRA® cortical screw - Dia. 3.5 mm x L. 42 mm - size 5 & 6


HINTEGRA® Total Ankle Prosthesis Associated Instruments for Primary Surgery



Generic base
(upper tray)



Trial tray



Talar base

HINTEGRA® Total Ankle Prosthesis Associated Instruments for Primary Surgery

Instruments for LOWER TRAY

Talar Base	
309300	Drilling guide talus - Size 0
309301	Drilling guide talus - Size 1
309302	Drilling guide talus - Size 2
309303	Drilling guide talus - Size 3
309304	Drilling guide talus - Size 4
309305	Drilling guide talus - Size 5
309306	Drilling guide talus - Size 6
309360	Tal. cut. guide-Right - Size 0
309361	Tal. cut. guide-Right - Size 1
309362	Tal. cut. guide-Right - Size 2
309363	Tal. cut. guide-Right - Size 3
309364	Tal. cut. guide-Right - Size 4
309365	Tal. cut. guide-Right - Size 5
309366	Tal. cut. guide-Right - Size 6
309370	Tal. cut. guide-Left - Size 0
309371	Tal. cut. guide-Left - Size 1
309372	Tal. cut. guide-Left - Size 2
309373	Tal. cut. guide-Left - Size 3
309374	Tal. cut. guide-Left - Size 4
309375	Tal. cut. guide-Left - Size 5
309376	Tal. cut. guide-Left - Size 6
309601	Screwdriver tibial screws
309608	Spacer h 12 mm
119609	Talar impactor
119641	Cutting guide tibial revision
309200	Reamer talar dia. 6 mm
309309	Peg drill talus - dia. 4.5 mm
309697	Holder trial talar implant
309699	Impactor trial talar implant

Trial Tray

HINTEGRA® Tibial Trial-Right	
119646	Tibial Trial-Right - Size 0
119690	Tibial Trial-Right - Size 1
119691	Tibial Trial-Right - Size 2
119692	Tibial Trial-Right - Size 3
119693	Tibial Trial-Right - Size 4
119694	Tibial Trial-Right - Size 5
119778	Tibial Trial-Right - Size 6
HINTEGRA® Tibial Trial-Left	
119647	Tibial Trial-Left - Size 0
119695	Tibial Trial-Left - Size 1
119696	Tibial Trial-Left - Size 2
119697	Tibial Trial-Left - Size 3
119698	Tibial Trial-Left - Size 4
119699	Tibial Trial-Left - Size 5
119779	Tibial Trial-Left - Size 6
HINTEGRA® Talar Trial-Right	
309680	Talar Trial-Right - Size 0
309681	Talar Trial-Right - Size 1
309682	Talar Trial-Right - Size 2
309683	Talar Trial-Right - Size 3
309684	Talar Trial-Right - Size 4
309685	Talar Trial-Right - Size 5
309686	Talar Trial-Right - Size 6
HINTEGRA® Talar Trial-Left	
309690	Talar Trial-Left - Size 0
309691	Talar Trial-Left - Size 1
309692	Talar Trial-Left - Size 2
309693	Talar Trial-Left - Size 3
309694	Talar Trial-Left - Size 4
309695	Talar Trial-Left - Size 5
309696	Talar Trial-Left - Size 6
HINTEGRA® Trial inlay	
119665	Trial inlay - thickness 5mm
119666	Trial inlay - thickness 6mm
119667	Trial inlay - thickness 7mm
119669	Trial inlay - thickness 9 mm

Sterilization container

309900	HINTEGRA® Container (Upper + Lower) including:		
309902	Stainless steel base upper, including:	309904	Stainless steel base lower, including:
309920	Generic base	309960	Talar base
309940	Pins insert	309980	Thermo insert trial implants
309942	Lid pins insert	309970	Lid trial implants
		996350	Lid
		119909	Silicone wedge

Instruments for UPPER TRAY

309625	Tibial v positioning
309630	V translation block
119664	Distractor Hintermann® big model
119662	Tibial - Talar Holder
119751	Tibial impactor
309604	Pins L. 30 mm
309605	Pins L. 70 mm
309606	Pins L. 45 mm
309607	Depth gauge tibial
159127	Drilling guide dia. 2.7mm
309615	Alignment rod tibial
309620	Connector v - tibial rod
309622	Saw blade Aesculap attachment
309623	Saw blade AO Synthes attachment
309624	Saw blade Stryker 6 attachment
309626	Saw blade Stryker attachment
309627	Saw blade Conmed Linvatec attachment
309635	Tib. cut. guide large
309636	Tib. cut. guide medium / standard
309637	Tib. cut. guide small
309645	Screwdriver Hex dia. 3.5 mm
309652	Imp / Extract pins
309656	Cutting block-talar cut / Medium - Standard
309657	Cutting block-talar cut - Small
309661	Holder PE inlays
309753	Cutting block proximal - tibial
309773	Cutting block distal - tibial
997034	Curved chisel 20mm r40mm
997035	Straight chisel 10mm
275560	Parallel pliers 7-1/4"
119611	Drill dia. 2.5 mm x L. 150 mm
115116	K-Wire dia. 1.6 x L. 150 mm
115225	K-Wire dia. 2.5 x L. 200 mm

HINTEGRA® Total Ankle Prosthesis – Flat Cut – Revision Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).



HINTEGRA® Total Ankle Prosthesis Flat Cut – Tibial Component – Revision Surgery

- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).
- Two thicknesses: 8 and 12 mm.



Reference Description

HINTEGRA® - revision tibial component - right - 8 mm

301222	HINTEGRA® revision tibial component-right-8 mm - size 2
301223	HINTEGRA® revision tibial component-right-8 mm - size 3
301224	HINTEGRA® revision tibial component-right-8 mm - size 4
301225	HINTEGRA® revision tibial component-right-8 mm - size 5

HINTEGRA® - revision tibial component - left - 8 mm

302222	HINTEGRA® revision tibial component-left-8 mm - size 2
302223	HINTEGRA® revision tibial component-left-8 mm - size 3
302224	HINTEGRA® revision tibial component-left-8 mm - size 4
302225	HINTEGRA® revision tibial component-left-8 mm - size 5

Reference Description

HINTEGRA® - revision tibial component - right - 12 mm

301232	HINTEGRA® revision tibial component-right-12 mm - size 2
301233	HINTEGRA® revision tibial component-right-12 mm - size 3
301234	HINTEGRA® revision tibial component-right-12 mm - size 4
301235	HINTEGRA® revision tibial component-right-12 mm - size 5

HINTEGRA® - revision tibial component - left - 12 mm

302232	HINTEGRA® revision tibial component-left-12 mm - size 2
302233	HINTEGRA® revision tibial component-left-12 mm - size 3
302234	HINTEGRA® revision tibial component-left-12 mm - size 4
302235	HINTEGRA® revision tibial component-left-12 mm - size 5

HINTEGRA® Total Ankle Prosthesis Flat Cut – Talar Component – Revision Surgery / Flat Talus

- For flat talus and revision surgeries.
- Double coating: Porous Titanium & Hydroxyapatite (Ti and HAP), for optimal bone ingrowth (osteointegration).



Reference Description

HINTEGRA® - flat cut talar component - right

301121	HINTEGRA® revision talar component-right - size 1
301122	HINTEGRA® revision talar component-right - size 2
301123	HINTEGRA® revision talar component-right - size 3
301124	HINTEGRA® revision talar component-right - size 4
301125	HINTEGRA® revision talar component-right - size 5

Reference Description

HINTEGRA® - flat cut talar component - left

302121	HINTEGRA® revision talar component-left - size 1
302122	HINTEGRA® revision talar component-left - size 2
302123	HINTEGRA® revision talar component-left - size 3
302124	HINTEGRA® revision talar component-left - size 4
302125	HINTEGRA® revision talar component-left - size 5

HINTEGRA® Total Ankle Prosthesis – Cortical screw

- Stainless steel: M30NW.
- Sterile.



Reference Description

303312S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 12 mm - talus
303316S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 16 mm - talus
303320S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 20 mm - talus
303334S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 34 mm - tibia - size 2
303338S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 38 mm - tibia - size 3
303340S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 40 mm - tibia - size 4
303342S	HINTEGRA® cortical screw - dia. 3.5 mm x l. 42 mm - tibia - size 5 & 6

HINTEGRA® Total Ankle Prosthesis – PE Inlay – Revision Surgery (See page 19 for product references.)

Size of the PE inlay has to be the same as the talar component.




HINTEGRA® Total Ankle Prosthesis – Flat Cut Instruments for Revision Surgery

Instruments

119790	FC Talar Trial for Drilling-Right - Size 1
119791	FC Talar Trial for Drilling-Right - Size 2
119792	FC Talar Trial for Drilling-Right - Size 3
119793	FC Talar Trial for Drilling-Right - Size 4
119794	FC Talar Trial for Drilling-Right - Size 5
119657	Drill dia. 4.5 mm x L. 22 mm
119659	FC Talar trial for drilling Impactor
119795	FC Talar Trial for Drilling-Left - Size 1
119796	FC Talar Trial for Drilling-Left - Size 2
119797	FC Talar Trial for Drilling-Left - Size 3
119798	FC Talar Trial for Drilling-Left - Size 4
119799	FC Talar Trial for Drilling-Left - Size 5
119631	FC Spacer 4 mm for revision Tibial Component 8 mm
119632	FC Spacer 8 mm for revision Tibial Component 12 mm
119757	FC spacer
119760	FC Talar Cutting Guide-Right - Size 1
119761	FC Talar Cutting Guide-Right - Size 2
119762	FC Talar Cutting Guide-Right - Size 3
119763	FC Talar Cutting Guide-Right - Size 4
119764	FC Talar Cutting Guide-Right - Size 5
119765	FC Talar Cutting Guide-Left - Size 1
119766	FC Talar Cutting Guide-Left - Size 2
119767	FC Talar Cutting Guide-Left - Size 3
119768	FC Talar Cutting Guide-Left - Size 4
119769	FC Talar Cutting Guide-Left - Size 5
309 708	Tibial depth gauge for tibial component 8 mm
309712	Tibial depth gauge for tibial component 12 mm

Sterilization container

119965	HINTEGRA® FC revision container includes:
119978	Container pocket
119966	HINTEGRA® FC container base
119967	Lid for HINTEGRA® FC container base



HINTEGRA® Sensitive Total Ankle Prosthesis – Primary Surgery

The HINTEGRA® Sensitive Total Ankle Prosthesis is also a three component design specially made for allergic patients.

- Materials:
 - > Talar & Tibial components:
 - » Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
 - > PE inlay:
 - » Polyethylene (UHMWPE).



HINTEGRA® Sensitive Total Ankle Prosthesis – Tibial Component – Primary Surgery

- Anatomical shape for minimal contact stress to the bone.
- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 6 right and 6 left components.

Reference	Description
HINTEGRA® sensitive - tibial component - right	
351251	HINTEGRA® sensitive - tibia-right - size 1
351252	HINTEGRA® sensitive - tibia-right - size 2
351253	HINTEGRA® sensitive - tibia-right - size 3
351254	HINTEGRA® sensitive - tibia-right - size 4
351255	HINTEGRA® sensitive - tibia-right - size 5
351256	HINTEGRA® sensitive - tibia-right - size 6
HINTEGRA® sensitive - tibial component - left	
352251	HINTEGRA® sensitive - tibia-left - size 1
352252	HINTEGRA® sensitive - tibia-left - size 2
352253	HINTEGRA® sensitive - tibia-left - size 3
352254	HINTEGRA® sensitive - tibia-left - size 4
352255	HINTEGRA® sensitive - tibia-left - size 5
352256	HINTEGRA® sensitive - tibia-left - size 6

HINTEGRA® sensitive - tibial component - right

351251 HINTEGRA® sensitive - tibia-right - size 1

351252 HINTEGRA® sensitive - tibia-right - size 2

351253 HINTEGRA® sensitive - tibia-right - size 3

351254 HINTEGRA® sensitive - tibia-right - size 4

351255 HINTEGRA® sensitive - tibia-right - size 5

351256 HINTEGRA® sensitive - tibia-right - size 6

HINTEGRA® sensitive - tibial component - left

352251 HINTEGRA® sensitive - tibia-left - size 1

352252 HINTEGRA® sensitive - tibia-left - size 2

352253 HINTEGRA® sensitive - tibia-left - size 3

352254 HINTEGRA® sensitive - tibia-left - size 4

352255 HINTEGRA® sensitive - tibia-left - size 5

352256 HINTEGRA® sensitive - tibia-left - size 6



HINTEGRA® Sensitive Total Ankle Prosthesis – Talar component with Pegs – Primary Surgery

- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- Anatomical conical shape.
- 5 sizes and 5 left components.

Reference	Description
HINTEGRA® - talar component with pegs - right	
351152	HINTEGRA® sensitive - talus -right - size 2
351153	HINTEGRA® sensitive - talus -right - size 3
351154	HINTEGRA® sensitive - talus -right - size 4
351155	HINTEGRA® sensitive - talus -right - size 5
351156	HINTEGRA® sensitive - talus -right - size 6
HINTEGRA® - talar component with pegs - left	
352152	HINTEGRA® sensitive - talus-left - size 2
352153	HINTEGRA® sensitive - talus-left - size 3
352154	HINTEGRA® sensitive - talus-left - size 4
352155	HINTEGRA® sensitive - talus-left - size 5
352156	HINTEGRA® sensitive - talus-left - size 6



HINTEGRA® Sensitive Total Ankle Prosthesis – PE Inlay (See page 45 for product references.)

- HINTEGRA® PE inlays are “sensitive compatible” since metallic rods are made of titanium Alloy (TA6V).
- Size of the PE inlay has to be the same as the talar component.



HINTEGRA® Sensitive Total Ankle Prosthesis – (Optional) Tibial cortical screw

- Titanium alloy (TA6V).
- Optionally to be put only with tibial components.

Reference	Description
303530S	Cortical Screw Ti - Dia. 3.5 mm x L. 30 mm - Tibia - Size 1
303534S	Cortical Screw Ti - Dia. 3.5 mm x L. 34 mm - Tibia - Size 2
303538S	Cortical Screw Ti - Dia. 3.5 mm x L. 38 mm - Tibia - Size 3
303540S	Cortical Screw Ti - Dia. 3.5 mm x L. 40 mm - Tibia - Size 4
303542S	Cortical Screw Ti - Dia. 3.5 mm x L. 42 mm - Tibia - Size 5 & 6



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Revision Surgery

The HINTEGRA® Sensitive Total Ankle Prosthesis is also a three component design specially made for allergic patients.

- Materials:
 - > Talar & Tibial components:
 - » Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
 - > PE inlay:
 - » Polyethylene (UHMWPE).



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Tibial Component – Revision Surgery

- Anatomical shape for minimal contact stress to the bone.
- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 4 right and 4 left components.

Reference	Description	Reference	Description
HINTEGRA® - Revision Tibial Component - Right - 8 mm		HINTEGRA® - Revision Tibial Component - Right - 12 mm	
351522	HINTEGRA® sensitive FC - tibia-right 8 mm - size 2	351532	HINTEGRA® sensitive FC - tibia-right 12 mm - size 2
351523	HINTEGRA® sensitive FC - tibia-right 8 mm - size 3	351533	HINTEGRA® sensitive FC - tibia-right 12 mm - size 3
351524	HINTEGRA® sensitive FC - tibia-right 8 mm - size 4	351534	HINTEGRA® sensitive FC - tibia-right 12 mm - size 4
351525	HINTEGRA® sensitive FC - tibia-right 8 mm - size 5	351535	HINTEGRA® sensitive FC - tibia-right 12 mm - size 5
HINTEGRA® - Revision Tibial Component - Left - 8 mm		HINTEGRA® - Revision Tibial Component - Left - 12 mm	
352522	HINTEGRA® sensitive FC - tibia-left 8 mm - size 2	352532	HINTEGRA® sensitive FC - tibia-left 12 mm - size 2
352523	HINTEGRA® sensitive FC - tibia-left 8 mm - size 3	352533	HINTEGRA® sensitive FC - tibia-left 12 mm - size 3
352524	HINTEGRA® sensitive FC - tibia-left 8 mm - size 4	352534	HINTEGRA® sensitive FC - tibia-left 12 mm - size 4
352525	HINTEGRA® sensitive FC - tibia-left 8 mm - size 5	352535	HINTEGRA® sensitive FC - tibia-left 12 mm - size 5



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut Talar component with Pegs – Revision Surgery / Flat talus

- Cobalt Chromium alloy (CoCr) with Titanium Nitride (TiN) coating.
- 4 right and 4 left components.

Reference	Description
HINTEGRA® - Talar Component with Pegs - Right	
351512	HINTEGRA® sensitive FC - talus-right - size 2
351513	HINTEGRA® sensitive FC - talus-right - size 3
351514	HINTEGRA® sensitive FC - talus-right - size 4
351515	HINTEGRA® sensitive FC - talus-right - size 5
HINTEGRA® - Talar Component with Pegs - Left	
352512	HINTEGRA® sensitive FC - talus-left - size 2
352513	HINTEGRA® sensitive FC - talus-left - size 3
352514	HINTEGRA® sensitive FC - talus-left - size 4
352515	HINTEGRA® sensitive FC - talus-left - size 5



HINTEGRA® Sensitive Total Ankle Prosthesis – PE Inlay (See page 45 for product references.)

- HINTEGRA® PE inlays are “sensitive compatible” since metallic rods are made of titanium Alloy (TA6V).
- Size of the PE inlay has to be the same as the talar component.



HINTEGRA® Sensitive Total Ankle Prosthesis – Flat Cut – Cortical screw

Titanium alloy (TA6V).

Reference	Description
303512S	Cortical Screw Ti - dia.3.5 x L. 12 mm - Talus
303516S	Cortical Screw Ti - dia.3.5 x L. 16 mm - Talus
303520S	Cortical Screw Ti - dia.3.5 x L. 20 mm - Talus
303534S	Cortical Screw Ti - dia. 3.5 x L. 34 mm - Tibia - Size 2
303538S	Cortical Screw Ti - dia. 3.5 x L. 38 mm - Tibia - Size 3
303540S	Cortical Screw Ti - dia. 3.5 x L. 40 mm - Tibia - Size 4
303542S	Cortical Screw Ti - dia. 3.5 x L. 42 mm - Tibia - Size 5



PANTA® / PANTA® XL Arthrodesis Nail

The PANTA® Nail is intended for use in tibiototalcalcaneal arthrodesis and treatment of trauma of the hindfoot and distal tibia.

- Cannulated Titanium Nail.
- Radio-transparent instrumentation system to apply compression.
- Multi-axial fixation:
 - > Oblique antero-posterior calcaneal screws.
 - > Medio-lateral tibial screws.
- 2 kind of screws:
 - > Oblique antero-posterior calcaneal screws.
 - » fully threaded
 - » partially threaded
- 5 mm diameter Titanium screw.
- Screw lengths: 20 to 110 mm.
- Material of the nail, the screws and the end cap: Titanium alloy TA6V.

PANTA® Nails (sterile)

Reference Dia. - length

500050	Dia. 10 mm - L. 150 mm
500080	Dia. 10 mm - L. 180 mm
500150	Dia. 11 mm - L. 150 mm
500180	Dia. 11 mm - L. 180 mm
500250	Dia. 12 mm - L. 150 mm
500280	Dia. 12 mm - L. 180 mm
500350	Dia. 13 mm - L. 150 mm
500380	Dia. 13 mm - L. 180 mm

PANTA® XL Nails (sterile)

To be used only with Partially threaded Screws

Reference Dia. - length

510111	Dia. 11mm - L 210 mm
510141	Dia. 11mm - L 240 mm
510211	Dia. 12 mm - L 210 mm
510241	Dia. 12 mm - L 240 mm
510311	Dia. 13mm - L 210 mm
510341	Dia. 13mm - L 240 mm

End cap (sterile)

Reference Description

500001	End cap
510005	Locking end cap



PANTA® Arthrodesis Nail – Fully threaded screws dia. 5 mm (not for use with PANTA® XL Arthrodesis Nail)

Reference	Length
501020	L. 20 mm
501022	L. 22.5 mm
501025	L. 25 mm
501027	L. 27.5 mm
501030	L. 30 mm
501035	L. 35 mm
501040	L. 40 mm
501045	L. 45 mm
501050	L. 50 mm
501055	L. 55 mm
501060	L. 60 mm
501065	L. 65 mm
501070	L. 70 mm
501075	L. 75 mm
501080	L. 80 mm
501085	L. 85 mm
501090	L. 90 mm
501095	L. 95 mm
501100	L. 100 mm
501105	L. 105 mm
501110	L. 110 mm


PANTA® / PANTA® XL Arthrodesis Nail – Partially threaded lag screws dia. 5 mm

Reference	Length
511020	L. 20 mm
511022	L. 22.5 mm
511025	L. 25 mm
511027	L. 27.5 mm
511030	L. 30 mm
511035	L. 35 mm
511040	L. 40 mm
511045	L. 45 mm
511050	L. 50 mm
511055	L. 55 mm
511060	L. 60 mm
511065	L. 65 mm
511070	L. 70 mm
511075	L. 75 mm
511080	L. 80 mm
511085	L. 85 mm
511090	L. 90 mm
511095	L. 95 mm
511100	L. 100 mm
511105	L. 105 mm
511110	L. 110 mm



PANTA® Arthrodesis Nail – Instruments

Reference	Description
519003	Drill dia. 4.3 mm for Calcaneum Screws
519004	Drill dia. 4.3 mm for Tibial Screws
519005	Drill dia. 5.0 mm – AO Attachment
519006	Long drill for calcaneus screws dia. 4.3 mm
519007	Cannulated Drill dia. 7.0 mm
519009	Cannulated Drill dia. 9.0 mm
519011	Reamer dia. 11 mm
519012	Reamer dia. 12 mm
519013	Reamer dia. 13 mm
519014	Reamer dia. 10.5 mm
519015	Reamer dia. 11.5 mm
519016	Reamer dia. 12.5 mm
519017	Reamer dia. 13.5 mm
519021	T Handle AO Attachment
519028	Internal Protecting Sleeve dia. 3.2 mm
519029	Central Internal Protecting Sleeve dia. 9.0 mm
519030	External Protecting Sleeve
519032	Pin dia. 3.2 mm x L. 400 mm
519040	Trocard Awl
519110	Support Device
519120	Nail Fixation Axis
519121	Toothed Wheel
519131	Threaded Axis
519130	Compression Device
519133	Teflon Ring
519135	Compression Wheel
519160	Depth Gauge
519175	Compression Rod
519179	Short Drill Guide dia. 4.3 mm L. 70 mm
519180	Long Drill Guide dia. 4.3 mm L. 120 mm
519181	Drill Guide dia. 5.0 mm
519184	Short Drill Guide dia. 7.0 mm L. 70 mm
519185	Drill Guide dia. 7.0 mm L. 120 mm
519290	Hexagonal screwdriver tip dia. 3.5 mm
519295	Hexagonal screwdriver dia. 3.5 mm
539015	Tap dia. 6.5 mm

PANTA® – Sterilization Container

519920	PANTA® sterilization container including:
519910	PANTA® base
519911	PANTA® Insert
519920	Lid
119909	Silicone wedge

PANTA® – Optional Instruments

519010	Reamer dia. 10 mm
519020	Quick Coupling AO
519002	Short drill for calcaneus screws used with medium drill guide dia. 4.3 mm
519008	Long drill for calcaneus screws used with medium drill guide dia. 4.3 mm
519034(S)	Guide wire dia. 3.2 mm, L. 600 mm (sterile upon request)
519178	Medium drill guide dia. 4.3 mm
519183	Medium soft tissue protector dia. 7 mm




PANTA® XL – Associated Instruments (optional)

539020	Conical reamer dia. 10.5 mm
539021	Conical reamer dia. 11.5 mm
539022	Conical reamer dia. 12.5 mm
539023	Conical reamer dia. 13.5 mm
519912	Reamer holder



(S) means available in sterile and non sterile

PANTA[®] Arthrodesis Nail - Associated Instruments



PANTA® Removal Set

Reference	Description
519209	Sliding hammer
519210	Handle extractor
309645	Screwdriver Hex dia.3.5 mm
119552	External protection sleeve dia.6.5 mm
Sterilization container	
519950	PANTA® Removal sterilization container
119909	Silicone wedge
519951	Basis
996100	Lid



ADVANSYS® TTC Fusion Plate

The TTC Plates are intended for arthrodesis of the ankle joint and calcanea, osteotomies, fusions and replantations of small bones in the ankle and distal tibia.

- 4 sizes of implants
 - > 4 right & 4 left plates (in 6,7,8,9 holes configurations).
- Screw dia.4.5 mm for fixation of the tibia and talus – Length 14 mm to 65 mm.
- Screw dia.6.5 mm for fixation of the calcaneus – Length 20 mm to 65 mm.
- Lock-screws for 4.5 mm dia. screws.
- Lock-screws for 6.5 mm dia. screws.
- Material:
 - > Plate: Stainless Steel 316LVM.
 - > Screw: Titanium TA6V.



Advansys® TTC Fusion Plate

Reference	Description
181001S	ADVANSYS® TTCP Right-6 holes
181002S	ADVANSYS® TTCP Right-8 holes
181003S	ADVANSYS® TTCP Right-7 holes
181004S	ADVANSYS® TTCP Right-9 holes
181011S	ADVANSYS® TTCP Left-6 holes
181012S	ADVANSYS® TTCP Left-8 holes
181013S	ADVANSYS® TTCP Left-7 holes
181014S	ADVANSYS® TTCP Left-9 holes



Surfix® Locking screw
Sterile Titanium dia. 4.5 mm + lock-screw

Reference	Length
285414s	L. 14 mm
285416s	L. 16 mm
285418s	L. 18 mm
285420s	L. 20 mm
285422s	L. 22 mm
285424s	L. 24 mm
285426s	L. 26 mm
285428s	L. 28 mm
285430s	L. 30 mm
285435s	L. 35 mm
285440s	L. 40 mm
285445s	L. 45 mm
285450s	L. 50 mm
285455s	L. 55 mm
285460s	L. 60 mm
285465s	L. 65 mm
185400s	SURFIX® Lock Screw dia. 4.5 mm



Surfix® Locking screw
Sterile Titanium dia. 6.5 mm + lock-screw

Reference	Length
285620s	L. 20 mm
285625s	L. 25 mm
285630s	L. 30 mm
285635s	L. 35 mm
285640s	L. 40 mm
285645s	L. 45 mm
285650s	L. 50 mm
285655s	L. 55 mm
285660s	L. 60 mm
285665s	L. 65 mm
185600s	SURFIX® lock Screw dia. 6.5 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

ADVANSYS® TTC Fusion Plate Instruments



Reference	Description
219845	Screwdriver hexa 2.5 x L. 191 mm
219865	Screwdriver hexa dia.3.5 x L. 230 mm
219445	AO Screwdriver tip hexa dia.2.5 x L. 76 mm
219465	AO Screwdriver tip hexa dia.3.5 x L. 76 mm
219645	Drilling guide dia.3.0 mm
219665	Drilling guide dia.3.5 mm
219545	AO Drill dia.3.0 x L. 190 mm
219565	Drill dia.3.5 x L. 200 mm
219900	Compression forceps L. 260 mm
219765	Bending forceps L. 193 - hole dia.6.5 mm
219345	Depth gauge hole dia.4.5 mm
219365	Depth gauge hole dia.6.5 mm
115116	K-Wire dia.1.6 x L. 150 mm
115425	K-Wire dia.2.5 x L. 200 mm
169001	TTCP trial implant 6 holes
169002	TTCP trial implant 8 holes
169003	TTCP trial implant 7 holes
169004	TTCP trial implant 9 holes

Sterilization container

169101	ADVANSYS® Hindfoot sterilization container including:
169111	Base
996200	Lid
169104	Cylinder
169105	TTCP trial implant module
119909	Silicone wedge

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm





TIBIAXYS® Ankle Arthrodesis Plating System

TIBIAXYS® Distal Tibia and Fibula Osteotomy Plating Systems

The TIBIAXYS® Plates are intended for fixation of arthrodesis osteotomies and fractures of ankle joint, distal tibia and fibula.

- Each plate is anatomically contoured and available in 1 size.
- 2 systems of 2 anterior plates for ankle fusion (green left and blue right systems, including 1 antero-medial and 1 antero-lateral plates).
- 3 different plates for distal tibia osteotomies (medial plate, left lateral / right medial anterior plate, left medial/ right lateral anterior plate).
- 2 different plates for fibula osteotomies (4 and 6 holes).
- dia. 3.5 mm Surfix® or Surfix®-Alpha screws and lock-screws.
- dia. 4.0 mm cortical screws (only for 1 hole in each plate for the ankle fusion).
- Material:
 - > Plates: Titanium alloy TA6V.
 - > Surfix® and Surfix®-Alpha screws and lock-screws: Titanium alloy TA6V.
 - > Cortical screws: Titanium alloy TA6V.



Reference	Description
150010S	TIBIAXYS® Arthrodesis Plate / Left antero-medial
150020S	TIBIAXYS® Arthrodesis Plate / Left antero-lateral
150110S	TIBIAXYS® Arthrodesis Plate / right antero-medial
150120S	TIBIAXYS® Arthrodesis Plate / right antero-lateral

Reference	Description
150030S	TIBIAXYS® Osteotomy Plate / Lateral-right / Medial-left
150130S	TIBIAXYS® Osteotomy Plate / Medial-right / Lateral-left
150040S	TIBIAXYS® Osteotomy Plate / Medial

Reference	Description
150514S	TIBIAXYS® Fibula Osteotomy Plate 4 holes
150516S	TIBIAXYS® Fibula Osteotomy Plate 6 holes

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

TIBIAXYS® Ankle Arthrodesis Plating System

TIBIAXYS® Distal Tibia and Fibula Osteotomy Plating Systems



Surfix®-Alpha variable angle locking screws – Sterile Titanium
dia. 3.5 mm + lock-screw

Reference	Length
295310s	L. 10 mm
295312s	L. 12 mm
295314s	L. 14 mm
295316s	L. 16 mm
295318s	L. 18 mm
295320s	L. 20 mm
295322s	L. 22 mm
295324s	L. 24 mm
295326s	L. 26 mm
295328s	L. 28 mm
295330s	L. 30 mm
295332s	L. 32 mm
295334s	L. 34 mm
295336s	L. 36 mm
295338s	L. 38 mm
295340s	L. 40 mm
295344s	L. 44 mm
295348s	L. 48 mm
295350s	L. 50 mm
195300s	Lock-screw dia. 3.5 mm

Surfix® Locking screws
Sterile Titanium
dia. 3.5 mm + lock-screw

Reference	Length
285310s	L. 10 mm
285312s	L. 12 mm
285314s	L. 14 mm
285316s	L. 16 mm
285318s	L. 18 mm
285320s	L. 20 mm
285322s	L. 22 mm
285324s	L. 24 mm
285326s	L. 26 mm
285328s	L. 28 mm
285330s	L. 30 mm
285332s	L. 32 mm
285334s	L. 34 mm
285336s	L. 36 mm
285338s	L. 38 mm
285340s	L. 40 mm
285344s	L. 44 mm
285348s	L. 48 mm
285350s	L. 50 mm
185300s	Lock-screw dia. 3.5 mm

Sterile Titanium cortical screw
dia. 4.0 mm

Reference	Length
150240s	L. 40 mm
150242s	L. 42 mm
150246s	L. 46 mm
150250s	L. 50 mm
150255s	L. 55 mm
150260s	L. 60 mm
150265s	L. 65 mm
150270s	L. 70 mm
150275s	L. 75 mm
150280s	L. 80 mm
150285s	L. 85 mm
150290s	L. 90 mm
150200s	L. 100 mm

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm

TIBIAXYS® Instruments

Reference	Description
115116	K-Wire dia.1.6 x L. 150 mm
115225	K-Wire dia.2.5 x L. 200 mm
159103	Wedge thickness 3 mm
159106	Wedge thickness 6 mm
159109	Wedge thickness 9 mm
159130	Cortical screw drilling guide dia.3.0 mm
159400	Depth gauge for dia.4.0 mm screws
159635	Compression forceps guide
159740	Compression forceps screw dia.4.0 x L. 40 mm
159755	Compression forceps screw dia.4.0 x L. 55 mm
159760	Compression forceps screw dia.4.0 x L. 60 mm
219035	Surfix® Alpha drilling guide dia.2.7 mm
219135	Screwdriver hexalobular T10
219335	Depth gauge for dia. 3.5 mm screws
219435	AO Screwdriver tip hexa 2.0 x L. 76 mm
219445	AO Screwdriver tip hexa dia.2.5 x L. 76 mm
219535	AO drill dia.2.7 x L. 125 mm
219545	AO Drill dia.3.0 x L. 190 mm
219635	Drilling guide dia.2.7 mm
219835	Screwdriver hexa dia. 2.0 x L. 180 mm
219845	Screwdriver hexa dia. 2.5 x L. 191 mm
219960	Compression forceps L. 260 mm
159010	Trial plate Left Antero-Medial
159020	Trial plate Left Antero-Lateral
159110	Trial plate Right Antero-Medial (optional)
159120	Trial plate Right Antero-Lateral (optional)

Sterilization container

159991	Tibixys® sterilization container including
159960	Base
996200	Lid
119909	Silicone wedge

Trial container


159010	Left medial anterior trial plate
159020	Left lateral anterior trial plate
159110	Right medial anterior trial plate
159020	Right lateral anterior trial plate

Container

159071	Base
309942	Lid

Optional instrument

Reference	Description
MA3511	Surfix® screw extractor 3.5 mm




ACHILLON® Minimal Invasive Achilles Tendon Suture System – Single Use

The ACHILLON® System is a method to treat acute Achilles tendon ruptures.

- Needle driver.
- Suture size: USP (3.5 metric).
- Materials:
 - > Polycarbonate.
 - > Stainless steel 316LVM.

Reference	Description
119700	ACHILLON® + 2 Needles and 1 Needle Driver



Hintermann® Distractors

Range of 4 spreading forceps.

For each of the 4 models, two pins (2.5 mm for the large model, and 1.6 or 2.5 mm diameter for the small one) are slid into the appropriate holes of the spreading forceps, and are then introduced at appropriate localisation into each of the two bone fragments to be outspread. The spreader is then distracted, and locked when the distance between the two fragments is reached.

The flexible positioning of these spreading forceps enables them not to interfere with the main procedure. These spreaders are available in two different sizes (15 cm long, or 21 cm long), and with closed or outspread arms, making them fitting many possible indication.

Hintermann® Small Distractor – 15 cm

Reference	Description
119653	Hintermann® Distractor Small open handle
119654	Hintermann® Distractor Small closed handle

Hintermann® Large Distractor – 21 cm

Reference	Description
119664	Hintermann® Distractor Big model
119651	Hintermann® Distractor Large closed handle

Hintermann® Distractor open handle



Hintermann® Distractor closed handle



Osteoguides

- Suggested use:
 - > Scarf osteotomy.
 - > Chevron osteotomy.
 - > Basal osteotomy.
 - > Arthrodesis.
 - > Wedge osteotomies (Akin...).
- Details:
 - > Hole for K-Wire: Diameter 1.0 to 1.6 mm.
 - > Knurled handle for optimal grip.
 - > Blade thickness maxi 0.8 mm.
 - > Positioning tip.

Reference	Description
229011	Longitudinal
229012	Transversal



K-FIX® – K-Wire End Protector

The K-Fix® pin protector is a simple system to protect patients or hospital staff from K-Wires sharp ends.

- Delivered sterile.

Reference	Length
119050	K-FIX® for K-WIRE - Dia. 0.8 mm to 2.0 mm (5 units per pack)



Weil Clamp

Reference	Length
119106	Weil Clamp

General Optional Instruments

Reference	Description
998001	Scalpel Handle N° 3
998002	Scalpel Handle N° 4
998003	Beyer Bone Rongeur Forceps
998004	Backhaus Towel Forceps
998005	Beckmann Mini Retractor
998006	Museux Forceps (14 cm)
998007	Lange- Hohmann Retractor
998008	Liston Bone Cutting Forceps
998010	Ombredanne Square Points
998011	Wolkman Retractor
998012	Wolkman Curette
998013	Lambotte Raspatory 12 mm
998014	Guillaume Gouge 8 mm
998015	Mac Even Hammer
998016	Redon Needle
998017	Finochietto Needle Holder 13 cm
998018	Pauwels Osteotomes 15 mm
998019	Farabeuf Retractor 12 cm
998020	Dissection Forceps (GR) 14 cm
998021	Dissection Forceps Curved 14 cm
998022	Metzenbaum Scissors
998023	Straight Scissors
998024	Mayo Scissors
998025	Leriche Forceps
998026	Verbrugge Bone Holding Forceps
998027	Poirier Gouge (14 mm)
998031	Metatarsal Head Extractor
MA3511	Surfix® screw 3.5 mm extractor
119651	Hintermann® Distractor / Large Closed Jaws
119653	Hintermann® Distractor / Small Open Jaws
119654	Hintermann® Distractor / Small Closed Jaws
119106	Weil clamp
998028	Inge Retractor
998029	Banaleck Forceps Right
998030	Banaleck Forceps Left
119909	Silicone wedge
119989	Instrument Tray



MA3511

Surfix® screw extractor 3.5 mm

Post-Operative Shoes

Reference	Description
Flat Post Operative Shoe	
999001	Small
999002	Medium
999003	Large
999004	X-Large



Reference	Description
Wegged Hell Post Operative Shoe	
999100	Size 4 to 6.5
999101	Size 7 to 8.5
999102	Size 9 to 10.5







Integra® Upper Extremities Solutions

Pyrocarbon PIP Implant

The Pyrocarbon PIP is a 2-component total joint replacement made from biocompatible Pyrocarbon material.

PIP is indicated for use in arthroplasty of the proximal interphalangeal (PIP) joint when the patient:

- Has soft tissue and bone that can provide adequate stabilization and fixation under high-demand loading conditions after reconstruction;
- Needs a revision of a failed PIP prosthesis, or has pain, limited motion, or joint subluxation/dislocation secondary to damage or destruction of the articular cartilage.



PIP implants

References	Description
PIP-200-10D-WW	Proximal Interphalangeal Size 10 Distal
PIP-200-10P-WW	Proximal Interphalangeal Size 10 Proximal
PIP-200-20D-WW	Proximal Interphalangeal Size 20 Distal
PIP-200-20P-WW	Proximal Interphalangeal Size 20 Proximal
PIP-200-30D-WW	Proximal Interphalangeal Size 30 Distal
PIP-200-30P-WW	Proximal Interphalangeal Size 30 Proximal
PIP-200-40D-WW	Proximal Interphalangeal Size 40 Distal
PIP-200-40P-WW	Proximal Interphalangeal Size 40 Proximal



Optional Instrument Set Components

References	Description
OSG-200-10PV	PIP Volar approach osteotomy guide, Size 10
OSG-200-20PV	PIP Volar approach osteotomy guide, Size 20
OSG-200-30PV	PIP Volar approach osteotomy guide, Size 30
OSG-200-40PV	PIP Volar approach osteotomy guide, Size 40



PIP Instrument Set

References	Description
CSA-200-02	Instrument Base
CSA-000-01	Instrument Lid
CSA-200-03	Instrument Insert
ALG-100-00	Alignment Guide, E/M
AWL-100-01	Awl, Starter
AWL-200-00	PIP Alignment Awl
BRH-200-10D	PIP, Cutting Broach, Distal, Size 10
BRH-200-10P	PIP, Cutting Broach, Proximal, Size 10
BRH-200-20D	PIP, Cutting Broach, Distal, Size 20
BRH-200-20P	PIP, Cutting Broach, Proximal, Size 20
BRH-200-30D	PIP, Cutting Broach, Distal, Size 30
BRH-200-30P	PIP, Cutting Broach, Proximal, Size 30
BRH-200-40D	PIP, Cutting Broach, Distal, Size 40
BRH-200-40P	PIP, Cutting Broach, Proximal, Size 40
EXT-200-00	Extractor
IMP-200-00D	PIP, Impactor, Distal
IMP-200-00P	PIP, Impactor, Proximal
OSG-200-00	PIP Vertical Cut Guide
OSG-200-10P	PIP, Proximal Cut Guide, Size 10
OSG-200-20P	PIP, Proximal Cut Guide, Size 20
OSG-200-30P	PIP, Proximal Cut Guide, Size 30
OSG-200-40P	PIP, Proximal Cut Guide, Size 40
TMP-200-01	PIP Distal Template
TRL-202-10D	PIP Distal Ti Trial Size 10
TRL-202-10P	PIP Proximal Ti Trial Size 10
TRL-202-20D	PIP Distal Ti Trial Size 20
TRL-202-20P	PIP Proximal Ti Trial Size 20
TRL-202-30D	PIP Distal Ti Trial Size 30
TRL-202-30P	PIP Proximal Ti Trial Size 30
TRL-202-40D	PIP Distal Ti Trial Size 40
TRL-202-40P	PIP Proximal Ti Trial Size 40
INS-200-01	Complete Instrument Set


By Request it is possible to order a set with Volar saw guides that have to be used when performing the volar approach and additionally ordered

Silicon PIP Implant

Silicone PIP Finger Joint Prosthesis is intended for cementless replacement of the proximal interphalangeal (PIP) joint, where disabled by rheumatoid, degenerative, or traumatic arthritis.

SPIP implants

References	Description
SPIP-520-0-WW	Silicone PIP, Size 0
SPIP-520-1-WW	Silicone PIP, Size 1
SPIP-520-2-WW	Silicone PIP, Size 2
SPIP-520-3-WW	Silicone PIP, Size 3
SPIP-520-4-WW	Silicone PIP, Size 4
SPIP-520-5-WW	Silicone PIP, Size 5



SPIP Instrument Set

References	Description
CSA-520-05	Instrument Base
CSA-000-04	Instrument Lid
CSA-520-07	Instrument tray Base
CSA-520-06	Instrument tray Lid
ALG-100-00	Alignment Guide, E/M
AWL-520-01	Starter Awl
AWL-200-00	PIP Alignment Awl
TRL-520-0	Silicone PIP Sizing Trial Size 0
TRL-520-1	Silicone PIP Sizing Trial Size 1
TRL-520-2	Silicone PIP Sizing Trial Size 2
TRL-520-3	Silicone PIP Sizing Trial Size 3
TRL-520-4	Silicone PIP Sizing Trial Size 4
TRL-520-5	Silicone PIP Sizing Trial Size 5
OSG-442-00	THB Vertical Cut Guide
BRH-525-0P	Silicone PIP Proximal Broach Size 0
BRH-525-1P	Silicone PIP Proximal Broach Size 1
BRH-525-2P	Silicone PIP Proximal Broach Size 2
BRH-525-3P	Silicone PIP Proximal Broach Size 3
BRH-525-4P	Silicone PIP Proximal Broach Size 4
BRH-525-5P	Silicone PIP Proximal Broach Size 5
BRH-525-0D	Silicone PIP Distal Broach Size 0
BRH-525-1D	Silicone PIP Distal Broach Size 1
BRH-525-2D	Silicone PIP Distal Broach Size 2
BRH-525-3D	Silicone PIP Distal Broach Size 3
BRH-525-4D	Silicone PIP Distal Broach Size 4
BRH-525-5D	Silicone PIP Distal Broach Size 5
INS-520-01	Complete set

Pyrocarbon MCP Implant


The MCP set contains:

- 2 Component anatomical total joint replacement made out of biocompatible Pyrocarbon.
- MCP Silicon spacer implant silicone spacer.

MCP is indicated for use as a total joint replacement of index, long, ring, and small finger metacarpophalangeal (MCP) joints that exhibit symptoms of pain, limited motion, or inadequate bony alignment (i.e., subluxation/dislocation) secondary to articular destruction or degenerative disease related to rheumatoid arthritis, systemic lupus erythematosus, osteoarthritis, or post-traumatic arthritis where soft tissue reconstruction can provide adequate stabilization.

MCP implants

References	Description
MCP-110-05D-WW	Metacarpophalangeal Joint – Size 05, Distal
MCP-110-05P-WW	Metacarpophalangeal Joint – Size 05, Proximal
MCP-100-10D-WW	Metacarpophalangeal Joint – Size 10, Distal
MCP-100-10P-WW	Metacarpophalangeal Joint – Size 10, Proximal
MCP-100-20D-WW	Metacarpophalangeal Joint – Size 20, Distal
MCP-100-20P-WW	Metacarpophalangeal Joint – Size 20, Proximal
MCP-100-30D-WW	Metacarpophalangeal Joint – Size 30, Distal
MCP-100-30P-WW	Metacarpophalangeal Joint – Size 30, Proximal
MCP-100-40D-WW	Metacarpophalangeal Joint – Size 40, Distal
MCP-100-40P-WW	Metacarpophalangeal Joint – Size 40, Proximal
MCP-100-50D-WW	Metacarpophalangeal Joint – Size 50, Distal
MCP-100-50P-WW	Metacarpophalangeal Joint – Size 50, Proximal




Instrument Set, MCP and SMCP 05-40

References	Description
CSA-100-05	Instrument Base
CSA-000-01	Instrument Lid
CSA-100-06	Instrument Tray
ALG-100-00	Alignment Guide, E/M
AWL-200-00	Alignment Awl
AWL-100-01	Awl, Starter
BRH-100-05D	MCP Distal Broach Sz. 05
BRH-100-05P	MCP Proximal Broach Sz. 05
BRH-100-10D	MCP, Broach, Distal, Sz. 10
BRH-100-10P	MCP, Broach, Proximal, Sz. 10
BRH-100-20D	MCP, Broach, Distal, Sz. 20
BRH-100-20P	MCP, Broach, Proximal, Sz. 20
BRH-100-30D	MCP, Broach, Distal, Sz. 30
BRH-100-30P	MCP, Broach, Proximal, Sz. 30
BRH-100-40D	MCP, Broach, Distal, Sz. 40
BRH-100-40P	MCP, Broach, Proximal, Sz. 40
EXT-200-00	Extractor
EXT-100-01	MCP, Extractor, Implant
IMP-100-00D	MCP, Impactor, Distal
IMP-100-00P	MCP, Impactor, Proximal
OSG-100-00D	MCP, Osteotomy Guide, Distal
OSG-100-00P	MCP, Osteotomy Guide, Proximal
TRL-100-05D	MCP Distal Trial Sz. 05
TRL-100-05P	MCP Proximal Trial Sz. 05
TRL-100-10D	MCP, Trial, Distal, Sz. 10
TRL-100-10P	MCP, Trial, Proximal, Sz. 10
TRL-100-20D	MCP, Trial, Distal, Sz. 20
TRL-100-20P	MCP, Trial, Proximal, Sz. 20
TRL-100-30D	MCP, Trial, Distal, Sz. 30
TRL-100-30P	MCP, Trial, Proximal, Sz. 30
TRL-100-40D	MCP, Trial, Distal, Sz. 40
TRL-100-40P	MCP, Trial, Proximal, Sz. 40
TRL-500-05	SMCP Trial Sz. 05
TRL-500-10	SMCP Trial Sz 10
TRL-500-20	SMCP Trial Sz 20
TRL-500-30	SMCP Trial Sz 30
TRL-500-40	SMCP Trial Sz 40
INS-100-02	Complete set

Optional Instrument Set Components

References	Description
BRH-100-50D	MCP Distal Broach Size 50
BRH-100-50P	MCP Proximal Broach Size 50
TRL-100-50D	MCP Distal Trial Size 50
TRL-100-50P	MCP Proximal Trial Size 50
TRL-500-50	SMCP Trial, Sz 50



Silicon MCP Implant

Silicone MCP Finger Joint Prosthesis is intended for cementless replacement of the metacarpophalangeal (MCP) joint, where disabled by rheumatoid, degenerative, or traumatic arthritis.

MCP System is used for proper placement. See MCP System for instrument set.

- Instrumentation for implant placement is the same set as for the Pyrocarbon MCP implant.
- Surgical technique and sequence is identical to MCP pyrocarbon implant (NB Possible soft tissue repair).

SMCP implants

References	Description
SMCP-500-05-WW	Silicone MCP, Size 05
SMCP-500-10-WW	Silicone MCP, Size 10
SMCP-500-20-WW	Silicone MCP, Size 20
SMCP-500-30-WW	Silicone MCP, Size 30
SMCP-500-40-WW	Silicone MCP, Size 40
SMCP-500-50-WW	Silicone MCP, Size 50



Instrument Set, MCP and SMCP 05-40

References	Description
CSA-100-05	Instrument Base
CSA-000-01	Instrument Lid
CSA-100-06	Instrument Tray
ALG-100-00	Alignment Guide, E/M
AWL-200-00	Alignment Awl
AWL-100-01	Awl, Starter
BRH-100-05D	MCP Distal Broach Sz. 05
BRH-100-05P	MCP Proximal Broach Sz. 05
BRH-100-10D	MCP, Broach, Distal, Sz. 10
BRH-100-10P	MCP, Broach, Proximal, Sz. 10
BRH-100-20D	MCP, Broach, Distal, Sz. 20
BRH-100-20P	MCP, Broach, Proximal, Sz. 20
BRH-100-30D	MCP, Broach, Distal, Sz. 30
BRH-100-30P	MCP, Broach, Proximal, Sz. 30
BRH-100-40D	MCP, Broach, Distal, Sz. 40
BRH-100-40P	MCP, Broach, Proximal, Sz. 40
EXT-200-00	Extractor
EXT-100-01	MCP, Extractor, Implant
IMP-100-00D	MCP, Impactor, Distal
IMP-100-00P	MCP, Impactor, Proximal
OSG-100-00D	MCP, Osteotomy Guide, Distal
OSG-100-00P	MCP, Osteotomy Guide, Proximal
TRL-100-05D	MCP Distal Trial Sz. 05
TRL-100-05P	MCP Proximal Trial Sz. 05
TRL-100-10D	MCP, Trial, Distal, Sz. 10
TRL-100-10P	MCP, Trial, Proximal, Sz. 10
TRL-100-20D	MCP, Trial, Distal, Sz. 20
TRL-100-20P	MCP, Trial, Proximal, Sz. 20
TRL-100-30D	MCP, Trial, Distal, Sz. 30
TRL-100-30P	MCP, Trial, Proximal, Sz. 30
TRL-100-40D	MCP, Trial, Distal, Sz. 40
TRL-100-40P	MCP, Trial, Proximal, Sz. 40
TRL-500-05	SMCP Trial Sz. 05
TRL-500-10	SMCP Trial Sz 10
TRL-500-20	SMCP Trial Sz 20
TRL-500-30	SMCP Trial Sz 30
TRL-500-40	SMCP Trial Sz 40
INS-100-02	Complete set

Optional Instrument Set Components

References	Description
BRH-100-50D	MCP Distal Broach Size 50
BRH-100-50P	MCP Proximal Broach Size 50
TRL-100-50D	MCP Distal Trial Size 50
TRL-100-50P	MCP Proximal Trial Size 50
TRL-500-50	SMCP Trial, Sz 50

The Thumb Carpometacarpal (CMC) Joint

Integra® CMC Implants

PyroDisk™

The PyroDisk™ CMC implant is intended to replace the joint between the first metacarpal and the trapezium in cases of rheumatoid arthritis, traumatic arthritis, osteoarthritis or post fracture deformation or bone loss which present as either a painful, unstable thumb, or a thumb with limited range of motion.



PyroDisk™ implants

References	Description
PYD-420-145-WW	PYD Comp. Size 145
PYD-420-165-WW	PYD Comp. Size 165
PYD-420-167-WW	PYD Comp. Size 167
PYD-420-168-WW	PYD Comp. Size 168
PYD-420-187-WW	PYD Comp. Size 187
PYD-420-189-WW	PYD Comp. Size 189

PyroDisk™ instrument set

References	Description
CSA-420-02	Instrument Base
CSA-000-05	Instrument Lid
AWL-420-00	Starter Awl
GOU-420-00	Gouge
SHP-425-14	PyroDisk Shaper Size 14

SHP-425-16	PyroDisk Shaper Size 16
SHP-425-18	PyroDisk Shaper Size 18
TRL-420-145	PYD Trial Size 145
TRL-420-165	PYD Trial Size 165
TRL-420-167	PYD Trial Size 167
TRL-420-168	PYD Trial Size 168
TRL-420-187	PYD Trial Size 187
TRL-420-189	PYD Trial Size 189
INS-420-01	Complete set



NuGrip™


NuGrip™ CMC Implant is intended to replace the proximal end of the first metacarpal in cases of rheumatoid arthritis, traumatic arthritis, osteoarthritis, post fracture deformity or bone loss which presents as either a painful, unstable thumb, or a thumb with limited range of motion.

NuGrip™ implants

References	Description
NUG-443-10S	NuGrip™ Implant Size 10-S
NUG-443-10M	NuGrip™ Implant Size 10-M
NUG-443-20S	NuGrip™ Implant Size 20-S
NUG-443-20M	NuGrip™ Implant Size 20-M
NUG-443-30M	NuGrip™ Implant Size 30-M
NUG-443-30L	NuGrip™ Implant Size 30-L
NUG-443-40M	NuGrip™ Implant Size 40-M
NUG-443-40L	NuGrip™ Implant Size 40-L

NuGrip™ Disposables

References	Description
CMC-DIS-HAL	CMC Disposable Pack for Hall Micro Systems
CMC-DIS-STR	CMC Disposable Pack for Stryker Command II and TPS Micro Drill
CUT-443-S	Cannulated Spherical Cutter Size S
CUT-443-M	Cannulated Spherical Cutter Size M
CUT-443-L	Cannulated Spherical Cutter Size L



NuGrip™ Instruments

References	Description
CSA-443-02	Instrument Base
CSA-000-04	Instrument Lid
CSA-443-03	Instrument Tray
CSA-443-04	Instrument Tray Lid
ALG-100-00	Alignment Guide
AWL-100-01	Starter Awl
AWL-200-00	Alignment Awl
BRH-442-10	Cutting Broach Size 10
BRH-442-20	Cutting Broach Size 20
BRH-442-30	Cutting Broach Size 30
BRH-442-40	Cutting Broach Size 40
EXT-200-00	Trial Extractor
IMP-442-00	NuGrip™ Impactor
OSG-442-00	Vertical Cut Guide
SHP-443-S	Finish Shaper Size S
SHP-443-M	Finish Shaper Size M
SHP-443-L	Finish Shaper Size L
TMP-443-S	Size Template Size S
TMP-443-M	Size Template Size M
TMP-443-L	Size Template Size L
TRL-443-10S	Metal NuGrip™ Trial Size 10S
TRL-443-10M	Metal NuGrip™ Trial Size 10M
TRL-443-20S	Metal NuGrip™ Trial Size 20S
TRL-443-20M	Metal NuGrip™ Trial Size 20M
TRL-443-30M	Metal NuGrip™ Trial Size 30M
TRL-443-30L	Metal NuGrip™ Trial Size 30L
TRL-443-40M	Metal NuGrip™ Trial Size 40M
TRL-443-40L	Metal NuGrip™ Trial Size 40L
EXT-100-01	Implant Extractor
RET-443-00	Hohman Retractor
INS-443-00	Complete set

CMC Saddle™ Implant

Saddle™ CMC implant is intended to resurface the proximal end of the first metacarpal in cases of rheumatoid arthritis, traumatic arthritis, osteo-arthritis, or post-fracture deformation or bone loss which present as either a painful, unstable thumb, or a thumb with limited range of motion. This is an uncemented prosthesis designed for press fit use only.



Saddle™ CMC implant

References	Description
CMC-410-10	CMC Distal Hemi, Size 10
CMC-410-20	CMC Distal Hemi, Size 20
CMC-410-30	CMC Distal Hemi, Size 30
CMC-410-40	CMC Distal Hemi, Size 40

Instrument Set, Saddle™ CMC

References	Description
CSA-410-02	Instrument Base
CSA-000-01	Instrument Lid
ALG-100-00	Alignment Guide, E/M
AWL-100-01	Awl, Starter
AWL-200-00	Alignment Awl
BRH-410-10	CMC Broach Size 10
BRH-410-20	CMC Broach Size 20
BRH-410-30	CMC Broach Size 30
BRH-410-40	CMC Broach Size 40
EXT-200-00	Extractor
IMP-410-00	CMC Hemi Impactor
OSG-200-00	Vertical Cut Guide
TMP-410-01	CMC Sizing Template
TRL-410-10M	CMC Hemi Trial Metal Size 10
TRL-410-20M	CMC Hemi Trial Metal Size 20
TRL-410-30M	CMC Hemi Trial Metal Size 30
TRL-410-40M	CMC Hemi Trial Metal Size 40
INS-410-00	Complete set




PCS Spacer*

PyroSphere is intended to replace the joint between the first metacarpal and the trapezium in cases of rheumatoid arthritis, traumatic arthritis, osteoarthritis or post-fracture deformation or bone loss which present as either a painful, unstable thumb, or a thumb with limited range of motion.

PyroSphere is also intended for use in the fourth/fifth tarsometatarsal (TMT) joint involvement where degenerative or post-traumatic arthritis presents:

- decreased motion.
- arthritic changes and/or subluxation.
- unstable, stiff, or painful joints.
- degenerative joint disease of the midfoot.
- associated with gout or pseudogout.

References	Description
PCS-430-10-WW	PyroSphere Size 10
PCS-430-20-WW	PyroSphere Size 20
PCS-430-30-WW	PyroSphere Size 30
PCS-430-40-WW	PyroSphere Size 40
PCS-430-50-WW	PyroSphere Size 50



Instrument Set - PyroSphere™

References	Description
CSA-430-01 CSA-000-03	Instrument Base and Lid
EXT-430-00	Implant Extractor
TMP-400-10	Implant sizer instrument Size 10
TMP-400-20	Implant sizer instrument Size 20
TMP-400-30	Implant sizer instrument Size 30
TMP-400-40	Implant sizer instrument Size 40
TMP-400-50	Implant sizer instrument Size 50
INS-430-00	Instrument Set - PyroSphere

Options Items for Instrument Set

References	Description
TRL-435-10	Pyrosphere Trial Size 10
TRL-435-20	Pyrosphere Trial Size 20
TRL-435-30	Pyrosphere Trial Size 30
TRL-435-40	Pyrosphere Trial Size 40
TRL-435-50	Pyrosphere Trial Size 50
BRH-445-10	Finishing Shapers Size 10
BRH-445-20	Finishing Shapers Size 20
BRH-445-30	Finishing Shapers Size 30
BRH-445-40	Finishing Shapers Size 40
BRH-445-50	Finishing Shapers Size 50

* Available in 5 different sizes

Safeguard™ Mini-Open Carpal Tunnel Release System

Indications for use:

- The Safeguard™ System is used for the treatment of carpal tunnel syndrome in patients who fail to respond to conservative treatment.

Instruments

References	Description
080001	SafeGuard Guide

Disposable

References	Description
080003	SafeGuard Knife (Individual)
080005	SafeGuard 5-pack knives



Kompressor™ Compression Screw system

The use of the Kompressor compression screw is generally indicated for the reduction and fixation of fractures appropriate for the size of the device. It is indicated for use in the internal fixation of fractures, fusion, and revision. It is also indicated for reconstructive procedures where reduction and fixation of bone fragments is required (e.g. osteotomies).

Kompressor™ Screws system: 4.0 mm screws

References	Description
210414	Kompressor™ screw 14 mm length
210416	Kompressor™ screw 16 mm length
210418	Kompressor™ screw 18 mm length
210420	Kompressor™ screw 20 mm length
210422	Kompressor™ screw 22 mm length
210424	Kompressor™ screw 24 mm length
210426	Kompressor™ screw 26 mm length
210428	Kompressor™ screw 28 mm length
210430	Kompressor™ screw 30 mm length

Kompressor™ Screws system: 4.0 mm screws

References	Description
210432	Kompressor™ screw 32 mm length
210434	Kompressor™ screw 34 mm length
210310	Mini-Kompressor™ screw, 10 mm length
210312	Mini-Kompressor™ screw, 12 mm length
210314	Mini-Kompressor™ screw, 14 mm length
210316	Mini-Kompressor™ screw, 16 mm length
210318	Mini-Kompressor™ screw, 18 mm length
210320	Mini-Kompressor™ screw, 20 mm length
210322	Mini-Kompressor™ screw, 22 mm length
210324	Mini-Kompressor™ screw, 24 mm length
210326	Mini-Kompressor™ screw, 26 mm length



Disposables (Kompressor™)

References	Description
211622	2.2 mm Cannulated Drill Bit
211628	2.8 mm Cannulated Drill Bit
211409	0.9 mm K-wires

Disposables (Mini Kompressor™)

211625	Drill Bit, 2.5 mm, Cannulated
211635	3.5 mm Cannulated Drill Bit
211411	1.1 mm K-wires

Instruments

References	Description
211100	Handle, Drill Guide Assembly
211200	Guide Sleeve Assembly
211300	Center Punch
211501	K-Wire Depth Gauge, STD.
211511	K-Wire Depth Gauge, Mini
211600	3.5 mm Drill Guide
211601	2.5 mm Drill Guide
211610	2.8 mm Drill Guide
211611	2.2 mm Drill Guide
211700	Handle, AO couple
211725	Standard Placement Driver
211720	Mini Placement Driver
211801	Standard Solid Leading Driver
211811	Mini Solid Leading Driver
211900	Standard Trailing Driver
211910	Mini Trailing Driver

Sterilization container

215001	Sterilization Tray
--------	--------------------

Hand Fractures and Deformity Repair Set

It is a set to be used for fractures.

Only one instrumentation module in order to set up 4 screws:


- Bold® Compression Screw:
 - > 2.5 mm
 - > 3.0 mm
 - > 3.7 mm
- Qwix® Fixation Screw:
 - > 3.0 mm

And two associated optional modules :

- Implant module: for non sterile screws.
- Generic module: for other personal instrumentation and support for 2 or 3 floors for the chosen modules.

Support

References	Description
229996	Support - 2 floors
229990	Support - 3 floors
229970	Forefoot II- lid- generic



Upper Screw Module

Instrument module

229002	Screwdriver hexa 2.0 - cannulated
229005	Screwdriver- AO attachment -hexa 2.0 - cannulated
229003	Screwdriver- T7 star - dia. 0.9 mm - cannulated
229004	Screwdriver- AO attachment - T7 star - dia. 0.9 mm - cannulated
997311	Newdeal - handle - AO straight
229101	Ruler
997301	Depth Gauge
115070S	K-Wire dia. 1.0 L. 70 mm - Sterile
115100S	K-Wire dia. 1.0 L. 100 mm - Sterile
115008S	K-Wire dia. 0.8 L. 70 mm - Sterile
119023S	2 in 1 cannulated drill dia. 2.2/3.0 x L. 32 mm
119025S	2 in 1 cannulated drill dia. 2.2/3.0 x L. 22 mm
119027S	2 in 1 cannulated drill dia. 2.2/3.0 x L. 12 mm
229227S	2 in 1 cannulated drill dia. 1.9 x L. 12 mm
229237S	2 in 1 cannulated drill dia. 1.9 x L. 22 mm
229961	Forefoot II - tray - screw module
229992	Forefoot II - mat - screw module



Upper Screw Module

Sterilization Container

229965	Forefoot II - tray - non sterile implants
229975	Forefoot II - rack - non sterile implants



Upper Generic Module

Sterilization Container

229966	Forefoot II - tray - generic
229976	Forefoot II - mat - generic



Other Possible Drills

Reference	Description
159027S	2in1 Drill - dia. 2.2 mm x L12 mm - AO - Cannulated
159025S	2in1 Drill - dia. 2.2 mm x L22 mm - AO - Cannulated
159023S	2in1 Drill - dia. 2.2 mm x L32 mm - AO - Cannulated
119227S	2in1 Drill - dia. 2.2 mm x L12 mm - Non Cannulated
119225S	2in1 Drill - dia. 2.2 mm x L22 mm - Non Cannulated
119223S	2in1 Drill - dia. 2.2 mm x L32 mm - Non Cannulated
119027S	2in1 Drill - dia. 2.2 mm x L12 mm - Cannulated
119025S	2in1 Drill - dia. 2.2 mm x L22 mm - Cannulated
119023S	2in1 Drill - dia. 2.2 mm x L32 mm - Cannulated
159227S	2in1 Drill - dia. 2.2 mm x L12 mm - AO - Non Cannulated
159225S	2in1 Drill - dia. 2.2 mm x L22 mm - AO - Non Cannulated
159223S	2in1 Drill - dia. 2.2 mm x L32 mm - AO - Non Cannulated
229027S	2in1 Drill - dia. 1.9 mm x L12 mm - AO - Cannulated
229037S	2in1 Drill - dia. 1.9 mm x L22 mm - AO - Cannulated
229227S	2in1 Drill - dia. 1.9 mm x L12 mm - Cannulated
229237S	2in1 Drill - dia. 1.9 mm x L22 mm - Cannulated
229028S	2in1 Drill - dia. 1.9 mm x L12 mm - AO - Non Cannulated
229038S	2in1 Drill - dia. 1.9 mm x L22 mm - AO - Non Cannulated
229228S	2in1 Drill - dia. 1.9 mm x L12 mm - Non Cannulated
229238S	2in1 Drill - dia. 1.9 mm x L22 mm - Non Cannulated



QWIX® Fixation Screw, dia. 3.0 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction:

- Self-drilling and self-tapping.
- Totally intra-osseous fixation.
- K-Wire guided drilling and insertion.
- Material: Titanium alloy TA6V.

Qwix® dia. 3.0 mm

References	Description
111312(S)	L. 12 mm
111314(S)	L. 14 mm
111316(S)	L. 16 mm
111318(S)	L. 18 mm
111320(S)	L. 20 mm
111322(S)	L. 22 mm
111324(S)	L. 24 mm
111326(S)	L. 26 mm
111328(S)	L. 28 mm
111330(S)	L. 30 mm
111332(S)	L. 32 mm
111334(S)	L. 34 mm


Bold® Compression Screw dia. 2.5 mm / 3.0 mm / 3.7 mm (for associated instrumentation, see p. 8)

For fixation of bone fractures or for bone reconstruction:

- Totally intra-osseous implants.
- Large range of sizes:
 - > 10 to 30 mm for Bold® 2.5 by 2 mm increment.
 - > 10 to 34 mm for Bold® 3.0 by 2 mm increment.
 - > 14 to 34 mm for Bold® 3.7 by 2 mm increment.
- K-Wire guided drilling and insertion.
- Self tapping.
- True & controlled compression.
- Dual threads.
- Material: Titanium alloy TA6V.



Bold® dia 2.5 mm

Reference	Length
200010 (S)	L. 10 mm
200012 (S)	L. 12 mm
200014 (S)	L. 14 mm
200016 (S)	L. 16 mm
200018 (S)	L. 18 mm
200020 (S)	L. 20 mm
200022 (S)	L. 22 mm
200024 (S)	L. 24 mm
200026 (S)	L. 26 mm
200028 (S)	L. 28 mm
200030 (S)	L. 30 mm



Bold® dia 3.0 mm

Reference	Length
111010 (S)	L. 10 mm
111012 (S)	L. 12 mm
111014 (S)	L. 14 mm
111016 (S)	L. 16 mm
111018 (S)	L. 18 mm
111020 (S)	L. 20 mm
111022 (S)	L. 22 mm
111024 (S)	L. 24 mm
111026 (S)	L. 26 mm
111028 (S)	L. 28 mm
111030 (S)	L. 30 mm
111032 (S)	L. 32 mm
111034 (S)	L. 34 mm



Bold® dia 3.7 mm

Reference	Length
131014 (S)	L. 14 mm
131016 (S)	L. 16 mm
131018 (S)	L. 18 mm
131020 (S)	L. 20 mm
131022 (S)	L. 22 mm
131024 (S)	L. 24 mm
131026 (S)	L. 26 mm
131028 (S)	L. 28 mm
131030 (S)	L. 30 mm
131032 (S)	L. 32 mm
131034 (S)	L. 34 mm

UNI-CLIP® Notched & Standard Staples (for associated instrumentation, see p. 9)

For fixation of bone fractures or for bone reconstruction.

- Range sizes:
 - > Interaxis: 11, 12, 13, 15 and 20 mm.
 - > Leg lengths: 12, 13, 14, 15, 16, 17 and 20 mm.
- Material: Stainless steel 316L.



Standard Staples

Reference	Description
113113	Uni-Clip® interaxis 11.0 x L. 13 mm
113114	Uni-Clip® interaxis 11.0 x L. 14 mm
113115	Uni-Clip® interaxis 11.0 x L. 15 mm
113116	Uni-Clip® interaxis 11.0 x L. 16 mm
113117	Uni-Clip® interaxis 11.0 x L. 17 mm
113213	Uni-Clip® interaxis 12.0 x L. 13 mm
113214	Uni-Clip® interaxis 12.0 x L. 14 mm
113215	Uni-Clip® interaxis 12.0 x L. 15 mm
113216	Uni-Clip® interaxis 12.0 x L. 16 mm
113217	Uni-Clip® interaxis 12.0 x L. 17 mm
113313	Uni-Clip® interaxis 13.0 x L. 13 mm
113314	Uni-Clip® interaxis 13.0 x L. 14 mm
113315	Uni-Clip® interaxis 13.0 x L. 15 mm
113316	Uni-Clip® interaxis 13.0 x L. 16 mm
113317	Uni-Clip® interaxis 13.0 x L. 17 mm
113512	Uni-Clip® interaxis 15.0 x L. 12 mm
113820	Uni-Clip® interaxis 20.0 x L. 20 mm

Notched Staples

Reference	Description
213113(S)	Uni-Clip® interaxis 11.0 x L. 13 mm
213114(S)	Uni-Clip® interaxis 11.0 x L. 14 mm
213115(S)	Uni-Clip® interaxis 11.0 x L. 15 mm
213116(S)	Uni-Clip® interaxis 11.0 x L. 16 mm
213117(S)	Uni-Clip® interaxis 11.0 x L. 17 mm
213213(S)	Uni-Clip® interaxis 12.0 x L. 13 mm
213214(S)	Uni-Clip® interaxis 12.0 x L. 14 mm
213215(S)	Uni-Clip® interaxis 12.0 x L. 15 mm
213216(S)	Uni-Clip® interaxis 12.0 x L. 16 mm
213217(S)	Uni-Clip® interaxis 12.0 x L. 17 mm
213313(S)	Uni-Clip® interaxis 13.0 x L. 13 mm
213314(S)	Uni-Clip® interaxis 13.0 x L. 14 mm
213315(S)	Uni-Clip® interaxis 13.0 x L. 15 mm
213316(S)	Uni-Clip® interaxis 13.0 x L. 16 mm
213317(S)	Uni-Clip® interaxis 13.0 x L. 17 mm
213512(S)	Uni-Clip® interaxis 15.0 x L. 12 mm
213820(S)	Uni-Clip® interaxis 20.0 x L. 20 mm

SPIN[®] SNAP-OFF[®] Screw, dia. 2.0 mm (for associated instrumentation, see p. 8)

For bone fixation.

- Break-off barrel.
- Self-drilling and self-tapping.
- Material: Titanium alloy TA6V.

Reference	Description
112011(S)	Spin [®] dia. 2.0 x L. 11 mm
112012(S)	Spin [®] dia. 2.0 x L. 12 mm
112013(S)	Spin [®] dia. 2.0 x L. 13 mm
112014(S)	Spin [®] dia. 2.0 x L. 14 mm



Universal2™ Total Wrist Implant System

The Universal2™ Total Wrist System is indicated for use in patients suffering from pain and/or loss of function due to:

- Traumatic Arthritis
- Osteoarthritis
- Rheumatoid Arthritis
- SLAC Wrist

The Universal2™ Total Wrist may also be indicated for use in the revision of a failed implant, in younger patients with articular joint disease, or in difficult situations where clinical experience indicates that other reconstructive efforts are likely to achieve unsatisfactory results. Use should be limited to those conditions where preoperative exam indicates sufficient bone stock for the support of the implants.



References	Description
Radial Implants	
263100R	Radial Implant, Extra Small, Right
263100L	Radial Implant, Extra Small, Left
263200R	Radial Implant, Small, Right
263200L	Radial Implant, Small, Left
263300R	Radial Implant, Medium, Right
263300L	Radial Implant, Medium, Left
263400R	Radial Implant, Large, Right
263400L	Radial Implant, Large, Left
Carpal Implants	
261101	Carpal Implant, Extra Small
261201	Carpal Implant, Small
261301	Carpal Implant, Medium
261401	Carpal Implant, Large

References	Description
Poly Implants	
262100	Carpal Poly Implant, Extra Small Standard
262200	Carpal Poly Implant, Small, Standard
262300	Carpal Poly Implant, Medium Standard
262400	Carpal Poly Implant, Large, Standard
262101	Carpal Poly Implant, Extra Small +1
262201	Carpal Poly Implant, Small +1
262301	Carpal Poly Implant, Medium +1
262401	Carpal Poly Implant, Large +1
262102	Carpal Poly Implant, Extra Small +2
262202	Carpal Poly Implant, Small +2
262302	Carpal Poly Implant, Medium +2
262402	Carpal Poly Implant, Large +2

References	Description
Non Sterile Screws	
264515	4.5 mm Bone Screw, 15 mm non sterile
264520	4.5 mm Bone Screw, 20 mm non sterile
264525	4.5 mm Bone Screw, 25 mm non sterile
264530	4.5 mm Bone Screw, 30 mm non sterile
264535	4.5 mm Bone Screw, 35 mm non sterile
Sterile Screws	
264515S	4.5 mm Bone Screw, 15 mm –sterile
264520S	4.5 mm Bone Screw, 20 mm –sterile
264525S	4.5 mm Bone Screw, 25 mm –sterile
264530S	4.5 mm Bone Screw, 30 mm –sterile
264535S	4.5 mm Bone Screw, 35 mm –sterile

Universal2™ Total Wrist Implant System – Instruments



References	Description
Disposables	
268111	1.1 mm K-Wire 100 mm
268114	1.4 mm K-Wire 150 mm
268232	Radial Guide Rod, 150 mm
268233	Radial Guide Rod, 200 mm
268125	Drill Bit, 2.5 mm, Solid
211625	Drill Bit, 2.5 mm, Cannulated
268136	Drill Bit, 3.5 mm, Cannulated
268145	Drill Bit, 4.5 mm, Cannulated
Instruments	
260515	Trial screw, 4.5 mm Solid Cancellous 15 mm
260520	Trial screw, 4.5 mm Solid Cancellous 20 mm
260525	Trial screw, 4.5 mm Solid Cancellous 25 mm
260530	Trial screw, 4.5 mm Solid Cancellous 30 mm
260535	Trial screw, 4.5 mm Solid Cancellous 35 mm
265100	Trial, Carpal X-Small
265200	Trial, Carpal Small
265300	Trial, Carpal Medium
265400	Trial, Carpal Large
266100	Trial, Carpal Poly, X-Small

References	Description
266101	Trial Carpal Poly, X-Small +1
266102	Trial Carpal Poly, X-small +2
266200	Trial, Carpal Poly, Small
266201	Trial, Carpal Poly, Small +1
266202	Trial Carpal Poly, Small +2
266300	Trial, Carpal Poly, Medium
266301	Trial, Carpal Poly, Medium +1
266302	Trial Carpal Poly, Medium +2
266400	Trial, Carpal Poly, Large
266401	Trial Carpal Poly, Large +1
266402	Trial Carpal Poly, Large +2
267100LT	Trial, Radial X-Small, Left
267100RT	Trial, Radial X-Small, Right
267200LT	Trial, Radial Small, Left
267200RT	Trial, Radial Small, Right
267300LT	Trial, Radial Medium, Left
267300RT	Trial, Radial Medium, Right
267400LT	Trial, Radial Large, Left
267400RT	Trial, Radial Large, Right
268100	Modular Drill Guide
268120	Guide Plate Extra Small
268121	Guide Plate, Small
268122	Guide Plate, Medium
268123	Guide Plate Large
268220	Bone Awl

References	Description
268230LF	Radial Cutting Guide, Left
268230RT	Radial Cutting Guide, Right
268231	Radial Guide Bar
268325	Carpal Guide Bar, Minimal
268335	Carpal Guide Bar, Standard
268340	Carpal Cutting Guide
268410L	Broach, Radial X-Small, Left
268410R	Broach, Radial X-Small, Right
268420L	Broach, Radial Small, Left
268420R	Broach, Radial Small, Right
268430L	Broach, Radial Medium, Left
268430R	Broach, Radial Medium, Right
268440L	Broach, Radial Large, Left
268440R	Broach, Radial Large, Right
268510	Broach Handle
268610	T handle
268640	Carpal Poly Impactor
268650	Radial Impactor
268660	Poly Impactor
268700	AO Quick Couple Driver
268706	6mm Counter Sink
268725	Shaft, Hex Driver 2.5 mm
Sterilization container	
260100	Sterilization Tray
268104	Modular drill guide (K-Wire) sleeve

Integra® Total Wrist Fusion System ITWFS

The Integra® Total Wrist Fusion System is indicated for use in patients with:

- Post-traumatic arthritis of the joints of the wrist
- Rheumatoid wrist deformities requiring restoration
- Complex carpal instability
- Post-septic arthritis of the wrist
- Severe unremitting wrist pain related to motion
- Brachial plexus nerve palsies
- Tumor resection
- Spastic deformities
- Pain and/or loss of function due to osteoarthritis
- Revision of failed partial wrist fusions




ITWFS implants

References	Description
303100S	Straight plate
303101S	Standard bend plate
303102S	Short bend plate
303210CS	2.7 mm cortical screw 10 mm long
303212CS	2.7 mm cortical screw 12 mm long
303214CS	2.7 mm cortical screw 14 mm long
303216CS	2.7 mm cortical screw 16 mm long
303218CS	2.7 mm cortical screw 18 mm long
303220CS	2.7 mm cortical screw 20 mm long
303222CS	2.7 mm cortical screw 22 mm long
303224CS	2.7 mm cortical screw 24 mm long
303312CS	3.5 mm cortical screw 12 mm long
303314CS	3.5 mm cortical screw 14 mm long
303316CS	3.5 mm cortical screw 16 mm long
303318CS	3.5 mm cortical screw 18 mm long
303320CS	3.5 mm cortical screw 20 mm long
303322CS	3.5 mm cortical screw 22 mm long
303324CS	3.5 mm cortical screw 24 mm long
303326CS	3.5 mm cortical screw 26 mm long
303328CS	3.5 mm cortical screw 28 mm long
286210SND	2.7 mm surfix screw 10 mm long and 2.7 mm surfix washer
286212SND	2.7 mm surfix screw 12 mm long and 2.7 mm surfix washer
286214SND	2.7 mm surfix screw 14 mm long and 2.7 mm surfix washer
286216SND	2.7 mm surfix screw 16 mm long and 2.7 mm surfix washer
286218SND	2.7 mm surfix screw 18 mm long and 2.7 mm surfix washer
286220SND	2.7 mm surfix screw 20 mm long and 2.7 mm surfix washer

286222SND	2.7 mm surfix screw 22 mm long and 2.7 mm surfix washer
286224SND	2.7 mm surfix screw 24 mm long and 2.7 mm surfix washer
186200SND	2.7 mm surfix locking screw
286312SND	3.5 mm surfix screw 12 mm long and 3.5 mm surfix washer
286314SND	3.5 mm surfix screw 14 mm long and 3.5 mm surfix washer
286316SND	3.5 mm surfix screw 16 mm long and 3.5 mm surfix washer
286318SND	3.5 mm surfix screw 18 mm long and 3.5 mm surfix washer
286320SND	3.5 mm surfix screw 20 mm long and 3.5 mm surfix washer
286322SND	3.5 mm surfix screw 22 mm long and 3.5 mm surfix washer
286324SND	3.5 mm surfix screw 24 mm long and 3.5 mm surfix washer
286326SND	3.5 mm surfix screw 26 mm long and 3.5 mm surfix washer
286328SND	3.5 mm surfix screw 28 mm long and 3.5 mm surfix washer
186300SND	3.5 mm surfix locking screw

ITWFS instrumentation

References	Description
303400	2.0 mm drill bit
303401	2.5 mm drill bit
303402	2.7 mm drill bit
303403	Double-ended drill guide for 2.0 mm drill bit
303404	Double-ended drill guide for 2.5 mm drill bit
303405	Threaded drill guide for 2.0 mm drill bit
303406	Threaded drill guide for 2.7 mm drill bit
303407	Depth gauge
303408	Self-retaining torx driver
303409	2.5 mm self-retaining hex driver
303410	2.0 mm self-retaining hex driver
303411	Plate bender - (2.7 mm hole)
303412	Plate bender - (3.5 mm hole)
303450	Straight plate trial
303451	Standard bend plate trial
303452	Short bend plate trial
600731	Screw forceps



ITWFS set

References	Description
303510	Instrument case complete - EMEA
including:	
303502	Instrument case top level - US/EMEA
303504	Instrument case lid - US/EMEA

Spider™ – Mini-Spider™ / Limited Wrist Fusion System

SPIDER / MINI-SPIDER is indicated for use in patients suffering pain and/or loss of function due to:

- Osteoarthritis.
- Post-traumatic arthritis.
- Fractures.
- Revision of failed partial wrist fusions.
- Carpal Instability.
- Rheumatoid arthritis.

IMPLANTS Stainless steel

References	Description
070005	Spider™ Limited Wrist Fusion Plate Steel
070006	Mini Spider™ Limited Wrist Fusion Plate Steel
072408	2.4 mm Bone Screw, 8 mm Steel
072410	2.4 mm Bone Screw, 10 mm Steel
072412	2.4 mm Bone Screw, 12 mm Steel
072414	2.4 mm Bone Screw, 14 mm Steel
072810	2.8 mm Bone Screw, 10 mm Steel
072814	2.8 mm Bone Screw, 14 mm Steel
072818	2.8 mm Bone Screw, 18 mm Steel

IMPLANTS Titanium

References	Description
060005	Spider™ Limited Wrist Fusion Plate Titanium
060006	Mini Spider™ Limited Wrist Fusion Plate Titanium
062806	2.8 mm Titanium Bone Screw, 6 mm
062808	2.8 mm Titanium Bone Screw, 8 mm



062810	2.8 mm Titanium Bone Screw, 10 mm
062812	2.8 mm Titanium Bone Screw, 12 mm
062814	2.8 mm Titanium Bone Screw, 14 mm
062816	2.8 mm Titanium Bone Screw, 16 mm
062818	2.8 mm Titanium Bone Screw, 18 mm

Spider™ – Mini Spider™ / Limited Wrist Fusion System – Instruments


References	Description
Disposables	
600709	1.5 mm Drill Bit
600704	2.0 mm Drill Bit
600720	Std. Spider Rasp
600721	Mini Spider Rasp
Instruments	
600706	Adjustable Depth Gauge
600715	Spider Plate Holder
600716	Spider Rasp Cleaning Brush
600718	T-Handle
600719	Mini Spider Plate Holder
600722	1.5/2.0 mm Drill Guide
600724	2.5mm Screw Driver Bit QC
600726	Screw Driver Holding Sleeve
600730	4-1/2 Screw Driver Handle QC
600731	Screw Forceps
600740	Tap 4-1/2, 3.5 mm dia.
Sterilization container	
075020	Sterilization Case

PyroCarbon Lunate Implant

PyroCarbon Lunate is intended for replacement of the lunate bone in the proximal carpal row of the wrist in the presence of:

- Avascular necrosis (Kienboch's disease)
- Localized osteoarthritic changes
- Long-standing dislocations

The implant is available in 5 sizes for use in left or right applications.

Lunate Implants

References	Description
LUN-710-01-WW	PyroCarbon Lunate, Size 01
LUN-710-02-WW	PyroCarbon Lunate, Size 02
LUN-710-03-WW	PyroCarbon Lunate, Size 03
LUN-710-04-WW	PyroCarbon Lunate, Size 04
LUN-710-05-WW	PyroCarbon Lunate, Size 05



Lunate Instrument Set

References	Description
CSA-170-02	Instrument Base
CSA-710-01	Instrument Lid
SZT-715-01S	Sizing Template, Size 01
SZT-715-02S	Sizing Template, Size 02
SZT-715-03S	Sizing Template, Size 03
SZT-715-04S	Sizing Template, Size 04
SZT-715-05S	Sizing Template, Size 05
TRL-715-01T	Lunate Trial, Size 01
TRL-715-02T	Lunate Trial, Size 02
TRL-715-03T	Lunate Trial, Size 03
TRL-715-04T	Lunate Trial, Size 04
TRL-715-05T	Lunate Trial, Size 05

First Choice™ System (DRUJ)

The First Choice™ system (DRUJ) contains a Total Ulna Head replacement as well as a Partial Ulna head replacement in one system.

The instrumentation set provides the same instruments for both implant options, but there are additional instruments present that are needed for the Partial Ulna head replacement.

Within one set, 2 trays are available:

- Bottom tray : general instruments (reamers, awl etc...).
- Top tray: Trails partial Ulna head replacement and the dedicated cutting guide that is needed for the partial Ulna Head replacement surgical technique.

The same instrument set is needed to be used with both replacements.

There is a range sizes of Partial Distal Ulnar Radial head replacement available. The partial ulna head replacement is a monoblock implant

- 3 Stem sizes: 4.5 mm, 4.5 mm and 5.5 mm
- All these stem size are available with the following head sizes: -14.5, 16.0- 17.5 – 19.0

The Ascension First Choice™ Partial Ulnar Head Replacement

The Ascension® First Choice® is intended for replacement of the distal radioulnar joint.

Replacement of the distal ulnar head for rheumatoid, degenerative, or post-traumatic arthritis presenting the following:

- Pain and weakness of the wrist joint not improved by conservative treatment
- Instability of the ulnar head with x-ray evidence of dorsal subluxation and erosive changes
- Failed ulnar head resection



First Choice™ DRUJ Modular Ulnar Head



DRUJ implant

Description

14.5 mm head, 4.5 mm standard stem

16.0 mm head, 4.5 mm standard stem

17.5 mm head, 4.5 mm standard stem

19.0 mm head, 4.5 mm standard stem

14.5 mm head, 5.5 mm standard stem

16.0 mm head, 5.5 mm standard stem

17.5 mm head, 5.5 mm standard stem

19.0 mm head, 5.5 mm standard stem

14.5 mm head, 6.5 mm standard stem


16.0 mm head, 6.5 mm standard stem

17.5 mm head, 6.5 mm standard stem

19.0 mm head, 6.5 mm standard stem

First Choice DRUJ™ Instrument system

References	Description
DRUJ-610-1445	14.5 mm head, 4.5 mm standard stem
DRUJ-610-1645	16.0 mm head, 4.5 mm standard stem
DRUJ-610-1745	17.5 mm head, 4.5 mm standard stem
DRUJ-610-1945	19 mm head, 4.5 mm standard stem
DRUJ-610-1455	14.5 mm head, 5.5 mm standard stem
DRUJ-610-1655	16.0 mm head, 5.5 mm standard stem
DRUJ-610-1755	17.5 mm head, 5.5 mm standard stem
DRUJ-610-1955	19 mm head, 5.5 mm standard stem
DRUJ-610-1465	14.5 mm head, 6.5 mm standard stem
DRUJ-610-1665	16.0 mm head, 6.5 mm standard stem
DRUJ-610-1765	17.5 mm head, 6.5 mm standard stem
DRUJ-610-1965	19.0 mm head, 6.5 mm standard stem



Note: To perform a Partial Ulnar Head procedure, both the MUH & Partial instrument sets are needed.

Partial Ulnar Head Implants

References	Description
INS-610-00	Instrument Set
OSG-610-00	DRUJ Osteotomy Guide
TRL-610-H145-S45	DRUJ Trial Size 14.5-4.5
TRL-610-H160-S45	DRUJ Trial Size 16.0-4.5
TRL-610-H175-S45	DRUJ Trial Size 17.5-4.5
TRL-610-H190-S45	DRUJ Trial Size 19.0-4.5
TRL-610-H145-S55	DRUJ Trial Size 14.5-5.5
TRL-610-H160-S55	DRUJ Trial Size 16.0-5.5
TRL-610-H175-S55	DRUJ Trial Size 17.5-5.5
TRL-610-H190-S55	DRUJ Trial Size 19.0-5.5
TRL-610-H145-S65	DRUJ Trial Size 14.5-6.5
TRL-610-H160-S65	DRUJ Trial Size 16.0-6.5
TRL-610-H175-S65	DRUJ Trial Size 17.5-6.5
TRL-610-H190-S65	DRUJ Trial Size 19.0-6.5

Modular Ulnar Head Instruments

References	Description
INS-600-00	Instrument Set
AWL-100-01	Starter Awl
TRL-600-H160	MUH Head Trial Size 16.0
TRL-600-H175	MUH Head Trial Size 17.5
TRL-600-H190	MUH Head Trial Size 19.0
TRL-600-S45	MUH Stem Trial Size 4.5
TRL-600-S55	MUH Stem Trial Size 5.5
TRL-600-S65	MUH Stem Trial Size 6.5
TRL-600-MED	MUH Trial Collar Medium
TRL-600-LNG	MUH Trial Collar Long
OSG-600-00	MUH Resection Guide
BRH-600-45	MUH Reamer Size 4.5
BRH-600-55	MUH Reamer Size 5.5
BRH-600-65	MUH Reamer Size 6.5
IMP-600-00	MUH Stem Impactor
IMP-600-01	MUH Assembly Pad
OSG-600-01	MUH Osteotomy Guide
IMP-300-00	Implant Impactor
EXT-200-00	Trial Extractor

Total Modular Ulnar Head Implants

References	Description
MUH600H160WW	MUH Head Sz. 16.0
MUH600H175WW	MUH Head Sz. 17.5
MUH600H190WW	MUH Head Sz. 19.0
MUH600S45LWW	MUH Stem Sz. 4.5 L
MUH600S45MWW	MUH Stem Sz. 4.5 M
MUH600S45SWW	MUH Stem Sz. 4.5 S
MUH600S55LWW	MUH Stem Sz. 5.5 L
MUH600S55MWW	MUH Stem Sz. 5.5 M
MUH600S55SWW	MUH Stem Sz. 5.5 S
MUH600S65LWW	MUH Stem Sz. 6.5 L
MUH600S65MWW	MUH Stem Sz. 6.5 M
MUH600S65SWW	MUH Stem Sz. 6.5 S

Modular Ulnar Head Instruments

References	Description
INS-600-00	Instrument Set
AWL-100-01	Starter Awl
TRL-600-H160	MUH Head Trial Size 16.0
TRL-600-H175	MUH Head Trial Size 17.5
TRL-600-H190	MUH Head Trial Size 19.0
TRL-600-S45	MUH Stem Trial Size 4.5
TRL-600-S55	MUH Stem Trial Size 5.5
TRL-600-S65	MUH Stem Trial Size 6.5
TRL-600-MED	MUH Trial Collar Medium
TRL-600-LNG	MUH Trial Collar Long
OSG-600-00	MUH Resection Guide
BRH-600-45	MUH Reamer Size 4.5
BRH-600-55	MUH Reamer Size 5.5
BRH-600-65	MUH Reamer Size 6.5
IMP-600-00	MUH Stem Impactor
IMP-600-01	MUH Assembly Pad
OSG-600-01	Osteotomy Guide
IMP-300-00	Implant Impactor

Wifix® Plate – Surfix® Volar Plate for Distal Radius

This plate is indicated for:

Traumatology:

- Articular or distal complex fractures of the radius.
- Anterior or posterior comminuted fractures with bone defect.

Reconstructive surgery:

- Shortening osteotomy for Kienbock's disease.
- Correction of distal radius malunion particularly anterior additional bone graft.
- Lunar radio-scapho arthrodesis.
- Radio-carpal arthrodesis.
 - > 3 sizes.
 - > Symmetrical and bendable plate.
 - > Anatomical.
 - > Low profil design.
 - > Oblong hole with Surfix® fixation.
 - > Materials:
 - » Plate: Stainless steel 316L
Standard ISO 5832-1 & ASTM F138&139
 - » Screws and lock-screws: Titanium TA6V
Standard ISO 5832-3 & ASTM F136

Contra- indication: Complex distal radial fractures.

References	Description
PDRA1	Volar plate for distal radius - Size 1 - L. 44 mm
PDRA2	Volar plate for distal radius - Size 2 - L. 70 mm
PDRA3	Volar plate for distal radius - Size 3 - L. 96 mm

WIFIX® plate – Associated Surfix® screw

Surfix® screw + lock-screw diam 3.5 mm

References	Description
V35 012	L. 12 mm
V35 014	L. 14 mm
V35 016	L. 16 mm
V35 018	L. 18 mm
V35 020	L. 20 mm
V35 022	L. 22 mm
V35 024	L. 24 mm
V35 026	L. 26 mm



WIFIX® – Associated Instruments

References	Description
MA3501	Plates holder
MA5303	Plates bender
MA3504	Hexagonal screwdriver 2 mm
MA3505	Length gauge
MA3506	Drill guide diam. 2.5 mm
MA3508	Drill AO diam. 2.5 mm
FPDRA1	Trial PDRA1
FPDRA2	Trial PDRA2
FPDRA3	Trial PDRA3





Elbow Pathologies Solutions

Modular Radial Head (MRH) and Carbon Modular Radial Head (CMRH) are intended for:

- Replacement of the radial head for degenerative or post-traumatic disabilities presenting pain, crepitation, and decreased motion at the radio-humeral and/or proximal radio-ulnar joint with joint destruction or subluxation visible on x-ray resistance to conservative treatment.
- Primary replacement after fracture of the radial head.
- Symptomatic sequelae after radial head resection.
- Revision following failed radial head arthroplasty.

RadFx system allows surgeons to have all the tools necessary to repair or replace a radial head fracture in one instrument set.



Pyrcarbon Radial Head Implants (CMRH)

Modular Radial Head (MRH) and Carbon Modular Radial Head (CMRH) are intended for:

- Replacement of the radial head for degenerative or post-traumatic disabilities presenting pain, crepitation, and decreased motion at the radio-humeral and/or proximal radio-ulnar joint with:
 - > joint destruction or subluxation visible on x-ray
 - > resistance to conservative treatment
- Primary replacement after fracture of the radial head.
- Symptomatic sequelae after radial head resection.
- Revision following failed radial head arthroplasty.

RadFx™ system allows surgeons to have all the tools necessary to repair or replace a radial head fracture in one instrument set.

It offers patients a radial head replacement.

The CMRH is an anatomically designed, semi-constrained, single-use, modular (hemi) radial head replacement. It has a proximal articular Pyrocarbon head component and a distal intramedullary Titanium stem component.

Design Specific Features:

- 3 different Pyrocarbon head sizes (20, 22, 24).
- 4 different stem lengths (10, 20, 30, 40) and these 4 sizes are available in both standard or long collar size.
- 24 head/stem combinations possible.
- Anatomically curved, cement less, long stem option.
- Same Instrument set is used for both the CMRH and the MRH replacement.
- Colour coded broaches and trials.

CMRH implants

References	Description
CMRH-310-H20-WW	Carbon Radial Head size 20
CMRH-310-H22-WW	Carbon Radial Head size 22
CMRH-310-H24-WW	Carbon Radial Head size 24
CMRH-310-So1S-WW	CMRH Stem Standard Size 01
CMRH-310-So2S-WW	CMRH Stem Standard Size 02
CMRH-310-So3S-WW	CMRH Stem Standard Size 03
CMRH-310-So4S-WW	CMRH Stem Standard Size 04
CMRH-310-So1L-WW	CMRH Stem Long Size 01
CMRH-310-So2L-WW	CMRH Stem Long Size 02
CMRH-310-So3L-WW	CMRH Stem Long Size 03
CMRH-310-So4L-WW	CMRH Stem Long Size 04



Metal Radial Head

The Metal Radial Head (MRH) is an anatomically designed, semi-constrained, single-use, modular (hemi) radial head replacement. It has a proximal articular CoCr head component and a distal intramedullary titanium stem component.

- Sloping lateral surface is radial-ulnar notch and accommodates wrapping of annular ligament.
- Same instrument system is used for both the Pyrocarbon and the Metal proximal Radial Head replacement.

MRH Implants

References	Description
MRH-350-20L	Radial Head Size 20 Long
MRH-350-20S	Radial Head Size 20 Standard
MRH-350-22L	Radial Head Size 22 Long
MRH-350-22S	Radial Head Size 22 Standard
MRH-350-24L	Radial Head Size 24 Long
MRH-350-24S	Radial Head Size 24 Standard
MRH-350-01	Stem Size 01
MRH-350-02	Stem Size 02
MRH-350-03	Stem Size 03
MRH-350-04	Stem Size 04



Optional Instrument Set Components

References	Description
MRH-RMR-01	MRH Reamer Sz. 01
MRH-RMR-02	MRH Reamer Sz. 02
MRH-RMR-03	MRH Reamer Sz. 03
MRH-PLN-00	MRH Planar

Instrument Set, Metal Radial Head

References	Description
CSA-350-02	Instrument Base
CSA-000-01	Instrument Lid
AWL-100-01	Awl, Starter
IMP-300-00	Implant Impactor
IMP-350-00	Impactor, Stem, MRH
IMP-350-01	Back Table Assembly Pad, MRH
BRH-300-20	Radial Head Broach Size 20
BRH-300-22	Radial Head Broach Size 22
BRH-300-24	Radial Head Broach Size 24
BRH-350-04	Broach, MRH, Sz. 04
EXT-350-00	Stem Holder, MRH
OSG-350-01	Osteotomy Cut Guide, MRH
TRL-351-20L	MRH Trial, Ti, Sz. 20L
TRL-351-20S	MRH Trial, Ti, Sz. 20S
TRL-351-22L	MRH Trial, Ti, Sz. 22L
TRL-351-22S	MRH Trial, Ti, Sz. 22S
TRL-351-24L	MRH Trial, Ti, Sz. 24L
TRL-351-24S	MRH Trial, Ti, Sz. 24S
INS-350-00	Complete Instrument Set

RADFx™ Radial Head Fracture System

One set containing 2 trays :

- Bottom tray- CHRH or MRH set.
- Top tray- Compression screw set and proximal radial head fracture plates.
 - > Compression Screws.
 - » 1.5 mini compression
 - » 2.0 mm high compression
 - > Large and Small plates.
- Low profile titanium 1.2mm locking plates can accept locking and non-locking screws.

References	Description
28.25.010	Plate locking screw, 10 mm
28.25.012	Plate locking screw, 12 mm
28.25.014	Plate locking screw, 14 mm
28.25.016	Plate locking screw, 16 mm
28.25.018	Plate locking screw, 18 mm
28.25.020	Plate locking screw, 20 mm
28.25.022	Plate locking screw, 22 mm
28.25.024	Plate locking screw, 24 mm
28.25.026	Plate locking screw, 26 mm
28.25.110	Plate lag screw, 10 mm
28.25.112	Plate lag screw, 12 mm
28.25.114	Plate lag screw, 14 mm
28.25.116	Plate lag screw, 16 mm
28.25.118	Plate lag screw, 18 mm
28.25.120	Plate lag screw, 20 mm
28.25.122	Plate lag screw, 22 mm
28.25.124	Plate lag screw, 24 mm
28.25.126	Plate lag screw, 26 mm
501 009 240	RH Plate, Short
501 009 241	RH Plate, Long


519 010 150	Hi comp Screw, cannulated, 10 mm
519 012 150	Hi comp Screw, cannulated, 12 mm
519 014 150	Hi comp Screw, cannulated, 14 mm
519 016 150	Hi comp Screw, cannulated, 16 mm
519 018 150	Hi comp Screw, cannulated, 18 mm
519 020 150	Hi comp Screw, cannulated, 20 mm
519 022 150	Hi comp Screw, cannulated, 22 mm
519 024 150	Hi comp Screw, cannulated, 24 mm
519 026 150	Hi comp Screw, cannulated, 26 mm
520 028 150	Hi comp Screw, cannulated, 28 mm
521 030 150	Hi comp Screw, cannulated, 30 mm
515 010 070	Mini comp Screw, cannulated, 10 mm
515 012 070	Mini comp Screw, cannulated, 12 mm
516 014 070	Mini comp Screw, cannulated, 14 mm
515 016 070	Mini comp Screw, cannulated, 16 mm
515 018 070	Mini comp Screw, cannulated, 18 mm
515 020 070	Mini comp Screw, cannulated, 20 mm
515 022 070	Mini comp Screw, cannulated, 22 mm
515 024 070	Mini comp Screw, cannulated, 24 mm
515 026 070	Mini comp Screw, cannulated, 26 mm
516 028 070	Mini comp Screw, cannulated, 28 mm
517 030 070	Mini comp Screw, cannulated, 30 mm

Disposables for RadFX™

References	Description
503 006320	Ascension Micro Drill Cannula Depth Gauge
503 006343	Ascension Tissue protector
502 015206	Drill bit for plate screws
605.090.105	K wire (for cannulated screws)
605 120 100	K wire (To pin bone)

Instrument Set, RadFX™

References	Description
CSA-850-00	Lid
CSA-850-01	Base
CSA-850-02	Top insert
CSA-850-03	Tag
CSA-860-01	Compression Screw Lid
CSA-860-02	Base Compr. Screw set
CSA-860-03	Screw tray w/lid
503 004197	Tweezers
503 004267	TX 6 driver, cannulated
503 004270	TX 8 driver, cannulated
503 004331	Cannula mallet plug (tap)
INS-860-00	Compression Screw Set
CSA-870-01	Plate Screw lid
CSA-870-02	Base Plate Set
CSA-870-03	Plate and screw tray w/lid
PBI-800-00	Plate Bending iron
503 004157	Drill guide handle for plate screws
503 004170	Screw on drill guide for plate screws
503 004197	Pickups
503 004262	Depth gage for plate screws
503 004268	TX 8 driver, non-cannulated



Katalyst™ Bipolar Radial Head System

KATALYST™ Bipolar Radial Head System is indicated in:

- Comminuted radial fracture with or without associated elbow instability.
- Comminuted radial fracture with associated rupture of the interosseous membrane; the Essex-Lopresti lesion.
- Post-traumatic arthritis involving the radiocapitellar or proximal radioulnar joint.
- Primary osteoarthritis involving the radiocapitellar joint or proximal radioulnar joint.
- Rheumatoid arthritis involving the radiocapitellar joint or proximal radioulnar joint.
- Revision of failed radial head prosthesis.

Heads

References	Description
221418	18 mm Head & Poly Assembly
221421	21 mm Head & Poly Assembly
221424	24 mm Head & Poly Assembly

Stems

References	Description
221665	6.5 mm Stem Assembly, 20 mm length
221675	7.5 mm Stem Assembly, 20 mm Length
221800	Stem Replacement Set Screw



Katalyst™ Bipolar Radial Head System – Instruments

Instruments

References	Description
224165	Stem Trial 6.5 mm
224175	Stem Trial 7.5 mm
224318	Head Trial 18 mm
224321	Head Trial 21 mm
224324	Head Trial 24 mm
225365	6.5 mm Broach
225375	7.5 mm Broach
225400	Spacer
225405	First Spacer
225600	Cutting Template
225900	Head Sizing Template
226000	Starter Awl
600717	2.5 mm Hex Driver
225205	Pliers

Sterilization container

References	Description
221001	Sterilization tray



EndoRelease™ Endoscopic Cubital Tunnel Release System

Indications for use:

- The EndoRelease™ System is intended for use in endoscopic surgical treatment of cubital tunnel syndrome by releasing the fascia around the ulnar nerve.

Disposables

References	Description
096000	Blade
096005	Set of 5 blades



Endorelease™ Endoscopic Cubital Tunnel Release System – Instruments

Instruments

References	Description
095070	Small Spatula cbtl tunnel
095080	Large Spatula cbtl tunnel
095090	Small Cannula
095160	Small obturator
095250	Large obturator



Titan™ Humeral Resurfacing

For complete product information, please see the Instructions for Use accompanying each product.



TITAN HRA Implants

References	Description
HRA-900-40/16-WW	Titan HRA Size 40/16
HRA-900-44/15-WW	Titan HRA Size 44/15
HRA-900-44/18-WW	Titan HRA Size 44/18
HRA-900-48/17-WW	Titan HRA Size 48/17
HRA-900-48/20-WW	Titan HRA Size 48/20
HRA-900-52/18-WW	Titan HRA Size 52/18
HRA-900-52/21-WW	Titan HRA Size 52/21
HRA-900-56/21-WW	Titan HRA Size 56/21

Disposables, TITAN HRA Resurfacing

References	Description
605 303230	Partial Threaded Kirschner
605 300230	Non-Threaded Kirschner

Instrument Set, TITAN HRA Resurfacing

References	Description
CSA-900-02	HRA Case
CSA-900-01	HRA Lid
CSA-900-03	HRA Top Tray
CSA-900-04	HRA Bottom Tray
ADP-900-01	HRA Adapter Size1
ADP-900-02	HRA Adapter Size2
ADP-900-03	HRA Adapter Size 3
BRH-900-00	HRA Starter Broach
BRH-900-01	HRA Broach Size 1
BRH-900-02	HRA Broach Size 2
BRH-900-03	HRA Broach Size 3
CTD-900-01	HRA Center Drill Size 1
CTD-900-02	HRA Center Drill Size 2
CTD-900-03	HRA Center Drill Size 3
IMP-900-01	HRA Impactor
RET-900-01	HRA Retractor
SHP-900-40/16	HRA Shaper 40/16
SHP-900-44/15	HRA Shaper 44/15
SHP-900-44/18	HRA Shaper 44/18
SHP-900-48/17	HRA Shaper 48/17
SHP-900-48/20	HRA Shaper 48/20
SHP-900-52/18	HRA Shaper 52/18
SHP-900-52/21	HRA Shaper 52/21
SHP-900-56/21	HRA Shaper 56/21
SZR-900-40-16	Sizer/Centralizer 40-16
SZR-900-44-15	Sizer/Centralizer 44-15
SZR-900-44-18	Sizer/Centralizer 44-18
SZR-900-48-17	Sizer/Centralizer 48-17
SZR-900-48-20	Sizer/Centralizer 48-20
SZR-900-52-18	Sizer/Centralizer 52-18
SZR-900-52-21	Sizer/Centralizer 52-21
SZR-900-56-21	Sizer/Centralizer 56-21
TRL-900-40/16	HRA Trail 40/16
TRL-900-44/15	HRA Trial 44/15
TRL-900-44/18	HRA Trial 44/18
TRL-900-48/17	HRA Trial 48/17
TRL-900-48/20	HRA Trial 48/20
TRL-900-52/18	HRA Trial 52/18
TRL-900-52/21	HRA Trial 52/21
TRL-900-56/21	HRA Trial 56/21
INS-900-00	Total
INS-900-00	Complete Instrument Set

Titan™ Modular Total Shoulder System

Total Shoulder Arthroplasty or Hemiarthroplasty is indicated for:

- Severely painful and/or disabled joint resulting from osteoarthritis, traumatic arthritis or rheumatoid arthritis.
- Fracture-dislocations of the proximal humerus where the articular surface is severely comminuted, separated from its blood supply or where the surgeon's experience indicates that alternative methods of treatments are unsatisfactory.
- Other difficult clinical problems where shoulder arthrodesis or resection arthroplasty are not acceptable (e.g. - revision of a failed primary component).

Shoulder Hemiarthroplasty is also indicated for:

- Ununited humeral head fractures.
- Avascular necrosis of the humeral head.
- Rotator cuff arthroplasty.
- Deformity and/or limited motion.

The humeral component is intended for cemented or uncemented use. The glenoid component is intended for cemented use only.



Glenoid tray: Base

References	Description
238500	Straight Drill Shaft
238000	Angled Drill Shaft
IMP-0920-064-001	Glenoid impactor
PUN-0920-062-001	Glenoid Keel Punch, XS/S
PUN-0920-062-002	Glenoid Keel Punch, M/L
PUN-0920-072-001	Glenoid, Peg Punch, XS/S
PUN-0920-072-002	Glenoid, Peg Punch, M/L

Glenoid tray: Insert

References	Description
SZR-0920-052-000	Glenoid Sizer, Extra Small
SZR-0920-052-001	Glenoid Sizer, Small
SZR-0920-052-002	Glenoid Sizer, Medium
SZR-0920-052-003	Glenoid Sizer, Large
RMR-0920-057-000	Glenoid Reamer, Extra Small
RMR-0920-057-001	Glenoid Reamer, Small
RMR-0920-057-002	Glenoid Reamer, Medium
RMR-0920-057-003	Glenoid Reamer, Large
TRL-0920-030-00P	Glenoid, Peg Trial, Extra Small
TRL-0920-030-01P	Glenoid, Peg Trial, Small
TRL-0920-030-02P	Glenoid, Peg Trial, Medium
TRL-0920-030-03P	Glenoid, Peg Trial, Large
TRL-0920-030-00K	Glenoid, Keel Trial, Extra Small
TRL-0920-030-01K	Glenoid, Keel Trial, Small
TRL-0920-030-02K	Glenoid, Keel Trial, Medium
TRL-0920-030-03K	Glenoid, Keel Trial, Large
GDW-0920-061-501	2.0 mm x 150 mm Guide Pin
DRG-0920-049-001	Center Glenoid Drill Guide, XS/S
DRG-0920-049-002	Center Glenoid Drill Guide, M/L

HDR-0920-063-001

Glenoid Holder

HDL-0920-075-001	Anti-rotation Peg Holder
PEG-0920-074-001	Anti-rotation Peg
DRL-0920-056-001	Solid Center Starter Drill
DRL-0920-056-002	Solid Peg Drill
BSH-0920-048-002	Center Starter Drill Bushing
DRL-0920-047-001	Cannulated Center Starter Drill
DRL-0920-047-002	Cannulated Peg Drill
BSH-0920-048-001	Pin Bushing

Humeral tray: Base

References	Description
HDL-0920-043-001	Stem Trial Handle
TRL-0920-020-08STD	Body Trial, 8 Standard
TRL-0920-020-10SML	Body Trial, 10 Small
TRL-0920-020-10STD	Body Trial, 10 Standard
TRL-0920-020-10LRG	Body Trial, 10 Large
TRL-0920-020-12STD	Body Trial, 12 Standard
TRL-0920-020-14SML	Body Trial, 14 Small
TRL-0920-020-14STD	Body Trial, 14 Standard
TRL-0920-020-14LRG	Body Trial, 14 Large
TRL-0950-020-08SML	Fracture Body Trial, 8 Small
TRL-0950-020-08STD	Fracture Body Trial, 8 Standard
TRL-0950-020-08LRG	Fracture Body Trial, 8 Large
BSW-0920-021-01	Body Screws
TRL-0920-025-06	Humeral Stem Trial, 6 mm
TRL-0920-025-07	Humeral Stem Trial, 7 mm
TRL-0920-025-08	Humeral Stem Trial, 8 mm
TRL-0920-025-09	Humeral Stem Trial, 9 mm
TRL-0920-025-10	Humeral Stem Trial, 10 mm
TRL-0920-025-11	Humeral Stem Trial, 11 mm

TRL-0920-025-12	Humeral Stem Trial, 12 mm
TRL-0920-025-13	Humeral Stem Trial, 13 mm
TRL-0920-025-14	Humeral Stem Trial, 14 mm
TRL-0920-025-15	Humeral Stem Trial, 15 mm
TRL-0920-025-16	Humeral Stem Trial, 16 mm
INS-0920-046-001	Humeral Trial Inserter/Extractor
INS-0950-046-001	Fracture Trial Inserter/Extractor
RMR-0920-050-001	Calcar Planer, Small
RMR-0920-050-002	Calcar Planer, Large
ADT-0920-065-001	Taper Adapter
SCR-0920-066-001	Taper Adapter Shaft
IMP-0920-079-501	Head Impactor
SCR-0920-060-002	3/32 Hexdriver
G207153	Torque Limiter
SET188-A001	Head Extractor
STD-0920-071-001	Implant Stand, Base Assembly
SEP-0920-068-001	Body Separator
INS-0920-045-001	Humeral Implant Inserter/Extractor
INS-0950-045-001	Fracture Implant Inserter/Extractor
INS-0920-071-SML	Stand Insert, Small
INS-0920-071-STD	Stand Insert, Standard
INS-0920-071-LRG	Stand Insert, Large
SIS-0920-054-001	Stem Impaction Stand
1 MP-0920-055-001	Stem Impactor
G107991_B	Slotted Mallet

Humeral tray: Insert

References	Description
SET189-A001	Head Cutting Guide
AWL-0920-042-001	Starter Awl
G108404_B	T-Handle w/Large AO Connect
SET189-D007	Version Rod
3020	Pin Puller
PIN-0920-051-001	Fixation Pin
TMP-0920-040-001R	Head Cutting Template, Right
TMP-0920-040-001L	Head Cutting Template, Left
ROD-0920-040-002	Head Cutting Template, Handle
GAU-0920-058-001	Head Cutting Depth, Gauge
HSG-0920-041-001	Head Sizing Gauge, 38-46
HSG-0920-041-002	Head Sizing Gauge, 48-52
CVR-0920-076S	Trial Osteotomy Cover, Small
CVR-0920-076L	Trial Osteotomy Cover, Large
CVR-0920-077S	Spiked Osteotomy Cover, Small
CVR-0920-077L	Spiked Osteotomy Cover, Large
G107992_B	Driver Handle w/Small AO Connect
SCR-0920-060-001	1/8 Hexdriver
HSP-0920-070-001E	Head Sizing Plate, 38 mm
HSP-0920-070-002E	Head Sizing Plate, 40 mm

HSP-0920-070-003	Head Sizing Plate, 42 mm
HSP-0920-070-004	Head Sizing Plate, 44 mm
HSP-0920-070-005	Head Sizing Plate, 46 mm
HSP-0920-070-006	Head Sizing Plate, 48 mm
HSP-0920-070-007	Head Sizing Plate, 50 mm
HSP-0920-070-008	Head Sizing Plate, 52 mm
TRL-0920-010-3814E	Head Trial, Eccentric, 38x14 mm
TRL-0920-010-4025E	Head Trial, Eccentric, 40x15 mm
TRL-0920-010-4216E	Head Trial, Eccentric, 42x16 mm
TRL-0920-010-4416E	Head Trial, Eccentric, 44x16 mm
TRL-0920-010-4419E	Head Trial, Eccentric, 44x19 mm
TRL-0920-010-4614E	Head Trial, Eccentric, 46x14 mm
TRL-0920-010-4617E	Head Trial, Eccentric, 46x17 mm
TRL-0920-010-4620E	Head Trial, Eccentric, 46x20 mm
TRL-0920-010-4815E	Head Trial, Eccentric, 48x15 mm
TRL-0920-010-4818E	Head Trial, Eccentric, 48x18 mm
TRL-0920-010-4821E	Head Trial, Eccentric, 48x21 mm
TRL-0920-010-5019E	Head Trial, Eccentric, 50x19 mm
TRL-0920-010-5022E	Head Trial, Eccentric, 50x22 mm
TRL-0920-010-5220E	Head Trial, Eccentric, 52x20 mm
TRL-0920-010-4216C	Head Trial, Concentric, 42x16 mm
TRL-0920-010-4416C	Head Trial, Concentric, 44x16 mm
TRL-0920-010-4419C	Head Trial, Concentric, 44x19 mm
TRL-0920-010-4614C	Head Trial, Concentric, 46x14 mm
TRL-0920-010-4617C	Head Trial, Concentric, 46x17 mm
TRL-0920-010-4620C	Head Trial, Concentric, 46x20 mm
TRL-0920-010-4815C	Head Trial, Concentric, 48x15 mm
TRL-0920-010-4818C	Head Trial, Concentric, 48x18 mm
TRL-0920-010-4821C	Head Trial, Concentric, 48x21 mm
TRL-0920-010-5019C	Head Trial, Concentric, 50x19 mm
TRL-0920-010-5022C	Head Trial, Concentric, 50x22 mm
TRL-0920-010-5220C	Head Trial, Concentric, 52x20 mm

Retractor set

References	Description
RET-920-11S	Darrach Retractor, Small
RET-920-11L	Darrach Retractor, Large
RET-920-20	Deltoid Retractor, Large
RET-920-31	Hohmann Retractor
RET-920-41S	Spiked Glenoid Retractor, Small
RET-920-41L	Spiked Glenoid Retractor, Large
RET-920-50S	Fukuda Retractor, Small
RET-920-50L	Fukuda Retractor, Large
RET-920-60	Kolbel Retractor Frame
RET-920-61S	Kolbel Retractor Blade, Small
RET-920-61L	Kolbel Retractor Blade, Large
RET-920-61XL	Kolbel Retractor Blade, Extra Large

CoCr Humeral Head Implants

References	Description
MHH-0920-010-4216C	Humeral Head, Cen, Size 42x16 mm
MHH-0920-010-4416C	Humeral Head, Cen, Size 44x16 mm
MHH-0920-010-4419C	Humeral Head, Cen, Size 44x19 mm
MHH-0920-010-4614C	Humeral Head, Cen, Size 46x14 mm
MHH-0920-010-4617C	Humeral Head, Cen, Size 46x17 mm
MHH-0920-010-4620C	Humeral Head, Cen, Size 46x20 mm
MHH-0920-010-4815C	Humeral Head, Cen, Size 48x15 mm
MHH-0920-010-4818C	Humeral Head, Cen, Size 48x18 mm
MHH-0920-010-4821C	Humeral Head, Cen, Size 48x21 mm
MHH-0920-010-5019C	Humeral Head, Cen, Size 50x19 mm
MHH-0920-010-5022C	Humeral Head, Cen, Size 50x22 mm
MHH-0920-010-5220C	Humeral Head, Cen, Size 52x20 mm
MHH-0920-010-3814E	Humeral Head, Ecc, Size 38x14 mm
MHH-0920-010-4015E	Humeral Head, Ecc, Size 40x15 mm
MHH-0920-010-4216E	Humeral Head, Ecc, Size 42x16 mm
MHH-0920-010-4416E	Humeral Head, Ecc, Size 44x16 mm
MHH-0920-010-4419E	Humeral Head, Ecc, Size 44x19 mm
MHH-0920-010-4614E	Humeral Head, Ecc, Size 46x14 mm
MHH-0920-010-4617E	Humeral Head, Ecc, Size 46x17 mm
MHH-0920-010-4620E	Humeral Head, Ecc, Size 46x20 mm
MHH-0920-010-4815E	Humeral Head, Ecc, Size 48x15 mm
MHH-0920-010-4818E	Humeral Head, Ecc, Size 48x18 mm
MHH-0920-010-4821E	Humeral Head, Ecc, Size 48x21 mm
MHH-0920-010-5019E	Humeral Head, Ecc, Size 50x19 mm
MHH-0920-010-5022E	Humeral Head, Ecc, Size 50x22 mm
MHH-0920-010-5220E	Humeral Head, Ecc, Size 52x20 mm

Primary Body Implants

References	Description
BOD-0920-020-08STD	Humeral Body, Size 08 Standard
BOD-0920-020-10SML	Humeral Body, Size 10 Small
BOD-0920-020-10STD	Humeral Body, Size 10 Standard
BOD-0920-020-10LRG	Humeral Body, Size 10 Large
BOD-0920-020-12STD	Humeral Body, Size 12 Standard
BOD-0920-020-14SML	Humeral Body, Size 14 Small
BOD-0920-020-14STD	Humeral Body, Size 14 Standard
BOD-0920-020-14LRG	Humeral Body, Size 14 Large
BOD-0950-020-08SML	Fracture Body, Size 08 Small
BOD-0950-020-08STD	Fracture Body, Size 08 Standard
BOD-0950-020-08LRG	Fracture Body, Size 08 Large
BSW-0920-021-01	Humeral Body Screw

Press Fit Stem Implants

References	Description
STEM-0920-025-06	Humeral Stem, 6 mm
STEM-0920-025-07	Humeral Stem, 7 mm
STEM-0920-025-08	Humeral Stem, 8 mm
STEM-0920-025-09	Humeral Stem, 9 mm
STEM-0920-025-10	Humeral Stem, 10 mm
STEM-0920-025-11	Humeral Stem, 11 mm
STEM-0920-025-12	Humeral Stem, 12 mm
STEM-0920-025-13	Humeral Stem, 13 mm
STEM-0920-025-14	Humeral Stem, 14 mm
STEM-0920-025-15	Humeral Stem, 15 mm
STEM-0920-025-16	Humeral Stem, 16 mm

Cemented Stem Implants

References	Description
STEM-0950-025-06	Cemented Stem, 6 mm
STEM-0950-025-08	Cemented Stem, 8 mm
STEM-0950-025-10	Cemented Stem, 10 mm
STEM-0950-025-12	Cemented Stem, 12 mm
STEM-0950-025-14	Cemented Stem, 14 mm

Keeled Glenoid Implants

References	Description
GLN-0920-030-00K	Glenoid, Keel, Extra Small
GLN-0920-030-01K	Glenoid, Keel, Small
GLN-0920-030-02K	Glenoid, Keel, Medium
GLN-0920-030-03K	Glenoid, Keel, Large

Pegged Glenoid Implants

References	Description
GLN-0920-030-00P	Glenoid, Peg, Extra Small
GLN-0920-030-01P	Glenoid, Peg, Small
GLN-0920-030-02P	Glenoid, Peg, Medium
GLN-0920-030-03P	Glenoid, Peg, Large

PyroTitan™ Humeral Resurfacing Arthroplasty

The PyroTITAN™ HRA System is intended for resurfacing of the humeral head due to non-inflammatory or inflammatory arthritis (i.e. rheumatoid arthritis, osteoarthritis and some cases of osteonecrosis), mild or moderate humeral head deformity and / or limited motion, post-traumatic arthritis, focal and large (Hill-Sachs) osteochondral defects, and patients with an intact or raparable rotator cuff.



PyroTITAN HRA Implants

References	Description
CHRA-910-38/14	Pyro TITAN Implant, Size 38/14
CHRA-910-41/15	Pyro TITAN Implant, Size 41/15
CHRA-910-41/18	Pyro TITAN Implant, Size 41/18
CHRA-910-44/16	Pyro TITAN Implant, Size 44/16
CHRA-910-44/19	Pyro TITAN Implant, Size 44/19
CHRA-910-47/17	Pyro TITAN Implant, Size 47/17
CHRA-910-47/20	Pyro TITAN Implant, Size 47/20
CHRA-910-50/18	Pyro TITAN Implant, Size 50/18
CHRA-910-50/21	Pyro TITAN Implant, Size 50/21
CHRA-910-53/19	Pyro TITAN Implant, Size 53/19
CHRA-910-53/22	Pyro TITAN Implant, Size 53/22
CHRA-910-56/21	Pyro TITAN Implant, Size 56/21

Instrument Set, PyroTITAN HRA

References	Description
CSA-00-007	Lid
CSA-910-02	Tray Inset
CSA-910-03	Tray Base
OSG-910-01	Sizer, 38/14, 41/15, 41/18, 44/16
OSG-910-02	Sizer, 44/19, 47/17, 47/20, 50/18
OSG-910-03	Sizer, 50/21, 53/19, 53/22, 56/21
SZR-910-38	Centralizer, Size 38
SZR-910-41	Centralizer, Size 41
SZR-910-44	Centralizer, Size 44
SZR-910-47	Centralizer, Size 47
SZR-910-50	Centralizer, Size 50
SZR-910-53	Centralizer, Size 53
SZR-910-56	Centralizer, Size 56
SHP-910-38-14	Shaper, Size 38/14
SHP-910-41-15	Shaper, Size 41/15
SHP-910-41-18	Shaper, Size 41/18
SHP-910-44-16	Shaper, Size 44/16
SHP-910-44-19	Shaper, Size 44/19
SHP-910-47-17	Shaper, Size 47/17
SHP-910-47-20	Shaper, Size 47/20
SHP-910-50-18	Shaper, Size 50/18
SHP-910-50-21	Shaper, Size 50/21
SHP-910-53-19	Shaper, Size 53/19
SHP-910-53-22	Shaper, Size 53/22
SHP-910-56-21	Shaper, Size 56/21
BRH-910-38	Broach, Size 38
BRH-910-41	Broach, Size 41

BRH-910-44	Broach, Size 44
BRH-910-47	Broach, Size 47
BRH-910-50	Broach, Size 50
BRH-910-53	Broach, Size 53
BRH-910-56	Broach, Size 56
TRL-910-38-14	Trial, Size 38/14
TRL-910-41-15	Trial, Size 41/15
TRL-910-41-18	Trial, Size 41/18
TRL-910-44-16	Trial, Size 44/16
TRL-910-44-19	Trial, Size 44/19
TRL-910-47-17	Trial, Size 47/17
TRL-910-47-20	Trial, Size 47/20
TRL-910-50-18	Trial, Size 50/18
TRL-910-50-21	Trial, Size 50/21
TRL-910-53-19	Trial, Size 53/19
TRL-910-53-22	Trial, Size 53/22
TRL-910-56-21	Trial, Size 56/21
IMP-910-00	Impactor
MAL-910-00	Mallet
EXT-910-00	Extractor
INS-910-01	Instrument Set, PyroTITAN HRA

Instrument Set, Shoulder Retractors (Coated)

References	Description
CSA-925-01	Retractor Base
CSA-000-07	Lid
RET-920-10S	Coated Darrach Retractor, Small
RET-920-10L	Coated Darrach Retractor, Large
RET-920-20	Coated Deltoid Retractor, Large
RET-920-30	Coated Hohmann Retractor
RET-920-40S	Coated Glenoid Retractor, Small
RET-920-40L	Coated Glenoid Retractor, Large
RET-920-50S	Coated Fukuda Retractor, Small
RET-920-50L	Coated Fukuda Retractor, Large
RET-920-60	Coated Kolbel Retractor, Frame
RET-920-60S	Coated Kolbel Retractor, Small
RET-920-60L	Coated Kolbel Retractor, Large
RET-920-60XL	Coated Kolbel Retractor, Extra Large
INS-920-02	Instrument Set, Coated Shoulder Retractors

PSH plate – Surfix[®] Proximal humeral plate

These plates are indicated for the fixation of proximal humeral fractures.

- 3 sizes.
- Symmetrical and bendable plates.
- Low profile design.
- Stabilisation of the lesser and greater tubercle.
- Materials:
 - > Plate: Stainless steel 316L
Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V
Standard ISO 5832-3 & ASTM F136

Indications

- The indications for surgery of the whole range of SURFIX[®] implants are described in the attached document Indications of these directions for use and in the leaflet.

Contraindications

- Acute infections; Absence of musculo-cutaneous covering; Bone deterioration preventing a good fixation of the screws in the bone.


References	Description
PSH	Humeral proximal plate Size 1
PSH2	Humeral proximal plate Size 2
PSH3	Humeral proximal plate Size 3




PSH Plate – Associated Surfix[®] screws

Surfix[®] Cancellous bone screw
+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



Surfix® Technology System




Surfix® Concept

Drilling:

Start with the placement of the Metacarpal Surfix® 2.7 mm Screw. Use the 2.0 mm Black drill through the black drill bit (It is possible to attach the drill bit to the plate).

For placement of the Surfix® metacarpal 3.5 mm screw to the radius. It is possible to prepare the holes with the 2.7 mm (blue) drill through blue drill bit.




Positioning the screw:

Using the Star screwdriver (303408 = Black), insert the screw into the prepared hole until the plate is at the desired position relative to the bone.


Using the 2.0 mm Hex screwdriver (303410 = Blue), insert the screw into the prepared hole until the plate is at the desired position relative to the bone.

The screw should be fully seated in the plate. Clean the threaded hole before and after introducing the screw. Maintain co-axiality between the screw and the threaded hole.




Measuring:

Measure the necessary screw length using the depth gauge (303407).



Positioning the lock-screw:

Assemble the lock-screw to the screwdriver. The lock-screw should be inserted after each screw, and before preparation and insertion of the subsequent screw. This prevents potential damage to the thread.




Chamfering:

Counter sink the drill hole with the Star screwdriver (303408 = Black).


Counter sink the drill hole with the 2.0 mm Hex screwdriver (303410 = Blue).

Ensure that the threaded hole is not damaged when performing the chamfering.



Locking:

Fully seat the lock-screw with the screwdriver. The lock-screw should be flush with the top of the plate when it is fully inserted.





Integra®
Large Bones
Solutions
Knee and pediatric
Osteotomy plates

POF Plate

Closing Femoral Osteotomy Plate

The POF plate is indicated for distal femoral varus osteotomy and lateral gonarthrosis on genu-valgum.

- Fixation by 3 Surfix® screws into the femoral condyle for rotatory stability.
- Plate angulation adapted to the anatomy of the tibia.
- Shapeable Peg.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136

Reference	Description
POF	Distal Femoral Varus Osteotomy Plate



POF Plate

Associated Surfix® Screws

Surfix® Cortical bone screw
+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm
CC 055	Length 55 mm
CC 060	Length 60 mm
CC 065	Length 65 mm



Surfix® Cancellous bone screw
+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



PTO Plate

Medial Opening Tibial Osteotomy Plate

ORTHOPAEDIC: The osteotomy plates (PTOD and PTOG) are used for fixation of opening wedge high tibial osteotomy and osteoarthritis of the medial tibio-femoral compartment with genu-varum.

TRAUMATOLOGY: These plates may be used for internal fixations of fractures of the medial tibial plateau.

- Parallel Epiphysal Screws
- Strengthened bridge.
- Adapted angulation and curve.
- Shapeable Plate.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136



References	Description
PTOD	Right tibial opening plate
PTOG	Left tibial opening plate

PTO plate

Associated Surfix® Screws

Surfix® Cortical bone screw

+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm
CC 055	Length 55 mm
CC 060	Length 60 mm
CC 065	Length 65 mm



Surfix® Cancellous bone screw

+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



POM Plate

Less Invasive Medial Opening Tibial Plate than PTO Plate

Fixation of opening wedge high tibial osteotomy $6^{\circ} < 12^{\circ}$.

Osteoarthritis of the medial tibiofemoral compartment with genu-varum.

- Pre-oriented screw.
- Bendable plate.
- Adapted angulation and curve.
- Strengthened bridge.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136



References	Description
POMD	Less invasive plate for opening wedge high tibial osteotomy right side - size 1
POMG	Less invasive plate for opening wedge high tibial osteotomy left side - size 1

POM Plate

Associated Surfix® Screws

Surfix® Cortical bone screw
+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm
CC 055	Length 55 mm
CC 060	Length 60 mm
CC 065	Length 65 mm



Surfix® Cancellous bone screw
+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



PTC Plate

Lateral Closing Tibial Osteotomy Plate

ORTHOPAEDICS: The osteotomy plates (PTCD and PTCG) are used for fixation of closing wedge high tibial osteotomy or for osteoarthritis of the medial tibia femoral compartment with genu-varum.

TRAUMATOLOGY: These plates may be used for internal fixations of fractures of the lateral tibial plateau.

- Pre-shaped to the tibia anatomy.
- Shapeable pegs.
- Reinforced plates.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136

References	Description
PTCD	Proximal tibial osteosynthesis plate right side
PTCG	Proximal tibial osteosynthesis plate left side



PTC Plate

Associated Surfix[®] Screws

Surfix[®] Cortical bone screw

+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm
CC 055	Length 55 mm
CC 060	Length 60 mm
CC 065	Length 65 mm



Surfix[®] Cancellous bone screw

+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



PTF Plate

Lateral Closing Tibial Osteotomy Plate

ORTHOPAEDICS: The osteotomy plates (PTFD and PTFG) are used for fixation of closing wedge high tibial osteotomy or for osteoarthritis of the medial tibia femoral compartment with genu-varum.

TRAUMATOLOGY: These plates may be used for internal fixations of fractures of the lateral tibial plateau.

- Adapted angulation and curve to the anatomy of the tibia.
- Epiphyseal screws in the same plane.
- Adapted Angulation and Curve.
- Shapeable Plate.
- Strengthened bridge.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136

References	Description
PTFD	Right plate for closing wedge HTO
PTFG	Left plate for closing wedge HTO



PTF Plate

Associated Surfix® Screws

Surfix® Cortical bone screw
+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm
CC 055	Length 55 mm
CC 060	Length 60 mm
CC 065	Length 65 mm



Surfix® Cancellous bone screw
+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



Osteotomy Plates

ATAO® Bone Substitute

The substitute ATAO® is indicated for filling a small sized bone defects in orthopaedic surgery.

The proportion of hydroxyapatite(HA) and tricalcium phosphate(TCP) allows an optimal resorption rate compatible with the physiology of the bone cells.

References	Description
WED006	Bone substitute - 06°
WED008	Bone substitute - 08°
WED010	Bone substitute - 10°
WED012	Bone substitute - 12°
WED014	Bone substitute - 14°
FIL008	Granules 08 cC
FIL003	Granules 03 cC
FIL005	Granules 05 cC




*Made in France by KASIOS for Newdeal S.A.S.

Osteotomy Plates

Associated Instruments

References	Description
MA003	Surfix® plate bender
MA004	Hexagonal screwdriver 3.5 mm
MA005	Length gauge
MA006	Drilling guide diam 3.5 mm
MA007	Drilling guide diam 4.5 mm
MA008	Drill AO diam 3.5 mm
MA009	Drill AO diam 4.5 mm
MA010	cutting closing guide
MA038	K-wire diam 1.8 mm, L. 150 mm
MA039	Bone substitute holder
MA040	Trial for bone substitute
MA041	K-wire guide diam 1.8 mm for drill guide diam 3.5 mm
MA042	Surfix® distractor + 2 measurers
MA043	Compressor
MA047	Cutting opening guide
FPOM	Trial for POM plates
FPTOPTF	Trial plate for pto and ptf plates
FPOF	Trial for POF plates
FPTCD	Trial plate for PTCD
FPTCG	Trial plate for PTCG
SET000	Basic Tray
SET100	Osteotomy tray
SET001	Lid of osteotomy set



PSFE Plates

Proximal Varus and Derotation Osteotomy Plate

The osteotomy plates (PPSFE, PSFE and PPSFE D/G) are used for fixation of proximal femoral varus and derotation osteotomy for children.

- Materials:
 - Plate: Stainless steel 316 LVM
Standard ISO 5832-1 & ASTM F138&139
 - Screws and lock-screws: Titanium TA6V
Standard ISO 5832-3 & ASTM F136

References	Description
PPSFE	Varus osteotomy plate size 1
PSFE	Varus osteotomy plate size 2
PSFE D	Reinforced varus osteotomy plate - right side
PSFE G	Reinforced varus osteotomy plate - left side



PPSFE



PSFE



PSFE D
PSFE G

PSFE Plates

PSFE Plates Associated Surfix® Screws

Surfix® Cortical bone screw
+ lock-screw diam 3.5 mm

References	Description
CS 016	Length 16 mm
CS 018	Length 18 mm
CS 020	Length 20 mm
CS 022	Length 22 mm
CS 024	Length 24 mm
CS 026	Length 26 mm
CS 028	Length 28 mm
CS 030	Length 30 mm
CS 032	Length 32 mm



Surfix® Cortical bone screw
+ lock-screw diam 4.5 mm

References	Description
CC 014	Length 14 mm
CC 016	Length 16 mm
CC 018	Length 18 mm
CC 020	Length 20 mm
CC 022	Length 22 mm
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm



PSFE Plates

PSFE Associated Surfix[®] Screws

Surfix[®] Cancellous bone screw
+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm



Surfix[®] Cannulated cancellous bone screw + lock-screw diam 6.5 mm

References	Description
VCA 020	Length 20 mm
VCA 025	Length 25 mm
VCA 030	Length 30 mm
VCA 035	Length 35 mm
VCA 040	Length 40 mm
VCA 045	Length 45 mm
VCA 050	Length 50 mm
VCA 055	Length 55 mm
VCA 060	Length 60 mm



PDTL Plate

Bendable Plates for Osteosynthesis

Fixation of epiphyseal, metaphyseal, diaphyseal fractures.

- Bendable plate.
- 4 sizes.
- Adapted angulation and curve.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136



References	Description
PD3TL52	Bendable plate 3 holes
PD4TL72	Bendable plate 4 holes
PD5TL92	Bendable plate 5 holes
PD7TL132	Bendable plate 7 holes

PDTL Reinforced Plate

Bendable and Reinforced Plates for Osteosynthesis

Fixation of epiphyseal, metaphyseal, diaphyseal fractures.

- Bendable plate.
- Reinforced bridge.
- 4 sizes.
- Adapted angulation and curve.
- Materials:
 - > Plate: Stainless steel 316 LVM Standard ISO 5832-1 & ASTM F138&139
 - > Screws and lock-screws: Titanium TA6V Standard ISO 5832-3 & ASTM F136



References	Description
PD2TL45	Bendable and reinforced plate 2 holes
PD4TL77	Bendable and reinforced plate 4 holes
PD4TL112	Bendable and reinforced plate 4 holes

PDTL Plates

Associated Surfix® Screws

Surfix® Cortical bone screw

+ lock-screw diam 3.5 mm

References	Description
CS 016	Length 16 mm
CS 018	Length 18 mm
CS 020	Length 20 mm
CS 022	Length 22 mm
CS 024	Length 24 mm
CS 026	Length 26 mm
CS 028	Length 28 mm
CS 030	Length 30 mm
CS 032	Length 32 mm



Surfix® Cortical bone screw

+ lock-screw diam 4.5 mm

References	Description
CC 024	Length 24 mm
CC 026	Length 26 mm
CC 028	Length 28 mm
CC 030	Length 30 mm
CC 032	Length 32 mm
CC 035	Length 35 mm
CC 040	Length 40 mm
CC 045	Length 45 mm
CC 050	Length 50 mm



Surfix® Cancellous bone screw


+ lock-screw diam 6.5 mm

References	Description
SP 020	Length 20 mm
SP 025	Length 25 mm
SP 030	Length 30 mm
SP 035	Length 35 mm
SP 040	Length 40 mm
SP 045	Length 45 mm
SP 050	Length 50 mm
SP 055	Length 55 mm
SP 060	Length 60 mm
SP 065	Length 65 mm
SP 070	Length 70 mm
SP 075	Length 75 mm
SP 080	Length 80 mm
SP 085	Length 85 mm
SP 090	Length 90 mm
SP 095	Length 95 mm
SP 100	Length 100 mm



Pediatric Set

Reference	Description
MA003	SURFIX® Plate bender
MA004	Hexagonal screwdriver 3.5 mm
MA005	Length gauge
MA006	Drilling guide dia. 3.5 mm
MA031	Drill guide dia. 2.7 mm
MA034	Hexagonal screwdriver 2.5 mm
MA036	Hexagonal canulated screwdriver dia. 3.5 mm
MA048	Varus Angle measurer
MA049	Ruler
MA050	Angle measurer / K-wire Guide
FPPSFE	Varus osteotomy trial plate size 1
FPSFE	Varus osteotomy trial plate size 2
FPSFED	Reinforced varus osteotomy trial plate - Right side
FPSFEG	Reinforced varus osteotomy trial plate - Left side
FPD3TL52	Bendable Trial 3 holes
FPD4TL72	Bendable Trial 4 holes
FPD5TL92	Bendable Trial 5 holes
FPD7TL132	Bendable Trial 7 holes
FPD2TL45	Bendable and reinforced trial 2 holes
FPD4TL77	Bendable and reinforced Trial 3 holes
FPD4TL112	Bendable and reinforced Trial 4 holes
MA008	Drill AO dia. 3.5 mm
MA032	Drill dia. 2.7 mm, L. 125 mm
MA037	Canulated Drill dia. 3.5 mm
MA038	K-wire dia. 1.8 mm, L. 150 mm
MA041	K-wire guide dia. 1.8 mm for drill guide dia. 3.5 mm
SET400	Tray for Pediatric
SET401	Lid for Pediatric





Integra[®] References

Symbols

05-0016.....	40	28.25.126.....	105	119 525.....	37
05-0017.....	40	028, CC.....	120, 121, 122, 123, 124, 126, 129	119 531.....	37
05-0106.....	40	028, CS.....	126, 129	119 532.....	37
05-0108.....	40	030, CC.....	120, 121, 122, 123, 124, 126, 129	119 536.....	37
05-0109.....	40	030, CS.....	126, 129	129 710.....	37
05-0110.....	40	030, SP.....	120, 121, 122, 123, 124, 127, 129	159 004 (S).....	11
05-0112.....	40	030, VCA.....	127	159 023 (S).....	10
05-0206.....	40	032, CC.....	120, 121, 122, 123, 124, 126, 129	159 025 (S).....	10
05-0208.....	40	032, CS.....	126, 129	159 027 (S).....	10
05-0209.....	40	035, CC.....	120, 121, 122, 123, 124, 126, 129	159 223 (S).....	10
05-0210.....	40	035, SP.....	120, 121, 122, 123, 124, 127, 129	159 225 (S).....	10
05-0212.....	40	035, VCA.....	127	159 227 (S).....	10
05-1002.....	40	040, CC.....	120, 121, 122, 123, 124, 126, 129	189 210.....	18
05-1006.....	40	040, SP.....	120, 121, 122, 123, 124, 127, 129	189 215.....	18
05-1008.....	40	040, VCA.....	127	189 220.....	18
05-1009.....	40	045, CC.....	120, 121, 122, 123, 124, 126, 129	195 100S.....	18
05-1010.....	40	045, SP.....	120, 121, 122, 123, 124, 127, 129	219 027.....	18
05-1012.....	40	045, VCA.....	127	219 127.....	18
05-1020.....	40	050, CC.....	120, 121, 122, 123, 124, 129	219 227.....	18
05-5001.....	40	050, SP.....	120, 121, 122, 123, 124, 127, 129	229 027 S.....	10
014, CC.....	126	050, VCA.....	127	229 028 S.....	10
016, CC.....	126	055, CC.....	120, 121, 122, 123, 124	229 037 S.....	10
016, CS.....	126, 129	055, SP.....	120, 121, 122, 123, 124, 127, 129	229 038 S.....	10
018, CC.....	126	055, VCA.....	127	229 227 S.....	10
018, CS.....	126, 129	060, CC.....	120, 121, 122, 123, 124	229 228 S.....	10
020, CC.....	126	060, SP.....	120, 121, 122, 123, 124, 127, 129	229 237 S.....	10
020, CS.....	126, 129	060, VCA.....	127	229 238 S.....	10
020, SP.....	120, 121, 122, 123, 124, 127, 129	065, CC.....	120, 121, 122, 123, 124	295 110S.....	18
020, VCA.....	127	065, SP.....	120, 121, 122, 123, 124, 129	295 112S.....	18
022, CC.....	126	070, SP.....	120, 121, 122, 123, 124, 129	295 114S.....	18
022, CS.....	126, 129	075, SP.....	120, 121, 122, 123, 124, 129	295 116S.....	18
024, CC.....	120, 121, 122, 123, 124, 126, 129	080, SP.....	120, 121, 122, 123, 124, 129	295 118S.....	18
024, CS.....	126, 129	085, SP.....	120, 121, 122, 123, 124, 129	295 120S.....	18
025, SP.....	120, 121, 122, 123, 124, 127, 129	090, SP.....	120, 122, 123, 124, 129	295 122S.....	18
025, VCA.....	127	095, SP.....	120, 121, 122, 123, 124, 129	295 124S.....	18
026, CC.....	120, 121, 122, 123, 124, 126, 129	100, SP.....	120, 121, 122, 123, 124, 129	295 126S.....	18
026, CS.....	126, 129	115 008 S.....	10	295 128S.....	18
28.25.010.....	105	115 070.....	18	295 130S.....	18
28.25.012.....	105	115 070 (S).....	10, 11	295 132S.....	18
28.25.014.....	105	115 100 (S).....	10, 11	295 134S.....	18
28.25.016.....	105	115 116.....	18	299 005.....	18
28.25.018.....	105	119 004 (S).....	11	299 020.....	18
28.25.020.....	105	119 006 (S).....	11	299 050.....	18
28.25.022.....	105	119 016 (S).....	11	309 708.....	49
28.25.024.....	105	119 023 (S).....	10	320 100S.....	18
28.25.026.....	105	119 025 (S).....	10	320 200S.....	18
28.25.110.....	105	119 027 (S).....	10	329 002.....	18
28.25.112.....	105	119 223 (S).....	10	329 500.....	18
28.25.114.....	105	119 225 (S).....	10	329 900.....	18
28.25.116.....	105	119 227 (S).....	10	329 901.....	18
28.25.118.....	105	119 501.....	37	329 902.....	18
28.25.120.....	105	119 506.....	37	501 009 240.....	105
28.25.122.....	105	119 510.....	37	501 009 241.....	105
28.25.124.....	105	119 516.....	37	502 015206.....	105
		119 522.....	37	503 004157.....	106

503 004170	106	080003	82	111332S	26
503 004197	106	080005	82	111334	14, 87
503 004262	106	095070	108	111334S	26
503 004267	106	095080	108	111424	34
503 004268	106	095090	108	111426	34
503 004270	106	095160	108	111428	34
503 004331	106	095250	108	111430	34
503 006320	105	096000	108	111432	34
503 006343	105	096005	108	111434	34
515 010 070	105	105426	36	111436	34
515 012 070	105	105428	36	111440	34
515 016 070	105	105430	36	111445	34
515 018 070	105	105432	36	111450	34
515 020 070	105	105434	36	111455	34
515 022 070	105	105436	36	111460	34
515 024 070	105	105438	36	111530S	35
515 026 070	105	105440	36	111535S	35
516 014 070	105	105445	36	111540S	35
516 028 070	105	105450	36	111545S	35
517 030 070	105	105455	36	111550S	35
519 010 150	105	105460	36	111555S	35
519 012 150	105	105640	36	111560S	35
519 014 150	105	105645	36	111565S	35
519 016 150	105	105650	36	111570S	35
519 018 150	105	105655	36	111575S	35
519 020 150	105	105660	36	111580S	35
519 022 150	105	105665	36	111740S	35
519 024 150	105	105670	36	111745S	35
519 026 150	105	105675	36	111750S	35
520 028 150	105	105680	36	111755S	35
521 030 150	105	105685	36	111760S	35
605.090.105	105	105690	36	111765S	35
605 120 100	105	111010	14, 87	111770S	35
605 300230	111	111012	14, 87	111775S	35
605 303230	111	111014	14, 87	111780S	35
060005	94	111016	14, 87	111785S	35
060006	94	111018	14, 87	111790S	35
062806	94	111020	14, 87	111795S	35
062808	94	111022	14, 87	111800S	35
062810	94	111024	14, 87	111805S	35
062812	94	111026	14, 87	111810S	35
062814	94	111028	14, 87	111815S	35
062816	94	111312	14, 87	111820S	35
062818	94	111314	14, 87	112011	15, 16, 89
070005	94	111316	14, 87	112012	15, 16, 89
070006	94	111318	14, 87	112013	15, 16, 89
072408	94	111320	14, 87	112014	15, 16, 89
072410	94	111322	14, 87	112111	15
072412	94	111324	14, 87	112111S	15
072414	94	111326	14, 87	112112	15
072810	94	111328	14, 87	112112S	15
072814	94	111328S	26	112113	15
072818	94	111330	14, 87	112113S	15
075020	95	111330S	26	112114(s)	15
080001	82	111332	14, 87	112114S	15

113113.....	88	117350.....	24	119641.....	47
113114.....	88	117440.....	24	119645.....	21
113115.....	88	117445.....	24	119646.....	47
113116.....	88	117450.....	24	119647.....	47
113117.....	88	119023S.....	85, 86	119651.....	65, 67
113213.....	88	119025S.....	85, 86	119653.....	65, 67
113214.....	88	119027S.....	85, 86	119654.....	65, 67
113215.....	88	119050.....	17, 66	119657.....	49
113216.....	88	119104.....	34	119659.....	49
113217.....	88	119105.....	10	119662.....	47
113313.....	88	119106.....	17, 66, 67	119664.....	47, 65
113314.....	88	119133.....	27	119665.....	47
113315.....	88	119135.....	27	119666.....	47
113316.....	88	119138.....	34	119667.....	47
113317.....	88	119139.....	34	119669.....	47
113512.....	88	119150.....	35	119690.....	47
113820.....	88	119151.....	35	119691.....	47
114438.....	34	119152.....	35	119692.....	47
115008S.....	85	119153.....	35	119693.....	47
115070.....	27	119154.....	35	119694.....	47
115070S.....	85	119155.....	35	119695.....	47
115100.....	25, 27, 42	119156.....	35	119696.....	47
115100S.....	85	119170.....	35	119697.....	47
115101.....	43	119171.....	35	119698.....	47
115116.....	27, 34, 38, 47, 63	119172.....	35	119699.....	47
115216.....	21, 25, 28, 60	119173.....	35	119700.....	64
115225.....	37, 38, 47, 63	119174.....	35	119751.....	47
115425.....	60	119175.....	35	119757.....	49
115516.....	35	119176.....	35	119760.....	49
117010.....	25	119177.....	35	119761.....	49
117012.....	25	119223S.....	86	119762.....	49
117014.....	25	119225S.....	86	119763.....	49
117016.....	25	119227S.....	86	119764.....	49
117018.....	25	119311.....	11	119765.....	49
117020.....	25	119510.....	37	119766.....	49
117022.....	25	119516.....	34	119767.....	49
117024.....	25	119531.....	38	119768.....	49
117026.....	25	119541.....	37	119769.....	49
117028.....	25	119546.....	37	119778.....	47
117030.....	25	119551.....	37	119779.....	47
117032.....	25	119552.....	37, 58	119790.....	49
117034.....	25	119555.....	37	119791.....	49
117110.....	25	119561.....	37, 38	119792.....	49
117112.....	25	119562.....	37	119793.....	49
117114.....	25	119566.....	37	119794.....	49
117116.....	25	119580.....	38	119795.....	49
117118.....	25	119582.....	38	119796.....	49
117150.....	24	119590.....	38	119797.....	49
117155.....	24	119592.....	38	119798.....	49
117160.....	24	119601.....	21	119799.....	49
117250.....	24	119609.....	47	119825.....	39
117255.....	24	119611.....	21, 47	119831.....	39
117260.....	24	119618.....	25	119850.....	39
117340.....	24	119631.....	49	119859.....	39
117345.....	24	119632.....	49	119860.....	39

119861	39	121805S	35	139216	21
119862	39	121810S	35	139311-313	21
119863	39	121815S	35	139321-323	21
119864	39	121820S	35	139331-333	21
119865	39	121950	35	139341-343	21
119867	39	121951	35	139351-353	21
119880	39	129710	25, 27, 28, 37	139940	21
119881	39	129714	25, 27, 28	139941	21
119882	39	129716	25, 27, 28	139942	21
11990919, 21, 25, 27, 29, 35, 37, 42, 43, 47, 56, 58, 60, 63, 67		129718	25, 27, 28	139943	21
119935	34	129720	25, 27, 28	139953	21
119936	34	129722	27, 28	141009	39
119937	34	129724	25, 27, 28	141010	39
119938	34	129726	25, 27, 28	141011	39
119949	37	129728	25, 27, 28	141012	39
119951	37	129730	25, 27, 28	141013	39
119953	37	129731	25	141014	39
119955	38	129732	25	141015	39
119956	37	129733	25	141017	39
119957	37	129734	25	150010S	61
119958	37	129735	25	150020S	61
119959	38	129736	25	150030S	61
119965	49	129752	27, 28	150040S	61
119966	49	129930	25	150110S	61
119967	49	129931	25	150120S	61
119974	37	129932	25	150130S	61
119975	38	129933	25	150200S	62
119976	37	129935	28	150240S	62
119977	38	129936	25	150242S	62
119978	49	129938	28	150246S	62
119989	67	131014S	14, 87	150250S	62
121530S	35	131016S	14, 87	150255S	62
121535S	35	131018S	14, 87	150260S	62
121540S	35	131020S	14, 87	150265S	62
121545S	35	131022S	14, 87	150270S	62
121550S	35	131024S	14, 87	150275S	62
121555S	35	131026S	14, 87	150280S	62
121560S	35	131028S	14, 87	150285S	62
121565S	35	131030S	14, 87	150290S	62
121570S	35	131032S	14, 87	150514S	61
121575S	35	131034S	14, 87	150516S	61
121580S	35	139100-106	21	159010	63
121740S	35	139107	21	159020	63
121745S	35	139111-116	21	159023S	86
121750S	35	139120	21	159025S	86
121755S	35	139121	21	159027S	27, 86
121760S	35	139121-126	21	159071	63
121765S	35	139128	21	159103	63
121770S	35	139129	21	159106	63
121775S	35	139140	21	159109	63
121780S	35	139150	21	159110	63
121785S	35	139201-205	21	159120	63
121790S	35	139210	21	159127	47
121795S	35	139211	21	159130	63
121800S	35	139211-215	21	159223S	86

159225S	86	186300SND	92	211200	83
159227S	86	195100S	26	211300	83
159400	63	195200S	26	211409	83
159635	63	200010(S)	14	211411	83
159740	63	200010S	87	211501	83
159755	63	200012S	14, 87	211511	83
159760	63	200014S	14, 87	211600	83
159960	63	200016S	14, 87	211601	83
159991	63	200018S	14, 87	211610	83
169001	60	200020S	14, 87	211611	83
169002	60	200022S	14, 87	211622	83
169003	60	200024S	14, 87	211625	83, 91
169004	60	200026S	14, 87	211628	83
169021	42	200028S	14, 87	211635	83
169022	42	200030S	14, 87	211700	83
169023	42	201010	22	211720	83
169031	42	201020	22	211725	83
169032	42	201030	22	211801	83
169033	42	201040	22	211811	83
169041	42	202010	22	211900	83
169042	42	202020	22	211910	83
169051	42	202030	22	213113	16, 88
169052	42	202040	22	213114	16, 88
169101	60	203010	22	213115	16, 88
169102	42	203020	22	213116	16, 88
169104	60	203030	22	213117	16, 88
169105	60	203040	22	213213	16, 88
169106	42, 43	203100	22	213214	16, 88
169107	42	203200	22	213215	16, 88
169111	60	203300	22	213216	16, 88
169112	42	203415	22	213217	16, 88
181001S	59	203500	22	213313	16, 88
181002S	59	203600	22	213314	16, 88
181003S	59	205000	22	213315	16, 88
181004S	59	210310	83	213316	16, 88
181011S	59	210312	83	213317	16, 88
181012S	59	210314	83	213512	16, 88
181013S	59	210316	83	213820	16, 88
181014S	59	210318	83	215001	83
181021S	41	210320	83	219027	27
181022S	41	210322	83	219035	63
181023S	41	210324	83	219127	27
181031S	41	210326	83	219135	63
181032S	41	210414	83	219227	27
181033S	41	210416	83	219335	42, 43, 63
181041S	41	210418	83	219345	60
181042S	41	210420	83	219365	60
181051S	41	210422	83	219435	42, 43, 63
181052S	41	210424	83	219445	60, 63
185100S	26	210426	83	219465	60
185200S	26	210428	83	219535	42, 43, 63
185300S	41, 43, 62	210430	83	219545	60, 63
185400S	59	210432	83	219565	60
185600S	59	210434	83	219627	27
186200SND	92	211100	83	219635	42, 43, 63

219645	60	229975	12, 85	262201	90
219665	60	229976	13, 85	262202	90
219735	42, 43	229990	8, 84	262300	90
219765	60	229992	85	262301	90
219835	42, 43, 63	229996	8, 84	262302	90
219845	60, 63	230001S	19	262400	90
219865	60	230002S	19	262401	90
219900	60	239000	19	262402	90
219960	63	239001	19	263100L	90
221001	107	239011	19	263100R	90
221418	107	239012	19	263200L	90
221421	107	239021	19	263200R	90
221424	107	239022	19	263300L	90
221665	107	239030	19	263300R	90
221675	107	239500	19	263400L	90
221800	107	250100	20	263400R	90
224165	107	250101	20	264515	90
224175	107	250102S	20	264515S	90
224318	107	250103	20	264520	90
224321	107	250104	20	264520S	90
224324	107	250105	20	264525	90
225205	107	250106	20	264525S	90
225365	107	250201	20	264530	90
225375	107	250202	20	264530S	90
225400	107	250203	20	264535	90
225405	107	250204	20	264535S	90
225600	107	250205	20	265100	91
225900	107	250311	20	265200	91
226000	107	250312	20	265300	91
229002	10, 85	250321	20	265400	91
229003	10, 85	250322	20	266100	91
229004	10, 85	250323	20	266101	91
229005	10, 85	250331	20	266102	91
229011	10, 17, 66	250332	20	266200	91
229012	10, 17, 66	250333	20	266201	91
229027S	86	250341	20	266202	91
229028S	86	250342	20	266300	91
229037S	86	250343	20	266301	91
229038S	86	250351	20	266302	91
229101	10, 11, 85	250352	20	266400	91
229102	11	250353	20	266401	91
229201	11	260100	91	266402	91
229202	11	260515	91	267100LT	91
229227S	85, 86	260520	91	267100RT	91
229228S	86	260525	91	267200LT	91
229237S	85, 86	260530	91	267200RT	91
229238S	86	260535	91	267300LT	91
229951	10	261101	90	267300RT	91
229952	11	261201	90	267400LT	91
229961	10, 85	261301	90	267400RT	91
229962	11	261401	90	268100	91
229965	12, 85	262100	90	268104	91
229966	13, 85	262101	90	268111	91
229970	8, 10, 11, 84	262102	90	268114	91
229971	10	262200	90	268120	91

268121.....	91	285124S	26	285650s	59
268122	91	285126S	26	285655s	59
268123	91	285210S	26	285660s	59
268125	91	285212S	26	285665s	59
268136	91	285214S	26	286210SND	92
268145	91	285216S	26	286212SND	92
268220.....	91	285218S	26	286214SND	92
268230LF.....	91	285220S	26	286216SND	92
268230RT	91	285222S	26	286218SND	92
268231.....	91	285224S	26	286220SND	92
268232	91	285226S	26	286222SND	92
268233	91	285310S	41, 62	286224SND	92
268325	91	285312S	41, 43, 62	286312SND	92
268335	91	285314S	41, 43, 62	286314SND	92
268340.....	91	285316S	41, 43, 62	286316SND	92
268410L.....	91	285318S	41, 43, 62	286318SND	92
268410R.....	91	285320S	41, 43, 62	286320SND	92
268420L	91	285322S	41, 43, 62	286322SND	92
268420R.....	91	285324S	41, 43, 62	286324SND	92
268430L	91	285326S	41, 43, 62	286326SND	92
268430R.....	91	285328S	41, 43, 62	286328SND	92
268440L	91	285330S	41, 43	290150S	26
268440R.....	91	285330S	62	290155S	26
268510	91	285332S	41, 43, 62	290160S	26
268610	91	285334S	41, 43, 62	290250S	26
268640	91	285336S	43	290255S	26
268650	91	285336S	62	290260S	26
268660	91	285338S	43	290340S	26
268700	91	285338S	62	290345S	26
268706	91	285340S	43	290350S	26
268725	91	285340S	62	290440S	26
275560.....	47	285344S	62	290445S	26
278002S	29	285348S	62	290450S	26
278005S	29	285350S	62	295110S	26
278007S	29	285414S	59	295112S	26
278008S	29	285416S	59	295114S	26
278009S	29	285418S	59	295116S	26
278010S	29	285420S	59	295118S	26
278100	29	285422S	59	295120S	26
278110	29	285424S	59	295122S	26
278119	29	285426S	59	295124S	26
278121	29	285428S	59	295126S	26
278122.....	29	285430S	59	295210S	26
278123.....	29	285435S	59	295212S	26
278124.....	29	285440S	59	295214S	26
278900	29	285445S	59	295216S	26
278901	29	285450S	59	295218S	26
278902.....	19, 29	285455S	59	295220S	26
285110S	26	285460S	59	295222S	26
285112S	26	285465S	59	295224S	26
285114S	26	285625S	59	295226S	26
285116S	26	285630S	59	295310S	62
285118S	26	285635S	59	295310S	62
285120S	26	285640S	59	295312S	62
285122S	26	285645S	59	295314S	62

295316S	62	300406	45	302222	48
295318S	62	300407	45	302223	48
295320S	62	300409	45	302224	48
295322S	62	300505	45	302225	48
295324S	62	300506	45	302232	48
295326S	62	300507	45	302233	48
295328S	62	300509	45	302234	48
295330S	62	300605	45	302235	48
295332S	62	300606	45	303100S	92
295334S	62	300607	45	303101S	92
295336S	62	300609	45	303102S	92
295338S	62	301110	45	303210CS	92
295340S	62	301111	45	303212CS	92
295344S	62	301112	45	303214CS	92
295348S	62	301113	45	303216CS	92
295350S	62	301114	45	303218CS	92
299005	27	301115	45	303220CS	92
299010	27	301116	45	303222CS	92
299020	27	301121	48	303224CS	92
299050	27	301122	48	303312CS	92
299070	27	301123	48	303312S	48
299080	27	301124	48	303314CS	92
299150	27	301125	48	303316CS	92
299155	27	301201	44	303316S	48
299160	27	301202	44	303318CS	92
299250	27	301203	44	303320CS	92
299255	27	301204	44	303320S	48
299260	27	301205	44	303322CS	92
299340	27	301206	44	303324CS	92
299345	27	301222	48	303326CS	92
299350	27	301223	48	303328CS	92
299440	27	301224	48	303328S	45
299445	27	301225	48	303330S	45
299450	27	301232	48	303334S	45, 48
299900	27	301233	48	303338S	45, 48
299901	27	301234	48	303340S	45, 48
299903	27	301235	48	303342S	45, 48
299904	27	302110	45	303400	93
300005	45	302111	45	303401	93
300006	45	302112	45	303402	93
300007	45	302113	45	303403	93
300009	45	302114	45	303404	93
300105	45	302115	45	303405	93
300106	45	302116	45	303406	93
300107	45	302121	48	303407	93
300109	45	302122	48	303408	93
300205	45	302123	48	303409	93
300206	45	302124	48	303410	93
300207	45	302125	48	303411	93
300209	45	302201	44	303412	93
300305	45	302202	44	303450	93
300306	45	302203	44	303451	93
300307	45	302204	44	303452	93
300309	45	302205	44	303502	93
300405	45	302206	44	303504	93

303510	93	309680	47	351513	53
303512S	53	309681	47	351514	53
303516S	53	309682	47	351515	53
303520S	53	309683	47	351522	52
303530S	51	309684	47	351523	52
303534S	51, 53	309685	47	351524	52
303538S	51, 53	309686	47	351525	52
303540S	51, 53	309690	47	351532	52
303542S	51, 53	309691	47	351533	52
309200	47	309692	47	351534	52
309300	47	309693	47	351535	52
309301	47	309694	47	352152	51
309302	47	309695	47	352153	51
309303	47	309696	47	352154	51
309304	47	309697	47	352155	51
309305	47	309699	47	352251	50
309306	47	309712	49	352252	50
309309	47	309753	47	352253	50
309360	47	309773	47	352254	50
309361	47	309900	47	352255	50
309362	47	309902	47	352256	50
309363	47	309904	47	352512	53
309364	47	309920	47	352513	53
309365	47	309940	47	352514	53
309366	47	309942	47	352515	53
309370	47	309960	47	352522	52
309371	47	309970	47	352523	52
309372	47	309980	47	352524	52
309373	47	330021S	43	352525	52
309374	47	330023S	43	352532	52
309375	47	330025S	43	352533	52
309376	47	330217S	43	352534	52
309601	47	330220S	43	352535	52
309604	47	330225S	43	500001	54
309605	47	330230S	43	500050	54
309606	47	330420S	43	500080	54
309607	47	330425S	43	500150	54
309608	47	330430S	43	500180	54
309615	47	339001	43	500250	54
309620	47	339003	43	500280	54
309622	47	339004	43	500350	54
309623	47	339005	43	500380	54
309624	47	339920	43	501020	55
309625	47	339950	43	501022	55
309626	47	351152	51	501025	55
309627	47	351153	51	501027	55
309630	47	351154	51	501030	55
309635	47	351155	51	501035	55
309636	47	351251	50	501040	55
309637	47	351252	50	501045	55
309645	47, 58	351253	50	501050	55
309652	47	351254	50	501055	55
309656	47	351255	50	501060	55
309657	47	351256	50	501065	55
309661	47	351512	53	501070	55

A

AWL-100-01	72, 104, 75, 77, 79, 99, 80, 99
AWL-200-00	72, 73, 75, 77, 79, 80
AWL-420-00	78
AWL-520-01	73

B

BOD	114
BOD-0920-020-08STD	114
BOD-0920-020-10LRG	114
BOD-0920-020-10SML	114
BOD-0920-020-10STD	114
BOD-0920-020-12STD	114
BOD-0920-020-14LRG	114
BOD-0920-020-14SML	114
BOD-0920-020-14STD	114
BOD-0950-020-08LRG	114
BOD-0950-020-08SML	114
Bold® Compression Screw	87
BRH-100-05D	75, 77
BRH-100-05P	75, 77
BRH-100-10D	75, 77
BRH-100-10P	75, 77
BRH-100-20D	75, 77
BRH-100-20P	75, 77
BRH-100-30D	75, 77
BRH-100-30P	75, 77
BRH-100-40D	75, 77
BRH-100-40P	75, 77
BRH-100-50D	75, 77
BRH-100-50P	75, 77
BRH-200-10D	72
BRH-200-10P	72
BRH-200-20D	72
BRH-200-20P	72
BRH-200-30D	72
BRH-200-30P	72
BRH-200-40D	72
BRH-200-40P	72
BRH-300-20	104
BRH-300-22	104
BRH-300-24	104
BRH-350-04	104
BRH-410-10	80
BRH-410-20	80
BRH-410-30	80
BRH-410-40	80
BRH-442-10	79
BRH-442-20	79
BRH-442-30	79
BRH-442-40	79
BRH-445-10	81
BRH-445-20	81
BRH-445-30	81
BRH-445-40	81

BRH-445-50	81
BRH-525-0D	73
BRH-525-0P	73
BRH-525-1D	73
BRH-525-1P	73
BRH-525-2D	73
BRH-525-2P	73
BRH-525-3D	73
BRH-525-3P	73
BRH-525-4D	73
BRH-525-4P	73
BRH-525-5D	73
BRH-525-5P	73
BRH-600-45	99
BRH-600-55	99
BRH-600-65	99
BRH-900-00	111
BRH-900-01	111
BRH-900-02	111
BRH-900-03	111
BRH-910-38	116
BRH-910-41	116
BRH-910-44	116
BRH-910-47	116
BRH-910-50	116
BRH-910-53	116
BRH-910-56	116
BSW-0920-021-01	114

C

CC 024	120, 121, 122, 123, 124, 129
CC 026	120, 121, 122, 123, 124, 129
CC 028	120, 121, 122, 123, 124, 129
CC 030	120, 121, 122, 123, 124, 129
CC 032	120, 121, 122, 123, 124, 129
CC 035	120, 121, 122, 123, 124, 129
CC 040	120, 121, 122, 123, 124, 129
CC 045	120, 121, 122, 123, 124, 129
CC 050	120, 121, 122, 123, 124, 129
CC 055	120, 121, 122, 123, 124
CC 060	120, 121, 122, 123, 124
CC 065	120, 121, 122, 123, 124
Cemented Stem Implants	114
CHRA-910-38/14	116
CHRA-910-41/15	116
CHRA-910-41/18	116
CHRA-910-44/16	116
CHRA-910-44/19	116
CHRA-910-47/17	116
CHRA-910-47/20	116
CHRA-910-50/18	116
CHRA-910-50/21	116
CHRA-910-53/19	116
CHRA-910-53/22	116

CHRA-910-56/21	116
CMC-410-10	80
CMC-410-20	80
CMC-410-30	80
CMC-410-40	80
CMC-DIS-HAL	79
CMC-DIS-STR	79
CMC Saddle™ implant	80
CMRH-310-H20-WW	103
CMRH-310-H22-WW	103
CMRH-310-H24-WW	103
CMRH-310-So1L-WW	103
CMRH-310-So1S-WW	103
CMRH-310-So2L-WW	103
CMRH-310-So2S-WW	103
CMRH-310-So3L-WW	103
CMRH-310-So3S-WW	103
CMRH-310-So4L-WW	103
CMRH-310-So4S-WW	103
CoCr Humeral Head Implants	114
CS 016	129
CS 018	129
CS 020	129
CS 022	129
CS 024	129
CS 026	129
CS 028	129
CS 030	129
CS 032	129
CSA-000-01	72, 104, 75, 77, 80
CSA-000-03	81
CSA-000-04	73, 79
CSA-000-05	78
CSA-000-07	116
CSA-100-05	75, 77
CSA-100-06	75, 77
CSA-170-02	96
CSA-200-02	72
CSA-200-03	72
CSA-350-02	104
CSA-410-02	80
CSA-420-02	78
CSA-430-01	81
CSA-443-02	79
CSA-443-03	79
CSA-443-04	79
CSA-520-05	73
CSA-520-06	73
CSA-520-07	73
CSA-710-01	96
CSA-850-00	106
CSA-850-01	106
CSA-850-02	106
CSA-850-03	106
CSA-860-01	106

CSA-860-02	106
CSA-860-03	106
CSA-870-01	106
CSA-870-02	106
CSA-870-03	106
CSA-900-01	111
CSA-900-02	111
CSA-900-03	111
CSA-900-04	111
CSA-910-02	116
CSA-910-03	116
CSA-925-01	116
CTD-900-01	111
CTD-900-02	111
CTD-900-03	111
CUT-443-L	79
CUT-443-M	79
CUT-443-S	79

D

DCG-890-00	23
DRL-890-00	23
DRUJ-610-1445	98
DRUJ-610-1455	98
DRUJ-610-1465	98
DRUJ-610-1645	98
DRUJ-610-1655	98
DRUJ-610-1665	98
DRUJ-610-1745	98
DRUJ-610-1755	98
DRUJ-610-1765	98
DRUJ-610-1945	98
DRUJ-610-1955	98
DRUJ-610-1965	98
DRUJ implant	98

E

Elbow pathologies solutions	102
ENDORELEASE™ ENDOSCOPIC CUBITAL TUNNEL RELEASE SYSTEM	108
EXT-100-01	75, 77, 79
EXT-200-00	72, 75, 77, 79, 80, 99
EXT-350-00	104
EXT-430-00	81
EXT-910-00	116

F

FIL003	125
FIL005	125
FIL008	125
FPD2TL45	130
FPD3TL52	130

FPD4TL72	130
FPD4TL77	130
FPD4TL112	130
FPD5TL92	130
FPD7TL132	130
FPDRA1	100
FPDRA2	100
FPDRA3	100
FPOF	125
FPOM	125
FPPSFE	130
FPSFE	130
FPSFED	130
FPSFEG	130
FPTCD	125
FPTCG	125
FPTOPTF	125

G

GDW-890-00	23
GLN-0920-030-00K	114
GLN-0920-030-00P	114
GLN-0920-030-01K	114
GLN-0920-030-01P	114
GLN-0920-030-02K	114
GLN-0920-030-02P	114
GLN-0920-030-03K	114
GLN-0920-030-03P	114
GOU-420-00	78

H

Hand Fractures and Deformity repair set	84
HINTEGRA® Total Ankle Prosthesis	44, 50, 52
HINTEGRA® Total Ankle Prosthesis – Flat Cut	49
Hintermann® distractors	65
HRA-900-40/16-WW	111
HRA-900-44/15-WW	111
HRA-900-44/18-WW	111
HRA-900-48/17-WW	111
HRA-900-48/20-WW	111
HRA-900-52/18-WW	111
HRA-900-52/21-WW	111
HRA-900-56/21-WW	111

I

I.CO.S.® Ideal COmpression Screw	36
I.CO.S.® Removal Set	38
IMP-100-00D	75, 77
IMP-100-00P	75, 77
IMP-200-00D	72
IMP-200-00P	72
IMP-300-00	104, 99, 99

IMP-350-00	104
IMP-350-01	104
IMP-410-00	80
IMP-442-00	79
IMP-600-00	99
IMP-600-01	99
IMP-890-00MT	23
IMP-890-00PP	23
IMP-900-01	111
IMP-910-00	116
INS-200-01	72
INS-350-00	104
INS-430-00	81
INS-600-00	99
INS-610-00	99
INS-860-00	106
INS-890-00	23
INS-900-00	111
INS-910-01	116
INS-920-02	116
Integra® CMC implants	78
ITWFS instrumentation	93
ITWFS Integra® Total Wrist Fusion System	92
ITWFS set	93

K

KALIX® II Flat Foot Implant	39
KATALYST™ TELESCOPING BIPOLAR RADIAL HEAD SYSTEM	107
Keeled Glenoid Implants	114
K-FIX® – K-wire end protector	66
KOMPRESSOR™ COMPRESSION SCREW SYSTEM	83

L

Large QWIX®	35
LUN-710-01-WW	96
LUN-710-02-WW	96
LUN-710-03-WW	96
LUN-710-04-WW	96
LUN-710-05-WW	96

M

MA003	125, 130
MA004	125, 130
MA005	125, 130
MA006	125, 130
MA007	125
MA008	125, 130
MA009	125
MA010	125
MA031	130
MA032	130

MA034.....	130
MA036.....	130
MA037.....	130
MA038.....	125, 130
MA039.....	125
MA040.....	125
MA041.....	125, 130
MA042.....	125
MA043.....	125
MA047.....	125
MA048.....	130
MA049.....	130
MA050.....	130
MA3501.....	100
MA3504.....	100
MA3505.....	100
MA3506.....	100
MA3508.....	100
MA3511.....	41, 42, 43, 59, 60, 61, 62, 67
MA5303.....	100
MAL-910-00.....	116
MCP-100-10D-WW.....	74
MCP-100-10P-WW.....	74
MCP-100-20D-WW.....	74
MCP-100-20P-WW.....	74
MCP-100-30D-WW.....	74
MCP-100-30P-WW.....	74
MCP-100-40D-WW.....	74
MCP-100-40P-WW.....	74
MCP-100-50D-WW.....	74
MCP-100-50P-WW.....	74
MCP-110-05D-WW.....	74
MCP-110-05P-WW.....	74
MCP implants.....	74
Metal radial head.....	104
MGT-890-10MT.....	23
MGT-890-10PPH.....	23
MGT-890-10PPT.....	23
MGT-890-20MT.....	23
MGT-890-20PPH.....	23
MGT-890-20PPT.....	23
MGT-890-30MT.....	23
MGT-890-30PPH.....	23
MGT-890-30PPT.....	23
MGT-890-40MT.....	23
MGT-890-40PPH.....	23
MGT-890-40PPT.....	23
MHH-0920-010-3814E.....	114
MHH-0920-010-4015E.....	114
MHH-0920-010-4216C.....	114
MHH-0920-010-4216E.....	114
MHH-0920-010-4416C.....	114
MHH-0920-010-4416E.....	114
MHH-0920-010-4419C.....	114
MHH-0920-010-4419E.....	114

MHH-0920-010-4614C.....	114
MHH-0920-010-4614E.....	114
MHH-0920-010-4617C.....	114
MHH-0920-010-4617E.....	114
MHH-0920-010-4620C.....	114
MHH-0920-010-4620E.....	114
MHH-0920-010-4815C.....	114
MHH-0920-010-4815E.....	114
MHH-0920-010-4818C.....	114
MHH-0920-010-4818E.....	114
MHH-0920-010-4821C.....	114
MHH-0920-010-4821E.....	114
MHH-0920-010-5019C.....	114
MHH-0920-010-5019E.....	114
MHH-0920-010-5022C.....	114
MHH-0920-010-5022E.....	114
MHH-0920-010-5220C.....	114
MHH-0920-010-5220E.....	114
MIS-890-00.....	23
MRH-350-01.....	104
MRH-350-02.....	104
MRH-350-03.....	104
MRH-350-04.....	104
MRH-350-20L.....	104
MRH-350-20S.....	104
MRH-350-22L.....	104
MRH-350-22S.....	104
MRH-350-24L.....	104
MRH-350-24S.....	104
MRH Implants.....	104
MRH-PLN-00.....	104
MRH-RMR-01.....	104
MRH-RMR-02.....	104
MRH-RMR-03.....	104
MSR-890-10/20.....	23
MSR-890-30/40.....	23

N

NUG-443-10M.....	79
NUG-443-10S.....	79
NUG-443-20M.....	79
NUG-443-20S.....	79
NUG-443-30L.....	79
NUG-443-30M.....	79
NUG-443-40L.....	79
NUG-443-40M.....	79
NuGrip™.....	79

O

OSG-100-00D.....	75, 77
OSG-100-00P.....	75, 77
OSG-200-00.....	72, 80
OSG-200-10P.....	72

OSG-200-10PV.....	72
OSG-200-20P.....	72
OSG-200-20PV.....	72
OSG-200-30P.....	72
OSG-200-30PV.....	72
OSG-200-40P.....	72
OSG-200-40PV.....	72
OSG-350-01.....	104
OSG-442-00.....	73, 79
OSG-600-00.....	99
OSG-600-01.....	99
OSG-610-00.....	99
OSG-910-01.....	116
OSG-910-02.....	116
OSG-910-03.....	116
Osteoguides.....	66

P

PANTA® Arthrodesis Nail.....	54
PANTA® Removal set.....	58
PBI-800-00.....	106
PCS-430-10-WW.....	81
PCS-430-20-WW.....	81
PCS-430-30-WW.....	81
PCS-430-40-WW.....	81
PCS-430-50-WW.....	81
PCSTM spacer*.....	81
PD2TL45.....	128
PD3TL52.....	128
PD4TL72.....	128
PD4TL77.....	128
PD4TL112.....	128
PD5TL92.....	128
PD7TL132.....	128
PDRA1.....	100
PDRA2.....	100
PDRA3.....	100
PDTL.....	128, 129
Pediatric set.....	130
Pegged Glenoid Implants.....	114
PIP-200-10D-WW.....	72
PIP-200-10P-WW.....	72
PIP-200-20D-WW.....	72
PIP-200-20P-WW.....	72
PIP-200-30D-WW.....	72
PIP-200-30P-WW.....	72
PIP-200-40D-WW.....	72
PIP-200-40P-WW.....	72
PIP implants.....	72
PIP Instrument Set.....	72
PIS-890-00.....	23
POF.....	120
POM.....	122
POMD.....	122

POMG.....	122
Post-Op Shoes	68
PPSFE.....	126
Press Fit Stem Implants.....	114
Primary Body Implants.....	114
PSFE	126, 127
PSFE D	126
PSFE G	126
PSH	117
PSH2	117
PSH3	117
PSH plate – Surfix® Proximal humeral plate	117
PSR-890-10/20.....	23
PSR-890-30/40	23
PTC.....	123
PTCD.....	123, 124
PTCG.....	123
PTF	124
PTFG	124
PTO	121
PTOD	121
PTOG	121
Pycarbon Radial Head implants (CMRH).....	103
PYD-420-145-WW	78
PYD-420-165-WW	78
PYD-420-167-WW.....	78
PYD-420-168-WW	78
PYD-420-187-WW	78
PYD-420-189-WW	78
PyroCarbon Lunate Implant.....	96
Pyrocarbon MCP implant	74
Pyrocarbon PIP implant	72
PyroDisk™	78
PyroDisk™ implants.....	78
PyroDisk™ instrument set	78
PyroTITAN HRA	116
PyroTITAN HRA Implants.....	116

Q

QWIX® Fixation Screw dia. 3.0 mm.....	87
QWIX® Fixation Screw dia.4.3 mm.....	34

R

RadFX™ implants	105
RET-443-00	79
RET-900-01	111
RET-920-10L.....	116
RET-920-10S.....	116
RET-920-20	116
RET-920-30	116
RET-920-40L.....	116
RET-920-40S.....	116
RET-920-50L.....	116

RET-920-50S	116
RET-920-60.....	116
RET-920-60L.....	116
RET-920-60S.....	116
RET-920-60XL	116

S

Safeguard™ Mini-open carpal tunnel release system	82
SET000	125
SET001	125
SET100	125
SET400.....	130
SET401.....	130
Shoulder Retractors (Coated)	116
SHP-425-14	78
SHP-425-16	78
SHP-425-18	78
SHP-443-L	79
SHP-443-M.....	79
SHP-443-S	79
SHP-900-40/16.....	111
SHP-900-44/15.....	111
SHP-900-44/18.....	111
SHP-900-48/17	111
SHP-900-48/20	111
SHP-900-52/18	111
SHP-900-52/21.....	111
SHP-900-56/21.....	111
SHP-910-38-14	116
SHP-910-41-15.....	116
SHP-910-41-18.....	116
SHP-910-44-16	116
SHP-910-44-19.....	116
SHP-910-47-17.....	116
SHP-910-47-20	116
SHP-910-50-18.....	116
SHP-910-50-21.....	116
SHP-910-53-19.....	116
SHP-910-53-22.....	116
SHP-910-56-21.....	116
Silicon MCP implant.....	76
Silicon PIP implant	73
SMCP-500-05-WW.....	76
SMCP-500-10-WW	76
SMCP-500-20-WW.....	76
SMCP-500-30-WW.....	76
SMCP-500-40-WW.....	76
SMCP-500-50-WW.....	76
SMCP implants	76
SP 020	120, 121, 122, 123, 124, 117
SP 025.....	120, 121, 122, 123, 124, 117
SP 030	120, 121, 122, 123, 124, 117
SP 035.....	120, 121, 122, 123, 124, 117
SP 040	120, 121, 122, 123, 124, 117

SP 045	120, 121, 122, 123, 124, 117
SP 050	120, 121, 122, 123, 124, 117
SP 055.....	120, 121, 122, 123, 124, 117
SP 060.....	120, 121, 122, 123, 124, 117
SP 065	120, 121, 122, 123, 124, 117
SP 070	120, 121, 122, 123, 124, 117
SP 075.....	120, 121, 122, 123, 124, 117
SP 080	120, 121, 122, 123, 124, 117
SP 085.....	120, 121, 122, 123, 124, 117
SP 090.....	120, 121, 122, 123, 124, 117
SP 095	120, 121, 122, 123, 124, 117
SP 100.....	120, 121, 122, 123, 124, 117
SPIDER™ – MINI-SPIDER™ / LIMITED WRIST FUSION SYS- TEM	94, 95
SPIN® Snap-off® screw – dia. 2.0 mm	89
SPIP-520-0-WW	73
SPIP-520-1-WW	73
SPIP-520-2-WW	73
SPIP-520-3-WW.....	73
SPIP-520-4-WW	73
SPIP-520-5-WW	73
SPIP implants	73
SPIP Instrument Set	73
STEM-0920-025-06.....	114
STEM-0920-025-07	114
STEM-0920-025-08.....	114
STEM-0920-025-09.....	114
STEM-0920-025-10.....	114
STEM-0920-025-11.....	114
STEM-0920-025-12.....	114
STEM-0920-025-13.....	114
STEM-0920-025-14	114
STEM-0920-025-15.....	114
STEM-0920-025-16.....	114
STEM-0950-025-06.....	114
STEM-0950-025-08.....	114
STEM-0950-025-10	114
STEM-0950-025-12	114
STEM-0950-025-14	114
Surfix® Concept	118
Surfix® technology	118
SZR-900-40-16	111
SZR-900-44-15	111
SZR-900-44-18	111
SZR-900-48-17.....	111
SZR-900-48-20	111
SZR-900-52-18	111
SZR-900-52-21	111
SZR-900-56-21	111
SZR-910-38.....	116
SZR-910-41.....	116
SZR-910-44	116
SZR-910-47.....	116
SZR-910-50	116
SZR-910-53.....	116

SZR-910-56	116
SZT-715-01S	96
SZT-715-02S	96
SZT-715-03S	96
SZT-715-04S	96
SZT-715-05S	96

T

The Thumb carpometacarpal (CMC) joint.....	78
TITAN HRA Implants	111
TITAN HRA Resurfacing	111
TITAN® Humeral Resurfacing Arthroplasty	110
TMP-200-01	72
TMP-400-10	81
TMP-400-20	81
TMP-400-30	81
TMP-400-40	81
TMP-400-50	81
TMP-410-01	80
TMP-443-L	79
TMP-443-M	79
TMP-443-S	79
TRL-100-05D	75, 77
TRL-100-05P	75, 77
TRL-100-10D	75, 77
TRL-100-10P	75, 77
TRL-100-20D	75, 77
TRL-100-20P	75, 77
TRL-100-30D	75, 77
TRL-100-30P	75, 77
TRL-100-40D	75, 77
TRL-100-40P	75, 77
TRL-100-50D	75, 77
TRL-100-50P	75, 77
TRL-202-10D	72
TRL-202-10P	72
TRL-202-20D	72
TRL-202-20P	72
TRL-202-30D	72
TRL-202-30P	72
TRL-202-40D	72
TRL-202-40P	72
TRL-351-20L	104
TRL-351-20S	104
TRL-351-22L	104
TRL-351-22S	104
TRL-351-24L	104
TRL-351-24S	104
TRL-410-10M	80
TRL-410-20M	80
TRL-410-30M	80
TRL-410-40M	80
TRL-420-14S	78
TRL-420-16S	78

TRL-420-167	78
TRL-420-168	78
TRL-420-187	78
TRL-420-189	78
TRL-435-10	81
TRL-435-20	81
TRL-435-30	81
TRL-435-40	81
TRL-435-50	81
TRL-443-10M	79
TRL-443-10S	79
TRL-443-20M	79
TRL-443-20S	79
TRL-443-30L	79
TRL-443-30M	79
TRL-443-40L	79
TRL-443-40M	79
TRL-500-05	75, 77
TRL-500-10	75, 77
TRL-500-20	75, 77
TRL-500-30	75, 77
TRL-500-40	75, 77
TRL-500-50D	75, 77
TRL-520-0	73
TRL-520-1	73
TRL-520-2	73
TRL-520-3	73
TRL-520-4	73
TRL-520-5	73
TRL-600-H160	99
TRL-600-H175	99
TRL-600-H190	99
TRL-600-LNG	99
TRL-600-MED	99
TRL-600-S45	99
TRL-600-S55	99
TRL-600-S65	99
TRL-610-H145-S45	99
TRL-610-H145-S55	99
TRL-610-H145-S65	99
TRL-610-H160-S45	99
TRL-610-H160-S55	99
TRL-610-H160-S65	99
TRL-610-H175-S45	99
TRL-610-H175-S55	99
TRL-610-H175-S65	99
TRL-610-H190-S45	99
TRL-610-H190-S55	99
TRL-610-H190-S65	99
TRL-715-01T	96
TRL-715-02T	96
TRL-715-03T	96
TRL-715-04T	96
TRL-715-05T	96
TRL-890-10MT	23

TRL-890-10PPH	23
TRL-890-10PPT	23
TRL-890-20MT	23
TRL-890-20PPH	23
TRL-890-20PPT	23
TRL-890-30MT	23
TRL-890-30PPH	23
TRL-890-30PPT	23
TRL-890-40MT	23
TRL-890-40PPH	23
TRL-890-40PPT	23
TRL-900-40/16	111
TRL-900-44/15	111
TRL-900-44/18	111
TRL-900-48/17	111
TRL-900-48/20	111
TRL-900-52/18	111
TRL-900-52/21	111
TRL-900-56/21	111
TRL-910-38-14	116
TRL-910-41-15	116
TRL-910-41-18	116
TRL-910-44-16	116
TRL-910-44-19	116
TRL-910-47-17	116
TRL-910-47-20	116
TRL-910-50-18	116
TRL-910-50-21	116
TRL-910-53-19	116
TRL-910-53-22	116
TRL-910-56-21	116

U

UNI-CLIP®	88
UNI-CP™	43
UNIVERSAL 2™ TOTAL WRIST IMPLANT SYSTEM	90, 91
Upper screw module	85

V

V35 012	100
V35 014	100
V35 016	100
V35 018	100
V35 020	100
V35 022	100
V35 024	100
V35 026	100

W

WED006	125
WED008	125
WED010	125

WEDo12.....	125
WEDo14.....	125
WIFIX [®] plate – Surfix [®] Volar plate for distal radius	100



250
240
230
220
210
200
190
180
170
160
150
140
130
120
110
100
90
80
70
60
50
40
30
20
10
MILLIMETERS





Distributed by

Integra LifeSciences Services (France) SAS
 Sales & Marketing EMEA
 Immeuble Séquoia 2 • 97 allée Alexandre Borodine
 Parc technologique de la Porte des Alpes
 69800 Saint Priest • FRANCE
 Phone: +33 (0)4 37 47 59 00 • Fax: +33 (0)4 37 47 59 99
 emea.info@integralife.com • integralife.eu

Customer Services
 International: +33 (0)4 37 47 59 50 • +33 (0)4 37 47 59 25 (Fax) • csmea@integralife.com
 United Kingdom: +44 (0)1 264 345 781 • +44 (0)1 264 363 782 (Fax) • csuk.ortho@integralife.com
 France: +33 (0)4 37 47 59 10 • +33 (0)4 37 47 59 29 (Fax) • cs-ortho@integralife.com
 Benelux: +32 (0)2 257 4130 • +32 (0)2 253 2466 (Fax) • custsvcbenelux@integralife.com
 Switzerland: +41 (0)2 27 21 23 30 • +41 (0)2 27 21 23 99 (Fax) • custsvcsuisse@integralife.com

©2013 Integra LifeSciences Corporation. All rights reserved. ILS 02-02-004-01-12 / ILS 02-02-006-01-12 / ILS 02-02-005-01-12 / ILS 08-02-009-02-12
 PRODUCTS FOR SALE IN EUROPE, MIDDLE-EAST AND AFRICA ONLY.

Availability of these products might vary from a given country or region to another, as a result of specific local regulatory approval or clearance requirements for sale in such country or region. Always refer to the appropriate instructions for use for complete clinical instructions. Non contractual document. The manufacturer reserves the right, without prior notice, to modify the products in order to improve their quality. WARNING: Applicable laws restrict these products to sale by or on the order of a physician. All text and diagrams contained in this catalog are intended as guidelines and are for illustrative purposes only. No information included in this catalog is intended to replace the Directions for Use packaged with each Integra product. Always refer to the appropriate Directions for Use for complete clinical instructions. For technical and additional clinical information on Integra products, please contact your Integra local sales representative. Nugrip, Safeguard, Saddle, Kompressor, Universal2, Spider, First Choice, RADFx, Catalyst, EndoRelease, Titan, PyroTitan, PyroDisc are trademarks of Integra LifeSciences Corporation or its subsidiaries. Qwix, Bold, Uni-Clip, Spin, Snap-Off, Wifix, Surfix, Integra and the Integra logo are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. All the references numbers mentioned on this document are CE marked according to European council directive 93/42/EEC on medical devices and its relatives, unless specifically identified as "NOT CE MARKED."