

MicroFrance® Duncavage Sphenoid Sinus Instrument Set

Currently, surgeons who operate within the sphenoid to remove posterior cells are forced to invert frontal sinus instruments in order to access the floor of the sphenoid sinus. All of the Duncavage instruments have a 25° bend superiorly and a 25° bend inferiorly at the distal tips, allowing the surgeon to access the sphenoid sinus quickly and accurately. They also feature reverse-action handles that reduce the amount of movement at the distal tips when actuated.

Ordering Information

Product	Description	Quantity
MCEN150*	$. {\sf MicroFrance}^{\tt @} {\sf Duncavage} {\sf Sphenoid} {\sf Sinus} {\sf Set} \dots \dots \dots$	1 set
Set Includes:		
MCEN151	. Sinus Vertical Cup Forceps	1
MCEN152	. Sinus Horizontal Cup Forceps	1
MCEN153	. Sinus Vertical Through-Cut Forceps	1
MCEN154	. Sinus Horizontal Through-Cut Forceps	1
MCEN155	. Sinus Blunt Probe	1
MCEN156	. Sinus Blunt Probe with Suction	1
MCEN157	. Sinus Freer Elevator	1
MCEN158	. Sinus Freer Elevator with Suction	1
3717002	. Instrument Tray, One-Level (52.0 cm L x 25.0 cm W x 2.8 cm	m D) 1

*Note: This set is only available in the US. If you are not in the US, please order the above items individually.



Medtronic ENT Medtronic USA, Inc.

6743 Southpoint Drive North Jacksonville, FL 32216 USA

www.MedtronicENT.com Toll free: (800) 874-5797 Fax: (800) 678-3995

International Telephone Numbers

Adriatic Regional Office 385-1-488-1120 Australia 1800-668-670 Baltic Regional Office 37-1-67560226 Belgium 32-2456-09-09 Canada 1800-217-1617 China 86-21-50800998 Czech Republic 420-2-9657-9580 France 33-470-679-800 Germany 49-2159-8149-209 Greece 30-210-67-79-099 Hong Kong 852-2919-1312 Hungary 36-30-5052987 India 91-22-26836733 Israel 972-9-972-4400 Italy 39-02-24137-324
Japan 81-6-4795-1506
Korea 82-2-3404-3600
Lebanon 961-1-370-670
Luxembourg 32-2456-09-09
Netherlands 31-45-566-8800
Poland 48-22-465-6942
Russian Federation 7-495-580-73 77
Singapore 65-6776-6255
South Africa 27-11-466-1820
Spain 34-91-625-05-40
Taiwan 886-2-2183-6000
UK 44-1923-205-166
USA 1-904-296-9600





This literature is intended for the exclusive use of physicians. Rx only.
© 2010, Medtronic, Inc. MicroFrance® is a registered trademark of Medtronic, Inc. 890512 07.10 2010-395 Rev. 2

890512_DuncavageSphenoid_R2.indd 2 7/6/10 9:12 AM